



ImageWorks
250 Clearbrook Road, Suite 240
Elmsford, NY 10523

EVA Sensor and EVAsoft Integration

EVA sensors acquire images directly into EVAsoft or almost any other imaging software. There are 2 ways to operate EVA sensors in imaging software:

1. Direct Driver Integration
2. TWAIN Driver Integration

Option 1: Direct Driver Integration. In this option the EVA Driver is directly incorporated into an imaging software and the EVA sensor will capture images directly into that software. This is how EVA sensors work in EVAsoft, and in:

- | | |
|--------------------------------|---------------------------------|
| ✓ AFP Digital | ✓ Lightyear (Apteryx) |
| ✓ Apteryx (XVLite and XVision) | ✓ LiveDDM |
| ✓ ADSTRA | ✓ Mediadent |
| ✓ CADI software/Synca | ✓ Suni (Apteryx) |
| ✓ DentalEye | ✓ TigerView |
| ✓ Mogo Imaging | ✓ XDR |
| ✓ GE Healthcare (Apteryx) | ✓ <i>Toothpics (Mac only)</i> |
| ✓ ImageFX | ✓ <i>RadioVision (Mac only)</i> |
| ✓ J. Morita (TigerView) | |

Option 2: TWAIN Driver Integration. In this option the EVA TWAIN Driver is incorporated into an imaging software and the EVA sensor will capture images into that imaging software. Examples of EVA sensors capturing into another imaging software using the EVA TWAIN Driver:

- | | |
|----------------------------------|-----------------------------|
| ✓ Dentrax Image | ✓ Dexis' Imaging Software |
| ✓ Eaglesoft Imaging | ✓ Schick's Imaging Software |
| ✓ PracticeWorks Imaging Software | ✓ Kodak's Imaging Software |
| ✓ Easy Dental Imaging Software | ✓ Gendex's Imaging Software |
| ✓ SoftDent Imaging Software | ✓ MANY OTHER SOFTWARE |

EVAsoft can bridge to all popular practice management software with one single left click of a mouse. This means that when you are in a patient's record or chart in a practice software you can automatically launch EVAsoft in that patient's record, ready to take x-rays in the Acquire Tab. There are 2 ways:

1. ProIntegrate Practice Software Bridge (must be purchased from Dent-X)
2. Built-in bridge, already provided by practice software

Option 1: EVAsoft Bridges to Practice Software via ProIntegrate:

ABELDent	D4W
ACE Dental	DataTeamDDS
AltaPoint	DESCO
Autopia Dental Works	dental.com
AVI Mark	OrthoTrack
BridgeIT	Dental-Exec / Perio-Exec
Computer Age Dentist (CADwin)	Dentisoft
Cornerstone Vet (IDEXX)	Dentrax
C.T.R.L.	Discus



ImageWorks
250 Clearbrook Road, Suite 240
Elmsford, NY 10523

Doc Aide	Open Dental
EagleSoft	Patient Base
Easy Dental	Power Practice Pro
Exact Dental System	PracticeWorks
First Pacific	PerioVision
Focus	Radiology Lab (RLMS)
Innova	SoftDent
Intellect Dental System (IDS)	Sorodex
Intellident	TDOCS
Julie	Time
Maxident	Visiodent
Maximus	WinDent
Mediadent	WinOMS
MeerMANAGER	XLDent
Oasis	

Option 2: The Bridge is already built into the Practice Software:

ACE	ExcelDent
Autopia Dental Works	Intellect Dental System (IDS)
AVImark	EZ2000
DentiMax	Mogo
Open Dental	Practice-Web
Genesis	

Disclaimer: This list is accurate to the best of our knowledge at the time of publication and is distributed as an informational service to our customers. However as each integration involves parties beyond ImageWorks control, ImageWorks cannot offer any warranty implied, expressed or written regarding this list. The names on this list are variously trademark and copyright by their respective owners. By inclusion in this list ImageWorks does not express a recommendation or endorsement and does not intend to imply any endorsement by any company named of ImageWorks or its products.