

ImageWorks

8 Westchester Plaza, Suite 112
Elmsford, NY 10523

Prolmage to EVAsoft Database Migration Tool

Part Number: **9992810105**

Serial Number: **UT-1055**

Dear Valued Customer:

Enclosed please find the USB key for the Prolmage to EVAsoft Database Migration Tool. Install this on your host/server or standalone computer after you have installed EVAsoft. This will enable EVAsoft to import your existing Prolmage Database and Images into EVAsoft's own Database and DICOM format image files. The USB key can only be used once on a single computer but can be used to import multiple Prolmage databases into a single EVAsoft database.

EVAsoft sequentially imports and resaves every image in the Prolmage Database when the import is run. The Prolmage database and images are left intact on the hard drive. It is therefore necessary that you ensure that you have sufficient empty space on your hard drive to run import tool. We recommend at least 1.5 times the size of your existing Prolmage Database and images are available. The time to run the import tool will vary by size of the original DB (in terms of number of images) and the speed of the local computer (and its hard drive). The table below is an estimate of the time to import a Prolmage Database into EVAsoft. Your actual times will vary. The migration time is significant for databases that were used in busy offices for a number of years so please plan accordingly.

No X-Rays in Prolmage DB	Estimated Time to to Import into EVAsoft	Approximate size of Prolmage Data on HD
1000	35 minutes	2 GB
5000	3.5 hours	10 GB
10000	7 hours	20 GB
20000	14 hours	40 GB
50000	35 hours	100 GB
100000	70 hours	200 GB

If you need assistance visit our website, email or give us a call.

www.imageworkscorporation.com
support@imageworkscorporation.com
1-914-592-6100

Thank you for choosing EVAsoft,

ImageWorks Team



8 Westchester Plaza, Suite 112

Elmsford, NY 10523

Phone: 914-592-6100

www.ImageWorksCorporation.com

Directions for EVAsoft – ProImage Database Migration

Install and configure EVAsoft prior to following these instructions.

This should be run from the EVAsoft Host.

- Step 1..... Close (exit) EVAsoft if it is running.
- Step 2..... Insert the Database Migration USB drive.
- Step 3..... Double-click on the **Setup.exe** file.
- Step 4..... Click on **Next** to begin the installation guide.
- Step 5..... Click on **Install** to begin the installation.
- Step 6..... Wait for the installation to complete, click **Finish**.
- Step 7..... Launch EVAsoft from the desktop shortcut.
- Step 8..... Register EVAsoft if still running in 30-Day Trial mode.
- Step 9..... In EVAsoft go to **Tools → Import ProImage Database**.
- Step 10..... Read the Warning message, verify that you will have enough free hard-drive space and enough time to complete the database migration and then click **Continue**.
- Step 11..... Select your ProImage database with the windows browsing dialog screen and click **Open**.
- Step 12..... Depending on the size of your ProImage Database and the number of images to be imported, this process can take many hours (approx 2-3 seconds per image). The status will be displayed in the lower left hand corner of EVAsoft while the import is running.
- Step. 13..... When the migration is complete a dialog will be displayed “The ProImage Database was successfully imported.”
- Step. 14..... Click on **OK**.