BD ClipVac™ Presurgical Hair Removal System Instructions For Use

REF 5500E / 5506A / 5506E / 5506U / 5506K

en English2

fr Français.....15

de Deutsch.....28

it Italiano41

es Español54

pt Português67

nl Nederlands...80

sv Svenska93

da Dansk.....106

fi Suomi.....119


no Norsk132

el Greek145

tr Türkçe.....158

pl Polski.....171

cz Český184

Intended Use

The intended use of the BD ClipVac™ system is to vacuum clipped hair and airborne contaminants generated from the clipping process. It should be used with the BD Surgical Clipper (REF 5513E). Only trained healthcare personnel should use BD ClipVac™ system and the BD Surgical Clippers. If the BD Surgical Clippers and BD ClipVac™ system are being used for the first time, see the training materials provided by BD.

How Supplied:

The BD ClipVac™ Vacuum Unit (REF 5500E) is a non-sterile reusable device. Included in this box:

1. One vacuum unit
2. IFU

The BD ClipVac™ Charging Adapter & Battery will be included in a separate box (REF 5506A/5506E/5506U/5506K) that contains:

1. One battery
2. One charging adapter (charger, power supply & power cord)

Note: The BD ClipVac™ Disposable Filter Assembly is packaged and sold separately (REF 5575).

Limitations on Reprocessing:

Repeated cleaning has minimal effect on the system. End of life is determined by wear and inadvertent damage due to use.


BD ClipVac™ Parts Identification:

REF 5500E: Vacuum Unit


REF 5505: Replacement Battery

REF 5506A/5506E/5506U/5506K: Charging Adapter & Battery

REF 5575: Disposable Filter Assembly


Step 1. Remove the battery from charging stand and review the LCD display ① to ensure battery has sufficient charge for the procedure prior to use. See Charging the Battery section.


Step 2. To insert battery in vacuum unit, retract battery latch ② and place in open slot with the battery tab opposite the battery latch.

Note: Battery will be properly seated when latch mechanism is slightly over the battery.


Step 3. Remove the disposable filter assembly from package. Uncoil tubing ③ by gripping near the nozzle ④ and filter body ⑤, gently pulling taut to remove any twist and straighten tubing.


Step 4. Place the filter body onto the vacuum port end of the vacuum unit.


Note: The filter body will rotate freely when attached.


Step 5. Place nozzle of disposable filter assembly over the clipper after the surgical blade has been installed on the clipper handle. Nozzle will click into place when properly seated to the clipper. Ensure the retention knob ⑥ on either side of the nozzle is properly seated prior to use.


Step 6. Turn ON the vacuum unit by flipping the switch ⑦.
Note: Ensure BD Surgical Clipper is turned ON prior to starting clipping procedure.


Step 7. While holding loose skin taut, clip hair using a slow steady forward motion avoiding downward pressure.
Note: The bottom side of the clipper blade should remain flat and gently rest on the surface of the skin while clipping. To ensure hair is properly collected by the vacuum, never tilt the blade edge into the skin.
Note: For optimal performance, use on dry hair.

Step 8. When clipping is complete, turn OFF unit and gently twist the filter body to remove disposable filter assembly away from the vacuum unit.

Step 9. Discard filter body into general waste and recycle nozzle and tubing according to local regulations. To detach the filter body from the tubing, pull with a gentle twist.

Step 10. After disinfecting and/or manual cleaning, place vacuum unit in the designated area for next use. If necessary, return battery to the charging adapter.

Note: The BD ClipVac™ Vacuum Unit may be placed on the bed next to the patient or secured to the bedrail.

Charging the Battery

Plug the charging cord into the charging adapter. A quick red and green LED flashing is an indication that charger is ON. Place the battery into the battery bay ensuring that the 5-way connector is fully seated. The LCD fuel gauge on the battery will begin flashing to indicate it is receiving a charge. The LED indicators on the front of the charging adapter will provide the charging status. Below are the LED indicators on the charging adapter:

- Green Flashing.....Battery is charging
- Green Solid.....Battery is fully charged
- Red Solid.....Battery failure, replace battery

Charging notes:

- Maximum time required for a charging adapter to charge an empty battery shall be 4 hours.
- Charge at room temperature of +41 °F (+5 °C) to +95 °F (+35 °C)
- Avoid charging in direct sunlight or near radiating heat.
- If not used for extended periods of time, charge battery every six (6) months to preserve battery life.
- The lithium-ion battery must be disposed of in accordance with local regulations or by using the RBRC 1-800 number listed on the battery to locate your local recycling point. Never incinerate the lithium battery. Do not dispose of the battery with general waste.

Disinfection Instructions:

The BD ClipVac™ Vacuum Unit should be thoroughly disinfected, after each use, using one of the recommended disinfectant solutions listed below:

- Isopropyl Alcohol/Quaternary Ammonia solution (e.g., validated disinfectant PDI Super Sani-Cloth®) containing the active ingredients:
 - 0.25% n-Alkyl (68% C₁₂, 32% C₁₄) dimethyl ethylbenzyl ammonium chloride
 - 0.25% n-Alkyl (60% C₁₄, 30% C₁₆, 5% C₁₂, 5% C₁₈) dimethyl benzyl ammonium chloride
 - 55.00% Isopropyl Alcohol
- Chlorinated bleach solution (e.g. validated disinfectant PDI Sani-Cloth® Bleach) containing the active ingredient:
 - 0.63% sodium hypochlorite (i.e., 6300 ppm sodium hypochlorite)

1. Thoroughly apply to all surfaces of the BD ClipVac™ Vacuum Unit and maintain a wet surface for at least:
 - a. 3 minutes for Isopropyl Alcohol/Quaternary Ammonia solutions or
 - b. 4 minutes for Chlorinated bleach solution.
2. Surfaces must remain visibly wet with the disinfectant for the complete timeframe indicated above.
3. Repeat the disinfection steps above if the unit dries prior to the allotted timeframe.
4. Allow the BD ClipVac™ Vacuum Unit to air dry completely prior to storage.

Fuel Gauge LCD display on top of the battery:

- 5 Bars = 80% to 100% charged
- 4 Bars = 60% to 79% charged
- 3 Bars = 40% to 59% charged
- 2 Bars = 20% to 39% charged
- 1 Bar = 0% to 19% charged

If the device comes in contact with bodily fluid, the following manual cleaning instructions can be used:

1. Using enzymatic/neutral pH detergent, prepare per the manufacturer's instructions utilizing tap water entailing a temperature range of 27°C to 44°C (81°F to 111°F).
2. Do not submerge or spray the device with cleaning solution.
3. Use a soft, clean cloth saturated with the cleaning solution (not dripping) to manually clean all external contaminated surfaces of the device for a minimum of 5 minutes. Be careful not to cause any liquid to drip into the charger or vacuum. In conjunction with this step, use an appropriately sized brush (standard soft bristle toothbrush and lumen brush) to remove all visible soil from the device; paying particular attention to cracks, crevices and other hard to clean areas. Note: It is recommended that enzymatic or neutral pH detergent solution should be changed when it becomes grossly contaminated (bloody and/or turbid).
4. Rinse thoroughly by using a soft, clean cloth saturated with lukewarm (27°C - 44°C) tap water (not dripping) for a minimum of 30 seconds. Repeat 2 additional times for a total of 3 total rinse cycles. Be careful not to cause any liquid to drip into the charger or vacuum.
5. Dry the device with a clean, lint-free towel.
6. Visually examine the device for cleanliness.
7. If visible soil remains, repeat cleaning procedure until the device is thoroughly clean.

Inspection/Maintenance

- The BD ClipVac™ Vacuum Unit requires no routine maintenance aside from cleaning system after each patient use. If a replacement is needed, contact your local BD representative or customer service. Service or attempted repair performed by unqualified personnel may result in a risk of injury, electric shock, fire, or permanent damage to the equipment and will void any warranty.
- Never use BD ClipVac™ Vacuum Unit or a battery that appears to be defective, damaged, or not working properly.

Limited Warranty:

- The BD ClipVac™ Presurgical Hair Removal System from BD (REFs 5500E, 5505, 5506A/5506E/5506U/5506K) is warranted against defects in material and workmanship under normal use for a period of two (2) years from the date of unit placement. If this product is determined to be defective with the warranty period, it will be repaired or replaced without charge by BD, if you contact your local representative or customer service representative. They will provide the proper return authorization. Please include a description of the defect or damage with the returned unit.
- This warranty does not cover the Disposable Filter Assembly (REF 5575) nor does it cover damage caused by misuse, negligence, accident, or abuse, alteration or modification of the unit by anyone other than BD or its authorized agent, or affixation of any accessory not manufactured or authorized by BD, or used in combination with other surgical clippers.
- BD is not liable in either tort or contract for any cost, expense, loss, or damage whether direct, consequential, special, collateral, or incidental arising out of the use, misuse, or inability to use this product. Except as expressly provided above, BD disclaims all warranties including, but not limited to the implied warranties of merchantability and fitness for particular purpose.

Disposal:

- Disposable assembly can be discarded into general waste or per hospital protocol.
- The nozzle and tubing may be recycled. To remove the filter body from the tubing, simply pull the tube from the filter body.
- The lithium-ion battery must be disposed of in accordance with local regulations or by using the RBRC 1-800 number listed on the battery to locate your local recycling point. Never incinerate the lithium battery. Do not dispose of the battery with general waste.

Caution:

- Use the BD ClipVac™ system as described in the Operating Instruction. Avoid direct patient contact with the nozzle. Direct patient contact should be with the bottom side of the clipper blade.
- Never use BD ClipVac™ Vacuum Unit, battery or charging adaptor that appears to be defective, damaged, or not working properly.
- Do not cover vacuum unit to prevent obstructing air flow as this may cause overheating of vacuum unit.
- Do not re-use disposable filter assembly, it is intended to be single use only.
- Clean the vacuum unit after each patient use, refer to cleaning and disinfecting instructions.
- Do not submerge or spray the vacuum unit, battery, or charging adapter with any liquid.
- Do not sterilize the vacuum unit, battery or charging adapter.
- Do not short the battery terminals.
- Avoid excessive physical shock or vibration to the battery.
- Do not place the battery and charging adapter in a location where it is subjected to direct sunlight and keep away from other external heat sources.
- Do not dispose the battery with general waste. Never incinerate the lithium ion battery.
- Do not replace charging adapter cord. Use only the battery and charging adapter supplied by BD.
- Do not attempt to modify or repair the vacuum unit, battery, or charging adapter as this can result in damage of the BD ClipVac™ system and will void the warranty.

Symbol Legend:

 Catalog Number

 Do not reuse

 Consult instructions for use

 Humidity limitation

 Batch Code

 Switch "On"
(Power)

 Direct current

 Temperature limit

 Manufacturer

 Switch "Off"
(Power)

 Type BF applied part

 Atmospheric pressure limitation

 Date of Manufacture

 Non-ionizing
electromagnetic radiation

 Authorized representative in
the European Community

Information for Users on Collection and Disposal of Old Equipment and used Batteries

These symbols on the products, packaging, and/or accompanying documents mean that used electrical and electronic products and batteries should not be mixed with general household waste. For proper treatment, recovery and recycling of old products and used batteries, please take them to applicable collection points, in accordance with your national legislation and the Directives 2012/19/EU and 2006/66/EC. By disposing of these products and batteries correctly, you will help to save valuable resources and prevent any potential negative effects on human health and the environment which could otherwise arise from inappropriate waste handling.

For more information about collection and recycling of old products and batteries, please contact your local municipality, your waste disposal service or the point of sale where you purchased the items. Penalties may be applicable for incorrect disposal of this waste, in accordance with national legislation.

Information on Disposal in other Countries outside the European Union

These symbols are only valid in the European Union. If you wish to discard these items, please contact your local authorities or dealer and ask for the correct method of disposal.

BD ClipVac™ Vacuum Unit Specifications:

BD ClipVac™ Vacuum Unit Mode of Operation	Continuous
BD ClipVac™ Disposable Filter Assembly Mode of Operation	Single Use
BD ClipVac™ Vacuum Unit	14.4VDC - 2.8A
BD ClipVac™ Vacuum Unit Weight	725g/25.6oz
Battery Type	Lithium ion with LCD Display
Battery Weight	234g/8.25 oz.
Full Charge Time	Approximately 3.5 hours
Continuous Run Time	Approximately 75 minutes on a full charge
Battery Life	Minimum 260 charges
Power Supply	100 – 240 VAC, 24V 2.71A DC
Ambient Operating Temperature, Humidity, and Atmospheric Pressure	10°C to 40°C (50°F to 104°F) 15% to 95% RH non-condensing 70 - 106kpa.
Storage Temperature, Humidity, and Atmospheric Pressure	-20°C to 50°C ≤ 80% RH -4°F to 122°F ≤ 80% RH 50 - 106kpa

Charging Adapter Specification:

Weight	235g/8 oz.
Height	58mm/2-1/4"
Length	180mm/7"
Width	92mm/3-1/2"
Mating Connector	5-blade standard battery connector
Communications Compliances	System Management Bus Rev. 1.0, Smart Battery Data Specification rev 1.0 and Smart Battery Charger specification Rev. 1.0 compliant
Power Supply	100-240VAC, 24V, 2.71A DC, 119x38x60mm, 500g/18oz, UL Listed, CE Compliant, CEC Compliant
Mains Cord	5506A - Type B 5506E - Type F 5506K - Type I 5506U - Type G

EMC Information of REF 5500

REF 5500E needs special precautions regarding EMC and needs to be used according to the following EMC information. REF 5500E complies with the EMC requirements of the IEC 60601-1-2.

The BD ClipVac™ system is intended for use in the electromagnetic environment specified in the following tables.

Warning:

REF 5500E should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, REF 5500E should be observed to verify normal operation in the configuration in which it will be used.

REF 5500E has been tested and found to comply with the limits for the medical devices to the IEC 60601-1-2. These limits are designed to provide reasonable protection against harmful interference. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to the other devices which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures;

- Reorient or relocate the receiving device.
- Increase the separation between the equipment.
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) are connected.
- Consult the manufacturer or field service for help.

Guidance and manufacturer's declaration — electromagnetic emissions

REF 5500E is intended for use in the electromagnetic environment specified below.
The customer or the user of REF 5500E/5506 series should insure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	REF 5500E uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	REF 5500E is suitable use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration — electromagnetic immunity

REF 5500E/5506 series is intended for use in the electromagnetic environment specified below.
The customer or the user of REF 5500E/5506 series should insure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 ¹ kV air ±8 ² kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV for differential mode ±2 kV common mode	±1 kV for differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<0% U_T (>100% dip in U_T) for 0.5 cycle (phase angle 0°,45°,90°,135°,180°,225°,270°) 0% U_T (100% dip in U_T) for 1 cycles 70% U_T (30% dip in U_T) for 25/30 cycles <0% U_T (>100% dip in U_T) for 250 cycle	<0% U_T (>100% dip in U_T) for 0.5 cycle (phase angle 0°,45°,90°,135°,180°,225°,270°) 0% U_T (100% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <0% U_T (>100% dip in U_T) for 250 cycle	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

¹REF 5500E - Vacuum unit and charger compliance level to ±15 kV air discharge.

²REF 5506A/5506E/5506U/5506K - Battery compliance level to ±8 kV air discharge. Battery fuel gauge LCD display may be damaged at ±15 kV air. But this is not a hazard to the patient.

Note: U_T is the a.c. mains voltage prior to application of the test level

Guidance and manufacturer's declaration — electromagnetic immunity

REF 5500E/5506 series is intended for use in the electromagnetic environment specified below.
The customer or the user of REF 5500E/5506 series should insure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms 150 kHz to 80 MHz (ISM and amateur radio bands)	3 Vrms (150 kHz to 80 MHz) 6 Vrms (150 kHz to 80 MHz)	Portable and mobile RF communications equipment should be used no closer to any part of REF 5500E, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strength from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol:
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m 80 MHz to 2.7 GHz	
Enclosure port immunity IEC 61000-4-3	27 V/m, 385MHz 28 V/m, 450MHz 9 V/m, 710 MHz to 780 MHz 28 V/m, 810 MHz to 2450 MHz 9 V/m, 5240 MHz to 5785 MHz	27 V/m, 385MHz 28 V/m, 450MHz 9 V/m, 710 MHz to 780 MHz 28 V/m, 810 MHz to 2450 MHz 9 V/m, 5240 MHz to 5785 MHz	$d = 6/\sqrt{fP}$ 380MHz to 5.8GHz

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^aField strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which REF 5500E is used exceeds the applicable RF compliance level above, REF 5500E should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating REF 5500E.

^bOver the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distance between portable and mobile RF communications equipment and REF 5500E.

REF 5500E/5506 series is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of REF 5500E/5506 series can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and REF 5500E/5506 series as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800MHz to 2.7GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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©2019 BD. BD, the BD Logo and BD ClipVac are
Becton Dickinson Pty Ltd.
Becton Dickinson Limited