Touch 'n Heat Model 5004

Instruction Guidelines

Warnings and Precautions

DANGER: Not for use in the presence of flammable anesthetics.

WARNING: Use of other accessories which are not authorized for use in connection with the Touch 'n Heat 5004 may cause malfunction and may compromise patient safety.

CAUTION: This device has been investigated with regard to safety from electrical shock and fire hazard as well as electromagnetic compatibility (EMC). The device has not been investigated for other physiological effects.

CAUTION: For use by qualified and trained personnel only.

Package Contents

P/N	Description
973-0212	Touch 'n Heat Unit (110V)
973-0213	Touch 'n Heat Unit (220V)
973-0217	Narrow Posterior Heat Carrier
973-0222	110 V Charge Adapter (for 110V unit only)
973-0223	220 V Charge Adapter (for 220V unit only)

Bottom Panel Symbols

	Attention! Consult the manual before use.
	CSA mark with "C/US" indicator for certified products.
MC 157498 (LR 113060)	CSA file number.
	Type BF Patient Applied Part

Back Panel Symbol

	Direct Current.
--	-----------------


Packaging/ Handling Symbols

	This way up.
	Fragile, handle with care.
	Keep away from rain.
	Temperature limits for storage and transportation.

Index:	page
English	1
Deutsch	4
Français	7
Italiano	10
Dansk	13
Español	16
Suomi	19
Nederlands	22
Svenska	25
Português	28
Norsk	31
Ελληνικά	34
日本語	37

Product Features


- The unit operates from a rechargeable battery.**
A sturdy, sealed lead-acid battery provides quick reliable power. At the maximum setting the battery will provide 30 to 45 minutes of heating power before being depleted.
- Recharge the battery with the charge adapter.**
The unit is delivered fully charged.
To recharge the battery, use only the plug-in adapter which comes with the Touch 'n Heat Model 5004. Remove the tip from the probe handle and place the metal part of the probe handle in the probe holder. Put the unit in the STORE mode. Plug the charge adapter into a grounded wall outlet and insert the jack at the end of the cord into the receptacle at the back of the unit. The yellow light will turn on.
Note: The Touch 'n Heat will not operate while charge adaptor is plugged into unit.
For longest battery life, keep the unit on charge at all times when not in use.
The battery is fully charged when the yellow CHRG light turns off. If the charge adapter is reinserted when the unit is fully charged, the CHRG light will flash once.
The yellow LOW light will light up when it is needed to fully recharge the battery.
- The front panel displays all power controls and indicators.**
The Touch 'n Heat Model 5004 panel is shown below.

- Power Control Dial
- Heat Light (Green)
- Charge Light (Yellow)
- Low Battery Light (Yellow)
- Selector Push Button

The **power control dial** (1) adjusts the amount of heat supplied to the tip. The heat increases linearly over a range of 0 (low power) to 10 (maximum power), so a setting of 5 equals 50% of the maximum power output. The green **Heat light** (2) comes on when the tip is heating. The **push button switch** (5) selects the **USE** or the **STORE** mode. (Button is in for USE mode, out for STORE mode.)


4. The Touch 'n Heat Model 5004 probe handle has special features.

The pinvise nut (1) provides a convenient holder for the different tips. The **contact spring** (2) on the pinvise nut allows you to activate the Touch 'n Heat from any direction while wearing gloves. When you want to heat a tip, depress the spring so that it contacts the **shaft** (3).

5. Using the selector button.

Whenever the unit is not being used, it should be put in the STORE mode with the selector button. The unit is deactivated in the STORE mode, preventing battery from draining. This will also prevent the tip from heating in case the touch switch is activated accidentally. This feature should always be used when changing tips.

- | | |
|---|---|
| 973-0214 Heat Carrier, Standard Anterior | 973-0219 Hot Pulp Test Tip |
| 973-0215 Heat Carrier, Standard Posterior | 973-0220 Plugger, Thick (1mm) |
| 973-0216 Heat Carrier, Narrow Anterior | 973-0221 Plugger, Thin (0.5 mm) |
| 973-0217 Heat Carrier, Narrow Posterior | 973-0224 Pin Vise Nut with Contact Spring |
| 973-0218 Touch 'n Heat Spoon Tip | |

Touch 'n Heat User Guidelines

1. Installing the tip.

Put the unit in the STORE mode using the panel switch. Loosen the pinvise nut a few turns and insert the tip as far as it will go. The thicker part of the tip should disappear into the pinvise. Hand tighten the nut. Adjust the intensity to your desired setting.


2. Tip life depends on temperature.

The tips have life expectancies which greatly depend on the operating temperature. At the maximum heat setting, the heat carrier will run between 1200 and 1400°F. The tip life at this temperature is short. If the average operating temperature is lower, the heat carrier's life expectancy dramatically increases. All tips will gradually lose their heating efficiency over a period of time.

3. Temperature versus power setting.

The temperature of a certain type of tip depends on its size and the applied heating power. The smaller the tip, the higher the temperature will be at a certain power. The following graph shows the temperature curves for the various tips when they are operating in open air.

Tip Temperature vs. Power Setting

4. Touch mode operation.

The touch mode has the advantage of making it possible for you to enter the tooth with a cool instrument and then rapidly heat the gutta-percha inside the tooth. Since the instrument remains hot, the gutta-percha releases without sticking. The recommended power setting for the touch mode is in the range of 6 to 10, depending on the technique used. On the other hand, if you want to remove gutta-percha, simply use a lower power setting or let the tip cool down in the canal so that the gutta-percha adheres to the tip.

Tips and Accessories

All tips are heated using a patented method with the heat generated close to the tip's point. The accessories available for the Touch 'n Heat include a flat spoon, two plugger sizes, narrow and standard anterior or posterior heat carriers plus a hot pulp test tip.

1. Using the pluggers.

The narrow plugger (0.5mm) is similar to the narrow posterior but has a flat end instead of being pointed. Many practitioners use this tip in the vertical condensation technique. They claim it is easier to remove gutta-percha with this tip than with the narrow posterior tip. The thick plugger (1.0mm) is used to prepare post space.

2. Using the Hot Pulp Test Tip (HPTT).

The HPTT is used to test for heat sensitivity. It is important to use some medium to transfer the heat from the pulp test tip to the tooth. A ball of gutta-percha may be softened and molded around the tip edges. Dip the gutta-percha ball in talcum powder or vasoline to prevent it from sticking to the tooth. Consult the temperature/power chart for your desired power setting. Place the tip in firm contact with the tooth, making sure a good contact is established. The black heat shielding tube on the HPTT minimizes your risk of burning the patient.

Frequently Asked Questions about the Touch 'n Heat

1. What size is the smallest tip?

The smallest tip we have are the narrow posterior and the narrow anterior tips. They both correspond to a #40 file. The straight portions of the posterior and the anterior tips are 18 mm and 21 mm, respectively. They are made from material of the same taper.

2. How hot can a tip get?

At maximum power, a tip in the air can reach 1200 to 1400°F in approximately 8 seconds. You can adjust the heat down to about 100°F by using the lowest power setting. A chart included in the instruction guidelines shows how the tip temperatures depend on the power setting.

3. Isn't this high temperature dangerous?

Like a dental drill, needle or scalpel, the unit should only be used by a trained professional. You should use the high heat setting to sear off gutta-percha in the lateral condensation technique where a cluster of cones are cut through. The high heat is also useful in preventing gutta-percha from sticking to the tip in the vertical condensation technique. The idea is to transfer a relatively large amount of heat in a very short time (1-2 seconds), so that the tip can be withdrawn while it is still hot. On the other hand, if you want to remove gutta-percha, you should use a lower heat setting; i.e., a cooler tip.

4. What is the normal power setting?

The maximum heat setting should be used only when searing off gutta-percha. For the vertical condensation technique, we recommend a power setting in the range of 6-8.

5. **How long will the tip last?**

It depends on how well you take care of the tip. We have had customers using the same tip for two years while other customers will use up a tip in two weeks. The life of a tip will be reduced by:

- a. running it too hot for long periods of time
- b. bending or putting too much mechanical force on it.

The tips should be used only to transfer heat—not to pack down gutta-percha. Our recommendation is that you use a cold plugger for that purpose. Tip life will also be greatly extended if you use the touch mode consistently so that the tip is heated only when needed.

6. **How does the touch switch work and why do I need a contact clip?**

A spring extending backward from the pinvis mount completes the circuit by making contact with the probe shaft.

7. **What kind of batteries are used?**

The Touch ‘n Heat uses sealed lead-acid batteries. Unlike the common nickel-cadmium batteries, you should never run down the lead-acid batteries completely. Total depletion will destroy them. Our Touch ‘n Heats have low battery warning lights to prevent this from happening. Always keep the unit on charge when not in use.

8. **How long do the batteries hold a charge?**

The Touch ‘n Heat battery will last 30–45 minutes and can be recharged in 8 hours. For your safety, the Touch ‘n Heat is designed not to function while being charged.

9. **When will the batteries have to be replaced?**

The battery manufacturer specifies an estimated battery life of 3 to 5 years when kept properly charged. When the batteries are depleted, they must be removed and replaced by SybronEndo since the batteries are soldered directly into the unit.

10. **How should I use the Touch ‘n Heat?**

Practitioners who ask this general question probably do not know enough about the warm gutta-percha technique. We usually recommend attending some courses on this subject. You should not start using the Touch ‘n Heat without having some course instruction and prior practice on extracted teeth.

Cleaning, Sterilization, and Maintenance

Sterilization

Pluggers and tips must be cleaned and sterilized before every use:

- 1. Clean with water and mild, non-abrasive, detergent such as dish washing liquid. Dry thoroughly.
- 2. Double wrap the pluggers and tips using a typical central service wrapping technique and place into autoclave unit avoiding contact with other instruments.
- 3. Operate sterilization cycle at 250°F for 30 minutes or 270°F for 25 minutes. Follow the instructions supplied with your autoclave unit for specific operating instructions.

Expected Life of Accessories

All accessories are reusable only for a limited number of procedures. This number is highly dependent on how and for what duration the instrument is used during each procedure and, therefore, will vary significantly from user to user. If the device function appears to be erratic, an accessory may have exceeded its useful life and should be replaced.

Cleaning

The Touch ‘n Heat 5004 enclosure and the handpiece should be wiped down with a cloth soaked in any of the common non-abrasive, mild, cleaning liquids containing alcohol. Do not submerge the device or let liquid enter the enclosure and handpiece.

Technical Specifications

Dimension/Weight

Casing of unit:	Weight:	1.4 lb. (0.64 kg)
	Dimensions:	Height: 2.2 inches /55 mm
		Width: 3.1 inches /79 mm
		Depth: 5.2 inches /132 mm

Classification: IEC 601-1 Internally Powered Equipment

Mode of Operation: IEC 601-1 Continuous Operation

Duty Cycle: 2 sec. ON, 5 sec. OFF for 5 mins.
@ 30 mins. Interval

Degree of safety of application in presence of a flammable anesthetic mixture:

Unit unsuitable for use in presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

Conditions for use: +10°C to +40°C
30% to 75% relative humidity

Storage and transport conditions: -10°C to +45°C
10% to 95% relative humidity

Atmospheric Pressure: 101.1 kPa

Origin: SybronEndo
1332 South Lone Hill Avenue
Glendora, CA 91740, USA

Model: Touch ‘n Heat 5004

Power Supply: rechargeable sealed lead-acid battery

Warranty Information

Returns and Warranty

The Touch ‘n Heat is designed for lasting durability and predictable results. The unit comes with a one-year warranty.

If the unit should fail to operate correctly, please contact your SybronEndo customer service representative at **1 800-346-3636**. Outside the United States, please call customer service at 714-516-7979 or contact your local dealer or supplier to coordinate the returns/repairs with SybronEndo.

Warranty repairs can only be performed by SybronEndo or by authorized agencies with original factory parts. Any unauthorized repairs will void the warranty.

- 1. Use the original packaging when shipping or storing your Touch ‘n Heat.
- 2. Many problems can be solved over the phone. Do not hesitate to contact us if you experience difficulties when using the Touch ‘n Heat.
- 3. For returns, call us for a Return Good Authorization (RGA) number.
- 4. Damages incurred during shipping due to improper packaging are not covered by the warranty. If the original box and/or foam packing are not available, please call SybronEndo for packing instructions.
- 5. Mark all outer packaging with the RGA number, your name, address, and phone number.
- 6. Shipping costs are not covered by the warranty.

Disclaimer

The responsibility of SybronEndo, as a manufacturer of electrical/medical devices, extends to the technical safety features of the device only if maintenance and repairs are carried out by SybronEndo or by authorized agencies with original factory parts.

For safety reasons, this product should be used with accessories manufactured and sold by SybronEndo. Any use of non-authorized accessories or not following any of the instructions for use is done so at the operator’s risk and voids the warranty.

SybronEndo does not assume any responsibility for poor results due to operator error or equipment malfunction.

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