en

fr

es

Directed Flow Disposable Impression Tray

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3M ESPE
Dental Products
St. Paul, MN 55144-1000 U.S.A.

3M ESPE Customer Care/MSDS Information:
U.S.A. 1-800-634-2249 and Canada 1-888-363-3685.
3M, ESPE, Penta and Position are trademarks of 3M Company
or 3M Deutschland GmbH. Used under license in Canada.
© 2012, 3M. All rights reserved.

en ENGLISH

Product Description

The disposable Impression Tray is made of plastic and can be used for all types of precision and preliminary impressions. Due to the special design of the tray, the pressure during impression taking is higher than with usual stock trays. This reduces the tendency to form fin-shaped voids caused by the material's characteristics of flow. The tray can be customized with tungsten carbide burs. Three different sizes of Trays for upper and lower jaw each are available.

 Instructions for Use should not be discarded for the duration of product use.

Indications

Suitable impression materials:

- Preliminary impression materials, e.g. Position Penta Quick
- Alginates
- Elastomeric precision impression materials

Types of impression:

- All types of precision and preliminary impressions
- Impressions for the fabrication of temporary restorations
- Opposing jaw impressions
- Impressions for the fabrication of splints

Preparatory Steps

Selection of a Tray of Suitable Size

- ▶ Select a tray and check the size inside the patient's mouth.
 - There must be sufficient space between the wall of the tray and the row of teeth to ensure that there is no permanent deformation introduced when the tray is removed from the row of teeth. As a rule of thumb, this space should be 3–4 times the depth of the undercut of the row of teeth viewed from occlusal.

- The tray can be autoclaved once at 121 °C/250 °F for 15 minutes, if the size does not fit exactly.
- ▶ Interfering areas, e.g. at the lip frenulum or dorsal end of the tray, can be removed with a tungsten carbide bur. Remove any fins or grinding burs.
- ▶ Clean the fitting tray and dry thoroughly with a stream of air.
- ▶ Block the distal end of the tray with wax or putty silicone, if required.

Adhesive

The Impression Tray is equipped with a self-adhesive fleece, making the use of an adhesive unnecessary.

Impression-taking

In the course of positioning common stock trays in the patient's mouth, fin-shaped voids caused by the material's flow characteristics are often introduced. These voids later become apparent in the model.

In contrast, the palatal and lingual areas of the Impression Trays are special-shaped such that the impression material flows towards the teeth when the tray is inserted into the patient's mouth. This increases the pressure during impression taking and reduces flow-related voids. In order to make use of this effect, though, the tray must be filled with impression material all the way up to the retention edge. The upper jaw tray is designed to contain an additional palatal groove serving to take up any excess impression material filled onto the tray. This groove prevents the excessive impression material from flowing into the pharynx, where it would cause a gag reflex.

- ▶ Fill the tray with impression material up to the retention edge.
- ▶ Then slowly insert the tray parallel to the tooth axes.
- ▶ Allow the impression material to set.
- ▶ Once the impression material is completely set, remove the tray from the mouth.

Disinfection of the Impression

- ▶ Place the impression in a standard disinfectant solution. Allow the disinfectant solution to act for the period of time recommended by the manufacturer.
- ▶ After disinfection rinse the impression under running tap water for approx. 15 sec.

Disposal

The Impression Tray is made of polypropylene and can be recycled in accordance with local regulations.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

Limitation of Liability

Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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