DISTRIBUTEURS / DISTRIBUTORS :

FR	Laboratoires Septodont 58, rue du Pont de Créteil 94107 Saint-Maur-des-Fossés Cedex
GB	SEPTODONT Ltd. Units R & S - Orchard Business Centre St Barnabas Close Allington, Maidstone, Kent ME 16 0JZ
ES	LABORATORIOS MAB DENTAL, S.A. C/Córcega 116 SP - 08029 - BARCELONA
DE	SEPTODONT GmbH Felix - Wankel - Str. 9 53859 Niederkassel - Mondorf
IT	GIOVANNI OGNA & FIGLI SpA Via Figini 41 20053 MUGGIO (MI)
BE	Distributeur/Verdeler BENELUX : SEPTODONT NV-SA Grondwetlaan 87 - Av. de la Constitution 87 Brussel 1083 Bruxelles

PL	Przedstawiciel w Polsce : Septodont-Polska sp.z o.o. ul. Kosiarzy 35 - 02-953 Warszawa.
SE	BIGMAN AB Järnväggsgatan 5 SE-172 75 Sundbyberg
US	Septodont, Inc. 245-C Quigley Blvd. New Castle, Delaware 19720
CA	Septodont au Canada Inc. CP 68 - Cambridge, ON N1R 5S9 CANADA
AU	Sponsor : Specialities SEPTODONT PTY Ltd. 14 Cliffbrook Crescent, Emu Plains, NSW 2750

Ultra
Safety Plus

Ultra
Safety Plus **XL**

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INSTRUCTIONS FOR USE - *For professional use only in the practice of dentistry.*

Step 1 : Remove the sterile injectable system from its protective wrapping by tearing it off. Gripping the barrel firmly, fully insert the anaesthetic cartridge into the open end of the injectable system.

Step 2 : Grip the plunger handle, putting the thumb behind the finger holder. Introduce the handle tip into the barrel of the injectable system, behind the cartridge.

Step 3 : Now slide the sheath protecting the needle backwards, towards the handle until it CLICKS (the click is made as the sheath hits the handle and locks the unit together).

NOTE : Failure to retract the sheath fully until you hear the click may result in the system disassembling during use.

Step 4 : All movements now are away from the needle.

Remove the needle cap and discard. The system is now ready to use.

Passive aspiration : At the base of the cartridge barrel of the injectable system you will see there is a small protuberance. It appears as a blob of glue holding the centred needle, the needle-end that penetrates the diaphragm of the cartridge when inserted. At the start of the injection, the diaphragm is pressed against the protuberance, a depression occurs and when released (injection stopped) the diaphragm moves back away from the protuberance; Aspiration occurs.

Active aspiration : Obtained by the silicone top of the plunger handle creating vacuum when, thumb in ring, the practitioner pulls back. The bung follows the plunger tip providing active aspiration, best observed when a minimum of 0.25 ml to 0.35 ml of solution (providing space) has been expelled or used from the cartridge.

Step 5 : When using only one cartridge.

NOTA : During multiple injection procedures using one cartridge the practitioner may safely retain the syringe for further use by movement of the sheath towards the needle until it reaches the holding position (A). Should you need to insert a second cartridge, follow from step 7 onwards.

Complete the procedure, now slide the protective sheath towards the locking position (B) which is the second notch at the end of the barrel. This has now locked the needle safely in the protective sheath.

Step 6 : Removing the plunger handle after use. After you have locked the injectable system in the (B) position, hold the barrel of the injectable system with one hand and using the other hand place a finger in the ring of the plunger handle and pull backwards until the plunger is fully retracted. After you have fully retracted the plunger, peel off the handle in one movement. Now the injectable system can be disposed of safely and the handle autoclaved.

Step 7 : Inserting a second cartridge.

NOTE : During procedures that require more than one cartridge, retract the sheath to the holding position (A) as a needle stick prevention device. Should the system be inadvertently fully locked into position (B), no attempt whatsoever should be made to unlock it, use a new injection system.

Hold the barrel of the injectable system with one hand and using the other hand place a finger in the ring of the plunger handle and pull backwards until the plunger is fully retracted. Now that you have fully retracted the plunger, peel off the handle in one movement.

Step 8 : Take hold of the plunger handle of the injection system and pull the grip handle back towards the ring. Now insert the tip of the plunger into the empty cartridge, which is inside the injectable system. Pull out the cartridge attached to the plunger by the silicone tip, remove cartridge from the plunger and dispose of safely. You are now ready to insert a fresh new cartridge and proceed as from step 1.

PRECAUTIONS FOR USE

1. Ensure protective sheath is attached to the handle with an audible click confirming proper locked on position.

2. Prior to injection remove and dispose of cap. **WARNING: DO NOT RECAP.**

3. WARNING: DO NOT BEND A NEEDLE PRIOR TO INJECTION.

4. If inadvertently when moving protective sheath to the holding position "position A", you miss and reach "position B" the locked position, discard the product and start again with a new piece. **WARNING: DO NOT TRY TO FORCE THE PROTECTIVE SHEATH OFF THE LOCK.**

5. When finished and protective sheath is moved to "position B" and locked by hearing an audible click, remove cartridge by using "replacement cartridge routine" (see step 7). **WARNING: DO NOT TRY TO BREAK OR TEST THE LOCK.** Dispose in separate glass container and the instrument then into a sharps container. The handle is now ready for autoclaving, please note autoclaving instructions.

6. Always place contaminated injection system locked "needle first" into sharps container.

7. Choose sharps container to suit consumption, minimum size 6 litres or more.

WARNING NEVER TAP INSTRUMENT WITH HAND INTO SEMI FULL SHARPS BOX.

8. Always have sharps container adjacent to where procedure is taking place to facilitate immediate and easy disposal of contaminated sharps.

CLEANING AND STERILISATION OF THE HANDLE

Caution: Clean and sterilise the handle before the first and after every use.

- Remove all traces of blood and saliva by brushing and washing the handle.

- Disinfectant solution with a (pH<10) may be used, in accordance with manufactures instructions for diluting and use.

- Sterilize in autoclave between 134° C and 138° C, for the duration of the machine cycle. These cycles may vary according to the regulations and/or the law in force in the country involved.

- Store the Handle in a sterile place.

- End of life is normally determined by wear and damage due to use.

USEFUL TIPS WHEN FIRST USING THE ULTRA SAFETY PLUS XL NEEDLE:



Its larger bore requires less pressure for injection. The excellent flexibility of the XL needle enables gentle insertion of the needle into the tissues.

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