POUCHES FOR STERILIZATION

Formula:

Polyamide+additives.

GENERAL DIRECTIONS

This product is approved for use in sterilization applications in normal steam at temperatures ranging from 250°F (121°C) - 270°F (132°C) for a maximum of 30 minutes and dry heat up to 374°F (190°C) for a maximum of 30 minutes. The melting point of the material is 419°F (215°C). Follow sterilizer manufacturer's directions for sterilization. (Note: change in color of indicator does not guarantee sterilization).

Directions for use:

A/ROLLS:

1. Add an additional 2" if sterilizer tape is used for sealing. (For perforated tubing, tear at perforation then proceed to step 3.)
2. Seal one end of the tubing (to form a bag) using a heat sealer or sterilizer tape. (the selector must be positioned on 3).
3. Place the instrument inside the bag and expel excess air by carefully pressing the bag outward toward the open end, then seal the end with a heat sealer or autoclave tape.
4. Place each sealed bag flat in the sterilizer. Sterilizer must be checked and calibrated periodically.
5. Allow pouches to cool in the sterilizer 5-10 minutes after cycle is completed. To facilitate cooling, open door slightly (follow sterilizer manufacturer's directions).
6. The bagged instruments are to be displayed on trays, side by side, without piling them up and with avoiding contact with the walls

INDICATOR TAPE SEALING

1. Cut the indicator tape twice the width of the tubing plus an additional 1" for the tab.
 2. Lay the tape down with the adhesive surface up.
 3. Place the tubing onto the tape 1/2" from the edge.
 4. Fold remaining indicator tape over the tubing.
 5. check that the colour coded indicator has changed.
- Pink-à Brown

B/SELF-SEALING POUCHES:

Unstick the paper tape protecting the adhesive strip at the extremity of the pouch.

Fold along the dotted line and stick the flap to seal.

1. Place the instrument inside the bag and expel excess air by carefully pressing the bag outward toward the open end.
2. Seal the pouch by folding flap over where it states, "Fold Here" and seal with the tape provided on pouch.
3. Place each sealed pouch flat in the sterilizer. Sterilizer must be checked and calibrated periodically.
4. Allow pouches to cool in the sterilizer 5-10 minutes after cycle is completed. To facilitate cooling, open door slightly (follow sterilizer manufacturer's directions).
5. "Tear at Notch" to open sterilized pouch.
6. The bagged instruments are to be displayed on trays, side by side, without piling them up and with avoiding contact with the walls

WARNING !

Nylon Sterile

Pouches/ROLLS are for single-use only.

Do not reuse.

TROUBLESHOOTING TIPS

Pouches bursts during sterilization

Cause:

- a. The air was not completely expelled before sealing the pouch.
- b. In steam sterilizers, check that air trap jet is clean.
- c. Improper cycle used for sterilization.
- d. Pouch was over filled.

Pouches stick together after the sterilization process

Cause:

- a. Improper sterilization temperature used.
 - b. Folded bags were placed in sterilizer.
- Pouches becomes brittle or cracked

Cause:

- a. Pouch was used more than once.
 - b. Improper temperature and/or time cycle used.
- Moisture in Pouches

Cause:

- a. Pouch was removed from the sterilizer before the cycle was completed.
- b. Improper drying cycle after steam sterilization was used.
- c. A small amount of condensation was present immediately following the drying cycle. This will dissipate within 5-10 minutes.

When ordering, mention either rolls (tubing) or self-sealing pouches.

N.B: In no way, our responsibility should be involved in case of misuse of the sterilizer or a non observance of the maximum temperature as indicated for the different pouches.

General Note: Sterilizers should be checked periodically for correct time, temperature and pressure.

Storage: Store Self Seal Pouches or Rolls at room temperature between 59°-86°F (15°-30°C). Protect from light.

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Hersteller: StePac L.A Ltd. - Tefen Industrial Park, 24959 Israel.