

# shade modification

COLOR MODIFICATION TINTS AND MASKING AGENTS

FARBCHARAKTERISIERUNGS - UND OPAQUERMATERIAL

AGENTES MASCARADORES E MODIFICADORES DE COR

TINTES MODIFICADORES DE COLOR Y AGENTES ENMASCARADORES

AGENTS MASQUANT ET MODIFICATEUR DE TEINTES

MODIFICATORE DI TINTE PER CARATTERIZZAZIONE E OPACIZZANTI

KLEUR MODIFICATIE TINTEN EN MASKERENDE KLEUREN

TINTS TIL FARVEMODIFICERING OG MASKERING

FARGEMODIFISERINGS - OG MASKERINGSSETT

FÄRGNYANS MODIFIERARE ELLER SOM TÄCKFÄRG

ΧΡΩΣΤΙΚΕΣ ΤΡΟΠΟΠΟΙΗΣΗΣ ΧΡΩΜΑΤΟΣ ΚΑΙ ΠΑΡΑΓΟΝΤΕΣ  
ΕΠΙΚΑΛΥΨΗΣ

色調修正マスキングキット

色調修正



SDI

# shade modification

SDI's Tints and Masking Agents simulate a life like, natural tooth

The natural tooth color is derived from the color saturated dentin in conjunction with light penetrating the translucent enamel and undergoing reflection, transmission and absorption

The SDI Shade Modification system simply and quickly achieves optimum esthetics with paint on, single component, light cured liquid Tints and Masking Agents




## Tooth Shades

The selection of the appropriate Masking Agent, Tint and restorative material is controlled by the multichromatic nature of the tooth which consists of four color zones:

1. Cervical zone - Red, Yellow, Orange or Brown
2. Body zone - Grey, Yellow, Universal
3. Interproximal zone - Differing between patients
4. Incisal zone - Grey, Blue

Darker colors in gingival areas become lighter towards the incisal edge.



## Low Viscosity

The ideal low viscosity allows for quick, easy and direct paint on application. The low viscosity resin eliminates the appearance of brush strokes and micron thin layers permit fine veneers.

## Direct Application

The SDI Shade Modification system's syringe and bent tip delivery system facilitates controlled direct extrusion for faster procedures.

## Masking Agents

SDI Masking Agents hide enamel anomalies and discolorations to achieve natural looking restorations.

## Tints

The chart is a resource guide for color choice. Mix shades using the layer technique to create customized color options.

|                            |                                                                                                                       |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Blue or Blue + Grey</b> | Incisal translucency<br>Reduce brightness<br>Cover fluorosis<br>Hide hypocalcification                                |
| <b>Blue</b>                | Cover tetracycline stains<br>Simulate incisal translucency                                                            |
| <b>Yellow</b>              | Cervical darkening<br>Decrease brightness or increase intensity<br>Cover bases                                        |
| <b>Brown + Yellow</b>      | Dark cervical areas<br>Gingival highlights                                                                            |
| <b>Brown</b>               | Cracks<br>Tobacco stains<br>Dark grooves<br>Check line<br>Decrease brightness or increase color intensity             |
| <b>Grey</b>                | Reduce intensity of other shades<br>Incisal translucency<br>Tobacco stains                                            |
| <b>Red</b>                 | Create gingival<br>High smile line<br>Neutralize green<br>To eliminate bluish / grey of teeth stained by tetracycline |
| <b>Blue + Red</b>          | Cover tetracycline stains                                                                                             |
| <b>Orange</b>              | Cervical darkening<br>Create cracks and dark stains<br>Decrease brightness or increase color intensity                |
| <b>White</b>               | Create opacity and masking<br>Mask dark background                                                                    |
| <b>Clear</b>               | Reduce intensity of tints<br>Reduce opacity of masking agent                                                          |
| <b>Green</b>               | To match porcelain with greenish cast                                                                                 |

## Multiple Uses

Nine Tints and four Masking Agents ensure life like restorations are possible for even the most difficult restorative situations.

## Indications

- Enamel • Dentin • Composite • Compomer • Porcelain
- Etched glass ionomer • Metal • Carbon posts


## INSTRUCTIONS

Clean and isolate tooth

- 1 Etch tooth surface with Super Etch 37% phosphoric acid for 20 seconds
- 2 Wash thoroughly with water
- 3 Remove excess water. Keep moist
- 4 Apply Stae dentin/enamel adhesive to saturate all internal surfaces, or other bonding agent according to manufacturer's instructions






- 5 Apply appropriate Masking Agent or Tint in layers less than 0.5mm



- 6 Light cure for 20 seconds. Additional Masking Agent or Tint layers may be applied
- 7 Cover with a thin translucent layer of Ice Incisal, Wave, or other restorative material, according to manufacturer's instructions





# shade modification

## Shade Modification kit

Reorder 7900000

9 x 1g Tints: Blue, Yellow, Clear, Red, White, Orange, Grey, Brown, Green

4 x 1g Masking Agents: Grey, Universal, White, Dark Yellow  
accessories

## Masking Agent 1g refill

Reorder

|         |             |
|---------|-------------|
| 7970202 | Universal   |
| 7970404 | Dark Yellow |
| 7975151 | Grey        |
| 7975252 | White       |

## Shade Tint 1g refill

Reorder

|         |        |
|---------|--------|
| 7981414 | Yellow |
| 7982424 | Grey   |
| 7983333 | Brown  |
| 7985252 | White  |
| 7985353 | Blue   |
| 7985454 | Clear  |
| 7985555 | Green  |
| 7985656 | Orange |
| 7985757 | Red    |

Made in Australia by SDI Limited

Bayswater, Victoria 3153

Australia 1 800 337 003

Austria 00800 022 55 734

Brazil 0800 770 1735

France 00800 022 55 734

Germany 0800 100 5759

Ireland 01 886 9577

Italy 800 780625

New Zealand 0800 734 034

Spain 00800 022 55 734

United Kingdom 00800 022 55 734

USA & Canada 1 800 228 5166

[www.sdi.com.au](http://www.sdi.com.au)



+D036M7900011+

M790001 B