Integra® Miltex® Dental Wax

EN – English	2
FR – Français	4
IT – Italiano	6
ES – Español	8
DE – Deutsch	10
NL – Nederlands	12

Integra® Miltex® Dental Wax

Rx ONLY

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a licensed medical professional.

Description

Integra® Miltex® Beauty Pink Wax and Beauty White Wax are soft impression waxes to record accurate impressions and bite registrations. Beauty Pink Wax may also be used for fabrication of occlusal rims for denture tooth set-up and try-in procedures.

Indications for Use

Integra Miltex Beauty Pink Wax and Beauty White Wax is indicated for use to record the relationship of the upper & lower teeth in dentulous patients. Beauty Pink Wax is also indicated for use in the fabrication of occlusal rims for denture tooth set-up and try-in procedures.

Contraindications

- None Known

Warnings

- Do not reuse. Dental waxes are designed for single use only.
- Used waxes should be considered a biohazard, use a responsible and safe method of disposal.

Precautions

- Read the Instructions for Use and keep them in a safe place.
- Use the product only in accordance with its intended use, see **Indications for Use**.
- Dental waxes should be handled and operated by personnel completely familiar with their use, assembly and disassembly.
- It is recommended to wear appropriate personal protection equipment (PPE) as promulgated by OSHA & AORN when handling, preparing and using dental waxes.
- Do not use a flame to soften the wax.
- Do not overheat or melt the wax. The wax's melting point is approximately 140°F (60°C).
- Store waxes in a cool or refrigerated location.

Adverse Reactions

- None known

Instructions for Use

Occlusal Rim Fabrication

1. Soften the wax sheet in hot water to ensure even heat distribution.
2. For occlusal rim fabrication: Fold the sheet across the short end and manipulate to create a horseshoe shape to attach to the baseplate.
3. Cool to room temperature.
4. The rims are now ready for patient try-in for jaw relation and denture tooth set-up.

Bite registration

1. Soften the wax sheet in hot water to ensure even heat distribution.
2. Fold the sheet across the short end to the desired width and shape for the bite registration required.

3. Place in the patient's mouth and close into occlusion.
4. Allow wax to cool before removing from the patient's mouth.

Storage Conditions

Store in a cool, well-ventilated or refrigerated place. Do not store near strong oxidizers.

Returned Goods Policy

Products must be returned in unopened packages with manufacturer's seals intact to be accepted for replacement or credit unless returned due to a complaint of product defect. Determination of a product defect will be made by Integra. Products will not be accepted for replacement if they have been in the possession of the customer for more than 90 days.

Repairs and Maintenance

Should your instruments require repair or maintenance, contact Integra for return authorization and address. Instruments returned to Integra for repair must have a statement testifying that each instrument has been thoroughly cleaned and sterilized. Failure to supply evidence of cleaning and disinfection will result in a cleaning charge and delayed processing of your instrument repair.

Product Information Disclosure

INTEGRA AND ITS SUBSIDIARIES ("INTEGRA") AND MANUFACTURER EXCLUDE ALL WARRANTIES, EXCEPT INTEGRA'S APPLICABLE STANDARD WARRANTY WHETHER EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. NEITHER INTEGRA NOR MANUFACTURER SHALL BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE, DIRECTLY OR INDIRECTLY ARISING FROM USE OF THIS PRODUCT. NEITHER INTEGRA NOR MANUFACTURER ASSUME NOR AUTHORIZE ANY PERSON TO ASSUME FOR THEM ANY OTHER OR ADDITIONAL LIABILITY OR RESPONSIBILITY IN CONNECTION WITH THESE PRODUCTS.

Symbols Used On Labeling

	Manufacturer ¹
	Catalog Number
	Lot Number
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or practitioner
	Caution! See Warnings and Precautions
	Consult Instructions for use
	Do Not Re-use
	Temperature Limitation
	Authorized Representative in the European Community
	The product meets the essential requirements of the Medical Device Directive 93/42/EEC

¹ Company responsible for a device marketed under its own name regardless of whether "manufactured for" or "manufactured by" the company.

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Miltex, Integra und das Integra-Logo sind eingetragene Handelsmarken der Integra LifeSciences Corporation oder deren Tochtergesellschaften in den Vereinigten Staaten und bzw. oder anderen Ländern. © 2017 Integra LifeSciences Corporation. Alle Rechte vorbehalten.
Integra York PA, Inc.
Integra LifeSciences Services (France) SAS