



Technical Information Clorox Healthcare® VersaSure™ Cleaner Disinfectant Wipes



OVERVIEW

Clorox Healthcare® VersaSure™ Wipes are engineered as a ready-to-use cleaner disinfectant system for healthcare facilities.

These EPA-registered wipes contain quaternary ammonium and other ingredients to kill TB and other healthcare-relevant pathogens in 2 minutes or less on a variety of hard, nonporous healthcare surfaces.

INGREDIENT INFORMATION

Active ingredient:



n-Alkyl Dimethyl Benzyl Ammonium Chloride – 0.25%

n-Alkyl Dimethyl Ethyl Benzyl Ammonium Chloride – 0.25%

AVAILABLE FORMATS

Clorox Healthcare® VersaSure™ Cleaner Disinfectant Wipes are available in multiple sizes and formats to accomplish different disinfecting jobs.

At a Glance Clorox Healthcare® VersaSure™ Cleaner Disinfectant Wipes	
EPA Registration Number	67619-37
Product Classification	Pre-moistened wipes
Dilution	Ready to use
Disinfectant Type	1-step Cleaner-disinfectant
Active Ingredient Chemistry	Quaternary ammonium
Odor	Fragranced
Shelf Life	Expiration date listed on each package: month, day, year
Toxicity Signal Word	Caution
HMIS Rating	2-2-0
PPE Required	Follow your facility's protocols

	Designed for	Item #	Quantity	Shelf life
Multi-Purpose Wipes (6.75"x8") 	<ul style="list-style-type: none">Daily disinfectionPatient room surfacesMedical equipment	31757	85 ct. canister, 6/case	24 months
		31760	30 ct. soft pack, 24/case	24 months
		31758	150 ct. carton, 6/case	24 months
Terminal Wipes (12"x12") 	<ul style="list-style-type: none">Discharge disinfectionLarge areas	31759	110 ct. bucket, 2/case	24 months
		31761	110 ct. refill, 2/case	24 months



Technical Information Clorox Healthcare® VersaSure™ Cleaner Disinfectant Wipes



DISINFECTION EFFICACY CLAIMS

Bacteria - 2 minutes	
<i>Bordetella pertussis</i>	Multidrug-resistant [MDR] <i>Enterococcus faecium</i>
<i>Burkholderia cepacia</i>	Multidrug-resistant <i>Klebsiella pneumoniae</i>
<i>Campylobacter jejuni</i>	<i>Mycobacterium bovis</i> -BCG (TB)
Carbapenem-resistant <i>Escherichia coli</i>	Penicillin-resistant <i>Streptococcus pneumoniae</i>
Carbapenem-resistant <i>Klebsiella pneumoniae</i>	<i>Pseudomonas aeruginosa</i> (Pseudomonas)
Community-associated methicillin-resistant <i>Staphylococcus aureus</i> (CA-MRSA)*	<i>Salmonella enterica</i> (Salmonella)
<i>Enterobacter aerogenes</i>	<i>Serratia marcescens</i>
<i>Escherichia coli</i> (E.coli)	<i>Shigella dysenteriae</i>
<i>Escherichia coli</i> O157:H7	<i>Staphylococcus aureus</i>
<i>Klebsiella pneumoniae</i>	<i>Staphylococcus epidermidis</i> (Coagulase-negative staphylococci)
<i>Legionella pneumophila</i>	<i>Streptococcus pyogenes</i>
Linezolid-resistant <i>Staphylococcus aureus</i> (LRSA)	Vancomycin-intermediate-resistant <i>Staphylococcus aureus</i> (VISA)
Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA)	Vancomycin-resistant <i>Enterococcus faecalis</i> (VRE)
Multidrug-resistant [MDR] <i>Acinetobacter baumannii</i>	Vancomycin-resistant <i>Staphylococcus aureus</i> (VRSA)


Viruses - 30 seconds
Avian influenza Type A Virus (H5N1)
Bovine Viral Diarrhea Virus
Cytomegalovirus
Duck Hepatitis B virus as surrogate for Human Hepatitis B
Herpes Simplex Virus Type 1
Herpes Simplex Virus Type 2
Human Immunodeficiency Virus Type 1
Influenza A (H1N1) Virus
Influenza A Virus
Measles Virus
MERS (Middle East Respiratory Syndrome)-Associated Coronavirus (MERS-CoV)
Mumps virus
Respiratory Syncytial Virus (RSV)
Viruses - 2 minute
Rotavirus
Fungi - 2 minute
<i>Candida albicans</i>



Technical Information Clorox Healthcare® VersaSure™ Cleaner Disinfectant Wipes



DISINFECTION DIRECTIONS FOR USE



- 1.** Don PPE, then REMOVE gross soil if visible.*
- 2.** WIPE surface with wipe until completely wet.
- 3.** To DISINFECT, allow surface to remain wet for contact time listed below:
Bacteria - 2 minutes
Viruses† - 30 seconds

Special Instructions for Using This Product to Clean and Decontaminate Against HIV-1, HBV and HCV on Surfaces/Objects Soiled with Blood/Body Fluids:

- **Personal Protection:** When handling items soiled with blood or body fluids, use disposable impervious gloves, gowns, masks and eye coverings.
- **Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces and other objects before applying this product.
- **Contact Time:** Allow surfaces to remain wet for 1 minute, let air dry. For all other organisms, see directions for contact times.
- **Disposal of Infectious Materials:** Use disposable impervious gloves, gowns, masks and eye coverings. Blood and other body fluids must be autoclaved and disposed of according to local regulations for infectious waste disposal.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that **(1)** is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or **(2)** contacts intact mucous membranes, but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semicritical medical devices prior to sterilization or high-level disinfection.

* This is a one-step product. Precleaning is required only if blood, vomit, fecal matter or other soil is visible, or in advance of disinfecting. Use product as directed. For unfamiliar surfaces, test in a small, inconspicuous area first. Consult medical equipment warranty information and equipment care and cleaning guides prior to selecting and using any surface disinfectant.

† Rotavirus — 2 minute contact time



Technical Information Clorox Healthcare® VersaSure™ Cleaner Disinfectant Wipes



PRODUCTS
& SAFETY

SAFETY INFORMATION

Reference the Safety Data Sheet (SDS) for complete product safety information.

EPA Label Statements

This chemical is a pesticide product registered by the Environmental Protection Agency and is subject to certain labeling requirements under federal pesticide law. These requirements differ from the classification criteria and hazard information required for safety data sheets and for workplace labels of non-pesticide chemicals. Following is the hazard information as required on the pesticide label:

CAUTION: Causes moderate eye irritation. Avoid contact with skin, eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Prolonged or frequently repeated skin contact may cause allergic reactions in some individuals. **FIRST AID: IF ON SKIN:** Take off contaminated clothing. Rinse skin immediately with plenty of water for 15–20 minutes. Call a poison control center or doctor for treatment advice. **IF IN EYES:** Hold eye open and rinse slowly and gently with water for 15–20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor or when going for treatment. You may also contact: 1-800-446-1014 for emergency medical treatment information. **PHYSICAL AND CHEMICAL HAZARDS:** Keep product and empty container away from heat and sources of ignition.

National Fire Protection Association (NFPA) Ratings

Health Hazard: 2

Flammability: 2

Instability: 0

Physical and Chemical Hazards: none

Hazardous Materials Identification System (HMIS) Ratings

Health Hazard: 2

Flammability: 2

Physical Hazard: 0

Personal Protection: X