



Limit uncertainty with precision craftsmanship and accurate administration of dental anesthetics

N-Tralig® Intraligamentary Syringe

The N-Tralig Intraligamentary syringe is designed for the delivery of anesthetic into the periodontal ligament. Intraligamentary injections provide controlled administration of anesthesia to one tooth at a time.¹ An optional 30° angled barrel hub adapter is available to help facilitate injections into more difficult to reach areas of the oral cavity.

- N-Tralig syringe is made of German stainless steel
- Precise delivery of 0.2cc of anesthetic per ratchet¹
- Positive ratchet action helps to protect the anesthetic cartridge against breakage

Reference	Description
76-55SSOB	N-Tralig Intraligamentary Syringe
76-55MSSOB	N-Tralig Intraligamentary Syringe, Metric
PN2SS	N-Tralig Cartridge Barrel Solid Wall, Stainless Steel
PN3SS	N-Tralig Cartridge Barrel Shielded Windows, Stainless Steel
PN4	Replacement Cap for N-Tralig Intraligamentary Syringe
76-56 ²	30° Dental Syringe Adapter
76-56M ²	30° Dental Syringe Adapter, Metric


¹ Refer to the N-Tralig intraligamentary syringe directions for use for important information concerning anesthetic delivery.

² 30° adapter hub is sold separately. Adapter is also compatible with 76-105 Miltex Articulating Barrel Dental Syringe.

Aspirating Syringes

Available in Cooke-Waite (CW) and "A" Type

Standard Reference	Metric Reference	Description
76-70	76-70M	Aspirating Syringe, 1.8cc, Cook-Waite Type (CW)
76-80	76-80M	Aspirating Syringe, 1.8cc, Type "A"
76-105	---	Articulating Barrel Syringe, 1.8cc
76-30	76-30M	Syringe, Aspirating (CW), Blue Silicone Grips
76-31	76-31M	Syringe, Aspirating (CW), Pink Silicone Grips
76-32	76-32M	Syringe, Aspirating (CW), Red Silicone Grips
76-33	76-33M	Syringe, Aspirating (CW), Purple Silicone Grips



Petite Aspirating Syringes

Oval thumb ring provides a reduced diameter compared to a round thumb ring

Flared ring handle base may help to ease pressure while dispensing

Shortened harpoon rod may help to reduce strain during aspiration

Standard Reference	Metric Reference	Description
77-70	---	Petite Aspirating Syringe, 1.8cc, Cook-Waite Type (CW)
76-40	76-40M	Syringe, Petite Aspirating (CW), Blue Silicone Grips
76-41	76-41M	Syringe, Petite Aspirating (CW), Pink Silicone Grips
76-42	76-42M	Syringe, Petite Aspirating (CW), Red Silicone Grips
76-43	76-43M	Syringe, Petite Aspirating (CW), Purple Silicone Grips



77-70

Self-Aspirating Syringes

Blunt harpoon rod design helps facilitate aspiration by the means of counter pressure

Standard Reference	Metric Reference	Description
76-90	76-90M	Lightweight, Self-Aspirating Syringe, 1.8 cc
76-100	---	Self-Aspirating Syringe, 1.8cc



76-100

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

integralife.com/contact

