



Replacement Bulbs

Photo Curing Bulbs

Item#	Bulb Type/ANSI	Volts/Watts	Equipment Model
QUA-175	JCR/M12V75W	12V/75W	Apollo CP (Lares); Aristocrat I (Healthco); Coltolux II, 3, 4 & 75 (Coltene); Compolite 660 (Motloid); Cromalux 100 (Mega Physik); Degulux (Degussa); Elipar 2500 (3M); Executor (Pro Den); Heliolux GT & DLX (Ivoclar); Henry Schein Economy Unit & CU-1000; Hilux 200 & 350 (First Medica); Litex 660 & 680 (Dentamerica); LTI-5000 (Lasac); Megalux (Megadental); Nou-Lite (Nouvag); Omega 750 (Ryan Ind.); Op II (Almore); Optilux 400, 401 & 403 (Kerr); Ortholux & Ortholux XT (Unitek); Penn Med 300 & 310; PowerLite 75 & 100 (Penn-Med); QHL75 (Caulk); Spectrum 800 (Caulk); Super-Dent 2001 & 2500 (Darby); VCL 100, 200 & 300 (Kerr); Visilux I & II (3M); Wite Lite 1500 (Wkyle); XL3000 (3M)
QUA-195	JCR/M12V100W	12V/100W	Ceri-Light (DenMat); EFOS Lite; Focus (Teledyne Getz); Prismetics Lite (Caulk); Profile VCS (SS White); Sunspot (Kinetic Instr); Translux CL (Kulzer)
QUA-205	JCR/M14V35W	14V/35W	Aristocrat II (Healthco); EFOS-35; Elcelux (Weil Dental); FARO P287; Heliolux I & II (Ivoclar); Max Lite (Caulk); Optilux 35, 50 & 150 (Kerr); Patterson TCL 390
QUA-220	EFN	12V/75W	Fotofil (J&J)
QUA-225		12V/80W	Maxima 480 & Maxima Pure Power (Schein); Optilux 500 & 501 (Kerr)
QUA-230		10V/52W	Arcus II (Litema); Hilux 100 (First Medica); Jetlite 4000 (J. Morita); Optilux 180 & 360 (Kerr); Spectrum (Caulk); XL 1500 (3M)
QUA-190	EFP	12V/100W	Opticure (NSK)
QUA-165	EFR	15V/150W	Coe Light; Denta Cam II (Fuji Optical); Heliomat H10 European (Ivoclar); Initiator (Solid State Systems); Pluraflex HL 150 (Litema); Superlight (Bosworth); Translux II, Translux EC & EC-S (Kulzer); VLC High Intensity (Healthco)
QUA-185	EJA	21V/150W	Coby Light; Dual Lux III (Sci-Can); Heliomat Single (Ivoclar); Marathon (Den-Mat); Max Plus (Caulk)
QUA-180	EJL	24V/200W	Prisma Lite (Caulk)
QUA-160	EJM	21V/150W	Command Light (Kerr)
QUA-200	EJV	21V/150W	200L Wand (Penn-Med); Fiberlite 170 & 176 (Dolan-Jenner); Marathon II (Den-Mat); Sunlite II (Kinetic Instr); Super-Dent 7500 (Darby); Triad TCUII (Dentsply); Visar II (Den-Mat)
QUA-215	ESD	120V/150W	200L Original Wand (Penn-Med); Marathon (Den-Mat); Polylite 500, 1000 & 2000 (Pro-Den); Super-Dent 5000 & 8000 (Darby); Witelite Wand (Phasealloy)



QUA-165



QUA-185

Operatory Bulbs

Item#	Bulb Type/ANSI	Volts/Watts	Equipment Model
QUA-120	Q150T4/CL/25	25V/150W	Pelton & Crane
QUA-125	Q75T3/CL/DC-28V	28V/75W	Star Lite Model-J (Ritter)
QUA-100	01007 (H3)	12V/55W	1300 Light (Summit Dental Systems); Aristocrat & Celebrity (Healthco); Belmont; Celux (Dentamerica); Chayes Virginia; Danserau; Performer (A-dec); SDS Fantastic (Summit Dental Systems); The Light (Fiad)
QUA-110	DZE/FDS	24V/150W	Daystar S-90B (Porter); Forest Dental 9000; Henry Schein Starbrite & Crusader; Knight; Koolbeam Light 45 (FARO); LuxStar (Marus); Sunlight 70, S-90 & 2000 (FARO); Super Star Lite Model-K (Ritter); Vacudent
QUA-180	EJL	24V/200W	Prisma Lite (Caulk); Rolux (Castle/MDT)
QUA-115	FCS	24V/150W	AXCS Light (DentalEZ); Delight (Planmeca); Eurodent; J. Morita; SDS 1350 (Summit Dental Systems); Siemens A1/E1 (Sirona); Sunlight 70 Original (FARO)
QUA-130	501-074	24V/60W	Belmont Clesta
QUA-131	Philips 14623	17V/95W	A-dec 6300; Midmark Dental Light



QUA-120

Fiber Optic Handpiece & Intraoral Imaging System Bulbs

Item#	Bulb Type/ANSI	Volts/Watts	Equipment Model
QUA-220	EFN	12V/75W	Fotofil (J&J); Logicam II (Synca); Perspective Dental Imaging (Dentsply); Spectra Vu (Vision Tools); Stylus 240 & 2000 (Cygnus); Trophy Intraoral Camera (Practice Works)
QUA-190	EFP	12V/100W	AcuCam; KS102 (Panasonic); Optica (KaVo); Optica Fiber Optic Unit (NSK); Prism (Pro Dentec)
QUA-165	EFR	15V/150W	DentaCam II, 50 & 100 (Fuji Optical); Vision Plus (Fuji Optical); Vista Cam (Air Techniques)
QUA-185	EJA	21V/150W	Oral Scan (Lester A. Dine)
QUA-210	EKE	21V/150W	Dual Lux II Illuminator (SciCan); Midmark 354 (Ritter); Ultraeye (Trojan)
QUA-135	—	3.0-3.3 DC/2.5W	KaVo LCM Illuminator Bulb #553-3881
QUA-170	ETJ	120V/250W	AcuCam; CygnaScope (Cygnus Tech); Dual Lux II (Sci-Can); Insight I & II (Midwest); Oral Vision; Polar II (KaVo); Quartz Fire SX (Lares); Tradition I (Midwest); Vicon DLS, LS & Polar II (KaVo); VLC Illumination (Healthco)



QUA-125

X-Ray Viewing Box Bulbs

Item#	Bulb Type/ANSI	Volts/Watts	Equipment Model
QUA-140	F8T5/D	—	Rinn View Boxes; Star Endodontic Viewer (Star X-Ray); R130 Cephalometric (Star Dental) (8 Watt Daylight - 5/8" Diameter Tube x 12" Length)
QUA-145	F14T8/D	—	DE100 & DE200 (Healthco); DE100, DE120, DE150, DE160, DE200 & SW100 (Star X-Ray); Executive & Table Top Viewer (Pelton & Crane); Henry Schein #94, #100, #200
QUA-155	F14T12/CW	—	Pelton & Crane Spirit & Executive Cabinet (14 Watt Cool White - 1 1/2" Diameter Tube x 15" Length)



QUA-110