



## VK Solutions and Tablets®

Powerful cleaning with minimal effort.

- Powerful ultrasonic energy quickly removes and breaks down hard to clean deposits like blood, bone, and food particles that chemical detergents can leave behind.
- Powerful, yet gentle on your instruments.
- Biodegradable and pH neutral\*.
- Non-toxic, non-corrosive formula contains no chlorine or harsh alkaline to damage your instruments\*.

\*Refers to VibraKleen E2 Ultrasonic Cleaning Tablets.

Giving a hand to oral health.



**KULZER**  
MITSUI CHEMICALS GROUP

# Your Complete Ultrasonic Cleaning System

Product	Features	Benefits	Uses	Packaging	Product Numbers
<b>VIBRAKLEEN E<sup>2</sup>® ULTRASONIC CLEANING SOLUTIONS</b>					
<b>VibraKleen E<sup>2</sup> Ultrasonic Cleaning Solutions</b>	<ul style="list-style-type: none"> <li>• Effervescent, enzymatic cleaning tablets</li> <li>• Non-corrosive and low-foaming</li> <li>• Biodegradable, non-toxic &amp; pH neutral</li> <li>• No chlorine or harsh alkaline</li> <li>• Drop directly into ultrasonic cleaner</li> <li>• Compact, convenient package of 20 tablets dilutes to 10 gallons</li> <li>• Clean mint scent</li> <li>• Effective as evacuation cleaner</li> <li>• Power of protease enzymes</li> </ul>	<ul style="list-style-type: none"> <li>• Breaks down hard-to-clean deposits in ultrasonic and evacuation systems</li> <li>• Easy on instruments</li> <li>• Environmentally friendly and safe to work with</li> <li>• Won't damage stainless steel equipment</li> <li>• No messy mixing of liquid concentrates</li> <li>• Saves on counter and storage space</li> <li>• Minimizes any odor from tank, pleasant to work with</li> <li>• One product for multiple uses</li> <li>• Effective, fast removal of blood, bone and food particles that chemical detergents leave behind</li> </ul>	<ul style="list-style-type: none"> <li>• For general purpose ultrasonic cleaning, evacuation system cleaning and pre-soak</li> </ul>	Box of 20 Box of 80  2 tablets = 1 gallon	<b>50037100</b> — Box of 20 (Case Lot 24) <b>50037101</b> — Clinic Pack— Box of 80 (Case Lot 12)
<b>DRI-CLAVE® VK ULTRASONIC CLEANING SOLUTIONS</b>					
<b>VK-1 General Purpose Cleaner</b>	<ul style="list-style-type: none"> <li>• Efficient cleansing time (2–10) minutes</li> <li>• 10:1 concentrated formulation</li> <li>• Color coded solution</li> <li>• Anticorrosive</li> </ul>	<ul style="list-style-type: none"> <li>• Protects your equipment</li> <li>• Avoid messy cleaning</li> <li>• Economical gallon container</li> <li>• Mistake-free identification</li> </ul>	<ul style="list-style-type: none"> <li>• Removes debris from instruments, burs, general equipment, ortho appliances, gold and porcelain</li> <li>• Use directly in ultrasonic tank</li> </ul>	Gallon— blue Liquid	<b>50036202*</b> <b>4000239*</b> (Canadian customers only)
<b>VK-3 Plaster and Stone Remover</b>	<ul style="list-style-type: none"> <li>• Efficient cleansing time (5–10) minutes</li> <li>• No diluting necessary—use at full strength</li> <li>• Color coded solution</li> </ul>	<ul style="list-style-type: none"> <li>• Protects your instruments</li> <li>• Avoid messy cleaning</li> <li>• Saves time</li> <li>• Economical gallon container</li> <li>• Mistake-free identification</li> </ul>	<ul style="list-style-type: none"> <li>• For removal of plaster and stone</li> <li>• Do not use with aluminum</li> <li>• Place in glass breaker inside ultrasonic cleaner</li> </ul>	Gallon— pink Liquid	<b>50036206*</b> <b>4000240*</b> (Canadian customers only)
<b>VK-4 Tartar and Stain Remover</b>	<ul style="list-style-type: none"> <li>• Efficient cleansing time (5–10) minutes</li> <li>• No diluting necessary—use at full strength</li> </ul>	<ul style="list-style-type: none"> <li>• Avoid messy cleaning</li> <li>• Ready-to-use formula</li> <li>• Saves time</li> <li>• Economical gallon container</li> </ul>	<ul style="list-style-type: none"> <li>• Removal of tartar and stains</li> <li>• 5–10 minutes only</li> <li>• Place in glass beaker inside ultrasonic cleaner</li> </ul>	Gallon— clear Liquid	<b>50036208*</b> <b>4000241*</b> (Canadian customers only)
<b>VK-6 Temporary Cement Remover</b>	<ul style="list-style-type: none"> <li>• Efficient cleansing time (5–10) minutes</li> <li>• No diluting necessary—use at full strength</li> <li>• Color coded solution</li> </ul>	<ul style="list-style-type: none"> <li>• Protects your instruments</li> <li>• Saves time</li> <li>• Economical gallon container</li> </ul>	<ul style="list-style-type: none"> <li>• Removal of temporary cement including cavit and eugenol-based temporary cement</li> <li>• 5–10 minutes only</li> <li>• Place in glass beaker inside ultrasonic cleaner</li> </ul>	Gallon— yellow Liquid	<b>50036212*</b> <b>4000242*</b> (Canadian customers only)
<b>VACUKLEEN E<sup>2</sup>® EVACUATION SYSTEM CLEANER</b>					
<b>VacuKleen E<sup>2</sup> Evacuation System Cleaner</b>	<ul style="list-style-type: none"> <li>• Highly effective powdered formula</li> <li>• Non-corrosive and low-foaming</li> <li>• Biodegradable &amp; non-toxic</li> <li>• Convenient, easy-mix packets</li> <li>• Power of protease enzymes</li> <li>• Fresh, clean, scent of lemon</li> <li>• Simple to use</li> </ul>	<ul style="list-style-type: none"> <li>• Removes hard-to-clean protein-based deposits that soapy solutions can leave behind</li> <li>• Environmentally friendly and safe to work with</li> <li>• Pre-measured packets—saves time and minimizes waste</li> <li>• Effective and quick dissolving of blood, tissue, saliva and encrusted residue in evacuation lines</li> <li>• Real lemon extracts gives a clean smell that eliminates evacuation system odor and makes the operatories smell fresh</li> <li>• Add powder packet to warm water to make a low-foaming solution</li> </ul>	<ul style="list-style-type: none"> <li>• Cleaning of your evacuation equipment</li> <li>• Only mix what you need so you eliminate waste and storage problems</li> </ul>	Box of 25 pouches Box of 60 pouches Box of 140 pouches	<b>50036100</b> — Unit Dose Pouch for Quarts Intro Pack— Box of 60 Pouches (Case Lot 24) <b>50036101</b> — Unit Dose Pouch for Quarts Clinic Pack— Box of 140 Pouches (Case Lot 12) <b>50036102</b> — Unit Dose Pouch for Gallons Clinic Pack— Box of 25 Pouches (Case Lot 24)

\* Each VK Solution Sold in Case Lots of Four (4) Gallons Only.