



ABRASIVES



DIAMONDS

SHOFU DENTAL CORPORATION

RETAIL PRICE GUIDE

Visit www.shofu.com | Call 800.827.4638




CEMENTS



PREVENTIVES



RESTORATIVES



CAD/CAM





Dear Valued Shofu Customer:

All of us at Shofu Dental Corporation want to thank you for your commitment to Shofu products. Your continued support enables Shofu Dental Corporation to increase our development of innovative products that add value for you and your staff.

Shofu has been known for our quality for more than 93 years in the dental industry. Our finishers and polishers, Brownies, Greenies, DuraWhites and porcelain remain top in the industry, and because of your commitment we have expanded and added new categories; including our state of the art composites with our Giomer technology. A complete list of our products is highlighted in this price guide or on our website www.shofu.com.

In presenting our products to the dental community, we offer very competitive pricing so that the value to you will be unsurpassed. If you compare Shofu products with those of our competitors, you will find the highest quality and lower prices.

For additional information, please contact your preferred Shofu dealer, our Customer Service Department at 1-800-82-SHOFU (800-827-4638) or your Territory Sales Manager.

Thank you again,

Randy Bailey
President



SHOFU DENTAL CORPORATION

TABLE OF CONTENTS

ABRASIVES

3/29



DIAMONDS/CARBIDES 43



CEMENTS

55



PREVENTIVES

65



RESTORATIVES

67



TECHNICAL DATA

81




SHOFU DENTAL CORPORATION

DEFINITIONS

Abbreviations used:

- **PN: PART NUMBER**
Number used to place an order
- **ISO#: INTERNATIONAL ORGANIZATION FOR STANDARDIZATION**
Numbers to identify dental rotary instruments – Nominal diameters and designation code
- **CA: CONTRA-ANGLE SHANK**
- **HP: HANDPIECE SHANK**
- **FG: FRICTION GRIP SHANK**



CUSTOMER SERVICE



Shofu Dental Corporation is located in San Marcos, California. Our hours of operation are Monday – Friday from 7:00 am to 5:00 pm P.S.T. You can reach us for questions or concerns by phone, fax or email.

1.800.827.4638 – phone

1.760.736.3276 – fax

customer-service@shofu.com – email

WEBSITE



Our website is designed to serve as a resource for dental professionals. It contains current product information and educational videos that can be found quickly with our easy-to-use navigation, which makes searching for what you need a breeze.

Visit us on the web: **www.shofu.com**

While every attempt has been made to ensure accuracy in this publication, Shofu Dental Corporation is not responsible for possible printing or other errors. Please confirm contents and pricing with your Sales or Customer Service Representative. Thank you!

ABRASIVES

STONES

DURA-GREEN & DURA-GREEN FINE STONES	4 - 5
NEW! DURA-GREEN DIA STONES (NEW FG SHAPES!)	6
DURA-WHITE STONES	7
BROWNIE, GREENIE & SUPERGREENIE	8

POLISHERS

CERAMASTER COARSE & CERAMASTER	9
CERAMISTÉ & CERAMISTÉ SOFT	10
COREMASTER COARSE & COREMASTER	11
COMPOMASTER COARSE & COMPOMASTER	12
COMPOSITE & COMPOSITE FINE	13
SUPER-SNAP FINISHING & POLISHING	14
SUPER-SNAP DISK REPLACEMENTS	15
SUPER-SNAP POLYSTRIPS	15
SUPER-SNAP X-TREME FINISHING & POLISHING	16
SUPER-SNAP SUPERBUFF DISK	17
SUPER-SNAP BUFF DISK & MINI BUFF DISK	17
SUPER-SNAP METAL & PLASTIC MANDRELS	18
ONEGLOSS & ONEGLOSS PS	19

LAB STONES & POLISHERS

ACRYPOINT	20
LAB SERIES STONES	21
HARDIE	22
UNMOUNTED WHEELS	23

POLISHING PASTES

DURA-POLISH & DURA-POLISH DIA	24
COMPOSITE POLISHING PASTE & DIRECTDIA POLISHING PASTE (INTRAORAL)	24

ACCESSORIES

METAL & PLASTIC MANDRELS	18
STERILIZABLE BUR BLOCKS	25
NEW! BURBUTLER UNIVERSAL SILICONE BUR BLOCK (NEW SIZES!)	26 - 27

ABRASIVES

DURA-GREEN® STONES



- Impregnated with silicon carbide
- Long lasting, fast cutting & always concentric in motion
- **Use Dura-Green Stones for:** composites, amalgams, precious alloys, non-precious alloys & porcelain
- **Use Dura-Green Fine Stones for:** composites & porcelain
- See page 23, Unmounted Finishing Wheels section, for Dura-Green Wheels

SHAPE CATEGORY	CN1	CN2	CN3	CY2	CY3	CY4	FL1	FL2	FL3	FL4	FL5
SHAPE											
ISO#	024	037	040	025	028	045	030	025	030	045	065
HP (12 PK)	PN 0002	-	-	0006	0007	0008	0009	0010	0011	0012	0013
HP FINE (12 PK)	PN -	-	-	-	-	-	-	-	011F	-	-
HP (72 PK)	PN -	-	-	-	-	-	-	HG10	-	-	-
CA (12 PK)	PN 0052	0053	0054	-	0057	-	0059	0060	0061	0062	-
FG (12 PK)	PN 0101	-	-	0103	-	-	-	0104	-	-	-

SHAPE CATEGORY	IC1	IC2	IC3	IC4	IC5	IC6	IC7	IC8	IC9	KN1	KN3
SHAPE											
ISO#	030	032	050	045	070	065	125	104	135	050	070
HP (12 PK)	PN 0014	0015	0016	0017	0018	0019	0020	0021	0022	0023	0025
HP FINE (12 PK)	PN -	-	-	-	-	-	-	-	022F	-	-
HP (72 PK)	PN HG14	-	HG16	-	HG18	-	HG20	-	HG22	-	-
CA (12 PK)	PN 0064	-	0066	-	0068	-	-	-	-	-	-
FG (12 PK)	PN -	0105	-	-	-	-	-	-	-	-	-

SHAPE CATEGORY	KN4	KN5	KN6	KN7	PC1	PC2	RD1	RD2	RD3	RD4	RE1
SHAPE											
ISO#	080	060	100	125	030	050	030	024	036	050	062
HP (12 PK)	PN 0026	0027	0028	0029	0030	0031	0032	0033	0034	0035	0036
HP FINE (12 PK)	PN -	-	-	-	030F	031F	-	-	-	-	-
HP (72 PK)	PN -	-	-	HG29	-	-	-	-	-	-	-
CA (12 PK)	PN -	-	-	-	0080	0081	0082	0083	0084	0085	0086
FG (12 PK)	PN -	-	-	-	-	-	0107	0108	-	-	-



SHAPE CATEGORY		RE2	TC1	TC2	TC3	TC4	WH1	WH2	WH3	WH4	WH5	WH6
SHAPE												
ISO#		062	020	025	030	040	040	050	065	060	065	130
HP (12 PK)	PN	0037	0039	0040	0041	0042	0043	0044	0045	0046	0047	0048
HP FINE (12 PK)	PN	-	039F	-	-	042F	-	-	045F	-	-	048F
HP (72 PK)	PN	-	-	-	-	HG42	-	-	HG45	-	-	HG48
CA (12 PK)	PN	-	0089	0090	-	0092	-	-	0095	0096	-	0098
FG (12 PK)	PN	-	0109	-	-	-	0110	-	-	-	-	-

PRICING - DURA-GREEN® STONES

PN

Product

Suggested Retail

See chartsHP Shank (12 per pack)\$27.10

See chartsHP Shank (1/2 Gross/One Box)\$138.35

See chartsCA & FG Shank (12 per pack)\$30.30

*See page 23, Unmounted Finishing Wheels section, for Dura-Green Wheels

ABRASIVES

NEW! DURA-GREEN® DIA STONES



- NEW FG SHAPES!
- Vitrified bonded diamond abrasive
- Improved cutting, better than carborundums
- Reduced chipping & contamination
- Use **Dura-Green DIA Stones** for: contouring zirconia, porcelain & ceramics

SHAPE CATEGORY	CY5	FL3	IC5	KN7	RD2	RD3	TC2	TC4	WH6
SHAPE									
ISO#	056	030	070	125	024	036	025	025	130
HP (2 PK*/3 PK)	PN 0156*	0158	0162	0154*	-	0159	0161	0155	0157*
CA (3 PK)	PN -	0163	-	-	-	0164	0165	-	-
FG (3 PK)	PN -	0168	-	-	0167	-	-	-	-

PRICING - DURA-GREEN® DIA STONES

PN	Product	Suggested Retail
See charts	HP, CA, FG Shank	\$54.10
0153	HP Shank Dura-Green DIA Assortment	\$54.10
CY5, KN7, TC4		
0166	CA Shank Dura-Green DIA Assortment	\$54.10
FL3, RD3, TC2		
0504	Shofu Adjuster	\$84.20



DURA-WHITE® STONES



- Micrograined aluminum oxide grit
- Vibration-free performance
- Long lasting
- **Use Dura-White Stones for:** fine finishing on enamel, composites, glass ionomer restoratives & porcelains

SHAPE CATEGORY		CN1	CY2	FL1	FL2	FL3	IC2	IC4	KN1	KN3	RD1	RD2
SHAPE												
ISO#		024	025	030	025	030	032	045	050	070	030	024
HP (12 PK)	PN	0201	-	0202	0203	0204	-	0205	-	0206	0207	0208
CA (12 PK)	PN	0221	-	0222	0223	0224	-	0225	-	-	0227	0228
FG (12 PK)	PN	0241	0243	0240	0244	-	0245	-	0246	-	0247	0248

SHAPE CATEGORY		RE1	TC1	TC4	WH1	WH4
SHAPE						
ISO#		062	020	040	040	060
HP (12 PK)	PN	0209	0210	0211	-	0212
CA (12 PK)	PN	0229	0230	-	-	0232
FG (12 PK)	PN	-	0249	-	0250	-

PRICING - DURA-WHITE® STONES

PN

Product

Suggested Retail

See charts HP, CA, FG Shank (12 per pack) \$30.30

ABRASIVES

BROWNIE®/GREENIE®/SUPERGREENIE POLISHERS



- Silicone polishers impregnated with high quality abrasives
- Use Brownies for pre-polishing, Greenies for polishing & Supergreenies for super-polishing
- Use Brownie/Greenie/Supergreenie Stones for: amalgams, precious alloys & non-precious alloys

SHAPE CATEGORY	Cup	Mini-Pt.	KN7	PC2	WH6	Floppie	Disk	Unm. Wheel
SHAPE								
ISO#	065	030	125	050	120	220	160	220
BROWNIE								
HP (12 PK)	PN -	-	0271	0272	0273	0458	-	0453
CA (12 PK)	PN 0401	0403	0281	0282	0283	-	0405	-
CA (72 PK)	PN H401	H403	-	-	-	-	-	-
FG (12 PK)	PN 0411	0413	-	-	-	-	-	-
FG (72 PK)	PN H411	H413	-	-	-	-	-	-
GREENIE								
HP (12 PK)	PN -	-	0274	0275	0276	0459	-	0454
CA (12 PK)	PN 0402	0404	0284	0285	0286	-	0406	-
CA (72 PK)	PN H402	H404	-	-	-	-	-	-
FG (12 PK)	PN 0412	0414	-	-	-	-	-	-
FG (72 PK)	PN H412	H414	-	-	-	-	-	-
SUPERGREENIE (yellow band)								
HP (12 PK)	PN -	-	0277	0278	0279	-	-	454A
CA (12 PK)	PN 402A	404B	0287	0288	0289	-	-	-
FG (12 PK)	PN 412A	414B	-	-	-	-	-	-

PRICING - BROWNIE®/GREENIE®/SUPERGREENIE POLISHERS

PN	Product	Suggested Retail
See chart.	HP, CA, FG Shank (12 per pack)	\$35.85
	KN7, PC2, WH6, Floppie, Disk	
See chart.	CA & FG Shank (12 per pack)	\$33.55
	Cup & Mini-Point	
See chart.	CA & FG Shank (1/2 Gross/OneBox)	\$165.45
	Cup & Mini-Point	
See chart.	Unmounted Wheels (12 per pack)	\$21.80



CERAMASTER®/CERAMASTER® COARSE POLISHERS



- Diamond impregnated silicone polishers
- Saves time by providing a simple, two-step system
- Use CeraMaster Coarse to polish to a satin finish & CeraMaster to super-polish
- Use CeraMaster for: porcelains, enamel & zirconia


SHAPE CATEGORY	Bullet	Knife	Cup	Mini-Pt.	KN7	WH6	Unm. Wheel	Kits
SHAPE								See Assort. Below
ISO#	040	065	090	030	125	120	220	
REGULAR (blue & white band)								
HP (3 PK)	PN 0126	-	-	0129	0127	0128	-	0130
CA (3 PK)	PN 0121	0122	0123	0124	-	-	-	0125
FG (3 PK)	PN -	-	-	0136	-	-	-	
COARSE (no band)								
HP (3 PK)	PN 126C	-	-	129C	127C	128C	132C (1 per pack)	130C
CA (3 PK)	PN 121C	122C	123C	124C	-	-	-	125C
FG (3 PK)	PN -	-	-	136C	-	-	-	

PRICING - CERAMASTER®/CERAMASTER® COARSE POLISHERS

PN	Product	Suggested Retail
0125	Assortment Package - CA Shank	\$67.50
	Bullet, Knife, Cup, Mini-Point	
0130	Assortment Package - HP Shank	\$67.50
	Bullet, Mini-Point, KN7, WH6	
125C	Coarse Assortment Package - CA Shank	\$67.50
	Bullet, Knife, Cup, Mini-Point	
130C	Coarse Assortment Package - HP Shank	\$67.50
	Bullet, Mini-Point, KN7, WH6	
See chart	Mini-Point FG Shank Polishing Points (3 per pack)	\$48.85
See chart	Polishing Points (3 per pack)	\$50.60
See chart	Unmounted Wheel (1 per pack)	\$53.70

ABRASIVES

CERAMISTÉ®/CERAMISTÉ® SOFT POLISHERS



- Silicon carbide impregnated polishers
- Three-step system
- Use Ceramisté Standard for pre-polishing, Ultra (yellow band) for polishing & Ultra II (white band) for super-polishing
- Use Ceramisté Soft PA (green band) & Ceramisté Soft PB (blue band) for smoothing rough surfaces & fragile porcelain margins
- Use **Ceramisté for:** porcelains & enamel


SHAPE CATEGORY	Bullet	Cup	Knife	Wheel	Midi-Pt.	KN7	PC2	WH6	Unm. Wheel
SHAPE									
ISO#	040	065	090	090	060	125	050	120	220
STANDARD (no band)									
HP (12 PK)	PN	-	-	-	-	0261	0262	0263	0451
CA (12 PK)	PN	0251	251A	0252	0253	253A	-	-	-
FG (12 PK)	PN	-	251B	-	-	253B	-	-	-
ULTRA (yellow band)									
HP (12 PK)	PN	-	-	-	-	0264	0265	0266	0452
CA (12 PK)	PN	0254	254A	0255	0256	256A	-	-	-
FG (12 PK)	PN	-	254B	-	-	256B	-	-	-
ULTRA II (white band)									
HP (12 PK)	PN	-	-	-	-	0267	0268	0269	452A
CA (12 PK)	PN	0257	257A	0258	0259	259A	-	-	-
FG (12 PK)	PN	-	257B	-	-	259B	-	-	-
SOFT PA (green band)									
HP (12 PK)	PN	-	-	-	-	-	-	0564	0595
SOFT PB (blue band)									
HP (12 PK)	PN	-	-	-	-	-	0563	0565	0596

PRICING - CERAMISTÉ®/CERAMISTÉ® SOFT POLISHERS

PN	Product	Suggested Retail
See chart.....	HP, CA, FG Shank (12 per pack)	\$38.10
See chart.....	Unmounted Wheel (12 per pack)	\$22.55
452A.....	HP Unmounted Wheel (12 per pack)	\$22.20
See chart.....	Soft PA/PB (12 per pack)	\$44.45
See chart.....	Soft PA/PB Unmounted Wheel (12 per pack)	\$31.75



COREMASTER®/COREMASTER® COARSE POLISHERS



- Diamond impregnated polishers eliminate microfractures
- Two grit sizes & three shapes
- Contouring & finishing
- No water cooling needed
- **Use CoreMaster for:** zirconia frameworks, alumina frameworks & porcelains


SHAPE CATEGORY	WH6
SHAPE	
ISO#	120
REGULAR (red grit)	
HP (2 PK)	PN 0140
COARSE (black grit)	
HP (2 PK)	PN 140C

PRICING - COREMASTER®/COREMASTER® COARSE POLISHERS

PN	Product	Suggested Retail
140C	HP Shank (2 per pack)	\$42.15
0140	HP Shank (2 per pack)	\$43.85

ABRASIVES

COMPOMASTER®/COMPOMASTER® COARSE POLISHERS



- Carefully balanced blend of diamond particles
- Produces the smoothest finish & luster
- A simple two-step system
- Use **CompoMaster for**: composite restorations

SHAPE CATEGORY	Bullet	Knife	Cup	Mini-Point	KN7	WH6	Unm. Wheel
SHAPE							
ISO#	040	090	065	030	125	120	220
REGULAR (white & yellow band)							
HP (3 PK)	PN 0116	-	-	0119	0117	0118	-
CA (3 PK)	PN 0111	0112	0113	0114	-	-	-
FG (3 PK)	PN -	-	-	0134	-	-	-
COARSE (no band)							
HP (3 PK)	PN 116C	-	-	119C	117C	118C	131C (1 per pack)
CA (3 PK)	PN 111C	112C	113C	114C	-	-	-
FG (3 PK)	PN -	-	-	134C	-	-	-

PRICING - COMPOMASTER®/COMPOMASTER® COARSE POLISHERS

PN	Product	Suggested Retail
0115	Assortment Package - CA Shank	\$67.50
	Bullet, Knife, Cup, Mini-Point	
0120	Assortment Package - HP Shank	\$67.50
	Bullet, Mini-Point, KN7, WH6	
115C	Coarse Assortment Package - CA Shank	\$67.50
	Bullet, Knife, Cup, Mini-Point	
120C	Coarse Assortment Package - HP Shank	\$67.50
	Bullet, Mini-Point, KN7, WH6	
See chart	Unmounted Wheel (1 per pack)	\$53.70
See chart	Mini-Point FG Shank Polishing Points (3 per pack)	\$48.85
See chart	Polishing Points (3 per pack)	\$50.60



COMPOSITE/COMPOSITE FINE POLISHERS



- Use on microfilled & hybrid composites
- Remove orthodontic adhesives
- Two-step system to pre-polish & polish
- Use **CompoSite for**: composites

SHAPE CATEGORY		Bullet	Knife	Wheel	Disk	Cup	Pear-Point	IC	Midi-Point
SHAPE									
ISO#		040	090	090	160	080	033	-	040
REGULAR (no band)									
CA (12 PK)	PN	0191	0192	0193	0194	0195	0196	0197	0198
FG (12 PK)	PN	-	-	-	-	-	196B	-	0199
FINE (white band)									
CA (12 PK)	PN	0291	0292	0293	0294	0295	0296	0297	0298
FG (12 PK)	PN	-	-	-	-	-	296B	-	0299

PRICING - COMPOMASTER®/COMPOMASTER® COARSE POLISHERS

PN **Product** **Suggested Retail**
 See chart CA & FG Shank (12 per pack) \$38.10

ABRASIVES

SUPER-SNAP® - FINISHING & POLISHING



- Elastic silicone shank mount for easy placing of the disk onto the mandrel
- Manufactured without a metal center
- Very flexible ultra-thin disks allow easy access to interproximal areas
- Two sizes – four grits (silicon carbide & aluminum oxide)
- Use **Super-Snap for:** composites

PRICING - SUPER-SNAP® - FINISHING & POLISHING

PN	Product	Suggested Retail
0500	Super-Snap Rainbow Technique Kit	\$130.25
180 assorted disks (100 Standard & 80 Mini), 40 Polystrips, 2 Dura-White Stones (CN1 & FL2), 1 CompoSite Fine Midi-Point, 4 CA metal mandrels arranged in a semicircular rainbow organizer tray		
0505	Super-Snap Mini Kit	\$37.10
48 assorted disks (32 Standard & 16 Mini), 1 CompoSite Fine Midi-Point, 1 CA metal mandrel		
0501	Super-Snap Singles	\$194.00
Mini & Standard - Coarse, Medium, Fine & Super Fine grits, Super-Snap plastic CA mandrel (8 disks per Single - 48 Singles per box)		
0502	Super-Snap Singles Standard Size Disk (12mm)	\$155.20
1 disk of Coarse, Medium, Fine & Super Fine grits, Super-Snap plastic CA mandrel (4 disks per Single - 48 Singles per box)		
0503	Super-Snap Singles Mini Size Disk (8mm)	\$155.20
1 disk of Coarse, Medium, Fine & Super Fine grits, Super-Snap plastic CA mandrel (4 disks per Single - 48 Singles per box)		



SUPER-SNAP® DISK REPLACEMENTS

- Black (Coarse): Two sizes – safe side down, disposable silicon carbide disk
- Violet (Medium): Two sizes – safe side down, double sided & safe side up, disposable silicon carbide disk
- Green (Fine): Two sizes – double sided, aluminum oxide grit
- Red (Super Fine): Two sizes – double sided, aluminum oxide grit

SHAPE CATEGORY		Contouring		Finishing				Polishing		Super-Polishing	
SHAPE			MINI 				MINI 		MINI 		MINI 
COLOR/GRIT		Black (Coarse)		Violet (Medium)				Green (Fine)		Red (Super Fine)	
REFILLS (50 PK)	PN	L506	L507	L508	L528	L509	L522	L501	L519	L502	L521

PRICING - SUPER-SNAP® DISK REPLACEMENTS

PN

Product

Suggested Retail

See chart Standard & Mini (50 per box) \$29.00

SUPER-SNAP® POLYSTRIPS

- Use Super-Snap Polystrips for: interproximal areas

SHAPE CATEGORY		Contouring	Finishing	Polishing	Super-Polishing
SHAPE					
COLOR/GRIT		Black/Violet (Course/Medium)		Green/Red (Fine/Super Fine)	
(100 PK)	PN	L525		L526	

PRICING - SUPER-SNAP® POLYSTRIPS

PN

Product

Suggested Retail

See chart Polystrips (100 per box) \$32.25

ABRASIVES

SUPER-SNAP X-TREME™ - FINISHING & POLISHING



Photo courtesy of
Sam Halabo, DDS

- Enhancement to the original green & red Super-Snap disks
- Increased flexibility, tactile feel & durability
- New 3D semispherical grit coating reduces clogging, secondary scratches & produces high gloss results – red disk only
- Two sizes – double-sided, disposable, aluminum oxide disks
- Long lasting with 40% greater gloss than leading competitor
- **Use Super-Snap X-Treme for:** composites

SHAPE CATEGORY		Polishing		Super-Polishing	
SHAPE					
			MINI		MINI
COLOR/GRIT		Green (fine)		Red (super fine)	
REFILLS (50 PK)	PN	L511	L531	L512	L532

PRICING - SUPER-SNAP X-TREME™ - FINISHING & POLISHING

PN	Product	Suggested Retail
0500XS	Super-Snap X-Treme Standard Ultra-Gloss Performance Kit (12mm)	\$188.75
50 black, 50 violet, 50 X-Treme green, 50 X-Treme red, 25 Super-Snap Buff Standard Disks, 4 Super-Snap CA mandrels, 3g DirectDia Polishing Paste Bonus Items: 2.2g Beautifil Flow Plus F00-A2 syringe, 6 compartment storage case		
0500XM	Super-Snap X-Treme Mini Ultra-Gloss Performance Kit (8mm)	\$188.75
50 black, 50 violet, 50 X-Treme green, 50 X-Treme violet, 25 Super-Snap Buff Mini Disks, 4 Super-Snap CA mandrels, 3g DirectDia Polishing Paste Bonus Items: 2.2g Beautifil Flow Plus F03-A2 syringe, 6 compartment storage case		
0506	Super-Snap X-Treme Standard Kit (12mm)	\$53.30
50 green & 50 red disks		
0507	Super-Snap X-Treme Mini Kit (8mm)	\$53.30
50 green & 50 red disks		
See chart.	Standard & Mini (50 per box)	\$28.05



SUPER-SNAP® SUPERBUFF® DISKS




- Paste-impregnated with aluminum oxide
- Disposable
- Eliminates the potential mess associated with polishing paste syringe systems
- **Use Super-Snap SuperBuff for:** composites

PRICING - SUPER-SNAP® SUPERBUFF® DISKS

PN	Product	Suggested Retail
0535	Super-Snap SuperBuff Set	\$32.25
	25 disks & 2 CA mandrels	
0536	Super-Snap SuperBuff Disk (50 refill disks)	\$48.40

SUPER-SNAP® BUFF DISK & MINI BUFF DISK




- Only a very light touch is needed to produce a high lustrous shine on porcelains & composites
- Flexible disks adapt to tooth's natural contours
- Best results are obtained in a combination with polishing pastes
- **Use Super-Snap Buff Disk & Mini Buff Disk for:** composites & porcelains

PRICING - SUPER-SNAP® BUFF DISK & MINI BUFF DISK

PN	Product	Suggested Retail
L523	Standard (25 per box & 1 metal mandrel)	\$29.00
L524	Mini (25 per box & 1 metal mandrel)	\$29.00

ABRASIVES

METAL & PLASTIC MANDRELS



- Durable
- Secure seating of the disk
- HP mandrels are stainless steel
- Super-Snap mandrels are available in HP, CA & FG shanks
- Super-Snap plastic mandrels are single-use in CA shank only
- Use **Super-Snap Metal & Plastic Mandrels for:** Super-Snap, Super-Snap X-Treme & SuperBuff products

PRICING - METAL & PLASTIC MANDRELS

PN	Product	Suggested Retail
5011	HP Mandrel (12 per pack)	\$29.95
0437	Super-Snap Metal - HP Shank (6 per box)	\$24.15
0438	Super-Snap Metal - FG Shank (6 per box)	\$24.15
0439	Super-Snap Metal - CA Shank (6 per box)	\$24.15
0440	Super-Snap Plastic - CA Shank (100 per box)	\$38.70



ONEGLOSS® & ONEGLOSS® PS



- OneGloss aluminum oxide polishers are disposable
- OneGloss complete kit with stainless steel CA mandrels
- OneGloss PS aluminum oxide polishers are mounted on sturdy plastic mandrels
- OneGloss PS are disposable & individually wrapped
- Available in three different shapes
- Durable enough to withstand polishing an entire arch
- **Use OneGloss & OneGloss PS for:** composites, temporary crowns, removal of orthodontic cements & removal of stains from enamel

SHAPE CATEGORY		Cup	Midi-Pt.	IC	Mini-Pt.	Assort. Kit	CA Mandrel
SHAPE						See Descriptions Below	
ISO#		060	060	100	030		
ONEGLOSS							
CA (50 PK)	PN	0181	0182	0183	0186 (12 pk)	0180	0435 (6 pk)
ONEGLOSS PS							
CA (50 PK)	PN	0176	0177	0178	-	0175	-

PRICING - ONEGLOSS® & ONEGLOSS® PS

PN

Product

Suggested Retail

0180 **OneGloss Set** **\$98.20**

20 of each Cup, Midi-Point & IC shapes, plus 3 CA mandrels

See chart **OneGloss Refill (50 per box)** **\$80.50**

Cup, Midi-Point or IC

0186 **OneGloss Mini-Point Refill (12 per pack)** **\$34.10**

0435 **OneGloss CA Mandrel (6 per box)** **\$21.30**

0175 **OneGloss PS Assortment Kit** **\$121.00**

15 Cups, 20 Midi-Points, 15 IC

See chart **OneGloss PS (50 per box)** **\$121.00**

Cup, Midi-Point or IC

0185 **OneGloss/SuperBuff Combo Set** **\$119.55**

20 each of OneGloss Cup, Midi-Point & IC shapes, 25 SuperBuffs disks, 5 CA mandrels

ABRASIVES

ACRYPOINT LAB POLISHERS

SHAPE CATEGORY	BP1	PC2
SHAPE		
ISO#	100	050
DARK GREY (course)		
HP (6 PK)	PN 0426	-
HP (12 PK)	PN -	0423
BROWN (medium)		
HP (6 PK)	PN 0427	-
HP (12 PK)	PN -	0424
LIGHT GREY (fine)		
HP (6 PK)	PN 0428	-
HP (12 PK)	PN -	0425


- Available in 3 silicon carbide grits
- Use Dark Grey (coarse) for adjustments, Brown (medium) for finishing & Light Grey (fine) for polishing
- Use AcryPoint for: acrylic resin denture material

PRICING - ACRYPOINT LAB POLISHERS

PN	Product	Suggested Retail
See chart	HP Shank (6 per pack)	\$44.45
BP1		
See chart	HP Shank (12 per pack)	\$44.45
PC2		



LAB SERIES STONES



- Aluminum oxide grit
- Free from carbon contaminants
- **Use Lab Series White & Pink for:** precious alloys
- **Use Lab Series Coral for:** non-precious alloys & model castings

SHAPE CATEGORY	No. 1	No. 2	No. 3	No. 4	No. 5	No. 6	Unm. Wheel
SHAPE							
ISO#	040	070	056	030	025	130	190
WHITE							
HP (12 PK)	PN	461A	462A	463A	464A	465A	466A
HP (72 PK)	PN	0461	0462	0463	0464	0465	0466
PINK							
HP (12 PK)	PN	471A	472A	473A	474A	475A	476A
HP (72 PK)	PN	0471	0472	0473	0474	0475	0476
CORAL							
HP (12 PK)	PN	481A	482A	483A	484A	485A	486A
HP (72 PK)	PN	0481	0482	0483	0484	0485	0486

PRICING - LAB SERIES STONES

PN	Product	Suggested Retail
See chart.....	White & Coral Unmounted Wheels (12 per pack)...	\$23.90
See chart.....	White, Pink, Coral Mounted (12 per pack).....	\$23.90
See chart.....	White, Pink, Coral Mounted (1/2 Gross)	\$92.70

ABRASIVES

HARDIE® POLISHERS

- Use **Hardie Polishers for:** hard & non-precious alloys

SHAPE CATEGORY	FL2	PC2	WH6	Floppie	BP1	Unm. Wheel #8
SHAPE						
ISO#	030	050	120	220	100	230
BLACK						
HP (12 pk)	PN	-	0362	-	0366	0368 (6 pk)
CA (12 pk)	PN	0390	-	-	-	-
BROWN						
HP (12 pk)	PN	0361	0363	0365	0367	0369 (6 pk)
CA (12 pk)	PN	0391	-	-	-	-

PRICING - HARDIE® POLISHERS

PN	Product	Suggested Retail
See chart.	HP & CA Shank (12 per pack)	\$34.75
	FL2, PC2, WH6, Floppie	
See chart.	HP Shank (6 per pack)	\$34.75
	BP1	
See chart.	Hardie Black Unmounted Wheel (12 per pack)	\$21.85



UNMOUNTED FINISHING WHEELS

- Shape & contour
- Use **Dura-Green Wheels (silicon carbide) & Chipless Wheels (silicon carbide)** for: porcelains & precious alloys
- Use **Cut-Off Disks (resin-bonded)**: non-precious alloys & model castings

SHAPE CATEGORY	No. 7	No. 7a	No. 8	Cut-Off Disk
SHAPE				
ISO#	190	190	190	250
THICKNESS	1.7mm	2.4mm	3mm	-
DURA-GREEN				
(12 PK)	PN 0150	PN 0151	PN 0152	-
CHIPLESS				
(12 PK)	PN 0455	PN 0456	PN 0457	-
UNM. WHEELS				
(12 PK)	PN -	PN -	PN -	0441 (100 pk)

PRICING - UNMOUNTED FINISHING WHEELS

PN	Product	Suggested Retail
See chart.....	Dura-Green Unmounted Wheel (12 per pack)	\$25.50
See chart.....	Chipless Unmounted Wheel (12 per pack)	\$25.80
See chart.....	Cut-Off Disk (100 per pack)	\$31.55

ABRASIVES

DURA-POLISH & DURA-POLISH DIA POLISHING PASTE



- Dura-Polish is impregnated with 73% by weight aluminum oxide particles
- **Use Dura-Polish for:** non-precious alloys, precious alloys, composites, hard resins & denture based resins
- Dura-Polish Dia is impregnated with 67% by weight diamond particles
- **Use Dura-Polish Dia for:** porcelains, composites & hard resins

PRICING - DURA-POLISH & DURA-POLISH DIA POLISHING PASTE

PN	Product	Suggested Retail
0552	Dura-Polish (20g)	\$88.60
0554	Dura-Polish Dia (5g)	\$88.60

COMPOSITE & DIRECTDIA POLISHING PASTE



- **CompoSite Polishing Paste** is an intraoral polishing cream with submicron aluminum oxide particles
- **Use CompoSite Polishing Paste for:** composites
- **DirectDia Polishing Paste** is a diamond polishing paste with 20% super fine diamond particles in a smooth, lime flavored intraoral paste
- DirectDia has optimal viscosity & excellent handling
- **Use DirectDia Polishing Paste for:** enamel or porcelains

PRICING - COMPOSITE & DIRECTDIA POLISHING PASTE

PN	Product	Suggested Retail
0551	CompoSite Polishing Paste (5g)	\$41.40
0558	DirectDia Polishing Paste (3g)	\$27.35



STERILIZABLE BUR BLOCKS



- Autoclavable for operator
- Four unique variations
- Hard-anodized aluminum
- Withstand most sterilizing techniques
- Holds 18 burs
- Upright for easy access
- Protective snap-lock lid & open sides for proper autoclaving
- **Use Sterilizable Bur Blocks for:** polishers, diamonds & carbides

PRICING - STERILIZABLE BUR BLOCKS

PN	Product	Suggested Retail
0960	Contra-Angle	\$46.00
0961	Friction Grip	\$46.00
0962	Handpiece	\$46.00

ABRASIVES

NEW! BURBUTLER™ UNIVERSAL SILICONE BUR BLOCK



- NEW SIZES!
- Ideal for all shanks including FG, CA & short shank – mix & match burs
- No spilling
- Will not magnetize or rust
- Long-lasting and durable, no moving parts or replacement plugs
- Easy to insert & remove burs
- Sterilizable
- **Fully autoclavable** – fast cooling
- Removable lid can be used as base during procedures to discard used rotary instruments
- Simple one hand lid removal
- Secure clear lid for easy viewing & safety
- 5, 10 & 25-hole blocks available in 5 colors
- 25-hole base has the option of a tall lid to include HP burs
- **Use BurButler for:** all burs, implant components, perio irrigation tips & other small items

PRICING - BURBUTLER™ 5-HOLE UNIVERSAL SILICONE BUR BLOCK

PN	Product	Suggested Retail
8205AOC	Amber Orange	\$34.25
8205APC	Amethyst Purple	\$34.25
8205DWC	Diamond White	\$34.25
8205RPC	Ruby Pink	\$34.25
8205SBC	Sapphire Blue	\$34.25
8205LID	Replacement Lid	\$11.40

PRICING - BURBUTLER™ 10-HOLE UNIVERSAL SILICONE BUR BLOCK

PN	Product	Suggested Retail
8210AOC	Amber Orange	\$39.45
8210APC	Amethyst Purple	\$39.45
8210DWC	Diamond White	\$39.45
8210RPC	Ruby Pink	\$39.45
8210SBC	Sapphire Blue	\$39.45
8210LID	Replacement Lid	\$13.45



PRICING - BURBUTLER™ 25-HOLE UNIVERSAL SILICONE BUR BLOCK

PN	Product	Suggested Retail
8225AOC	Amber Orange	\$55.05
8225APC	Amethyst Purple	\$55.05
8225DWC	Diamond White	\$55.05
8225RPC	Ruby Pink	\$55.05
8225SBC	Sapphire Blue	\$55.05
8225LID	Replacement Lid	\$27.00

PRICING - BURBUTLER™ 25-HOLE TALL UNIVERSAL SILICONE BUR BLOCK

PN	Product	Suggested Retail
8225TAOC	Amber Orange	\$68.60
8225TAPC	Amethyst Purple	\$68.60
8225TDWC	Diamond White	\$68.60
8225TRPC	Ruby Pink	\$68.60
8225TSBC	Sapphire Blue	\$68.60
8225TLID	Replacement Lid	\$30.10

ABRASIVES



SHOFU RAVES AND RECOMMENDATIONS

Super-Snap X-Treme Disks

"Most dentists are concerned about saving time in placing and curing composites, but much of the time-wasting involves the finishing and polishing.

To help facilitate the process, Super-Snap X-Treme has incorporated an innovative technology into its original green and red Super-Snap polishing disks.

Their 3-D semi-spherical surface coating prevents secondary surface scratches and clogging, as well as the heat that denatures the material, so dentists using the system need only the two kinds of heads on the mandrel to create the needed finishing and polishing." ~ *Dr. Fay Goldstep*

"Super Snap X-Treme is the only bur system for provisionals that I simply can't do without. Regardless of single or multiple units it's reliable, fast, precise and efficient!" ~ *Shannon Pace Brinker*

"The Shofu Super-Snap X-Treme Finishing and Polishing Kit is an excellent example of an anterior composite polishing system. This kit can serve as your one-stop anterior composite polishing solution. It is simple and intuitive, and it yields an excellent result. Great work, Shofu. Hard-hit double in the gap!" ~ *Dr. Joshua Austin*

ABRASIVES

FINISHING & POLISHING KITS

NEW! ATTACHMENT REMOVAL KIT	30
--	----

BURBUTLER™ BUR BLOCK KITS

NEW! ARTISTIC COMPOSITE KIT	30
NEW! PINK PREP KIT	31
NEW! COMPOSITE FINISHING KIT	32
NEW! COMPOSITE POLISHING KIT	32
NEW! CERAMISTÉ POLISHING KIT	33
NEW! METALS POLISHING KIT	33
NEW! ALL CERAMIC PREPARATION KIT	34
NEW! CADMASTER CAD/CAM KIT	34

STERILIZABLE BUR BLOCK KITS

PORCELAIN & ENAMEL ADJUSTMENT KIT	35
COMPOSITE TECHNIQUE KIT	35

CLASSIC PLASTIC KITS

PORCELAIN & ENAMEL ADJUSTMENT KIT	35
COMPOSITE TECHNIQUE KIT	35
COMPOSITE POLISHING KIT	36
COMPOSITE FINISHING KIT	36
ACRYLIC ADJUSTING & POLISHING KIT	37
AMALGAM POLISHING KIT	38
GOLD POLISHING KIT	38
PORCELAIN LAMINATE POLISHING KIT	39
PORCELAIN VENEER KIT	39
NP ALLOY ADJUSTMENT KIT	40

ABRASIVES QUICK GUIDES

CHAIRSIDE	41
LABORATORY	42

ABRASIVES

NEW! ATTACHMENT REMOVAL KIT



- A complete system for safe & easy removal of orthodontic attachment
- Inclusive of Shofu's award-winning finishers & polishers
- Selected instruments & the removal technique support minimally invasive dentistry
- BurButler is **fully autoclavable**
- **Use Attachment Removal Kit for:** removal of composite buttons

PRICING - ATTACHMENT REMOVAL KIT

PN	Product	Suggested Retail
ARK.	Attachment Removal Kit.	\$275.60
40 Super-Snap Singles Mini Size Disks (8mm), 3 Robot Carbide Finishers CTF, OneGloss PS (ten Cups, ten Midi-Points, ten IC), 50 Super-Snap SuperBuff Disks, DirectDia Polishing Paste, 10-hole BurButler with Lid Bonus Items: 2 (0.25g) Beautifil II A2 tips, 2 (0.21g) Beautifil Flow Plus F00-A2 tips, 2 (0.21g) Beautifil Flow Plus F03-A2 tips		

NEW! ARTISTIC COMPOSITE KIT



- Comprehensive finishing system inclusive of finest diamonds, carbides & stones
- Supported by an efficient artisan's finishing technique
- Excellent for anterior composite restorations
- Promotes minimally invasive dentistry
- BurButler is **fully autoclavable**
- **Use Artistic Composite Kit for:** finishing of composite restorations

PRICING - ARTISTIC COMPOSITE KIT

PN	Product	Suggested Retail
0980BB.	Artistic Composite Kit.	\$171.60
Robot FG Diamonds (840F, 840V, 835F, 793V), T&F Hybrid Point (0933), Robot Carbide FG Finisher (0718), Dura-Green Stone (FL2), Robot Carbide CTF-FG Finishers (0770, 0771, 0772), 10-Hole Sapphire Blue BurButler with Lid		



NEW! PINK PREP KIT



- Two grits: Regular & Super Fine
- Three shapes: inverted cone, round-tapered cylinder & flame
- Two Minimally Invasive Preparation Robot Diamond points
- Allows for minimally-invasive preparation of direct composite restorations
- BurButler is **fully autoclavable**
- Complements the Beautifil II PINK Kit (PN GMK1)
- **Use Pink Prep kit for:** composite restorations

PRICING - PINK PREP KIT

PN	Product	Suggested Retail
8205K1	PINK Prep Kit	\$69.95
Robot FG Diamonds (0948, 0949, 0855, 0856, 840V), 5-Hole Ruby Pink BurButler with Lid		

ABRASIVES

NEW! COMPOSITE FINISHING KIT



- NEW & IMPROVED!
- More shapes for a complete finishing system – 10 burs in each kit
- Dura-Green silicon carbide stones for adjusting
- Dura-White aluminum oxide stones for fine finishing
- Contra-Angle & Friction Grip
- Now in BurButler – **fully autoclavable**
- Use **Composite Finishing kit for:** composites

PRICING - COMPOSITE FINISHING KIT

PN	Product	Suggested Retail
0306BB	Composite Finishing CA Kit	\$91.75
	Dura-Green Stones (CN1, FL2, FL2, RD1, RD2), Dura-White Stones (CN1, FL1, FL2, RD1, RD2), 10-Hole BurButler with Lid	
0302BB	Composite Finishing FG Kit	\$91.75
	Dura-Green Stones (CN1, two FL2, RD1, RD2), Dura-White Stones (CN1, FL1, FL2, RD1, RD2), 10-Hole BurButler with Lid	

NEW! COMPOSITE POLISHING KIT



- NEW & IMPROVED!
- 10 shapes for a complete finishing & polishing system
- Dura-Green silicon carbide stones for adjusting
- Dura-White aluminum oxide stones for finishing
- CompoSite silicone points are impregnated with aluminum oxide for polishing
- Contra-Angle
- Now in BurButler – **fully autoclavable**
- Use **CompoSite Polishing kit for:** composites

PRICING - COMPOSITE POLISHING KIT

PN	Product	Suggested Retail
0310BB	CompoSite Polishing CA Kit	\$91.75
	Dura-Green Stones (CN1 & FL2), Dura-White Stones (CN1 & FL2), CompoSite Points (two Bullets, two Knives, two Cups), 10-Hole BurButler with Lid	



NEW! CERAMISTÉ POLISHING KIT



- 10 shapes for a complete polishing system
- Ceramisté Standard for pre-polishing
- Ceramisté Ultra for polishing
- Ceramisté Ultra II for super-polishing
- Contra-Angle & Friction Grip
- Available in BurButler – **fully autoclavable**
- **Use Ceramisté Polishing kit for:** porcelains & enamel

PRICING - CERAMISTÉ POLISHING KIT

PN	Product	Suggested Retail
0317BB	Ceramisté Polishing CA Kit	\$91.75
	Ceramisté Standard (two Cups & two Midi-Points), Ultra (two Cups & two Midi-Points), Ultra II (Cup & Midi-Point), 10-Hole BurButler with Lid	
0318BB	Ceramisté Polishing FG Kit	\$91.75
	Ceramisté Standard (two Cups & two Midi-Points), Ultra (two Cups & two Midi-Points), Ultra II (Cup & Midi-Point), 10-Hole BurButler with Lid	

NEW! METALS POLISHING KIT



- 10 shapes for a complete polishing system
- Brownie points for pre-polishing
- Greenie points for polishing
- Supergreenie points for super-polishing
- Contra-Angle
- Available in BurButler – **fully autoclavable**
- **Use Metals Polishing kit for:** amalgams, gold non-precious & precious alloys

PRICING - METALS POLISHING KIT

PN	Product	Suggested Retail
0304BB	Metals Polishing CA Kit	\$91.75
	Brownies (Mini-Point, Cup, WH6, PC2), Greenies (Mini-Point, Cup, WH6, PC2), Supergreenies (Mini-Point & Cup), 10-Hole BurButler with Lid	

ABRASIVES

NEW! ALL CERAMIC PREPARATION KIT



- Exceptional Robot Diamonds
- Two grits of Robot Diamond Points, Regular & Super Fine
- Use for anterior, premolar, molar preparations, inlay & laminate veneer
- Available in BurButler – **fully autoclavable**
- **Use All Ceramic Preparation Kit for:** ceramics

PRICING - ALL CERAMIC PREPARATION KIT

PN	Product	Suggested Retail
0955BB	ALL CERAMIC PREPARATION KIT	\$200.00
17 FG Robot Diamond Points PN: 0825, 0826, 0833, 0836, 0839, 0893, 0896, 0897, 0899, 808V, 825V, 826V, 833V, 836V, 839V, 893V, 899V, 25-Hole BurButler with Lid		

NEW! CADMASTER CAD/CAM Kit



- Dura-Green silicon carbide stones for contouring & finishing
- Dura-White aluminum oxide stones for finishing
- Dura-Green DIA vitrified bonded diamond abrasives for contouring & finishing
- CeraMaster Coarse to polish to a satin finish
- CeraMaster to super-polish
- ZilMaster for contouring, finishing, pre-polishing & polishing
- Available in BurButler – **fully autoclavable**
- **Use CADMaster CAD/CAM Kit for:** composites, porcelains, precious and non-precious alloys, amalgams, enamel, glass ionomer restorations, zirconia, ceramics & lithium disilicate

PRICING - CADMASTER CAD/CAM Kit

PN	Product	Suggested Retail
0370BB	CADMASTER CAD/CAM Kit	\$349.00
Dura-Green Stones (three FL2, WH6, IC9), Dura-White Stones (three FL2 & RD2), Dura-Green DIA Stones (TC4 & KN7), CeraMaster Points (two Bullets, two Mini-Points, Cup, WH6), CeraMaster Coarse Points (two Bullets, Mini-Point, Cup, WH6), ZilMaster Points (three WH6 - Fine, Medium, Coarse), 25-Hole BurButler with Tall Lid		



PORCELAIN (HP) & ENAMEL (CA) ADJUSTMENT KIT



- Dura-White micrograined aluminum oxide stones
- Ceramisté silicone polishers
- Handpiece & Contra-Angle kits
- Your choice of a sterilizable bur block or plastic kit
- **Use Porcelain & Enamel kit for:** porcelains & enamel

PRICING - PORCELAIN (HP) & ENAMEL (CA) ADJUSTMENT KIT

PN	Product	Suggested Retail
0301B	Sterilizable Bur Block HP Kit	\$95.40
0301	Classic Plastic HP Kit	\$52.45
0307	Classic Plastic CA Kit	\$52.45
Dura-White Stones (three), Ceramisté Points (nine)		

COMPOSITE TECHNIQUE KIT



- Dura-White Stones for finishing
- CompoSite silicone points are impregnated with aluminum oxide for polishing
- CompoSite Fine silicone points impregnated with zirconia for fine polishing
- Contra-Angle kit
- Your choice of a sterilizable bur block or plastic kit
- **Use CompoSite Technique kit for:** composites

PRICING - COMPOSITE TECHNIQUE KIT

PN	Product	Suggested Retail
0322B	Sterilizable Bur Block CA Kit	\$95.40
Dura-White Stones (two), CompoSite Points (six), CompoSite Fine Points (six)		
0322	Classic Plastic CA Kit	\$71.55
Dura-White Stones (two), CompoSite Points (eight), CompoSite Fine Points (eight)		

ABRASIVES

COMPOSITE POLISHING KIT



- Dura-Green silicon carbide stones for adjusting
- Dura-White aluminum oxide stones for finishing
- CompoSite silicone points are impregnated with aluminum oxide for polishing
- CompoSite Fine silicone points impregnated with zirconia for fine polishing
- Contra-Angle kit
- Your choice of a sterilizable bur block or plastic kit
- Use **CompoSite Polishing kit for:** composites

PRICING - COMPOSITE POLISHING KIT

PN	Product	Suggested Retail
0310	Classic Plastic CA Kit	\$52.45
Dura-Green Stones (CN1 & FL2), Dura-White Stones (CN1 & FL2), CompoSite Points (six), CompoSite Fine Points (two)		

COMPOSITE FINISHING KIT



- Dura-Green silicon carbide stones for adjusting
- Dura-White aluminum oxide stones for fine finishing
- Contra-Angle & Friction Grip kits
- Your choice of a sterilizable bur block or plastic kit
- Use **CompoSite Finishing kit for:** composites

PRICING - COMPOSITE FINISHING KIT

PN	Product	Suggested Retail
0306	Classic Plastic CA Kit	\$52.45
0302	Classic Plastic FG Kit	\$52.45
Dura-Green Stones (three CN1 & three FL2), Dura-White Stones (three CN1 & three FL2)		



ACRYLIC ADJUSTING & POLISHING KIT



- A complete chairside or lab system for adjusting & polishing acrylic materials
- Handpiece kit
- Carbide burs for smooth, fast adjustment
- AcryPoints are impregnated with high quality abrasives
- Your choice of a sterilizable bur block or plastic kit
- **Use Acrylic Adjusting & Polishing kit for:** acrylic resin denture material

PRICING - ACRYLIC ADJUSTING & POLISHING KIT

PN	Product	Suggested Retail
0321	Classic Plastic HP Kit	\$125.65
Cross-cut Carbide HP Cutters (two), assorted shapes of AcryPoint (ten)		

ABRASIVES

AMALGAM POLISHING KIT



- Provides fast, brilliant results without pumice or tin oxide
- Cups are for proximals & Mini-Points are for occlusals
- Brownies for pre-polishing
- Greenies for polishing
- Supergreenies for super-polishing
- Contra-Angle & Friction Grip kits
- Your choice of a sterilizable bur block or plastic kit
- **Use Amalgam Polishing kit for:** amalgams & non-precious alloys

PRICING - AMALGAM POLISHING KIT

PN	Product	Suggested Retail
0308	Classic Plastic CA Kit	\$52.45
0309	Classic Plastic FG Kit	\$52.45
Brownie (three Cups & three Mini-Points), Greenie (two Cups & two Mini-Points), Supergreenie (Cup & Mini-Point)		

GOLD POLISHING KIT



- Brownies for pre-polishing
- Greenies for polishing
- Supergreenies for super-polishing
- Handpiece & Contra-Angle kits
- Your choice of a sterilizable bur block or plastic kit
- **Use Gold Polishing kit for:** precious alloys

PRICING - GOLD POLISHING KIT

PN	Product	Suggested Retail
0303	Classic Plastic HP Kit	\$52.45
Brownies (two KN7, two PC2, two WH6), Greenies (KN7, PC2, WH6), Supergreenies (KN7, PC2, WH6)		



PORCELAIN LAMINATE POLISHING KIT



- Ceramisté Standard for pre-polishing
- Ceramisté Ultra for polishing
- Ceramisté Ultra II for super-polishing
- Contra-Angle & Friction Grip kits
- Your choice of a sterilizable bur block or plastic kit
- **Use Porcelain Laminate Polishing kit for:** porcelains & enamel

PRICING - PORCELAIN LAMINATE POLISHING KIT

PN	Product	Suggested Retail
0317	Classic Plastic CA Kit	\$56.50
0318	Classic Plastic FG Kit	\$56.50
Ceramisté Standard (four Cups & four Midi-Points), Ultra (two Cups & two Midi-Points), Ultra II (Cup & Midi-Point)		

PORCELAIN VENEER KIT



- Indicated for contouring, finishing & polishing of porcelain veneers
- Incorporates FG Robot- and Hybrid-Points & CA Dura-White and Ceramisté Points
- Use Robot- & Hybrid-Points to contour
- Use Dura-White to finish & pre-polish
- Use Ceramisté Points to polish & super-polish
- **Use Porcelain Veneer kit for:** porcelains

PRICING - PORCELAIN VENEER KIT

PN	Product	Suggested Retail
0300	Classic Plastic CA Kit	\$117.70
Regular-grit FG Diamond Point, Fine-grit FG Diamond Points (two), T&F Hybrid Points (three), CA Dura-White Stones (FL2 & RE1), CA Ceramisté Standard, Ultra & Ultra II (three Knives, three Wheels, two Cups)		

ABRASIVES

NP ALLOY ADJUSTMENT KIT



- Lab Series White Stones for adjusting
- Hardie Points for polishing
- Handpiece kit
- Available in a plastic kit
- Use NP Alloy Adjustment kit for: non-precious alloys

PRICING - NP ALLOY ADJUSTMENT KIT		
PN	Product	Suggested Retail
0311	Classic Plastic HP Kit	\$52.45
Lab Series White Stones (ten), Hardie Points (two PC2)		



ABRASIVE QUICK GUIDE - CHAIRSIDE

OPERATORY	CONTOURING/FINISHING	POLISHING	HIGH-LUSTER POLISHING
ENAMEL			
	PINK Prep Kit	OneGloss Dura-White CeraMaster Coarse Ceramisté Standard, Ultra	DirectDia Polishing Paste with Super-Snap Buff Disk Super-Snap SuperBuff Disk CeraMaster Ceramisté Ultra II
FILLINGS			
COMPOSITES	Dura-Green/Dura-Green Fine Dura-White OneGloss Super-Snap black/violet Robot Carbide 8/16 blades T&F Hybrid Points CompoSite Artistic Composite Kit	OneGloss Super-Snap green Super-Snap X-Treme green Robot Carbide CTF 30 blades CompoMaster Coarse CompoSite Fine	Super-Snap red Super-Snap X-Treme red CompoSite Polishing Paste with Super-Snap Buff Disk CompoMaster SuperBuff
AMALGAMS	Dura-Green	Brownie/Greenie	Supergreenie
GLASS-IONOMERS	Dura-White T&F Hybrid Points Robot Carbides		Ceramisté Ultra II
INLAYS, CROWNS, BRIDGES, FRAMEWORKS			
TITANIUM	Dura-Green Hardie Black	Brownie/Greenie Hardie Brown	Supergreenie
PRECIOUS ALLOYS	Dura-Green	Brownie/Greenie	Supergreenie
NON-PRECIOUS ALLOYS	Dura-Green Hardie Black	Brownie/Greenie Hardie Brown	Supergreenie
VENEERS, FACINGS, INLAYS, CROWNS			
CERAMICS (Zirconia, Lithium Disilicate, **IPS e.max®, **Bruxzir®, Alumina)	Dura-Green DIA Dura-Green Dura-White *CoreMaster Coarse	*CoreMaster CeraMaster Coarse	CeraMaster Ceramisté Ultra II DirectDia Polishing Paste
PORCELAINS	Dura-Green DIA Dura-Green/Dura-Green Fine Dura-White T&F Hybrid Points	CeraMaster Coarse Ceramisté Standard, Ultra Soft PA & PB	CeraMaster Ceramisté Ultra II DirectDia Polishing Paste
COMPOSITES	Dura-Green/Dura-White T&F Hybrid Points OneGloss	CompoMaster Coarse CompoSite/CompoSite Fine Ceramisté Soft PA & PB OneGloss	CompoMaster CompoSite Polishing Paste with Super-Snap Buff Disk CompoSite Fine
LAB COMPOSITES (**Lava™ Ultimate, Ceramage)	Dura-Green DIA Dura-Green	Ceramisté Soft PB CompoMaster Coarse *Dura-Polish	CompoMaster *Dura-Polish DIA
PARTIAL/FULL DENTURES			
DENTURE ACRYLICS	*Robot Carbide HP Cutters *AcryPoint, dark grey and brown	*AcryPoint, light grey	
MOUTHGUARDS/SPLINTS			
ACRYLICS (**Under Armour®)	Dura-Green DIA Dura-Green *AcryPoint, dark grey CoreMaster Coarse	*AcryPoint, brown *AcryPoint, light grey CoreMaster	

■ one-step procedure

■ one-step procedure and disposable

* Not for intraoral

**All trademarks belong to their manufacturers

■ two-step procedure

■ three-step procedure

■ multi-step procedure

ABRASIVES

ABRASIVE QUICK GUIDE - LABORATORY

LABORATORY	CONTOURING/FINISHING	POLISHING	HIGH-LUSTER POLISHING
VENEERS, FACINGS, INLAYS, CROWNS			
CERAMICS (Zirconia, Lithium Disilicate, **IPS e.max®, **Bruxzir®, Alumina)	Dura-Green DIA CoreMaster Coarse	CoreMaster Dura-Polish	Dura-Polish DIA
PORCELAINS	Dura-Green DIA Dura-Green/Dura-Green Fine Dura-White Robot Carbide CTF T&F Hybrid Points	CeraMaster Coarse Ceramisté Standard, Ultra Soft PA & PB Dura-Polish	CeraMaster Ceramisté Ultra II DirectDia Polishing Paste Dura-Polish DIA
LAB COMPOSITES (**Lava™ Ultimate, Ceramarge)	Dura-Green Dura-White OneGloss Robot Carbide CTF T&F Hybrid Points Robot Carbide HP Cutters	CompoMaster Coarse Ceramisté Soft PA & PB OneGloss Dura-Polish	CompoMaster CompoSite Paste Dura-Polish DIA
INLAYS, CROWNS, BRIDGES, FRAMEWORKS			
PRECIOUS ALLOYS	Lab Pink Lab White	Brownie/Greenie	Supergreenie Dura-Polish
NON-PRECIOUS ALLOYS	Robot Carbide Lab Coral Hardie Black	Hardie Brown	Dura-Polish
PARTIAL DENTURES			
MODEL CASTING	Robot Carbide Lab Coral Hardie Black	Hardie Brown	
PARTIAL/FULL DENTURES			
DENTURE ACRYLICS	Robot Carbide HP Cutters AcryPoint, dark grey and brown	AcryPoint, light grey	Dura-Polish
ORTHO APPLIANCES/SPLINTS/MOUTHGUARDS			
ACRYLICS (**Under Armour®)	AcryPoint, dark grey CoreMaster Coarse Robot Carbide HP Cutters	AcryPoint, brown AcryPoint, light grey CoreMaster	

■ one-step procedure

■ one-step procedure and disposable

** All trademarks belong to their manufacturers

■ two-step procedure

■ three-step procedure

DIAMONDS/CARBIDES

DIAMONDS

SHOFU ROBOT FG DIAMONDS	44 - 47
MINIMALLY INVASIVE DIAMOND TECHNIQUE	47
TPE DIAMOND POINTS, HYBRID POINTS & T&F HYBRID POINTS	47 - 49
GINGI CURETTAGE	49
ABRASIVE DRESSER	49

CARBIDES

SHOFU ROBOT CARBIDE FG FINISHERS	50
SHOFU ROBOT CARBIDE CTF-FG FINISHERS	51
SHOFU ROBOT CARBIDE HP CUTTERS	52

EQUIPMENT

ARTICULATORS	53
LAB AIR-Z	54

Official Partner



Minimally Invasive
Cosmetic Dentistry

DIAMONDS/CARBIDES

ROBOT® FG DIAMONDS



Grit:
 R=Regular
 C=Coarse
 SC=Super Coarse
 F=Fine
 SF=Super Fine

- Rated #1 in overall performance by a national research organization
- Economically priced
- Made by high technology robots in a computer-automated environment untouched by human hands
- Multi-coating technique with a special surface treatment for durability – maintains cutting efficiency with less clogging
- Smooth cutting & long lasting performance
- Grits diamonds are available in: Super Fine (SF), Fine (F), Standard (R), Coarse (C) & Super Coarse (SC)
- Use **Robot FG Diamonds for:** cavity & crown preparation

SHAPE CATEGORY	Straight Cylinder						Round End Straight Cylinder			
SHAPE										
ISO#	107/010	109/013	110/015	110/016	110/018	158/013	289/016	289/015	130/016	130/014
LENGTH	4.0	4.5	7.0	7.0	7.0	8.0	8.0	8.0	8.0	8.0
GRIT	R	R	R	R	R	SF	SC	C	R	F
PN (1 PK)	0800-1	0801-1	0803-1	0805-1	0807-1	808V-1	845X-1	845C-1	0846-1	846F-1

SHAPE CATEGORY	Flat End Tapered Cylinder									
SHAPE										
ISO#	170/012	170/017	170/017	170/022	171/016	172/019	172/021	173/017	173/015	174/020
LENGTH	4.0	4.4	5.0	5.0	7.0	9.0	9.0	10.0	10.0	12.0
GRIT	R	R	R	R	R	C	R	R	F	R
PN (1 PK)	0810-1	0812-1	0813-1	0814-1	0816-1	817C-1	0818-1	0820-1	820F-1	0822-1

SHAPE CATEGORY	Flat End Tapered Cylinder	Corner Round Tapered Cylinder								
SHAPE										
ISO#	173/022	196/013	544/018	544/018	544/021	544/021	545/018	545/018	553/019	553/019
LENGTH	10.0	5.0	4.0	4.0	5.0	5.0	7.0	7.0	10.0	10.0
GRIT	R	R	R	SF	R	SF	R	SF	R	SF
PN (1 PK)	0824-1	0830-1	0789-1	789V-1	0790-1	790V-1	0791-1	791V-1	0792-1	792V-1



SHAPE CATEGORY	Corner Round Tapered Cylinder						Needle Tapered Cylinder			
SHAPE										
ISO#	553/023	553/023	544/027	544/027	544/031	544/031	164/014	165/013	165/012	166/016
LENGTH	10.0	10.0	3.5	3.5	3.5	3.5	6.0	8.0	8.0	10.0
GRIT	R	SF	R	SF	R	SF	F	R	F	R
PN (1 PK)	0793-1	793V-1	0825-1	825V-1	0826-1	826V-1	847F-1	0848-1	848F-1	0850-1

SHAPE CATEGORY	Needle Tapered Cylinder			Round End Tapered Cylinder						
SHAPE										
ISO#	166/018	166/015	166/022	198/013	222/016	222/018	222/015	223/017	198/018	193/016
LENGTH	10.0	10.0	10.0	9.0	7.0	7.0	7.0	9.0	9.0	10.0
GRIT	SC	F	R	R	R	SC	F	R	SF	R
PN (1 PK)	850X-1	850F-1	0852-1	0831-1	0832-1	832X-1	832F-1	0833-1	833V-1	0834-1

SHAPE CATEGORY	Round End Tapered Cylinder									
SHAPE										
ISO#	223/018	223/019	223/018	223/016	198/021	198/021	223/023	223/023	198/023	198/023
LENGTH	9.0	9.0	9.0	9.0	9.0	9.0	9.0	11.5	9.0	9.0
GRIT	R	SC	C	F	R	SF	R	R	R	SF
PN (1 PK)	0835-1	835X-1	835C-1	835F-1	0836-1	836V-1	0837-1	0838-1	0839-1	839V-1

SHAPE CATEGORY	Flame									
SHAPE										
ISO#	249/014	249/016	249/015	249/013	249/012	243/009	249/017	249/016	243/012	249/023
LENGTH	7.5	7.5	7.5	7.5	7.5	4.0	7.5	7.5	6.0	9.0
GRIT	R	SC	C	F	SF	SF	R	F	SF	R
PN (1 PK)	0840-1	840X-1	840C-1	840F-1	840V-1	841V-1	0842-1	842F-1	843V-1	0844-1

DIAMONDS/CARBIDES

ROBOT® FG DIAMONDS

SHAPE CATEGORY	S Type Inverted Cone						Corn			
SHAPE										
ISO#	032/013	032/018	032/021	019/010	019/015	019/017	039/033	039/034	031/035	031/035
LENGTH	1.5	2.0	2.3	2.2	3.6	3.5	5.0	5.0	5.7	5.7
GRIT	R	R	R	R	R	R	R	SC	R	SF
PN (1 PK)	0855-1	0856-1	0858-1	0859-1	0860-1	0861-1	0898-1	898X-1	0899-1	899V-1

SHAPE CATEGORY	Inverted Cone					Pear				
SHAPE										
ISO#	012/014	012/016	013/025	225/020	226/024	237/015	237/018	237/025	237/033	237/034
LENGTH	1.4	1.6	1.7	4.4	7.0	4.5	4.5	6.5	5.0	5.0
GRIT	R	R	R	R	R	R	R	R	R	SC
PN (1 PK)	0864-1	0865-1	0866-1	0868-1	0870-1	0878-1	0879-1	0891-1	0892-1	892X-1

SHAPE CATEGORY	Wheel		Knife Edge Wheel		Round Edge Wheel				Egg	
SHAPE										
ISO#	042/044	041/052	304/028	304/050	068/034	068/046	068/047	068/046	277/029	277/029
LENGTH	1.4	0.8	0.8	1.2	1.3	1.8	1.9	1.8	4.9	4.9
GRIT	R	R	R	R	R	R	SC	C	R	SF
PN (1 PK)	0885-1	0886-1	0888-1	0890-1	0894-1	0895-1	895X-1	895C-1	0893-1	893V-1

SHAPE CATEGORY	Guide Groove Cutter		Bud							
SHAPE										
ISO#	522/026	552/036	243/017	243/014	257/022	243/022	243/023	243/022	243/020	243/019
LENGTH	3.6	4.0	3.0	3.4	3.5	5.0	5.0	5.0	5.0	5.0
GRIT	R	R	R	SF	R	R	SC	C	F	SF
PN (1 PK)	0896-1	0897-1	0880-1	880V-1	0881-1	0883-1	883X-1	883C-1	883F-1	883V-1



SHAPE CATEGORY	Round						
SHAPE							
ISO#	001/014	001/016	001/015	001/018	001/017	001/022	001/024
LENGTH	1.0	1.2	1.3	1.4	1.4	1.8	1.8
GRIT	R	R	SF	R	F	R	SC
PN (1 PK)	0872-1	0874-1	874V-1	0875-1	875F-1	0876-1	876X-1

SHAPE CATEGORY	Minimally Invasive Preparation					
SHAPE						
ISO#	001/009	001/011	001/013	237/008	237/010	237/012
LENGTH	0.7	0.9	1.1	2.0	2.0	2.0
GRIT	R	R	R	R	R	R
PN (1 PK)	0944-1	0945-1	0946-1	0947-1	0948-1	0949-1

PRICING - ROBOT® FG DIAMONDS (ONLY SOLD IN A PACKAGE OF 6)

PN	Product	Suggested Retail
See charts	Fine/Regular/Coarse (1 per pack)	\$8.85/6 pk \$53.10
See charts	Super Coarse (1 per pack)	\$8.85/6 pk \$53.10
See charts	Super Fine (1 per pack)	\$8.85/6 pk \$53.10
See charts	Minimally Invasive (1 per pack)	\$8.85/6 pk \$53.10
0955BB	All Ceramic Preparation Kit (see page 34)	\$200.00

TPE DIAMONDS POINTS

SHAPE CATEGORY	TPE									
SHAPE										
ISO#	150/010	150/010	150/012	150/012	150/012	150/014	150/014	150/014	150/016	150/016
LENGTH	9.0	9.0	9.0	9.0	6.0	9.0	9.0	6.0	9.0	9.0
GRIT	C	F	C	F	C	C	F	C	C	F
PN (1 PK)	680C-1	680F-1	682C-1	682F-1	682S-1	684C-1	684F-1	684S-1	686C-1	686F-1

DIAMONDS/CARBIDES

HYBRID POINTS® (DIAMONDS)

- Hybrid Points cut like a 6 fluted carbide bur

SHAPE CATEGORY	2	4	6	34	35	245	331	556/56	557/57	558/58
SHAPE										
ISO#	001/010	001/014	001/018	010/008	010/010	238/008	237/010	108/009	107/010	107/012
LENGTH	1.0	1.0	1.6	1.0	0.8	3.0	1.4	2.9	4.0	4.5
GRIT	R	R	R	R	R	R	R	R	R	R
PN (1 PK)	0901-1	0903-1	0905-1	0907-1	0909-1	0911-1	0913-1	0915-1	0917-1	0919-1

SHAPE CATEGORY	700/170	701/701
SHAPE		
ISO#	170/010	170/012
LENGTH	4.0	4.0
GRIT	R	R
PN (1 PK)	0921-1	0923-1

T&F HYBRID POINTS® (DIAMONDS)

- T&F Hybrid Points cut like a 12 fluted carbide bur

SHAPE CATEGORY	7006	7104	7106	7404	7406	7803	7901	7903	LT1	LT2
SHAPE										
ISO#	001/018	254/014	254/018	277/014	277/018	499/012	243/009	243/012	164/013	165/011
LENGTH	1.6	3.0	3.5	3.0	3.5	4.0	3.7	4.2	6.0	8.0
GRIT	SF	SF	SF	SF	SF	SF	SF	SF	SF	SF
PN (1 PK)	0925-1	0927-1	0929-1	0931-1	0933-1	0935-1	0937-1	0939-1	0941-1	0943-1



PRICING - TPE, HYBRID POINTS®, T&F HYBRID POINTS®

PN	Product	Suggested Retail
See charts	TPE - End Cutting (1 per pack)	\$18.45
See charts	TPE - Short End Cutting (1 per pack)	\$14.10
682S-1, 684S-1		
See charts	Hybrid Points (1 per pack)	\$8.85/6 pk \$53.10
0900	Hybrid Points Kit	\$158.70
12 Hybrid Points, 10 T&F Hybrid Points		
See charts	T&F Hybrid Points (1 per pack)	\$11.10/6 pk \$66.60
0950	T&F Hybrid Points Kit (all shapes).	\$104.80

GINGICURETTAGE® DIAMONDS

SHAPE CATEGORY	Gingicurettage										
SHAPE											
ISO#	012	014	017	012	016	017	019	021	013	015	021
PN (1 PK)	C061-1	C062-1	C063-1	C081-1	C082-1	C083-1	C084-1	C085-1	C101-1	C102-1	C103-1

PRICING - GINGICURETTAGE® DIAMONDS

PN	Product	Suggested Retail
See charts	Gingicurettage Diamonds (1 per pack).	\$16.60

ABRASIVE DRESSER



- Diamond coated grit
- Excellent tool for abrasive cleaning & shaping

PRICING - ABRASIVE DRESSER

PN	Product	Suggested Retail
0781	Abrasive Dresser.	\$68.45
1 Abrasive Dresser		

DIAMONDS/CARBIDES

ROBOT® CARBIDE FG FINISHERS



- High quality 10 & 12 blade tungsten carbide finishing head on a stainless steel shank
- Surface finish of all types of restorative materials, both anterior & posterior regions
- **Use Robot Carbide FG Finishers for:** all restorative materials

SHAPE CATEGORY	SH 41	SH 41	SH 46	SH 46	SH 379	SH 379	SH 379
SHAPE							
HEAD SIZE (MM)	014	018	014	018	014	018	023
HEAD LENGTH (MM)	-	-	3.5	3.8	3.1	3.5	4.2
# OF BLADES	12	12	10	10	10	10	10
MAXIMUM RPM	450,000	450,000	450,000	450,000	450,000	450,000	300,000
FG (3 pk)	PN 0703	0704	0709	0710	0717	0718	0719

SHAPE CATEGORY	SH 375R	SH 247	SH 246	SH 246	SH 246
SHAPE					
HEAD SIZE (MM)	014	009	009	010	012
HEAD LENGTH (MM)	8.0	3.2	3.6	3.6	3.6
# OF BLADES	10	10	10	10	10
MAXIMUM RPM	300,000	450,000	450,000	450,000	450,000
FG (3 pk)	PN 0727	0736	0741	0742	0743

PRICING - ROBOT® CARBIDE FG FINISHERS

PN	Product	Suggested Retail
See charts	FG Finishers (3 per pack)	\$23.50
0727	FG Finishers - SH 375R (3 per pack)	\$27.95



ROBOT® CARBIDE CTF-FG FINISHERS



- High quality tungsten carbide finishing head on a stainless steel shank
- Trim & finish subgingival cosmetic restorations without damaging the subgingival tissue
- Available in three blade types: 8, 16 (yellow-banded shank) & 30 (white-banded shank)
- 8 blades: rapid trimming & contouring; 16 blades: Fine Finishing (F); 30 blades: Ultra Fine Finishing (UF)
- **Use Robot Carbide CTF-FG Finishers for:** cosmetic restorations

SHAPE CATEGORY	SH 132	SH 132 F	SH 132 UF	SH 133	SH 133F	SH 133 UF
SHAPE						
HEAD SIZE (MM)	008	008	008	010	010	010
HEAD LENGTH (MM)	3.0	3.0	3.0	4.2	4.2	4.2
# OF BLADES	8	16	30	8	16	30
MAXIMUM RPM	300,000	300,000	300,000	450,000	450,000	450,000
FG (3 PK)	PN 0761-3	0762-3	0763-3	0764-3	0765-3	0766-3

SHAPE CATEGORY	SH 134	SH 134 F	SH 134 UF	SH 135	SH 135F	SH 135 UF
SHAPE						
HEAD SIZE (MM)	014	014	014	014	014	014
HEAD LENGTH (MM)	6.0	6.0	6.0	9.0	9.0	9.0
# OF BLADES	8	16	30	8	16	30
MAXIMUM RPM	450,000	450,000	450,000	300,000	300,000	300,000
FG (3 PK)	PN 0767-3	0768-3	0769-3	0770-3	0771-3	0772-3

PRICING - ROBOT® CARBIDE CTF-FG FINISHERS

PN	Product	Suggested Retail
See charts	CTF-FG Regular - 8 blades (3 per pack)	\$30.75
See charts	CTF-FG Fine - 16 blades (3 per pack)	\$35.15
See charts	CTF-FG Ultrafine - 30 blades (3 per pack)	\$39.55

DIAMONDS/CARBIDES

ROBOT® CARBIDE HP CUTTERS



- High quality tungsten carbide finishing head on a stainless steel shank
- Universal cross cut
- **Use Robot Carbide HP Cutters for:** predominantly base, noble and high noble metal alloys, chrome cobalt, model casts, plaster & acrylic

SHAPE CATEGORY	SH 71E	SH 71E	SH 72E	SH 73E	SH 73E	SH 77E	SH 77E	SH 78E	SH 79E	SH 79E
SHAPE										
HEAD SIZE (MM)	023	050	060	014	060	023	060	023	040	050
HEAD LENGTH (MM)	-	-	12.0	3.1	10.0	5.0	11.0	6.0	13.0	13.0
# OF BLADES	10/9	16/10	18/12	10/6	18/12	14/12	18/12	14/9	12/8	16/13
MAXIMUM RPM	50,000	50,000	50,000	50,000	50,000	50,000	50,000	50,000	50,000	50,000
HP (1/PK)	PN	0602	0603	0605	0606	0609	0610	0613	0614	0616
										0617

SHAPE CATEGORY	SH 251E	SH 351E	SH 137E	SH 261E	SH 295E	SH 129E	SH 33
SHAPE							
HEAD SIZE (MM)	060	060	016	023	023	023	009
HEAD LENGTH (MM)	14.0	11.0	4.0	13.0	15.0	8.0	4.2
# OF BLADES	18/12	18/12	10/6	10/9	14/12	14/12	6
MAXIMUM RPM	50,000	50,000	50,000	50,000	50,000	50,000	50,000
HP (1/PK)	PN	0620	0623	0626	0630	0631	0634
							0667 (6PK)

PRICING - ROBOT® CARBIDE HP CUTTERS

PN	Product	Suggested Retail
See charts	HP Cutters (1 per pack)	\$24.30
0602, 0606, 0610, 0614, 0626, 0630, 0631, 0634		
See charts	HP Cutters (1 per pack)	\$28.75
0603, 0609, 0613, 0605, 0616, 0617, 0620, 0623		
0667	HP Cutters (6 per pack)	\$16.95/6 pk \$101.70



ARTICULATORS




- Durable & light weight articulators to reproduce the natural movement of the patient's jaw on the basis of mandible/maxilla relationship
- Provide flexible & versatile working features that are only found in more complex & expensive instruments
- Competitively priced

PRICING - ARTICULATORS

PN	Product	Suggested Retail
5070	Handy IIM (without occlusal plate)	\$169.75
5071	Handy IIM Complete (with occlusal plate)	\$185.35
5072	Metal Mounting Plate	\$72.80
4 universal mounting plates for Handy IIM & ProArch		
5073	Incisal Rod	\$10.25
5074	Incisal Rod Lockscrew	\$5.85
5075	Incisal Guide Pin	\$5.85
5076	Occlusal Plate	\$22.90
5078	Spring	\$2.10
5079	Central Screw with 2 Knurled Knobs	\$18.45
4070	ProArch	\$243.60
4071	ProArch Complete (with occlusal plate)	\$273.15

EQUIPMENT

LAB AIR-Z HIGH SPEED TURBINE



- Robust, ultra-quiet, precision instrument designed for extended laboratory duty with speed of 360,000 rpm
- Technician can choose between Oil-Free & regular Lab Air-Z versions
- Proven technology for technicians who seek value for a fair price

PRICING - LAB AIR-Z HIGH SPEED TURBINE		
PN	Product	Suggested Retail
5239	Lab Air-Z System.	\$922.95
5339	Lab Air-Z Oil-Free System.	\$941.45
Includes hand piece, foot air controller, air filter/pressure regulator & mounting bracket		
5320	Lab Air-Z Handpiece Oil-Free.	\$664.55
5221	Lab Air-Z Turbine Cartridge	\$280.50
5321	Lab Air-Z Turbine Cartridge Oil-Free	\$302.65
5222	Metal Chuck	\$56.05
5223	Chuck Tool	\$21.40
5224	Cartridge Tool.	\$31.75
5230	Foot Aircontroller.	\$141.70
5234	Air Filter/Pressure Regulator	\$137.25

CEMENTS

BEAUTICEM SA	56
MONOCem SELF-ADHESIVE RESIN CEMENT	56
HY-BOND POLYCARBOXYLATE CEMENT & REFILLS	57
HY-BOND POLYCARBOXYLATE TEMPORARY CEMENT	57
HY-BOND ZINC OXIDE EUGENOL TEMPORARY CEMENT	58
HY-BOND ZINC PHOSPHATE CEMENT & REFILLS	58
GLASIONOMER BASE CEMENT & REFILLS	59
CX-PLUS GLASIONOMER & REFILLS (NOW IN CAPSULES!)	60
CORESHADE GLASIONOMER BASE CEMENT & REFILLS	61
IP TEMP CEMENT	61
RESICem UNIVERSAL RESIN CEMENT	62
CEMENT ACCESSORIES	63
CEMENT DURATION CHART	64

Official Partner



CEMENTS

BEAUTICEM SA



PRG
Technology
Giomar

- Self-etching, self-adhesive resin cement
- S-PRG (Surface Pre-reacted Glass Ionomer) – ability to recharge fluoride
- No primer required (except porcelain)
- High bond values across all substrates
- Auto-mixing & dual cure
- High levels of radiopacity
- Low film thickness (12µm)

PRICING - BEAUTICEM SA

PN	Product	Suggested Retail
3213	BeautiCem SA Clear	\$119.35
	5ml paste, 10 mixer tips & step card	
3214	BeautiCem SA Ivory	\$119.35
	5ml paste, 10 mixer tips & step card	

MONOCEM® SELF-ADHESIVE RESIN CEMENT



- Self-etch & self-adhesive
- Bonds to all substrates with high retention values
- Moisture tolerant with two second tack cure for easy clean-up
- Eliminates sensitivity & has high fluoride release
- 100% polymerization/dual cure: three minute auto cure, 20 second light cure
- Command set, curing begins at restoration's placement
- Very low film thickness (<12µm)

PRICING - MONOCEM® SELF-ADHESIVE RESIN CEMENT

PN	Product	Suggested Retail
3208A	MonoCem Translucent Kit	\$126.30
	2 (3.5g) Translucent, 15 applicator tips & 2 OneGloss PS Midi-Point polishers	
3209A	MonoCem Bleach White Kit	\$126.30
	2 (3.5g) Bleach White, 15 applicator tips & 2 OneGloss PS Midi-Point polishers	
3211	MonoCem Auto-mix Tips	\$12.75
	20 applicator tips	



HY-BOND® POLYCARBOXYLATE CEMENT (LUTING)



- Contains tannin fluoride, which helps seal dentinal tubules, reducing post-operative sensitivity
- Chemically bonds to tooth structure & metal appliances
- High adhesive strength
- Extremely low film thickness (19µm)
- High level of radiopacity

PRICING - HY-BOND® POLYCARBOXYLATE CEMENT

PN	Product	Suggested Retail
1160	HY-Bond Polycarboxylate Cement Kit.	\$136.80
	60g powder, 33ml liquid, mixing pad, spatula, dispenser & instructions	
1161	Polycarboxylate Powder 60g.	\$74.75
1162	Polycarboxylate Liquid 33ml	\$70.85

HY-BOND® POLYCARBOXYLATE TEMPORARY CEMENT



- 79% of independent evaluators said they would switch, "It's easier to clean-up & easier to remove"*
- Excellent for temporary crowns, bridgework, implant prostheses, long-term provisional & temporary cementation of permanent restorations
- Eugenol-free
- Helps prevent microleakage & improves the marginal seal
- Excellent bonding strength for provisional restorations
- Fluoride additive for protection of tooth structure & reduction of post-operative sensitivity
- Economically priced compared to other leading brands

**Results from a recent study by an internationally respected, independent research institute.*

PRICING - HY-BOND® POLYCARBOXYLATE TEMPORARY CEMENT

PN	Product	Suggested Retail
1173	HY-Bond Polycarboxylate Temporary Cement Kit. . .	\$198.95
	125g powder, 60ml liquid, mixing pad, spatula & instructions	
1176	HY-Bond Polycarboxylate Temporary Cement Kit. .	\$110.60
	60g powder, 30ml liquid, mixing pad, spatula & instructions	

CEMENTS

HY-BOND® ZINC OXIDE EUGENOL TEMPORARY CEMENT



- Tannin fluoride seals the dentinal tubules & reduces post-operative sensitivity
- Helps prevent the development of caries & the decalcification of dentin and enamel
- Excellent handling for cementing temporary crowns or filling materials
- High level of radiopacity
- Sets up promptly even in presence of moisture
- Easy to mix & clean-up

PRICING - HY-BOND® ZINC OXIDE EUGENOL TEMPORARY CEMENT

PN	Product	Suggested Retail
1180	Zinc Oxide Eugenol (ZOE) Temporary Cement Kit	\$53.30
	20g powder, 12ml liquid, dispenser & instructions	

HY-BOND® ZINC PHOSPHATE CEMENT



- Seals dentinal tubules & minimizes post-operative sensitivity
- Excellent for seating permanent prostheses
- Film thickness (25µm)
- High compressive strength (98 MPa)
- High level of radiopacity
- Apply directly to prepared dentin
- Neutralizes the harshness of phosphoric acid

PRICING - HY-BOND® ZINC PHOSPHATE CEMENT

PN	Product	Suggested Retail
1170	Zinc Phosphate Cement Kit	\$53.30
	60g powder, 25ml liquid, dispenser & instructions	
1171	Zinc Phosphate Powder 60g	\$38.10
1172	Zinc Phosphate Liquid 25ml	\$21.50



GLASIONOMER BASE CEMENT (DENTIN REPLACEMENT)



- Recommended as a supportive, protective base under stress-bearing amalgam, composite, porcelain & indirect composite restorations
- Biocompatible & chemically bonds to dentin and enamel
- High fluoride release reduces the occurrence of secondary caries
- Compressive strength (147 MPa)
- High level of radiopacity
- Low solubility
- Easy mixing & handling
- Pediatric temporary cusp replacement

PRICING - GLASIONOMER BASE CEMENT

PN	Product	Suggested Retail
1110	Glaslonomer Base Cement Kit	\$111.35
	15g yellow color powder, 5ml liquid, mixing pad, spatula, dispenser & instructions	
1111	Yellow Powder 15g	\$72.45
1113	White Powder 15g	\$72.45
1114	Base Cement Liquid 5ml	\$41.30

CEMENTS

CX-PLUS GLASIONOMER (TYPE I LUTING CEMENT)



- NOW IN CAPSULES!
- Indicated for cementing crowns, inlays, long-span bridges & orthodontic appliances
- Low film thickness (15µm) with high compressive strength (149 MPa)
- Bonds to both dentin & enamel to provide an excellent marginal seal to overcome post-operative sensitivity
- High fluoride release & high level of radiopacity
- Easy mixing, handling & resistant to water contamination at the initial stage
- 3 year shelf life

PRICING - CX-PLUS GLASIONOMER

PN	Product	Suggested Retail
1166	CX-Plus Cement Kit	\$131.30
	35g powder, 17ml liquid, mixing pad, spatula, dispenser & instructions	
1167	CX-Plus Powder 35g	\$85.90
1168	CX-Plus Liquid 17ml	\$52.45
1169	CX-Plus Triple Kit	\$267.35
	3 (35g) powder, 3 (17ml) liquid, mixing pad, spatula, dispenser & instructions	
3277	NEW! CX-Plus Capsules	\$263.85
	50 (0.53g) capsules & instructions	



CORESHADE GLASIONOMER (CORE BUILD-UP BASE CEMENT)



- Supportive and protective base under stress-bearing fillings & restorations
- High compressive strength (147 MPa)
- Excellent bonding to tooth structure, metal pins & posts
- Gray in color for easy identification during placement & finishing
- High fluoride release & high level of radiopacity
- Biocompatible

PRICING - CORESHADE GLASIONOMER

PN	Product	Suggested Retail
1115	CoreShade Base Cement Kit	\$159.10
	25g gray powder, 10ml liquid, mixing pad, dispenser & instructions	
1116	Gray Powder 25g	\$97.00
1117	CoreShade Liquid 10ml	\$69.25

IP TEMP CEMENT



- Bioactive, fluoride-releasing Giomer temporary luting cement for implants & provisional crowns
- Easily remove excess cement with water spray
- Very low film thickness (10µm)

PRICING - IP TEMP CEMENT

PN	Product	Suggested Retail
1197	IP Temp Cement 1-1 Set	\$77.10
	45g powder, 26ml liquid, measuring spoon, mixing pad & instructions	
1198	IP Temp Cement Powder	\$38.60
1199	IP Temp Cement Liquid	\$38.60

CEMENTS

RESICEM® UNIVERSAL RESIN CEMENT



- Designed specially for zirconia, alumina & Shofu HC Blocks – highest reported bond strengths of any resin cement
- Formulated adhesion & aesthetics
- Customize to be universal for all materials
- Self-conditioning
- High fluoride release
- Very low film thickness (9µm)

PRICING - RESICEM® UNIVERSAL RESIN CEMENT

PN	Product	Suggested Retail
3220	ResiCem Basic Kit	\$360.70
	2 (5ml) syringes, 3ml Primer A, 2.5ml Primer B, 20 mixing tips, 1 fine microtip handle, 50 fine microtip applicators & petit dish	
3221	ResiCem	\$141.90
	5ml syringe & 10 mixing tips	
3227	Primer A (3ml)	\$68.10
3228	Primer B (2.5ml)	\$90.80
3232	AZ Primer (5ml)	\$97.75
	for bonding to zirconia (zirconia primer)	
1740	ML Primer (5ml)	\$55.20
	for bonding to metal (metal primer)	
3234	Porcelain Primer (5ml)	\$87.30
	for bonding to porcelain & indirect composites	



CEMENT ACCESSORIES

PRICING - CEMENT ACCESSORIES

PN	Product	Suggested Retail
1163	Mixing Pad (6 per pack)	\$16.10
1165	Spatulas (6 per pack)	\$6.05
1139	Cocoa Butter (10g)	\$12.65



SHOFU RAVES AND RECOMMENDATIONS

MonoCem® Self-Adhesive Resin Cement

"The most important factor I look for in a resin cement is the elimination of post-operative sensitivity. When I cement all porcelain restorations, I absolutely do not want a patient calling me the next day and stating that their \$900 investment hurts them. Compressive strength, tensile strength, solubility and retentive strength are all nice numbers to have but they will not help my patient relations if my patient is hurting.

The self-etching removes the chance of my using chemicals which could produce sensitivity.

I trust the Shofu name and that the company will test their product before it gets to me and not let me be the test laboratory." ~ *Dr. Craig Schrager*

"MonoCem eliminates many steps to get the job done. It is a pure, self-etch and self-adhesive resin cement that has unlimited set time and will polymerized at 100%. It is backed by fast and professional results. I highly recommend it!" ~ *Dr. Marshall Oleck*

CEMENTS

CEMENT DURATION CHART

	MIXING	WORKING	SETTING
HY-BOND CEMENT			
POLYCARBOXYLATE	within 60 sec	3 min start of mixing	6 min end of mixing
POLYCARBOXYLATE TEMPORARY	within 60 sec until homogenous	-	5 min end of mixing
ZINC OXIDE EUGENOL TEMPORARY	within 60 sec	2 min end of mixing until crown is seated	5 min end of mixing
ZINC PHOSPHATE	within 60 sec	3 min/30 sec start of mixing	6 min/30 sec end of mixing
RESI GLASS	10 sec	3 min/30 sec start of mixing	7 min end of mixing
GLASIONOMER CEMENT			
BASE	not to exceed 30 sec	2 min/30 sec start of mixing	4 min end of mixing
CORE SHADE	not to exceed 30 sec	2 min/50 sec start of mixing	4 min end of mixing
CX-PLUS	30 sec	3 min/30 sec start of mixing	3 min/45 sec end of mixing
CX-PLUS CAPSULES	10 sec	3 min/10 sec start of mixing	4 min/10 sec end of mixing

PREVENTIVES

BEAUTISEALANT	66
PRG BARRIER COAT	66



Official Partner



PREVENTIVES

BEAUTISEALANT



PRG
Technology
Giomer

- Fluoride releasing pit & fissure sealant system
- Ideal handling & viscosity
- Radiopacity equal to dentin
- Unique bioactive properties of S-PRG fillers prevents plaque formation & build-up
- High shear bond strength without phosphoric acid etch

PRICING - BEAUTISEALANT

PN	Product	Suggested Retail
1798	BeautiSealant	\$99.50
	1.2g BeautiSealant paste, 3ml BeautiSealant primer, 15 needle tips, 50 microbrushes, 25 V-Dishes	
1799	BeautiSealant Paste 1.2g	\$39.10
1800	BeautiSealant Primer 3ml	\$55.90
1794	BeautiSealant Needle Tips.	\$19.55
	50 (27g) needle tips	

PRG BARRIER COAT



PRG
Technology
Giomer

- Quick & easy to apply
- HEMA, Acetone & Alcohol Free
- Fluoride release & rechargeability, through Giomer technology
- Anti-plaque effect
- Ideal for enhancing patient's intraoral health
- Indications other than hypersensitivity: Hard to brush areas (orthodontic brackets, crowded teeth) & high caries risk areas (exposed root surfaces, newly erupted molars, white spots)

PRICING - PRG BARRIER COAT

PN	Product	Suggested Retail
1786	PRG Barrier Coat	\$39.65
	8 Base Units 0.02ml each, 1 Active 0.09ml, 8 brushes	

RESTORATIVES

BEAUTIFIL II	68
NEW! BEAUTIFIL II LS (LOW SHRINK).	69
NEW! BEAUTIFIL II GINGIVA & ENAMEL SHADES	70
BEAUTIFIL FLOW	71
BEAUTIFIL FLOW PLUS (NOW IN TIPS!).	72 - 73
BEAUTIBOND	74
BEAUTIFIL FLOW PLUS/BEAUTIBOND – ZERO & LOW FLOW KITS.	74
BEAUTIFIL-BULK FLOWABLE.	75
BEAUTIFIL-BULK RESTORATIVE	75
BEAUTIFIL OPAQUER	76
FL-BOND II ADHESIVE SYSTEM	76
CERARESIN BOND	77
GLASIONOMER FX-II CAPSULES	77
GLASIONOMER TYPE II CEMENT RESTORATIVE & REFILLS	78
RESTORATIVES DURATION CHART.	78
CAD/CAM	
NEW! SHOFU BLOCK HC	79
CAMERA	
EYESPECIAL C-III	80

Official Partner



RESTORATIVES

BEAUTIFIL® II



PRG
Technology
Giomar



- 8-year University of Florida clinical study showed no failures, no post-operative sensitivity & high retention of color match and luster
- 13-year recall study: 96% of restorations with no secondary caries, in 66% retention rate (27 of 41)
- S-PRG (Surface Pre-reacted Glass Ionomer) – ability to release/recharge fluoride, strontium & four other ions
- Optimal for patients with a high caries index
- Tooth-like fluorescence & life-like aesthetics
- Eliminates shade variations associated with polymerization
- High level of radiopacity (3.4 Al:mm)
- Enamel-like surface hardness
- See page 30, for Artistic Composite Kit

PRICING - BEAUTIFIL® II

PN	Product	Suggested Retail
1775	Beautifil II 6 Color Set	\$264.10
6 (4.5g) syringes A2, A3, A3.5, B2, A3O, INC; Dura-White Stone, Super-Snap Singles, OneGloss Trial Set & paper pad		
1776	Beautifil II Tips 6 Color Set	\$493.00
120 (0.25g) tips 6 colors, 4.5g in A2, A3, A3.5, B2, A3O, INC; Dura-White Stone, Super-Snap Singles, OneGloss Trial Set & paper pad		
1777	Beautifil II Cosmetic Kit	\$249.60
3 (4.5g) syringes A2, A3, INC; Beautifil Opaquer-LO, Beautifil Flow F02-A2, Dura-White Stone, Super-Snap Singles, OneGloss Trial Set & paper pad		
Beautifil II Syringe Refill 4.5g		\$40.55
1401 - A1	1408 - B1	1415 - INC (incisal)
1402 - A2	1409 - B2	1416 - BW (bleach)
1403 - A3	1410 - B3	1419 - A2O
1404 - A3.5	1412 - C2	1420 - A3O
1405 - A4	1413 - C3	
Beautifil II Tips 0.25g x 20 tips		\$86.70
1752 - A1	1759 - B1	1770 - A2O
1753 - A2	1760 - B2	1771 - A3O
1754 - A3	1763 - C2	
1755 - A3.5	1766 - INC (incisal)	
1756 - A4	1767 - BW (bleach)	



NEW! BEAUTIFIL® II LS (LOW SHRINK)



- Low volumetric shrinkage (0.85%) & shrinkage stress (2.72 MPa)
- Great strength & high wear resistance
- Tooth-like aesthetics with natural fluorescence & chameleon effect
- Polishes in an instant
- Award-winning handling
- Sustained fluoride release/recharge with bioactive Giomer chemistry

Giomer
PRG
Technology

PRICING - BEAUTIFIL® II LS (LOW SHRINK)

PN

Product

Suggested Retail

BIILSK1 Beautifil II LS Intro Tips Kit \$149.00

20 (0.25g) tips A2 and A3.5

Bonus Items: 2 OneGloss PS Midi-Points

BIILSK2 Beautifil II LS A2 Bulk Tips Kit \$399.00

80 (0.25g) tips A2

Bonus Items: 10 (0.25) tips A2, 25-hole BurButler with Tall Lid

Beautifil II LS Syringe Refill 4g \$69.95

Y2270 - A1

Y2275 - B1

Y2280 - A2O

Y2271 - A2

Y2276 - B2

Y2281 - A3O

Y2272 - A3

Y2277 - B3

Y2282 - INC (incisal)

Y2273 - A3.5

Y2278 - C2

Y2283 - BW (bleach)

Y2274 - A4

Y2279 - C3

Beautifil II LS Tips 0.25g x 20 tips \$99.95

Y2284 - A1

Y2289 - B1

Y2294 - A2O

Y2285 - A2

Y2290 - B2

Y2295 - A3O

Y2286 - A3

Y2291 - B3

Y2296 - INC (incisal)

Y2287 - A3.5

Y2292 - C2

Y2297 - BW (bleach)

Y2288 - A4

Y2293 - C3

RESTORATIVES

NEW! BEAUTIFIL® II GINGIVA & ENAMEL SHADES



- Beautifil II Gingiva was created based on Shofu's highly acclaimed Ceramage line, a nano-ceramic indirect composite
- Beautifil II Gingiva aesthetically blends with different hues of the natural gingiva
- Beautifil II Enamel has optical characteristics of natural teeth
- Designed for artistic restoration of anterior regions
- S-PRG (Surface Pre-reacted Glass Ionomer) – ability to recharge fluoride
- High level of radiopacity (3.4 Al:mm)
- Complementary with the entire Beautifil line of restorative materials
- See page 31, for PINK Prep Kit

PRG
Technology
Giomer

PRICING - BEAUTIFIL® II GINGIVA & ENAMEL SHADES

PN	Product	Suggested Retail
----	---------	------------------

GMK1	Beautifil II PINK Kit	\$249.00
-----------------------	--	-----------------

2.5g Gingiva syringes in Light Pink, Dark Pink, Violet, Brown & Orange;

2g Beautifil Flow F02 Gingiva syringe

Bonus Items: Beautifil Flow Plus 2 (0.21g) tips F00-A2, Beautifil Flow Plus 2 (0.21g) tips F03-A2, Beautifil II 2 (0.25g) tips A2

GMK3	Beautifil II PINK & WHITE Kit	\$399.00
-----------------------	--	-----------------

2.5g Gingiva syringes in Light Pink, Dark Pink, Violet, Brown & Orange;

2g Beautifil Flow F02 Gingiva syringe; 2.5g Enamel syringes in Translucent, High-Value Translucent, Low-Value Translucent & Amber

Bonus Items: Beautifil Flow Plus 5 (0.21g) tips F00-A2, Beautifil Flow Plus 5 (0.21g) tips F03-A2, 10 (0.10ml) BeautiBond Unit Dose, 10 Super-Snap Singles Mini Size Disk (8mm)

Beautifil II Enamel 2.5g	\$37.50
---	----------------

Y2250 - T	Y2252 - LVT
------------------	--------------------

Y2251 - HVT	Y2253 - AM
--------------------	-------------------

Beautifil II Gingiva 2.5g	\$49.95
--	----------------

Y2256 - G-LP	Y2259 - G-Br
---------------------	---------------------

Y2257 - G-DP	Y2260 - G-Or
---------------------	---------------------

Y2258 - G-V	
--------------------	--



BEAUTIFIL® FLOW



PRE
Technology
Giomer



- Complements Beautifil II to further enhance its applications
- Two viscosities, Beautifil Flow F02 (low flow) & Beautifil Flow F10 (high flow), maintain the same aesthetics, physical properties, particle size & filler load
- Ideal as a base or liner to help prevent caries recurrence
- Root caries & enamel defects
- Deciduous dentition
- Repair of composite or ceramic veneers
- Operator can select viscosity depending on the clinical case
- Low Flow: effective to use where paste flows moderately & shape should be retained, such as cervical areas
- High Flow: effective to use where smooth flow is necessary such as uneven or narrow spaces & usage as liner

PRICING - BEAUTIFIL® FLOW

PN	Product	Suggested Retail
	Beautifil Flow F02 Low Flow 2g (5 needle tips)	\$34.35
1431 - A1	1434 - A3.5	1451 - G
1432 - A2	1435 - A4	
1433 - A3	1449 - A3O	
	Beautifil Flow F10 High Flow 2g (5 needle tips)	\$34.35
1461 - A1	1464 - A3.5	1482 - T
1462 - A2	1465 - A4	
1463 - A3	1479 - A3O	
1484	Beautifil Flow Needle Tips	\$22.35
50 needle tips		

RESTORATIVES

BEAUTIFIL FLOW PLUS®



PRG
Technology
Gomer



- NOW IN TIPS!
- A base, liner & final restorative all in one!
- Indicated in classes I, II, III, IV & V
- Physical properties rival hybrids
- Stackable & sculptable; stays put
- Self-polishing "leveling-effect"
- High radiopacity (2.3 AI:mm)
- S-PRG (Surface Pre-reacted Glass Ionomer) – ability to recharge fluoride
- Two viscosities: stackable F00 (zero flow) & self-leveling F03 (low flow)

PRICING - BEAUTIFIL FLOW PLUS®

PN	Product	Suggested Retail
----	---------	------------------

2000P	Beautifil Flow Plus Pedo Kit*	\$118.45
------------------------	--	-----------------

4 (2.2g) syringes F00 & F03 in A1, Bleach White, 2 Beautifil II tips, 2 BeautiBond, 6 assorted OneGloss PS, 2 Super-Snap Singles

2000S	Beautifil Flow Plus Standard Kit*	\$118.45
------------------------	--	-----------------

4 (2.2g) syringes F00 & F03 in A2, A3, 2 Beautifil II tips, 2 BeautiBond, 6 assorted OneGloss PS, 2 Super-Snap Singles

Beautifil Flow Plus F00 Syringe 2.2g	\$43.35
---	----------------

2000 - A0.5
2001 - A1
2002 - A2
2003 - A3
2004 - A3.5

2005 - A4
2006 - A0.5O
2007 - A1O
2008 - A2O
2062 - B1

2063 - B2
2064 - C2
2065 - D2
2010 - INC (incisal)
2011 - BW (bleach white)

Beautifil Flow Plus F03 Syringe 2.2g	\$43.35
---	----------------

2013 - A0.5
2014 - A1
2015 - A2
2016 - A3
2017 - A3.5
2018 - A4

2020 - A2O
2021 - A3O
2066 - B1
2067 - B2
2068 - C2
2069 - D2

2023 - INC (incisal)
2024 - BW (bleach white)
2026 - MI (milky white)
2027 - CV (cervical)

1484	Beautifil Flow Needle Tips	\$22.35
-----------------------	---	----------------

50 (20g) needle tips

*Limited time offer



PRICING - BEAUTIFIL FLOW PLUS®

PN	Product	Suggested Retail
NEW! Beautifil Flow Plus F00 Tips 0.21g x 20 \$93.35		
2200 - A0.5	2205 - A4	2207 - B2
2201 - A1	2211 - A0.5O	2208 - C2
2202 - A2	2212 - A1O	2209 - D2
2203 - A3	2213 - A2O	2215 - INC (incisal)
2204 - A3.5	2206 - B1	2216 - BW (bleach white)
NEW! Beautifil Flow Plus F03 Tips 0.21g x 20 \$93.35		
2218 - A0.5	2229 - A2O	2232 - INC (incisal)
2219 - A1	2230 - A3O	2233 - BW (bleach white)
2220 - A2	2224 - B1	2235 - MI (milky white)
2221 - A3	2225 - B2	2236 - CV (cervical)
2222 - A3.5	2226 - C2	
2223 - A4	2227 - D2	

RESTORATIVES

BEAUTIBOND®



- Available in bottles & unit doses
- Dual functional monomers providing equally high bond strength to both dentin & enamel
- Bond strengths comparable to leading multi-component brands
- Single coat for shorter working time
- Efficient 30 second application
- 5µm film thickness
- HEMA-free
- No mixing required
- Convenient unit dose

PRICING - BEAUTIBOND®

PN	Product	Suggested Retail
1781	BeautiBond (Bottle)	\$163.35
	6ml bottle, 50 microbrush tips, 25 V-dishes & instructions	
1782	BeautiBond (Unit Dose)	\$163.35
	50 (0.1ml) ampules, 50 microbrushes & instructions	

BEAUTIFIL FLOW PLUS®/BEAUTIBOND® – ZERO & LOW FLOW KITS



- **BeautiFlow Plus** a base, liner & final restorative all in one!
- Indicated in classes I, II, III, IV & V
- Two viscosities: stackable F00 (zero flow) & self-leveling F03 (low flow)
- **BeautiBond** has dual functional monomers providing equally high bond strength to both dentin & enamel
- Bond strengths comparable to leading multi-component brands
- 5µm film thickness

PRICING - BEAUTIFIL FLOW PLUS®/BEAUTIBOND® – ZERO & LOW FLOW KITS

PN	Product	Suggested Retail
2000BF00	Zero Flow Kit	\$318.85
	6 (2.2g) BeautiFlow Plus F00 syringes A1, A4, 2 - A2, 2 - A3; 50 (0.1ml) BeautiBond Unit Dose, 50 microbrushes	
2000BF03	Low Flow Kit	\$318.85
	6 (2.2g) BeautiFlow Plus F03 syringes A1, A4, 2 - A2, 2 - A3; 50 (0.1ml) BeautiBond Unit Dose, 50 microbrushes	



BEAUTIFIL-BULK FLOWABLE



Giomer
PRE
Technology



- Low shrinkage stress (2.06 MPa)
- 4mm depth of cure
- Ideal handling without the ooze
- Superior fluoride release & rechargability
- Self-leveling allows easy cavity wall adaptation
- High level of radiopacity (2.5 Al:mm)
- Base/liner in class I-II posterior restorations

PRICING - BEAUTIFIL-BULK FLOWABLE

PN

Product

Suggested Retail

Beautifil-Bulk Flowable Syringe 2.4g \$43.40

2030 - Universal

2031 - Dentin

Beautifil-Bulk Flowable Tips 0.23g x 20 tips \$93.35

2028 - Universal

2029 - Dentin

BEAUTIFIL-BULK RESTORATIVE



Giomer
PRE
Technology



- High fill ratio (87.0wt% 74.5vol%) low polymerization shrinkage (1.7%)
- Ideal balance of light diffusion & transmission ensures complete polymerization at 4mm depth of cure
- Outstanding Vickers hardness value (61 Hv)
- Optimum translucency creates aesthetic shades unaffected by surrounding intraoral color
- Sustained fluoride release & rechargability with Giomer technology
- High level of radiopacity (3.0 Al:mm)

PRICING - BEAUTIFIL-BULK RESTORATIVE

PN

Product

Suggested Retail

Beautifil-Bulk Restorative Syringe 4.5g \$40.55

2034 - Universal

2035 - A

Beautifil-Bulk Restorative Tips 0.25g x 20 tips \$86.70

2032 - Universal

2033 - A

RESTORATIVES

BEAUTIFIL® OPAQUER



- Mask discolored or stained teeth, metal color & pulp capping agent
- Light-curing, fluoride releasing restorative material
- Available in the two shades, Universal Opaque & Light Opaque

PRG
Technology
Giomar

PRICING - BEAUTIFIL® OPAQUER

PN	Product	Suggested Retail
	Beautiful Flow Opaquers 2g (5 needle tips)	\$40.55

1363 - UO (universal) 1364 - LO (light)

FL-BOND II ADHESIVE SYSTEM



- S-PRG (Surface Pre-Reacted Glass Ionomer) technology that fully reacts with alumina silica glass particles
- High levels of fluoride release & recharge
- Excellent bonding to tooth structure
- Bonds any composite to enamel & dentin for Class I-V & cervical erosions
- Self-etching saves time & reduces post-operative sensitivity
- High level of radiopacity
- Total procedure time: 35 seconds

PRG
Technology
Giomar

PRICING - FL-BOND II ADHESIVE SYSTEM

PN	Product	Suggested Retail
1306	FL-Bond II Complete Set	\$198.10
5ml bonding agent, 5ml primer, 25 red microbrushes, 25 yellow microbrushes, dish & instructions		
1307	Primer 5ml each	\$99.05
1308	Bonding Agent 5ml each	\$143.05
1309	Etchant 6ml each	\$48.40



CERARESIN BOND



- Porcelain & resin repair system
- Intra oral & lab repairs
- Bonds to porcelain, ceramic, resin, acrylic & metal (with metal primer)
- Resistant to moisture & thermal strains
- Quick & easy to use with a 3 year shelf life
- Compatible with all resins & composites

PRICING - CERARESIN BOND

PN	Product	Suggested Retail
1741	Cerarasin Bond	\$196.95
	6ml CRB Bond 1, 5mL CRB Bond 2, 25 yellow microbrushes, 25 pink microbrushes, 5 disposable dishes & instructions	
1740	ML Primer (5ml)	\$55.20
	for bonding to metal (metal primer)	
3232	AZ Primer (5ml)	\$97.75
	for bonding to zirconia (zirconia primer)	
Y0005	NEW! HC Primer (5ml)	\$110.00
	for bonding to Shofu Block HC	

GLASIONOMER FX-II CAPSULES



- No conditioner required
- Radiopaque & moisture tolerant
- When compared to leading competitors*
 - 39% greater fluoride release & recharge
 - Highest value in hardness & equal compressive strength
 - Significantly decreases acid erosion levels
- Ease of handling & packability
- Allows optimum shade matching & stability
- 3 year shelf life

**Shofu R&D test results*

PRICING - GLASIONOMER FX-II CAPSULES

PN	Product	Suggested Retail
	Glaslonomer FX-II 0.55g x 50 capsules	\$263.85
3270 - A2	3272 - A3.5	
3271 - A3	3273 - B2	

RESTORATIVES

GLASIONOMER TYPE II CEMENT RESTORATIVE & REFILLS



- Appropriate for restoring cervical erosions, Class III, Class V & pedodontic restorations
- Biocompatible with high fluoride release & recharge
- Available in 3 blendable shades: Universal, Yellow & Gray
- Chemically bonds to enamel & dentin
- High compressive strength (186 MPa)
- Coefficient of thermal expansion, similar to natural teeth & ensures marginal integrity
- Superb handling

PRICING - GLASIONOMER TYPE II CEMENT RESTORATIVE & REFILLS


PN	Product	Suggested Retail
1130	GlasIonomer Cement Restorative Type II	\$337.35
4-Color Set (15g each) Universal, Yellow, Gray & Brown powder, 16ml liquid, varnish, cocoa butter, 2 mixing pads, matrix strips, shade guide, 2 spatulas, 2 dispensers & instructions		
1131	Universal One-Color Set (for filling)	\$126.55
1132	Yellow One-Color Set (for filling)	
1133	Gray One-Color Set (for filling)	
15g color powder, 8ml liquid, varnish, cocoa butter, mixing pad, dispenser, spatula, matrix strips & instructions (each set)		
1141	Universal Color Powder 15g.	\$91.50
1142	Yellow Color Powder 15g	
1143	Gray Color Powder 15g	
1140	GlasIonomer Type II Liquid 8ml	\$41.30
1139	Cocoa Butter 10g	\$12.65
1148	Matrix Strips 20 pk.	\$7.85

GLASIONOMER TYPE II RESTORATIVE DURATION CHART

	MIXING	WORKING	SETTING
GLASIONOMER TYPE II	not to exceed 45 secs	2 min/50 secs start of mixing	3 min/15 secs end of mixing

CAD/CAM

NEW! SHOFU BLOCK HC - CEREC®



- Zirconia-filled, nano-ceramic CAD/CAM composite
- Indicated for full-contour anterior/posterior crowns, inlays/onlays & veneers
- High compressive (471 MPa) & flexural strength (191 MPa)
- High degree of milling accuracy with ideal marginal integrity
- Impervious to stains and color-retentive with self-polishing attributes

PN

PRICING - SHOFU BLOCK HC - CEREC®

Product

Suggested Retail

Block HC One-Layer 5 pcs \$139.95
(10mm x 12mm x 16mm)

2151C - W2-LT
2153C - A1-LT
2155C - A2-LT
2164C - A3-LT

2156C - A3.5-LT
2165C - B3-LT
2152C - A1-HT
2154C - A2-HT

2163C - A3-HT
2166C - 59 (incisal)
2167C - OC (occlusal)

Block HC Two-Layer 5 pcs. \$153.95
(10mm x 12mm x 16mm)

Z0011 - A2-2L
Z0012 - A3-2L

Z0013 - A3.5-2L

Block HC One-Layer 5 pcs \$158.75
(12mm x 14mm x 18mm)

2157C - W2-LT
2159C - A1-LT
2161C - A2-LT
2169C - A3-LT

2162C - A3.5-LT
2170C - B3-LT
2158C - A1-HT
2160C - A2-HT

2168C - A3-HT
2171C - 59 (incisal)
2172C - OC (occlusal)

Block HC Two-Layer 5 pcs. \$175.65
(12mm x 14mm x 18mm)

Z0008 - A2-2L
Z0009 - A3-2L

Z0010 - A3.5-2L

8200GC NEW! BlockButler \$199.00

CAD/CAM block storage, organizer & display. Holds 100 blocks, fits all brands*. Wall mount or store in drawers.

*With the exception of PlanMill or PrograMill blocks Trademarks belong to their manufactures

Contact for additional non-CEREC® selection

CAMERA

EYE SPECIAL C-III



- Designed for dentistry: 8 pre-set dental modes
- Ultra-lightweight ≈ 1 lb
- OSHA-compliant: water, chemical & scratch resistant
- HIPAA-compliant
- Wi-Fi capable

PRICING - EYESPECIAL C-III

PN	Product	Suggested Retail
E0001	EyeSpecial C-III	\$2,995.00
4 AA batteries, 4GB SDHC card, high performance Kenko 49mm macro lens, lens cap, computer connection cord, hand strap, flash calibration card, user manual & quick reference guide		
E0001CASE	EyeSpecial C-III Case	\$99.00
14"x 9"x 5.25" high density case with custom foam inserts		



SHOFU RAVES AND RECOMMENDATIONS

EyeSpecial C-III



EyeSpecial C-III Digital Dental Camera by Shofu Named 2014 Best New Product for Women by the **American Association of Women Dentists**. This product was judged as being the most beneficial to the practicing woman dentist and to improve her practice. Its ease of use and light weight, combined with the quality and enhanced features of other heavier intraoral digital cameras, was cited for its selection.

*Brinker, S. (2015, January). HIPAA compliance and digital photography with personal mobile devices. Dental Products Report, 76-80.

TECHNICAL DATA

RECOMMENDED SPEEDS	82 - 84
RECOMMENDED TYPES OF STERILIZATION.....	85
CERAMASTER CONVERSION CHART	85
DIAMOND CONVERSION CHART	87 - 89



TECHNICAL DATA

RECOMMENDED & MAXIMUM SPEEDS

	RECOMMENDED SPEEDS (RPM)	WATER SPRAY	MAXIMUM SPEED (RPM)
ABRASIVES			
AcryPoint	5,000 to 10,000	no	15,000
Brownie, Greenie, Supergreenie	5,000 to 7,000	yes	30,000
Brownie, Greenie, Supergreenie unmounted; Supergreenie PC2	5,000 to 7,000	yes	20,000
Brownie, Greenie, Supergreenie Floppie & Cup (FG)	5,000 to 7,000	no	15,000
CeraMaster Coarse & CeraMaster	10,000 to 15,000	no	20,000
Ceramisté	10,000 to 12,000	yes	20,000 (HP/CA) 15,000 (FG)
Ceramisté Soft	5,000 to 20,000	yes	20,000
Ceramisté Soft Unmounted	5,000 to 15,000	yes	15,000
Chipless Wheels	5,000 to 12,000	no	20,000
CompoMaster Coarse & CompoMaster	5,000 to 15,000	no	20,000 (HP/CA) 15,000 (Wheel)
CompoSite & CompoSite Fine	10,000 to 12,000	yes	15,000 (FG) 20,000 (CA)
CoreMaster Coarse & CoreMaster	10,000 to 15,000	no	20,000
CoreMaster unmounted	5,000 to 10,000	no	15,000
Dura-Green, Fine, DIA & Dura-White	5,000 to 20,000	no	30,000
Dura-Green unmounted wheels	5,000 to 15,000	no	30,000
Hardie Black & Hardie Brown	5,000 to 30,000	no	30,000
Hardie Black & Hardie Brown Floppie	5,000 to 15,000	no	15,000
Lab Series Stones	5,000 to 30,000	no	30,000
Lab Series Stones Unmounted	5,000 to 12,000	no	30,000
OneGloss	3,000 to 10,000	yes	15,000
OneGloss PS	3,000 to 5,000	yes	5,000
SuperBuff	8,000 to 12,000	moisten before using to generate paste	15,000
Super-Snap & Super-Snap X-Treme	10,000 to 12,000	no	15,000
Super-Snap Buff Disk	5,000 to 10,000	no	10,000



	RECOMMENDED SPEEDS (RPM)	WATER SPRAY	MAXIMUM SPEED (RPM)
SHOFU ROBOT FG POINT			
0897, 0893, 0897, 0899, 899V	-	yes	120,000
0822, 0826, 0885, 0886, 0890, 0892, 0894, 0895, 0898, 895C, 892X, 895X, 898X, 826V, 893V	-	yes	160,000
0790, 0792, 0793, 0800, 0803, 0810, 0814, 0818, 0820, 0824, 0825, 0830, 0831, 0833, 0834, 0835, 0836, 0837, 0838, 0839, 0840, 0844, 0846, 0848, 0850, 0852, 0858, 0866, 0870, 0876, 0878, 0880, 0881, 0883, 0888, 0891, 0896, 817C, 835C, 840C, 845C, 883C, 835X, 840X, 845X, 850X, 876X, 883X, 820F, 835F, 840F, 846F, 847F, 848F, 850F, 883F, 790V, 792V, 793V, 808V, 825V, 833V, 836V, 839V, 840V, 841V, 843V, 880V, 883V	-	yes	300,000
0789, 0791, 0801, 0805, 0807, 0812, 0813, 0816, 0832, 0842, 0855, 0856, 0859, 0860, 0861, 0864, 0865, 0868, 0872, 0874, 0875, 0879, 832X, 832F, 842F, 875F, 789V, 791V, 874V	-	yes	450,000
SHOFU ROBOT FG POINT MINIMALLY INVASIVE			
0944, 0945, 0946, 0947, 0948, 0949	-	yes	300,000
SHOFU ROBOT CARBIDE BURS			
0602, 0603, 0605, 0606, 0609, 0610, 0613, 0614, 0616, 0617, 0620, 0623, 0626, 0630, 0631, 0634, 0667	10-15,000	no	50,000
0719, 0727, 0761-3, 0762-3, 0763-3, 0770-3, 0771-3, 0772-3	-	yes	300,000
0703, 0704, 0709, 0710, 0717, 0718, 0736, 0741, 0742, 0743, 0764-3, 0765-3, 0766-3, 0767-3, 0768-3, 0769-3,	-	yes	450,000

TECHNICAL DATA

RECOMMENDED & MAXIMUM SPEEDS

	RECOMMENDED SPEEDS (RPM)	WATER SPRAY	MAXIMUM SPEED (RPM)
GINGI CURETTAGE			
C063	-	yes	450,000
C061, C062, C081, C082, C083, C084, C085	-	yes	300,000
C101, C102, C103	-	yes	160,000
T&F HYBRID POINTS			
0925, 0927, 0929, 0931, 0933, 0935, 0937, 0939	-	yes	450,000
0941, 0943	-	yes	300,000
HYBRID POINTS			
0901, 0903, 0905, 0907, 0909	-	yes	450,000
0911, 0913, 0915, 0917, 0919, 0921, 0923	-	yes	300,000
TPE DIAMONDS			
684C, 686C, 684F, 686F, 682S, 684S	-	yes	450,000
680C, 682C, 680F, 682F	-	yes	300,000



RECOMMENDED TYPES OF STERILIZATION

INSTRUMENTS	TYPES OF STERILIZATION				
	Dry Heat Oven	Steam Autoclave	Chemical Vapor	ETO (Ethylene Oxide)	Cold Sterilization
DIAMOND POINT					
Shofu Robot Points	√	√	√	√	√
Hybrid Points	√	√	√	√	√
TPE Points, GingiCurettage Points	√	√	√	√	√
CARBIDE BURS					
Shofu Carbide FG Finishers	√	√	√	√	√
Shofu Carbide CTF-FG Finishers	√	√	√	√	√
Shofu Carbide HP Cutters	√	√	√	√	√
STONES					
Dura-Green Stones, Fine & DIA	√	√	√	√	√
Dura-White Stones	√	√	√	√	√
Lab Series Stones	√	√	√	√	√
SILICONE POINTS					
Brownie, Greenie, Supergreenie	∅	√	√	√	√
Hardie, AcryPoint	∅	√	√	√	√
Ceramisté, Ceramisté Soft	∅	√	√	√	√
CeraMaster Coarse & CeraMaster	∅	√	√	√	√
CoreMaster Coarse & CoreMaster	∅	√	√	√	√
CompoSite, CompoSite Fine	∅	√	√	√	√
CompoMaster Coarse & CompoMaster	∅	√	√	√	√
OneGloss Mini Points	∅	√	√	√	√
BUR BLOCKS					
HP, CA, FG, CA/FG, Original CA/FG	√	√	∅	∅	√
BurButler	∅	√	∅	∅	√
DISPOSABLE FINISHING & POLISHING					
OneGloss PS & OneGloss	N/A	N/A	N/A	N/A	N/A
Super-Snap, Super-Snap X-Treme Disk	N/A	N/A	N/A	N/A	N/A
Super-Snap Buff Disk	N/A	N/A	N/A	N/A	N/A
Super-Snap SuperBuff Disk	N/A	N/A	N/A	N/A	N/A
Super-Snap Plastic Mandrel	N/A	N/A	N/A	N/A	N/A
LEGEND √ = RECOMMENDED ∅ = NOT RECOMMENDED N/A = DISPOSABLE Do NOT Use ultrasonic cleaner for stones, polishers & bur blocks					

TECHNICAL DATA

CERAMASTER® CONVERSION CHART

SHAPE	SHOFU CERAMASTER	KERR ROTARY	BRASSER DIALITE
BULLET			
Contra-Angle Handpiece	0121	P30032	N/A
Contra-Angle Handpiece	0126	P30044	H2D
Contra-Angle Handpiece	121C	P3032	N/A
Contra-Angle Handpiece	126C	P3044	H2DG
KNIFE			
Contra-Angle Handpiece	0122	P30036	W18DG
Contra-Angle Handpiece	N/A	N/A	N/A
Contra-Angle Handpiece	122C	P3036	W18D
Contra-Angle Handpiece	N/A	N/A	N/A
CUP			
Contra-Angle Handpiece	0123	P30035	W17D
Contra-Angle Handpiece	N/A	N/A	N/A
Contra-Angle Handpiece	123C	P3035	W17DG
Contra-Angle Handpiece	N/A	N/A	N/A
MINI-POINT			
Contra-Angle Handpiece	0124	P30033	W16D
Contra-Angle Handpiece	0129	N/A	N/A
Contra-Angle Handpiece	124C	P3033	W16DG
Contra-Angle Handpiece	129C	N/A	N/A
KN7			
Contra-Angle Handpiece	N/A	N/A	N/A
Contra-Angle Handpiece	0127	P30042	N/A
Contra-Angle Handpiece	N/A	N/A	N/A
Contra-Angle Handpiece	127C	P3042	N/A
WH6			
Contra-Angle Handpiece	N/A	N/A	N/A
Contra-Angle Handpiece	0128	P30043	R17D
Contra-Angle Handpiece	N/A	N/A	N/A
Contra-Angle Handpiece	128C	P3043	R17DG
UNMOUNTED WHEEL			
Handpiece	132C	P301	L26DBC



DIAMOND CONVERSION CHART

SHOFU	BRASSLER	KERR ROTARY	PREMIER	SS WHITE*
STRAIGHT CYLINDER				
0800-1	6836-012	835-010	514.3	835-008
0801-1	6836-014	835-012	514.4	835-012
0803-1	6837-014	836-014	514.7	837-014
0805-1	6837-016	837-014	515.7	837-016
0807-1	6837-016	837-016	515.8	N/A
808V-1	837KR-014	837KR-014	586.8	507
STRAIGHT CYLINDER (ROUND END)				
845X-1	5885-014	878-012	248.8	878K-014
845C-1	5885-014	878-012	248.8	878K-014
0846-1	6885-014	885-014	250.9	886-016
846F-1	6885-014	885-014	250.9	886-016
TAPERED CYLINDER (FLAT END)				
0810-1	N/A	845-014	722.6	846-014
0812-1	6846-012	845-016	722.8	846-016
0813-1	6846-016	846-016	724.8	846-016
0814-1	N/A	845-014	701.7	N/A
0816-1	6847-012	847-016	703.8	847-016
817C	5848-016	847-018	701.9	847-018
0818-1	6848-023	848-018	703.8	N/A
0820-1	6848-016	848-016	700.8	848-016
820F-1	6848-016	848-016	700.8	848-016
0822-1	N/A	848-018	798.10	848-018
0824-1	6848-018	848-018	727.10	N/A
0830-1	N/A	849-010	767.7	855-014
TAPERED CYLINDER (CORNER ROUND)				
0789-1	846KR-016	846KR	585.5KR	N/A
789V-1	846KR-016	846KR	585.5KR	N/A
0790-1	846KR-016	846KR	575.7KR	N/A
790V-1	846KR-016	846KR	575.7KR	N/A
0791-1	847KR-018	847KR-018	703.8KR	N/A
791V-1	847KR-018	847KR-018	703.8KR	N/A
0792-1	850KR-014	850KR-014	722.8KR	N/A
792V-1	850KR-014	850KR-014	722.8KR	N/A
0793-1	850KR-018	850KR-018	588.10KR	N/A
793V-1	850KR-018	850KR-018	588.10KR	N/A
0825-1	N/A	845KR-018	708.KR	N/A
825V-1	N/A	845KR-018	708.4KR	N/A
0826-1	N/A	845KR-025	708.4KR	N/A
826V-1	N/A	845KR-025	708.4KR	N/A
GUIDE GROOVE CUTTER				
0896-1	834-016	834-016	DC0.3	N/A
0897-1	834A-031	834A-031	DC0.5	N/A

TECHNICAL DATA

DIAMOND CONVERSION CHART

SHOFU	BRASSELER	KERR ROTARY	PREMIER	SS WHITE*
TAPERED CYLINDER (NEEDLE)				
847F-1	8852-012	858-014	207.7	858-014
0848-1	6858-014	859-014	608.9	858-014
848F-1	6858-014	859-014	608.9	858-014
0850-1	6859-018	859-016	207.10	859-016
850X-1	6859-018	859-016	207.10	859-016
850F-1	6859-018	859-016	207.10	859-016
0852-1	6859-021	859-021	209.10	N/A
TAPERED CYLINDER (ROUND END)				
0831-1	6850-014	850-014	781.8	N/A
0832-1	6855-012	855-014	767.9	856-016
832X-1	6855-012	855-014	767.9	856-016
832F-1	6855-012	855-014	767.9	856-016
0833-1	856-014	856-014	767.8/782.8	755 (1DT)
833V-1	856-014	856-014	767.8/782.8	755 (1DT)
0834-1	6850-014	850-014	777.8	854-016
0835-1	6850-016	850-018	770.9	856-018
835X-1	6850-016	850-018	770.9	856-018
835C-1	6850-016	850-018	770.9	856-018
835F-1	6850-016	850-018	770.9	856-018
0836-1	856-016	856-016	770.8/767.8	775/776
836V-1	8850-016	856-018	784.10	N/A
0837-1	6850-023	856-021	784.10	N/A
0838-1	6850-018	856-021	799.11	854-021
0839-1	856L-018	856L-018	780.9/770.8	765 (1 1/8DT)
839V-1	856L-018	856L-018	780.9/770.8	765 (1 1/8DT)
FLAME				
0840-1	6862-014	862-010	261.8	N/A
840X-1	6862-014	862-010	261.8	N/A
840C-1	6862-014	862-010	261.8	N/A
840F-1	6862-014	862-010	261.8	N/A
840V-1	6862-014	862-010	261.8	N/A
841V-1	6860-010	889-009	260.3	N/A
0842-1	6862-016	863-016	260.8	862-016
842F-1	6862-016	863-016	260.8	862-016
843V-1	8860-012	860-012	260.6.5	N/A
0844-1	6888-012	888-012	263.8	N/A
INVERTED CONE S-TYPE				
0855-1	6813-014	813-012	1210	N/A
0856-1	6813-016	813-014	1220	N/A
0858-1	6813-016	813-016	1220	N/A
0859-1	806-009	806-010	N/A	N/A
0860-1	806-016	806-012	N/A	N/A
0861-1	806-018	806-014	1250	806-016
EGG				
0893-1	379-023	379-023	N/A	N/A
893V-1	38010-024	379-023	287.4	N/A



SHOFU	BRASSLER	KERR ROTARY	PREMIER	SS WHITE*
CORN				
0898-1	6811-033	811-033	234	N/A
898X-1	6811-033	811-033	234	N/A
0899-1	N/A	N/A	N/A	N/A
899V-1	N/A	N/A	N/A	N/A
INVERTED CONE				
0864-1	805-016	805-010	310.1	805-016
0865-1	805-018	805-012	315.1.75	N/A
0866-1	805-023	805-018	324.1	N/A
0868-1	807-016	807-016	N/A	N/A
0870-1	807-018	807-018	N/A	N/A
PEAR				
0878-1	6830L-014	830L-014	362.3	N/A
0879-1	6830L-016	830L-016	364.5	N/A
0891-1	N/A	N/A	365.4	N/A
0892-1	35010-031	830-016	287.4	N/A
892X-1	35010-031	830-016	287.4	N/A
WHEEL				
0885-1	820-060	818-040	863	820-042
0886-1	818-050	818-035	N/A	N/A
WHEEL (KNIFE EDGE)				
0888-1	825-023	N/A	N/A	N/A
0890-1	825-050	N/A	N/A	N/A
WHEEL (ROUND END)				
0894-1	909-040	909-037	860	909-035
0895-1	6909-040	909-040	862	909-042
895X-1	6909-040	909-040	862	909-042
895C-1	6909-040	909-040	862	909-042
BUD				
0880-1	8368-016	368-016	292.3	N/A
880V-1	8368-016	368-016	292.3	N/A
0881-1	6368-016	368-018	283.4	868-024
0883-1	6368-023	368-023	285.5	868-024
883X-1	6368-023	368-023	285.5	868-024
883C-1	6368-023	368-023	285.5	868-024
883F-1	6368-023	368-023	285.5	868-024
883V-1	6368-023	368-023	285.5	868-024
ROUND				
0872-1	6801-014	801-014	115	801-014
0874-1	6801-016	801-016	120	N/A
874V-1	6801-016	801-016	120	N/A
0875-1	6801-018	801-018	125	801-018
875F-1	6801-018	801-018	125	801-018
0876-1	6801-023	801-021	130	801-021
876X-1	6801-023	801-023	130	801-021

*These conversions are drawn solely from SS White's multi-use Black Diamond & Revelation Diamond product lines. All trademarks belong to their manufacturers.

[illegible]

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.



SHOFU INC. 11 Kamitakamatsu-cho, Fukuino, Higashiyama-ku, Kyoto 605-0983, Japan (www.shofu.co.jp)
SHOFU INC. Singapore Branch 10 Science Park Road, #03-12 The Alpha, Science Park II, Singapore 117684 (www.shofu.com.sg)
SHOFU DENTAL TRADING (SHANGHAI) CO., LTD. No. 645 Jiye Road, Sheshan Industrial Park, Songjiang, Shanghai 201602, China (www.shofu.com.cn)
SHOFU DENTAL CORPORATION 1225 Stone Drive, San Marcos, CA 92078-4059, USA (www.shofu.com)
SHOFU UK Riverside House, River Lawn Road, Tonbridge, Kent, TN9 1EP, UK (www.shofu.co.uk)
SHOFU DENTAL GmbH   Am Brühl 17, 40878 Ratingen, Germany (www.shofu.de)

SHOFU is a registered trademark of SHOFU INC. All other trademarks and registered trademarks are the property of their respective holders. SHOFU INC. reserves the right to change specifications without notice.