


SIMPLY AESTHETIC CAD/CAM SOLUTIONS IMPROVE YOUR WORKFLOW



www.shofu.com | 800.827.4638

SIMPLY AESTHETIC

SHOFU BLOCK & DISK HC



CERAMIC-BASED RESTORATIVE

Shofu Block & Disk HC

A new hybrid ceramic restorative for CAD/CAM, HC Block/Disk features excellent physical and mechanical attributes with multiple shades and translucencies. The optical properties of HC Blocks/Disks resemble the translucency and fluorescence of natural teeth, permitting the material to seamlessly integrate with the surrounding tissue. The proprietary nano-architecture of HC Block/Disk demonstrates high flexural strength and wear resistance allowing for a fast milling and broad applications, including permanent and transitional anterior/posterior restorations for inlays/onlays, full-contour crowns, and implant-supported cases.

CAD/CAM STAIN KIT




CUSTOMIZE AESTHETICS

CAD/CAM Stain Kit

Universal and simple to use, LiteArt Stain Set for CAD/CAM & Hybrid Ceramics is design to improve the aesthetics of CAD/CAM restorations fabricated chairside. Incorporating a novel technology with the ability to mimic the natural dentition's optical properties, the assortment of stains featured in the LiteArt kit will help clinicians achieve life-like restorations, chairside. The proprietary chemical composition of the LiteArt supports a reliable adhesion to a variety of substrates, including resin-based direct and indirect materials. The low viscosity of the LiteArt stains allows for a precise application and prevents product's wasting. All stains can be layered alone or blended together to effortlessly replicate the various characteristics of natural dentition.

CAD/CAM SOLUTIONS

BLOCKBUTLER



CAD/CAM BLOCK ORGANIZER

BlockButler

Efficient storage and organization of CAD/CAM blocks is possible with the BlockButler. This 100-hole silicone bur holder can house blocks of all sizes and it can be mounted on a wall or placed in a drawer for immediate access to the blocks' inventory.

CADMASTER CAD/CAM KIT



CONTOUR, FINISH, AND POLISH

CADMaster CAD/CAM Kit

With the BurButler CADMaster Kit, clinicians will be able to contour, finish, and pre-polish CEREC restorations fabricated from modern CAD/CAM restorative materials, including Zirconia, lithium silicate and disilicate. Combining in one BurButler 25 diamond-filled abrasives in all shank-, head-, and grit-sizes, this universal abrasive system is a must-have for a CDA/CAM-focused dental practice of any size.

ALL CERAMIC PREPARATION KIT



QUICKLY PREP TEETH

All Ceramic Preparation Kit

The BurButler All Ceramic Preparation Kit features a selection of Robot Diamonds in FG shank, in 17 shapes and head sizes. These comprehensive system includes everything-abrasive a doctor will ever need to quickly prep a cavity for ceramic restorations fabricated from lithium silicate or disilicate, Zirconia, porcelain, leucite-reinforced ceramics, hybrid ceramics, CAD/CAM composites, and PFM.

IMPROVE YOUR WORKFLOW

MONOCEM & RESICEM CEMENT



STRONG SEAMLESS BOND

MonoCem

Formulated to simplify adhesive cementation to direct and indirect resin materials, MonoCem demonstrates excellent bond strength without the numerous intermediary steps. A moisture-tolerant, self-etch, self-adhesive resin cement with an unlimited setting time, MonoCem bonds to all substrates with high retention rates.

ResiCem

A dual-cure, high-strength resin cement formulated for long-term bonding to zirconia, ResiCem can be expanded using porcelain and metal primers to become a universal adhesive system for a wide range of clinical applications.

BEAUTICEM SA



BIOACTIVE SEAMLESS BOND

BeutiCem SA

A self-etch, self-adhesive dual-cure resin cement with only 12 μ m film thickness, BeutiCem SA demonstrates reliable adhesion across a variety of materials, including composites, alloys, and selected ceramic substrates. Proprietary Giomer technology incorporate into BeutiCem SA offers bioactive benefits of fluoride release/recharge, antiplaque attributes and acid neutralization.

CERARESin BOND



RELIABLE ADHESION

CeraResin Bond

Simple and intuitive to use, CeraResin Bond facilitates reliable adhesion to porcelain, ceramic, resin and acrylic materials. Designed specifically for the reparation of crowns and bridges, CeraResin Bond overcomes the existing clinical difficulties associated with the chairside repair processes.

RELIABLE & EASY

HC PRIMER



SURFACE PRE-TREATMENT

HC Primer

HC Primer provides an ideal surface pre-treatment for a reliable adhesion of modern hybrid-ceramic materials to a tooth abutment using a self-adhesive resin cement. Thanks to its innovative chemical composition, HC Primer thoroughly penetrates the matrix of hybrid ceramics effectively preparing the surface of a restoration for the subsequent cementation.

AZ PRIMER



SIMPLIFIES ADHESION

AZ Primer

Inclusive of proprietary monomers adapted for high-performance bonding, AZ Primer simplifies adhesion to Zirconia, providing the surface of a restoration with a reliable chemical retention.

ML PRIMER



ENHANCES BONDING

ML Primer

A single-bottle metal primer, ML Primer enhances the bonding of resin materials to a variety of semi-precious, precious and non-precious alloys and metals, including gold, titanium, palladium, chromium, stainless steel, and nickel.

ORDER INFORMATION

BLOCK HC - CEREC

10mm x 12mm x 16mm

SHADE	PN
ONE-LAYER 5 PCS.	
W2-LT	2151C
A1-LT	2153C
A2-LT	2155C
A3-LT	2164C
A3.5-LT	2156C
B3-LT	2165C
A1-HT	2152C
A2-HT	2154C
A3-HT	2163C
59 (Incisal)	2166C
OC (Occlusal)	2167C
TWO-LAYER 5 PCS.	
A2-2L	Z0011
A3-2L	Z0012
A3.5-2L	Z0013

BLOCK HC - CEREC

12mm x 14mm x 18mm

SHADE	PN
ONE-LAYER 5 PCS.	
W2-LT	2157C
A1-LT	2159C
A2-LT	2161C
A3-LT	2169C
A3.5-LT	2162C
B3-LT	2170C
A1-HT	2158C
A2-HT	2160C
A3-HT	2168C
59 (Incisal)	2171C
OC (Occlusal)	2172C
TWO-LAYER 5 PCS.	
A2-2L	Z0008
A3-2L	Z0009
A3.5-2L	Z0010

DISK HC - CEREC

98mm x 14mm

SHADE	PN
ONE-LAYER	
W2-LT	2180S
A1-LT	2181S
A2-LT	2182S
A3-LT	2183S
A3.5-LT	2184S
B3-LT	2185S
A1-HT	2186S
A2-HT	2187S
A3-HT	2188S
59 (Incisal)	2189S
OC (Occlusal)	2190S
TWO-LAYER 5 PCS.	
A2-2L	Z0014
A3-2L	Z0015
A3.5-2L	Z0016

*Visit www.shofu.com for the complete Shofu HC product line.

BLOCK HC - CERAMILL

12mm x 14mm x 18mm

SHADE	PN
ONE-LAYER 5 PCS.	
W2-LT	2157A
A1-LT	2159A
A2-LT	2161A
A3-LT	2169A
A3.5-LT	2162A
B3-LT	2170A
A1-HT	2158A
A2-HT	2160A
A3-HT	2168A
59 (Incisal)	2171A
OC (Occlusal)	2172A
TWO-LAYER 5 PCS.	
A2-2L	Z0014
A3-2L	Z0015
A3.5-2L	Z0016

CAD/CAM Stain Kit | PN 2150-SS

BlockButler | PN 8200GC

CADMaster CAD/CAM Kit | PN 0370BB

All Ceramic Preparation Kit | PN 0955BB

MonoCem Translucent Kit | PN 3208A

MonoCem Bleach White Kit | PN 3209A

ResiCem Basic Kit | PN 3220

BeautiCem SA Clear | PN 3213

BeautiCem SA Ivory | PN 3214

CRB - Ceraresin Bond | PN 1741

HC Primer | PN Y0005

AZ Primer | PN 3232

ML Primer | PN 1740



SHOFU INC. 11 Kamitakamatsu-cho, Fukuine, Higashiyama-ku, Kyoto 605-0983, Japan (www.shofu.co.jp)
SHOFU INC. Singapore Branch 10 Science Park Road, #03-12 The Alpha, Science Park II, Singapore 117684 (www.shofu.com.sg)
SHOFU DENTAL TRADING (SHANGHAI) CO., LTD. No. 645 Jiye Road, Sheshan Industrial Park, Songjiang, Shanghai 201602, China (www.shofu.com.cn)
SHOFU DENTAL CORPORATION 1225 Stone Drive, San Marcos, CA 92078-4059, USA (www.shofu.com)
SHOFU UK Riverside House, River Lawn Road, Tonbridge, Kent, TN9 1EP, UK (www.shofu.co.uk)
SHOFU DENTAL GmbH Am Brüll 17, 40878 Ratingen, Germany (www.shofu.de)

SN319-0518