



SAFETY DATA SHEET

page 1/5
Glaslonomer Cement CX-Plus
POWDER
Printing date: August 2, 2017

SECTION 1. Identification of the substance or mixture and of the supplier

- 1.1 Product identifier
Trade Name:
Glaslonomer Cement CX-Plus "POWDER"
- 1.2 Relevant identified uses of the substance or mixture and uses advised against
Relevant identified uses: Dental material
Uses advised against: No further data
- 1.3 Details of the supplier of the safety data sheet
Company/Undertaking identification
Manufacturer's Name: **SHOFU INC.**
Address: 11 Kamitakamatsu-cho, Fukuine, Higashiyama-ku, Kyoto 605-0983, JAPAN
Phone: +81-75-561-1112
Fax: +81-75-275-4795
Section in Charge: Quality Assurance Section
- 1.4 Emergency Telephone Number
+81-75-561-1112

SECTION 2. Hazards identification

- 2.1 GHS Classification
It is not possible to classify it.
- 2.2 Label elements
SYMBOL Void
SIGNAL WORD Void
HAZARD STATEMENTS
Void
- PRECAUTIONARY STATEMENTS
[Prevention]
Obtain special instruction before use.
Do not handle until all safety precautions have been read and understood.
[Response]
Get medical advice/attention if you feel unwell.
[Storage]
Store in a cool and dark area.
[Disposal]
Dispose of contents and container in accordance with regulation.
- 2.3 Other hazards
Results of PBT and vPvB assessment
PBT: Not applicable.
vPvB: Not applicable.

SECTION 3. Composition/information on ingredients

- 3.1 Chemical characterization: Mixtures
- 3.2 Description: Mixture of substances listed below with nonhazardous additions.



SAFETY DATA SHEET

page 2/5
Glaslonomer Cement CX-Plus
POWDER
Printing date: August 2, 2017

- 3.3 Dangerous components: Void
Fluoroaluminosilicate Glass [Cas No.65997-18-4] 90-95 %
Others
- 3.4 Additional information: For the wording of the listed risk phrases refer to section 2

SECTION 4. First-aid measures

- 4.1 Description of first aid measures
Eye contact: Rinse cautiously with water for several minutes. Remove contact lenses, if present. and easy to do. If eye irritation persists, get medical advice/attention.
Skin contact: Wash immediately with soap and plenty of water. If on skin, skin irritation, get medical advice/attention.
Ingestion: Rinse mouth and seek medical advice if necessary.
Inhalation: Remove victim to fresh air and keep at rest in a position comfortable for breathing. If symptom concerning breath goes out, call a POISON CENTER or doctor.
- 4.2 Most important symptoms and effects, both acute and delayed
No further relevant information available.
- 4.3 Indication of any immediate medical attention and special treatment needed
No further relevant information available.

SECTION 5. Fire-fighting measures

- 5.1 Extinguishing Media:
This product is not flammable.
- 5.2 Special hazards arising from the substance or mixture:
No further relevant information available.
- 5.3 Advice for firefighters:
No special measures required.

SECTION 6. Accidental release measures

- 6.1 Personal precautions, protective equipment and emergency procedures:
Avoid contact with eyes and skin.
- 6.2 Environmental Precautions:
Send to approved treatment/disposal company or dispose according to local, state and federal regulations.
- 6.3 Methods and material for containment and cleaning Up:
Wipe up and discard in a suitable container.
- 6.4 Reference to other section:
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

SECTION 7. Handling and storage

- 7.1 Precautions for safe handling:
Avoid inhaling and contact with eyes and skin.
Handle in a well ventilated area.
Handling of the product exclusively restrict to dental specialists, according to the instruction for use.



SAFETY DATA SHEET

page 3/5
Glaslonomer Cement CX-Plus
POWDER
Printing date: August 2, 2017

- 7.2 Conditions for safe storage, including any incompatibilities:
Store in a cool and dry conditions with lid tightly closed.
- 7.3 Specific end use(s):
No further relevant information available.

SECTION 8. Exposure controls/personal protection

- 8.1 Control parameters:
Ingredients with limit values that require monitoring at the workplace:
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- 8.2 Exposure controls:
Respiratory Protection:
Dust mask
- Skin Protection: Hand Protection
The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.
Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/the chemical mixture.
Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation.
- Material of gloves
The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.
 - Penetration time of glove material
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
 - For the permanent contact of a maximum of 15 minutes gloves made of the following materials are suitable:
Butyl rubber, BR
Nitrile rubber, NBR
- Eye Protection: Safety goggles

SECTION 9. Physical and chemical properties

- 9.1 Information on basic physical and chemical properties
- | | |
|------------------------------|-------------------------|
| Appearance/Odor/Colour: | White odourless powder. |
| Odour threshold | Not determined. |
| pH | Not determined. |
| Melting point/freezing point | Not determined. |
| Boiling Point: | Not determined. |
| Flash point: | Not determined. |
| Evaporation rate | Not determined. |
| Flammability (solid, gas) | Not applicable. |



SAFETY DATA SHEET

page 4/5
Glaslonomer Cement CX-Plus
POWDER
Printing date: August 2, 2017

Upper/lower flammability or explosive limits	Not determined.
Vapour pressure	Not determined.
Vapour density	Not determined.
Relative Density:	≈ 2.5 (water=1)
Solubility: water solubility	Insoluble.
Partition coefficient: n-octanol/water	Not determined.
Auto-ignition temperature	Not determined.
Decomposition temperature	Not determined.
Viscosity	Not determined.
Explosive properties	Not applicable.
Oxidising properties	Not applicable.

9.2 Other information

No further relevant information available.

SECTION 10. Stability and reactivity

10.1 Reactivity:

No further relevant information available.

10.2 Chemical stability:

Stable under normal temperatures and pressures.

10.3 Possibility of hazardous reactions:

No dangerous reactions known.

10.4 Condition to Avoid:

Avoid direct sunlight and high temperature.

10.5 Incompatible materials:

No further relevant information available.

10.6 Hazardous Decomposition Products:

None under normal conditions of storage and use.

SECTION 11. Toxicological information

11.1 Information on toxicological effects:

Acute toxicity:	Based on available data, the classification criteria are not met.
Skin corrosion/irritation:	Based on available data, the classification criteria are not met.
Eye damage/irritation:	Based on available data, the classification criteria are not met.
Sensitization to the respiratory tract:	Based on available data, the classification criteria are not met.
Skin sensitization:	Based on available data, the classification criteria are not met.
Germ cell mutagenicity/Genotoxicity:	Based on available data, the classification criteria are not met.
Carcinogenicity:	Based on available data, the classification criteria are not met.
Reproductive toxicity:	Based on available data, the classification criteria are not met.
Effects on or via lactation:	Lack of data.
Specific target organ toxicity (single exposure):	Based on available data, the classification criteria are not met.
Specific target organ toxicity (repeated exposure):	Based on available data, the classification criteria are not met.
Aspiration hazard:	Based on available data, the classification criteria are not met.



SAFETY DATA SHEET

page 5/5
Glaslonomer Cement CX-Plus
POWDER
Printing date: August 2, 2017

SECTION 12. Ecological information

- 12.1 Toxicity:
No further relevant information available.
- 12.2 Persistence and degradability:
No further relevant information available.
- 12.3 Bioaccumulative potential:
No further relevant information available.
- 12.4 Mobility in soil:
No further relevant information available.
- 12.5 Results of PBT and vPvB assessment:
Not applicable.
- 12.6 Other adverse effects:
No further relevant information available.

SECTION 13. Disposal considerations

- 13.1 Waste treatment methods:
Dispose of contents/container to in accordance with local/regional/national/international regulations.

SECTION 14. Transport information

- 14.1 UN number: Void
- 14.2 UN proper shipping name: Void
- 14.3 Transport hazard class(es): Void
- 14.4 Packing group: Void
- 14.5 Environmental hazards: No
- 14.6 Special precautions for user: Not applicable.
- 14.7 Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code:
Not applicable.

SECTION 15. Regulatory informati

Follow all regulations in your country.

SECTION 16. Other information

This product is intended for use by dental professionals. (instrument/material)

Abbreviations and acronyms:

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative