

SAFETY DATA SHEET

1. Identification of the substance or mixture and of the supplier*1.1 Product Name:***SILOCONE (POINT & WHEEL) P-type / Ceramiste (POINTS & WHEELS)***1.2 Company /Undertaking identification*Manufacturer's Name: **SHOFU INC.**

Address: 11 Kamitakamatsu-cho, Fukuine, Higashiyama-ku, Kyoto 605-0983 ,JAPAN

Phone: +81-75-561-1112

Fax: +81-75-561-2273

Section in Charge: Quality Assurance Section

Date Prepared: 2014.9.29

1.3 Emergency Telephone Number

+81-75-561-1112

2. Hazardous identification

GHS CLASSIFICATION

It is not possible to classify it.

LABEL ELEMENTS

SYMBOL	—
SIGNAL WORD	—
HAZARD STATEMENTS	—

PRECAUTIONARY STATEMENTS

[Prevention]

Do not handle until all safety precautions have been read and understood.

[Response]

If you feel un-well, get medical advice/attention.

[Storage]

Store in a dry place. Take measures to avoid mechanical and chemical damage.

[Disposal]

Dispose of contents and container in accordance with regulation.

3. Composition/information on ingredients

Working part :	Synthetic rubber
	Abrasive grain (SiC)
	Magnesium oxide (MgO)
	Others
Shank part :	Stainless steel shank

No hazardous material is present in this grinding stone.

4. First-aid measures

- 4.1 Eye contact :* Wash immediately plenty of water and seek medical advice, since the grinding dust cause scratched harm to the cornea if contact.
- 4.2 Skin contact :* N/A
- 4.3 Ingestion :* N/A
- 4.4 Inhalation:* Remove from exposure .(if inhaled the dust from grinding)

SAFETY DATA SHEET

5.Fire-fighting measures

- 5.1 Suitable Extinguishing Media:* CO₂, Dry chemical, Foam
- 5.2 Extinguishing Media Should not be used :* N/A
- 5.3 Special Exposure Hazards Arising:* In the case of fire, irritation gas and fume may emit.
- 5.4 Special Protective Equipment:* Wear fire protective cloth, safety goggles and dust mask.

6.Accidental release measures

- 6.1 Personal Precautions:* N/A
- 6.2 Methods for Cleaning Up:* Pick up and collect in a suitable container.

7.Handling and storage

- 7.1 Handling:* Use in accordance with instruction Regarding rotating speed.
- 7.2 Storage:* Store in a dry place. Take measures to avoid mechanical and chemical damage.

8.Exposure controls/personal protection

- | | Chemical name | CAS No. | ACGIH TLV |
|------------------------------------|---------------------------------------|---------|-----------|
| <i>8.1Exposure limits:</i> | N/A | | |
| <i>8.2 Respiratory Protection:</i> | Dust respirator and/or local exhaust. | | |
| <i>8.3 Hand Protection:</i> | N/A | | |
| <i>8.4 Eye Protection:</i> | Safety goggles | | |
| <i>8.5 Skin Protection:</i> | General protective clothing | | |

9.Physical and chemical properties

- 9.1 Appearance/Odor :* Dark gray rubber bonded grinding and polishing stone (points/wheels).
- 9.2 Boiling Point/Boiling Range:* N/D
- 9.3 Melting Point/Melting Range:* N/D
- 9.4 Flash point:* N/D
- 9.5 Auto flammability:* N/D
- 9.6 Explosive Properties:* N/D
- 9.7 Oxidizing Properties:* N/D
- 9.8 Vapour Pressure:* N/D
- 9.9 Relative Density:* N/D
- 9.10 Solubility: water solubility* Non-soluble

10.Stability and reactivity

- Stable in a normal temperature.
- 10.1 Condition to Avoid:* Avoid direct sunlight and high temperature.
- 10.2 Materials to Avoid:* N/A
- 10.3 Hazardous Decomposition Products:* Will not occur

11.Toxicological information

No information

SAFETY DATA SHEET**12.Ecological information**

No information

13.Disposal considerations

Dispose of contents/container to in accordance with local/regional/national/international regulations (to be specified).

14.Transport information

Airline restriction information ; It doesn't correspond.
UN No. ; It doesn't correspond.
Sea restriction information ; It doesn't correspond.
UN No. ; It doesn't correspond.

15.Regulatory information

Follow all regulations in your country.

16.Other information

This product is intended for use by dental professionals. (instrument/material)