

SAFETY DATA SHEET

1. Identification of the substance or mixture and of the supplier*1.1 Product Name:***ACRY POINT / BIG SILICONE POINT***1.2 Company /Undertaking identification*Manufacturer's Name: **SHOFU INC.**

Address: 11 Kamitakamatsu-cho, Fukuine, Higashiyama-ku, Kyoto 605-0983 ,JAPAN

Phone: +81-75-561-1112

Fax: +81-75-561-2273

Section in Charge: Quality Assurance Section

Date Prepared: 2014.9.29

1.3 Emergency Telephone Number

+81-75-561-1112

2. Hazardous identification**GHS CLASSIFICATION**

It is not possible to classify it.

LABEL ELEMENTS

SYMBOL	—
SIGNAL WORD	—
HAZARD STATEMENTS	—

PRECAUTIONARY STATEMENTS

[Prevention]

Do not handle until all safety precautions have been read and understood.

[Response]

If you feel un-well, get medical advice/attention.

[Storage]

Store in a dry place. Take measures to avoid mechanical and chemical damage.

[Disposal]

Dispose of contents and container in accordance with regulation.

3. Composition/information on ingredients

Working part :	Synthetic rubber
	Abrasive grain (SiC)
	Silicon oxide (SiO ₂)
	Others
Shank part :	Stainless steel shank

No hazardous material is present in this grinding stone.

4. First-aid measures

- 4.1 Eye contact :* Wash immediately plenty of water and seek medical advice, since the grinding dust cause scratched harm to the cornea if contact.
- 4.2 Skin contact :* N/A
- 4.3 Ingestion :* N/A
- 4.4 Inhalation:* Remove from exposure .(if inhaled the dust from grinding)

SAFETY DATA SHEET

5.Fire-fighting measures

- | | |
|---|---|
| <i>5.1 Suitable Extinguishing Media:</i> | CO ₂ , Dry chemical, Foam |
| <i>5.2 Extinguishing Media Should not be used :</i> | N/A |
| <i>5.3 Special Exposure Hazards Arising:</i> | In the case of fire, irritation gas and fume may emit. |
| <i>5.4 Special Protective Equipment:</i> | Wear fire protective cloth, safety goggles and dust mask. |

6.Accidental release measures

- | | |
|-------------------------------------|--|
| <i>6.1 Personal Precautions:</i> | N/A |
| <i>6.2 Methods for Cleaning Up:</i> | Pick up and collect in a suitable container. |

7.Handling and storage

- | | |
|----------------------|--|
| <i>7.1 Handling:</i> | Use in accordance with instruction Regarding rotating speed. |
| <i>7.2 Storage:</i> | Store in a dry place. Take measures to avoid mechanical and chemical damage. |

8.Exposure controls/personal protection

- | <i>8.1Exposure limits:</i> | Chemical name | CAS No. | ACGIH TLV |
|------------------------------------|---------------------------------------|---------|-----------|
| | N/A | | |
| <i>8.2 Respiratory Protection:</i> | Dust respirator and/or local exhaust. | | |
| <i>8.3 Hand Protection:</i> | N/A | | |
| <i>8.4 Eye Protection:</i> | Safety goggles | | |
| <i>8.5 Skin Protection:</i> | General protective clothing | | |

9.Physical and chemical properties

- | | |
|---|--|
| <i>9.1 Appearance/Odor :</i> | Gray, brown and white colored silicone rubber bonded grinding and polishing stone(points). |
| <i>9.2 Boiling Point/Boiling Range:</i> | N/D |
| <i>9.3 Melting Point/Melting Range:</i> | N/D |
| <i>9.4 Flash point:</i> | N/D |
| <i>9.5 Auto flammability:</i> | N/D |
| <i>9.6 Explosive Properties:</i> | N/D |
| <i>9.7 Oxidizing Properties:</i> | N/D |
| <i>9.8 Vapour Pressure:</i> | N/D |
| <i>9.9 Relative Density:</i> | N/D |
| <i>9.10 Solubility:</i> <i>water solubility</i> | Non-soluble |

10.Stability and reactivity

- | | |
|---|---|
| | Stable in a normal temperature. |
| <i>10.1 Condition to Avoid:</i> | Avoid direct sunlight and high temperature. |
| <i>10.2 Materials to Avoid:</i> | N/A |
| <i>10.3 Hazardous Decomposition Products:</i> | Will not occur |

11.Toxicological information

No information

SAFETY DATA SHEET

12.Ecological information

No information

13.Disposal considerations

Dispose of contents/container to in accordance with local/regional/national/international regulations (to be specified).

14.Transport information

Airline restriction information ; It doesn't correspond.

UN No. ; It doesn't correspond.

Sea restriction information ; It doesn't correspond.

UN No. ; It doesn't correspond.

15.Regulatory information

Follow all regulations in your country.

16.Other information

This product is intended for use by dental professionals. (instrument/material)