

Super-Snap[®]

FINISHING AND POLISHING



Super-Snap & SuperBuff



Super-Snap Disks
flexible and ultra-thin



Super-Snap Polystrips
great for interproximal areas



SuperBuff Disks
final polish on all composites

EASY PROCEDURES

Super-Snap is designed for easy, fast and safe contouring, finishing and polishing of all microfilled and hybrid composites. Super-Snap Polystrips are color coded with matching grit size for easy identification and offer the same results as the Super-Snap disks for interproximal areas. Use SuperBuff disks for the final polish on all composites. They are paste-impregnated disks that, when exposed to water, release a high quality polishing paste.

Super-Snap

- Elastic silicone shank mount for easy placement of the disk onto the mandrel
- Manufactured without a metal center
- Very flexible, ultra-thin disks allow easy access to interproximal areas
- Two sizes – four grits (silicon carbide and aluminum oxide)

SuperBuff

- Paste-impregnated with aluminum oxide
- Disposable
- Eliminates the potential mess associated with polishing paste syringe systems

Visit www.shofu.com or call 800.827.4638

SHAPE CATEGORY		Contouring		Finishing				Polishing		Super-Polishing											
Shape																					
Color/Grit		Black (coarse)		Violet (medium)				Green (fine)		Red (super fine)											
Refills (50 pk)	PN	L506	L507	L508	L528	L509	L522	L501	L519	L502	L521										

Order Information

Super-Snap Kits

- Part Number 0500 – Rainbow Technique Kit 100 standard disks, 80 mini disks, 40 polystrips, 2 Dura-White Stones (CN1 and FL2), CompoSite fine midi-point, 4 CA metal mandrels
- Part Number 0505 – Mini Kit 32 standard disks, 16 mini disks, CompoSite fine mini-point, CA metal mandrel

Super-Snap Singles – Individual Patient Use

- Part Number 0501 – Singles 2 disks of coarse, medium, fine and super fine grits, Super-Snap plastic CA mandrel

- Part Number 0502 – Standard 1 disk of coarse, medium, fine, super fine grits, CA mandrel
- Part Number 0503 – Mini 1 disk of coarse, medium, fine, super fine grits, CA mandrel
- See chart above for individual shapes

Super-Snap Buff and Super-Snap SuperBuff Kits


- Part Number L523 – Super-Snap Buff Disk 25 disks
- Part Number L524 – Super-Snap Buff Mini Disk 25 disks
- Part Number 0535 – SuperBuff Set 25 disks and 2 CA mandrels
- Part Number 0536 – SuperBuff Disk 50 disks (refills)

Super-Snap Polystrips

- Part Number L525 – Polystrips 100 black (coarse) and violet (medium)
- Part Number L526 – Polystrips 100 green (fine) and red (super fine)

Mandrels

- Part Number 0439 – CA Metal
- Part Number 0440 – CA Plastic



Shofu Dental Corporation
1225 Stone Drive
San Marcos, CA
92078-4059, USA

SN22-0814

Finishing and Polishing