

## SAFETY DATA SHEET

**1. Identification of the substance or mixture and of the supplier***1.1 Product Name:***Composite Polishing Paste***1.2 Company /Undertaking identification*Manufacturer's Name: **SHOFU INC.**

Address: 11 Kamitakamatsu-cho, Fukunine, Higashiyama-ku, Kyoto 605-0983 ,JAPAN

Phone: +81-75-561-1112

Fax: +81-75-561-2273

Section in Charge: Quality Assurance Section

Date Prepared: 2014.9.29

*1.3 Emergency Telephone Number*

+81-75-561-1112

**2. Hazardous identification**

## GHS CLASSIFICATION

HEALTH HAZARDS

ACUTE TOXICITY-ORAL  
EYE DAMAGE/IRRITATIONNot Classified  
Category 2B

The thing without mention is out of a classification object. Or cannot classify it.

## LABEL ELEMENTS

SYMBOL

—

SIGNAL WORD

Warning

HAZARD STATEMENTS

Causes eye irritation

## PRECAUTIONARY STATEMENTS

[Prevention]

Do not handle until all safety precautions have been read and understood.  
Wash hands thoroughly after handling.

[Response]

If in eyes, rinse cautiously with water for several minutes.  
Remove contact lenses, if present and easy to do.  
If eye irritation persists, get medical advice/attention.

[Storage]

Store away from sunlight and flame.

[Disposal]

Dispose of contents and container in accordance with regulation.

**3. Composition/information on ingredients***Composition**%*

Alumina

≒ 50

Glycerin

≒ 50

Benzalkonium Chloride

Peppermint oil

**4. First-aid measures***4.1 Eye contact :* If in eyes, rinse cautiously with water for several minutes.  
Remove contact lenses, if present and easy to do.

## SAFETY DATA SHEET

- If eye irritation persists, get medical advice/attention.
- 4.2 Skin contact :** If on skin, wash with plenty of soap and water. If on skin, skin irritation, get medical advice/attention.
- 4.3 Ingestion :** Rinse mouth. If you feel un-well, get medical advice/attention.
- 4.4 Inhalation:** If you feel un-well, get medical advice/attention.

**5.Fire-fighting measures**

- 5.1 Suitable Extinguishing Media:** CO2, dry chemical, foam
- 5.2 Extinguishing Media Should not be used :** N/A
- 5.3 Special Exposure Hazards Arising:** Irritating gases and fumes may emit.
- 5.4 Special Protective Equipment:** N/A

**6.Accidental release measures**

- 6.1 Personal Precautions:** N/A
- 6.2 Methods for Cleaning Up:** Wipe up and discard in a suitable container.

**7.Handling and storage**

- 7.1 Handling:** Keep away from naked flames.
- 7.2 Storage:** Keep away from direct sunlight or heat.

**8.Exposure controls/personal protection**

- | 8.1Exposure limits:                | Chemical name               | CAS No. | ACGIH TLV |
|------------------------------------|-----------------------------|---------|-----------|
|                                    | N/A                         |         |           |
| <b>8.2 Respiratory Protection:</b> | Not required.               |         |           |
| <b>8.3 Hand Protection:</b>        | Latex or rubber gloves      |         |           |
| <b>8.4 Eye Protection:</b>         | Safety goggles              |         |           |
| <b>8.5 Skin Protection:</b>        | General protective clothing |         |           |

**9.Physical and chemical properties**

- 9.1 Appearance/Odor :** White color rough paste
- 9.2 Boiling Point/Boiling Range:** N/D
- 9.3 Melting Point/Melting Range:** N/D
- 9.4 Flash point:** N/D
- 9.5 Auto flammability:** N/D
- 9.6 Explosive Properties:** N/D
- 9.7 Oxidizing Properties:** N/D
- 9.8 Vapour Pressure:** N/D
- 9.9 Relative Density:** N/D
- 9.10 Solubility:** water solubility Slightly soluble

**10.Stability and reactivity**

- Stable in a normal temperature.
- 10.1 Condition to Avoid:** Avoid exposure to direct sun light on the container.
- 10.2 Materials to Avoid:** Strong oxidizing materials
- 10.3 Hazardous Decomposition Products:** Will not occur

## SAFETY DATA SHEET

**11.Toxicological information**

No information

**12.Ecological information**

No information

**13.Disposal considerations**

Dispose of contents/container to in accordance with local/regional/national/international regulations (to be specified).

**14.Transport information**

Airline restriction information ; It doesn't correspond.  
UN No. ; It doesn't correspond.  
Sea restriction information ; It doesn't correspond.  
UN No. ; It doesn't correspond.

**15.Regulatory information**

Follow all regulations in your country.

**16.Other information**

This product is intended for use by dental professionals. (instrument/material)