



# BlockButler™

CAD/CAM Block  
Storage, Organizer, Display

Visit [www.shofu.com](http://www.shofu.com) | Call **800.827.4638**

PN 8200GC | LOT:

Distributed by: Shofu Dental Corporation, San Marcos CA, USA



## **BLOCKBUTLER™**

The only patented CAD/CAM dedicated block storage, organization and display system.

Effortlessly select desired blocks and quickly assess how many blocks are remaining or need to be replenished.

BlockButler may be wall mounted or stored in standard drawers.

## **MATERIAL:**

Medical grade silicone

## **CLEANING:**

Clean BlockButler with warm soapy water. Surface disinfectants may also be used.

## **ORGANIZATION:**


Organize the blocks by preference. For instance, by manufacturer, material, shade, translucency, and/or size. Each row or column will hold 10 blocks. Customize by labeling the rows/columns with a permanent marker or adhesive labels.

## **GRAY COLOR:**

The 18% neutral gray color gives the ideal neurological visual perception contrast for indicating the correct color CAD/CAM block for the patient.

## **WALL MOUNTED:**

There are mounting holes on the back side. Place two round head screws in a wall or a backboard.



## **DRAWER STORAGE:**

The BlockButler can be stored flat in a drawer or commercial toolbox.

## **DIMENSIONS:**

260mm x 240mm x 10mm