

SAFETY DATA SHEET

1. Identification of the substance or mixture and of the supplier*1.1 Product Name:*

**Dura-White Stones /WHITE POINTS /Lab Series Stones "white" (points & wheels)
/Ceramo-Metal (points & wheels)**

1.2 Company /Undertaking identification

Manufacturer's Name: **SHOFU INC.**

Address: 11 Kamitakamatsu-cho, Fukuine, Higashiyama-ku, Kyoto 605-0983 ,JAPAN

Phone: +81-75-561-1112

Fax: +81-75-561-2273

Section in Charge: Quality Assurance Section

Date Prepared: 2014.9.22

1.3 Emergency Telephone Number

+81-75-561-1112

2. Hazardous identification

GHS CLASSIFICATION

It is not possible to classify it.

LABEL ELEMENTS

SYMBOL	—
SIGNAL WORD	—
HAZARD STATEMENTS	—

PRECAUTIONARY STATEMENTS

[Prevention]

Do not handle until all safety precautions have been read and understood.

[Response]

If you feel un-well, get medical advice/attention.

[Storage]

Store in a dry place. Take measures to avoid mechanical and chemical damage.

[Disposal]

Dispose of contents and container in accordance with regulation.

3. Composition/information on ingredients

Working part :	Alumina grit(Al_2O_3)
	Others
Shank part :	Stainless steel shank

No hazardous material is present in this grinding stone.

4. First-aid measures

<i>4.1 Eye contact :</i>	Wash immediately plenty of water and seek medical advice, since the grinding dust cause scratched harm to the cornea if contact.
<i>4.2 Skin contact :</i>	N/A
<i>4.3 Ingestion :</i>	N/A
<i>4.4 Inhalation:</i>	Remove from exposure .(if inhaled the dust from grinding)

SAFETY DATA SHEET

5.Fire-fighting measures

This product is non-flammable.

<i>5.1 Suitable Extinguishing Media:</i>	N/A
<i>5.2 Extinguishing Media Should not be used :</i>	N/A
<i>5.3 Special Exposure Hazards Arising:</i>	N/A
<i>5.4 Special Protective Equipment:</i>	N/A

6.Accidental release measures

<i>6.1 Personal Precautions:</i>	N/A
<i>6.2 Methods for Cleaning Up:</i>	Pick up and collect in a suitable container.

7.Handling and storage

<i>7.1 Handling:</i>	Use in accordance with instruction Regarding rotating speed.
<i>7.2 Storage:</i>	Store in a dry place. Take measures to avoid mechanical and chemical damage.

8.Exposure controls/personal protection

<i>8.1Exposure limits:</i>	Chemical name	CAS No.	ACGIH TLV
	N/A		
<i>8.2 Respiratory Protection:</i>	Dust respirator and/or local exhaust.		
<i>8.3 Hand Protection:</i>	N/A		
<i>8.4 Eye Protection:</i>	Safety goggles		
<i>8.5 Skin Protection:</i>	General protective clothing		

9.Physical and chemical properties

<i>9.1 Appearance/Odor :</i>	White odorless Grinding stone(abrasive)
<i>9.2 Boiling Point/Boiling Range:</i>	N/D
<i>9.3 Melting Point/Melting Range:</i>	N/D
<i>9.4 Flash point:</i>	N/D
<i>9.5 Auto flammability:</i>	N/D
<i>9.6 Explosive Properties:</i>	N/D
<i>9.7 Oxidizing Properties:</i>	N/D
<i>9.8 Vapour Pressure:</i>	N/D
<i>9.9 Relative Density:</i>	N/D
<i>9.10 Solubility: water solubility</i>	Non-soluble

10.Stability and reactivity

	Stable in a normal temperature.
<i>10.1 Condition to Avoid:</i>	N/A
<i>10.2 Materials to Avoid:</i>	N/A
<i>10.3 Hazardous Decomposition Products:</i>	Will not occur

11.Toxicological information

No information

12.Ecological information

No information

SAFETY DATA SHEET

13. Disposal considerations

Dispose of contents/container to in accordance with local/regional/national/international regulations (to be specified).

14. Transport information

Airline restriction information ; It doesn't correspond.

UN No. ; It doesn't correspond.

Sea restriction information ; It doesn't correspond.

UN No. ; It doesn't correspond.

15. Regulatory information

Follow all regulations in your country.

16. Other information

This product is intended for use by dental professionals. (instrument/material)