

## SAFETY DATA SHEET

**1. Identification of the substance or mixture and of the supplier***1.1 Product Name:***Pink Points(points & wheels)****/ Dura-Pink Stones (points & wheels) / Lab Series Stones "pink" (points & wheels)***1.2 Company /Undertaking identification*Manufacturer's Name: **SHOFU INC.**

Address: 11 Kamitakamatsu-cho, Fukuine, Higashiyama-ku, Kyoto 605-0983 ,JAPAN

Phone: +81-75-561-1112

Fax: +81-75-561-2273

Section in Charge: Quality Assurance Section

Date Prepared: 2014.9.30

*1.3 Emergency Telephone Number*

+81-75-561-1112

**2. Hazardous identification**

## GHS CLASSIFICATION

It is not possible to classify it.

## LABEL ELEMENTS

SYMBOL —

SIGNAL WORD —

HAZARD STATEMENTS —

## PRECAUTIONARY STATEMENTS

[Prevention]

Do not handle until all safety precautions have been read and understood.

[Response]

If you feel un-well, get medical advice/attention.

[Storage]

Store in a dry place. Take measures to avoid mechanical and chemical damage.

[Disposal]

Dispose of contents and container in accordance with regulation.

**3. Composition/information on ingredients**

Working part :	Alumina grit( $\text{Al}_2\text{O}_3$ )
	Bonding glass phase
	Others
Shank part :	Stainless steel shank

No hazardous material is present in this grinding stone.

**4. First-aid measures**

*4.1 Eye contact :* Wash immediately plenty of water and seek medical advice, since the grinding dust cause scratched harm to the cornea if contact.

*4.2 Skin contact :* N/A

*4.3 Ingestion :* N/A

*4.4 Inhalation:* Remove from exposure .(if inhaled the dust from grinding)

## SAFETY DATA SHEET

**5.Fire-fighting measures**

**This product is non-flammable.**

- |   |     |
|---|-----|
| <i>5.1 Suitable Extinguishing Media:</i>            | N/A |
| <i>5.2 Extinguishing Media Should not be used :</i> | N/A |
| <i>5.3 Special Exposure Hazards Arising:</i>        | N/A |
| <i>5.4 Special Protective Equipment:</i>            | N/A |

**6.Accidental release measures**

- |                                     |  |
|-------------------------------------|--|
| <i>6.1 Personal Precautions:</i>    | N/A  |
| <i>6.2 Methods for Cleaning Up:</i> | Pick up and collect in a suitable container. |

**7.Handling and storage**

- |                      |  |
|----------------------|--|
| <i>7.1 Handling:</i> | Use in accordance with instruction Regarding rotating speed.                 |
| <i>7.2 Storage:</i>  | Store in a dry place. Take measures to avoid mechanical and chemical damage. |

**8.Exposure controls/personal protection**

- | <i>8.1Exposure limits:</i>         | Chemical name                         | CAS No. | ACGIH TLV |
|------------------------------------|---------------------------------------|---------|-----------|
|                                    | N/A                                   |         |           |
| <i>8.2 Respiratory Protection:</i> | Dust respirator and/or local exhaust. |         |           |
| <i>8.3 Hand Protection:</i>        | N/A                                   |         |           |
| <i>8.4 Eye Protection:</i>         | Safety goggles                        |         |           |
| <i>8.5 Skin Protection:</i>        | General protective clothing           |         |           |

**9.Physical and chemical properties**

- |  |  |
|--|--|
| <i>9.1 Appearance/Odor :</i>             | Pink odorless Grinding stone(abrasive) |
| <i>9.2 Boiling Point/Boiling Range:</i>  | N/D                                    |
| <i>9.3 Melting Point/Melting Range:</i>  | N/D                                    |
| <i>9.4 Flash point:</i>                  | N/D                                    |
| <i>9.5 Auto flammability:</i>            | N/D                                    |
| <i>9.6 Explosive Properties:</i>         | N/D                                    |
| <i>9.7 Oxidizing Properties:</i>         | N/D                                    |
| <i>9.8 Vapour Pressure:</i>              | N/D                                    |
| <i>9.9 Relative Density:</i>             | N/D                                    |
| <i>9.10 Solubility: water solubility</i> | Non-soluble                            |

**10.Stability and reactivity**

- |   |                                 |
|---|---------------------------------|
|   | Stable in a normal temperature. |
| <i>10.1 Condition to Avoid:</i>               | N/A                             |
| <i>10.2 Materials to Avoid:</i>               | N/A                             |
| <i>10.3 Hazardous Decomposition Products:</i> | Will not occur                  |

**11.Toxicological information**

No information

**12.Ecological information**

No information

# SAFETY DATA SHEET

---

**13.Disposal considerations**

---

Dispose of contents/container to in accordance with local/regional/national/international regulations (to be specified).

---

**14.Transport information**

---

Airline restriction information ; It doesn't correspond.

UN No. ; It doesn't correspond.

Sea restriction information ; It doesn't correspond.

UN No. ; It doesn't correspond.

---

**15.Regulatory information**

---

Follow all regulations in your country.

---

**16.Other information**

---

This product is intended for use by dental professionals. (instrument/material)