



**Works Better.
Lasts Longer.
Costs Less.**



IMPRESSION



- Work 1:00 minutes
- Intra Oral set 1:00 minute
- Total work/set time 2:00 minutes
- Durometer 60-65 (Shore A)
- Color: Clear

“Works Better”

Offers greater detail while allowing light to pass through for curing of dual cure or light-cured temporary materials.

“Lasts Longer”

Great visibility for voids or bubbles makes ClearBite an excellent choice as a matrix for quick and easy self-cure temporary materials.

“Costs Less”

Clear Bite is fast setting, dimensionally stable and flexible saving you time and money.

Defend ClearBite

Great visibility of any voids or bubbles makes ClearBite the perfect choice as a matrix for quick and easy self cure temporary materials. The clear material is also excellent to "cure through" when used with dual cure or light-cured temporary materials. The fast setting, unflavored, medium viscosity bite registration material is also dimensionally stable and flexible.

Packaging:	Item No:	Compares to*:
2 cartridges (50 ml. ea.) 6 pink mixing tips	#BR-9200	Clear Bite™ BLU-BITE® CLEAR

*Legal disclaimers: Clear Bite is a registered trademark of DenMat®, BLU-BITE® CLEAR is a registered trademark of Henry Schein®.

Use of these trademarks does not imply any affiliation with or endorsement by these companies. Mydent International is not responsible for any innacuracy or typographical errors in this document.

Request FREE Samples at: www.defend.com/freesamples

Mydent
International

80 Suffolk Court, Hauppauge, NY 11788
800-275-0020 • FAX 631-434-7750
www.defend.com • EMAIL: sales@defend.com



Follow us on Facebook.com/MydentDefend
Follow us on Twitter.com/Mydent_Defend