



TheSharpChoice®

NMSDC Certified Diversity Supplier

DATASHEET

## RELI® Blunt Fill Needles

Reli® Blunt Fill Needles eliminate the use of a sharp when reconstituting medication, while maintaining compatibility with most syringes available in the market.

**Note: Reli® Blunt Fill Needles are not for skin injection.**

- Angled bevel reduces coring of vial stopper
- 18G ensures rapid filling and smooth flow rate suitable for viscous medications
- 1.5" needle aids penetration of different septum barriers
- Eliminates the use of a sharp, reducing the risk of a needle-stick injury
- Compatible with most syringes available in the market
- Sterile, single-use, PVC-free, not made with natural rubber latex

BRAND (COMPANY NAME)	Angled Bevel	18G	1.5"	Red Shield
Reli® Blunt Fill Needle (MYCO Medical)	✓	✓	✓	✓
BD™ Blunt Fill Needle (BD)	✓	✓	✓	✓
JELCO® Blunt Fill Needle (Smiths Medical)	✓	✓	✓	

## Reli® Blunt Fill Needles (Sterile)

MYCO CODE	GAUGE	LENGTH	UNIT OF MEASURE
BFN18G101	18G	1"	100 ea/box
BFN18G151	18G	1.5"	100 ea/box
BFN20G101	20G	1"	100 ea/box

## Reli® Blunt Fill Needles with Filter (Sterile)

MYCO CODE	GAUGE	LENGTH	UNIT OF MEASURE
F-BFN18G151	18G	1.5"	100 ea/box



Angled bevel



Reli® Blunt Fill Needles  
with Filter (Sterile)



Reli® Blunt Fill Needles (Sterile)

As a leading NMSDC Diversity Certified Supplier, MYCO Medical is committed to “Bringing Reliability, Innovation and Value” to the Healthcare industry with an emphasis on “physician preference” surgical devices and disposables.

### To order, learn more or receive complimentary evaluation samples:

Please call MYCO Medical at (800) 454-6926 or (919) 460-2535 or email us at [info@mycomedical.com](mailto:info@mycomedical.com)

See the complete line of MYCO innovative disposable and reusable medical devices at  
[www.mycomedical.com](http://www.mycomedical.com)



2015 Production Drive  
Apex, NC 27539 USA

[info@mycomedical.com](mailto:info@mycomedical.com)  
[www.mycomedical.com](http://www.mycomedical.com)

Some medical devices are not licensed for sale in Canada.

© 2017 MYCO Medical, Inc. All rights reserved.  
LIT-00325 Rev 2 Eff 12.14.2017