



**ADU-20B**  
**ADU-22**

**Mobile  
Delivery  
Systems**



**ADVANCING DENTAL TECHNOLOGY**

# Mobile Dental Delivery Systems

- Dual handpiece control
- Water on/off and flow adjustment for highspeed handpiece
- 3-way air/water syringe
- Disc-type foot control
- High volume vacuum
- Durable casters for easy mobility between rooms
- Compatible with Midwest 4-hole or Borden 2/3-hole pneumatic handpieces. Specify connector type and tubing style (*coiled or straight*) when ordering
- Includes Ritter size stainless steel tray: 13.75" x 9.875" (34 x 25 cm)
- Quick connect to air source — no junction box required!
- Manufactured in the USA

## Premium

### ADU-20B *(pictured on opposite side)*

- Improved chassis design
- Automatic dual handpiece selection with individual air and water adjustments
- 2 Liter water reservoir
- Wet-dry foot control
- Saliva ejector suction
- Handpiece air pressure gauge
- Height adjustable: 23" to 32" (60 to 81 cm)

## Options

### Third Handpiece Connection

Water Connection for Scaler . . . . . AWS-1C

### Handpieces

4-hole highspeed . . . . . AHP-96  
4-hole highspeed . . . . . AHP-55  
4-hole highspeed E-style adapter AHP-50  
1:1 straight w/ Doriot post . . . . AHP-101  
Doriot style - Latch type burs . . . AHP-53

## Economical

### ADU-22

- 1 Liter water reservoir
- Height adjustable: 23" to 38" (59 to 97 cm)

## Roll From Room to Room



## Fiber Optic Upgrade

- Compatible with any standard fiber optic handpiece
  - Easy to replace halogen lamp provides brilliant, cool light for intra-oral use.
  - Halogen lamp, 3.5 VAC, 720 mA
  - ISO standard fiber optic handpiece terminal with stainless steel adapter and nut (*Midwest 5-hole*)
  - Air electric switch activated
- 5-hole . . . . . AA-19A-04  
6-pin . . . . . AA-19A-06