

Canada's #1 manufacturer of disposable dental and medical products



ABOUT MEDICOM



Corporate Headquarters in Montreal

Medicom[®] manufactures and distributes high quality infection control and patient care products worldwide. It is present on five continents through sales, logistics and manufacturing sites.

Founded by Ronald Reuben in 1988, Medicom[®] is privately owned and its product line spans dental, veterinary, industrial and medical single-use supplies, as well as a comprehensive line of oral health care products.

Leading Medicom[®] brand names include Denti-Care[®], Dry-Back[®], Duraflor[®] Halo, Safe+Mask[®], Safe Gauze[®] and SafeTouch[®] as well as AMD Ritmed medical products.

In a continued effort to "go green", as of March 2012 Medicom will be making its product catalog available to you in a paperless format only. We hope you share our enthusiasm for protecting the environment. Catalogs will be available for viewing and download on www.medicom.ca.



TABLE OF CONTENTS

SECTIONS	PAC	GES
	MASKS	. 4
	NON-WOVEN & COTTON PRODUCTS	. 8
	STERILIZATION PRODUCTS	12
	DISINFECTANTS	14
	PREVENTIVES	16
P	PAPER & PLASTIC PRODUCTS	22
	APPAREL	26
	OTHER	28
	MASKS ADDITIONAL INFORMATION	31





MEDICOM[®] Safe+Mask[®]

Different dental procedures involve different levels of spray, spatter and airborne particles. To ensure adequate protection and peace of mind, it is important to select the ideal face mask with proper bacterial and particle filtration and fluid resistance. While face masks may visibly look the same, their level of protection can differ significantly. To help distinguish the right mask for the right procedure, Medicom[®] Safe+Mask[®] has adopted ASTM 2100 standard to offer a complete range of procedural face masks for any situation.

Now you or your patient can finally choose the ideal face mask, easy as 1-2-3.





High Barrier mask ideal for oral exams of high-risk patients and procedures involving high levels of spray, spatter, moisture and airborne particles. Perfect for ultrasonic scalers, high-speed aerosols, air polishers and oral surgeries.

Safe+Mask

a * * • Premice

- BFE and PFE filtration efficiency > 98% @ 0.1micron assures a high level of protection against potential sources of infection.
- Shingle pleat design (will not lint, tear or shred), prevents the pooling of fluids and cross contamination, while dye free non-woven inner layer provides additional protection.
- Soft and strong, ample-length earloops so there is less pulling, and adjustable nosepiece for comfortable fit.

Item No.	Description	Unit/Pkg	Pkg/Case
2042	SafeMask Premier Elite - Blue	50	10
2046	SafeMask Premier Elite - Pink	50	10





High Barrier mask for oral exams of highrisk patients with an integrated optically clear, distortion free, wrap-around shield for full-face protection.

- Eve shield has anti-fog treatment to eliminate fogging and promote clear visibility. Mask and shield feature soft and strong earloops for custom fit and comfort.
- BFE and PFE filtration efficiency > 98% @ 0.1micron assures a high level of protection against potential sources of infection.
- Shingle pleat design (will not lint, tear or shred), prevents the pooling of fluids and cross contamination, while dye free non-woven inner layer provides additional protection.

Item No. Description Unit/Pkg Pkg/Case 2025 SafeMask ProShield 25 Δ W/Visor-Teal



Earloop Mask



Moderate Barrier mask for oral exams of known higher-risk patients and procedures involving moderate levels of spray, spatter and airborne particles such as endo or periodontic procedures, bleaching and lab work with exposure to blood or contaminated materials.

- Lower differential pressure offers excellent breathability and comfort with high levels of protection.
- Bacterial and particle filtration efficiency of 98% and more assures greater protection against potential sources of infection.
- Soft and strong, ample-length earloops so there is less pulling, and adjustable nosepiece for comfortable fit

Item No.	Description	Unit/Pkg	Pkg/Case
2040	SafeMask Premier Plus - Blue	50	10

Premier Earloop Mask



Low Barrier mask for oral exams of low-risk patients and procedures involving lowest levels of spray and spatter such as fluoride treatments, minor orthodontics and basic lab work.

- A daily favorite of most dental professionals designed for everyday protection and comfort.
- Non-woven layers and shingle pleat design prevent fluid pooling and cross contamination. Fluid resistant to help reduce potential exposure of the wearer to blood and body fluids.
- Soft and strong, ample-length earloops so there is less pulling, and adjustable nosepiece for comfortable fit.

Item No.	Description	Unit/Pkg	Pkg/Case
2011	SafeMask Premier Earloop - Lavender	50	10
2015	SafeMask Premier Earloop - Blue	50	10
2016	SafeMask Premier Earloop - Pink	50	10
2017	SafeMask Premier Earloop - Yellow	50	10
2018	SafeMask Premier Earloop - Teal	50	10





Low Barrier mask with comfortable tie-on laces. Suitable for procedures with low levels of spray and splatter such as fluoride treatments, minor orthodontics and basic lab work.

- Construction and materials providing adequate protection and comfort.
- Shingle pleat design prevents fluid pooling and cross contamination. Non-woven inner layer for superior protection.
- Adjustable nosepiece and tie-on for comfortable fit. Extra long laces fit with a variety of facial shapes and sizes.

Item No.	Description	Unit/Pkg	Pkg/Case
2000	SafeMask Tie On - Blue	50	6





Utility mask made of latex free, porous material which helps protect the wearer from body fluids and contamination caused by exhaled micro-organisms. Molded contour design increases comfort and allows for excellent breathability.

- Adjustable, aluminum nosepiece and stretchable headband allow for a custom, close fit.
- Materials and construction are soft and non-irritating.
- Comfortable to wear and good alternative to shingle pleat masks.

ltem No.	Description	Unit/Pkg	Pkg/Case
2020	SafeMask W/Alum Nose Guard - Blu	e 50	12
2021	SafeMask W/Alum Nose Guard - Pinl	< 50	12








NIOSH* approved duckbill and cone respirators feature a 95% filter efficiency level and provide protection against particulate aerosols free of oil.

- Complies with OSHA Bloodborne Pathogens Standard when properly fitted to face and worn with proper protective eyewear.
- Maximized air chamber, lightweight construction and soft inner lining contribute to user comfort and breathability.
- Adjustable nosepiece for custom fit. Secure seal against most face sizes and shapes.



* National Institute for Occupational Safety and Health



Eyewear for use when levels of fluids and splashes warrant extra protection. Wraparound, anti-fog lens provides added peripheral protection and optimal, clear visibility.

- Single-use unit is preassembled for timesaving and convenience, and allows for attachment of mask earloops for total protection.
- Lightweight, comfortable plastic frame adjusts, allowing for custom fit.
- Ergonomic design for all-day wearability.

Item No.	Description	Unit/Pkg	Pkg/Case
2900	SafeMask Protective Eyewear	10	6





Primary protection for simple, low-risk procedures involving very light dust, moisture, spray and spatter. Ideal for patient appointments and exams.

- Inner and outer layer constructed for basic fluid resistance, with ultrasonically sealed, soft, stretchy, latex free earloops for superior comfort.
- Adjustable nosepiece for custom and comfortable fit and easy breathing.
- Economic alternative for simple procedures.

Item No. 2115M

2116M

2117M

2118M





Designed with sensitive skin in mind. Its soft, lightweight and highly breathable inner layer is gentle, non-irritating and remains cool on the skin for ultimate ease and comfort.

- Cellulose-blend inner layer is white, extra soft and dye free which leaves skin feeling clean and dry. It does not lint, tear or shred.
- Meets ASTM 2100 standard for everyday, low barrier protection.
- Adjustable nose piece helps customize the fit. Earloops ties are soft, stretchy, latex free and sonically sealed.

Description	Unit/Pkg	Pkg/Case	Item No.	Description
Distech Earloop Mask - Blue	50	10	2080	SafeMask Sofskin E
Distech Earloop Mask - Pink	50	10	2081	SafeMask Sofskin E
Distech Earloop Mask - Yellow	50	10	2082	SafeMask Sofskin E
Distech Earloop Mask - Green	50	10	2083	SafeMask Sofskin Ea

ltem No.	Description	Unit/Pkg	Pkg/Case
2080	SafeMask Sofskin Earloop - White	50	10
2081	SafeMask Sofskin Earloop - Blue	50	10
2082	SafeMask Sofskin Earloop - Pink	50	10
2083	SafeMask Sofskin Earloop - Lavender	50	10





Creates a tight seal around the face with the help of the aluminum nose and chin pieces. This reduces gapping at the bottom and sides, provides good face coverage and better protects the wearer.

- Aluminum nose and chin pieces create a comfortable personalized fit around the circumference of the face every time.
- Meets ASTM 2100 performance standards for everyday, low and high barrier protection.
- Soft non-woven material. Earloops ties are soft, stretchy, latex free and sonically sealed.

ltem No.	Description	Unit/Pkg	Pkg/Case
2070	SafeMask TailorMade Low Barrier - Blue	50	10
2071	SafeMask TailorMade Low Barrier - Pink	50	10
2072	SafeMask TailorMade High Barrier - Blue	50	10







MEDICOM[®] US CATALOG 2012

MEDICOM[®] Safe-Gauze[®]

Clearly the best choice for absorption

Medicom[®] SafeGauze[®] Non-Woven sponges are the perfect combination of performance and price. SafeGauze[®] sponges are more absorbent than traditional cotton filled sponges, are nonlinting, maintain their shape during use and easily release from incisions improving the healing process. SafeGauze[®] sponges are costeffective and economical since 4 plies of nonwoven can substitute 8 plies of cotton gauze. They are made with extra soft, non-abrasive, premium quality non-woven material and have no chemical whiteners added.

See for yourself why SafeGauze® has been a trusted North American market leader for 20 years!



- Non-woven sponges minimize the number used per procedure, saving time and money, when compared with cotton-filled sponges.
- Extra soft, non-adherent, polyester/viscose blend for increased patient comfort.
- Latex free Non-sterile
- Made in the USA.

Item No.	Description	Unit/Pkg	Pkg/Case
4522	SafeGauze 2x2 N/W Sponges 4 Ply	200	20
4522-5	SafeGauze 2x2 N/W Sponges 4 Ply	200	25
4544	SafeGauze 4x4 N/W Sponges 4 Ply	200	10



Lyocell[®]. Available in 2x2 and 4x4 sizes.

• 100% biodegradable natural fibres using

• Luxuriously soft, highly absorbent and

SafeGauze Green 2x2 N/W Sponge 200

SafeGauze Green 4x4 N/W Sponge 200

Unit/Pkg Pkg/Case

20

10

Item No.

2102

2104

chlorine-free practices.

• Latex free. Non-sterile.

• Made in the USA.

non-linting.

Item No. Description

4722

4744







Increased patient comfort due to the soft, viscose/poly blend, non-woven material. An excellent and economical alternative to cotton gauze sponges.

Unit/Pkg

200

200

Pkg/Case

20

10

• Multipurpose, non-linting, and highly absorbent material saves time

2x2 N/W Gauze Sponges 4 Ply

4x4 N/W Gauze Sponges 4 Ply

- No chemical additives to avoid allergic reactions.
- Latex free. Non-sterile.

Description



Peace of mind with sterile, heat-sealed and secure envelope for added safety from cross contamination where additional caution is warranted.

- Absorbent, soft, non-linting and breathable material for increased patient comfort.
- Convenient, per-use, individual peel-open design.
- 2 sponges per envelope.

Item No.	Description	Env./Pkg	Pkg/Case
2240	2x2 N/W Sterile Gauze 4 Ply	100	30
9210	4x4 Sterile Drain Sponges 6 Ply	50	12





10 MEDICOM













Increased convenience with easy-to-use, individual, peel-open package design.

- Highly absorbent, soft sponges for optimal protection and patient comfort.
- Sterile, heat-sealed envelope for superior safety.
- 2 sponges per envelope.

Applicators are available in 3" and 6" lengths and in sterile or non-sterile versions. Sterile applicators are sealed in autoclavable bags of 100.

- Highly absorbent, extra-soft cotton tips contribute to patient comfort.
- High-quality construction for superior safety and dependability.
- Latex free, recyclable and degradable.

Item No.	Description	Sponges/Pkg	Pkg/Case	Item No.	Description	Unit/Pkg	Pkg/Case
4800	4x4 Sterile Gauze Double 8 Ply	100	12	803	Non Sterile Cotton Tip Applicator 3"	1000	10
				806	Non Sterile Cotton Tip Applicator 6"	1000	10
				804	Sterile Cotton Tip Applicator 3"	200	10
				807	Sterile Cotton Tip Applicator 6"	200	10

Flexible roll format holds its shape and performs well for a number of dental procedures. Suitable for sensitive patients due to easy removal and minimal irritation.

- Soft and pliable, Medium, #2 rolls made from 100% cotton.
- Highly absorbent and maintain their shape during use.
- Latex and chlorine free.

Item No.	Description	Unit/Pkg	Pkg/Case
4554	Medicom N/S Cotton Dental Rolls #2 24000 pcs/case	2000	12



Flexible roll format holds its shape and performs well for a number of dental procedures. Suitable for sensitive patients due to easy removal and minimal irritation.

- Soft and pliable, Medium, #2 rolls made from 100% cotton.
- Highly absorbent and maintain their shape during use.
- Latex and chlorine free.

tem No.	Description	Unit/Pkg	Pkg/Case
3554	Distech N/S Cotton Dental Rolls #2	2000	12



STERILIZATION PRODUCTS

SAFE-SEAL

To assure a high quality product, effective sterilization and compliance with ISO Standards, Safe-Seal Duet was put through a series of tests and validations. Safe-Seal Duet now contains an internal multi-variable indicator and an external process indicator for steam sterilization that responds to all variables of time, temperature and saturated steam. The product line meets CDC guidelines for effective sterilization and complies with both ISO 11140 and ISO 11607.





Assured sterilization with an internal multi-variable indicator and an external process indicator for steam sterilization. Meets CDC guidelines for effective sterilization and responds to all variables: time, temperature and saturated steam.

- Secure seal with extra wide chevron seal protects against instrument breaching while high quality porous medical-grade paper meets medical standard.
- Perforated fold facilitates easy and accurate folding of adhesive strip to form airtight seal. Thumb notch enables easy opening of pouch.
- Lead free, high-quality inks compatible with a large number of commonly used disinfectants and assure easy confirmation of sterilization process.

ltem No.	Description	Unit/Pkg	Pkg/Case
88000	SafeSeal Duet 2 3/4" x 9"	200	6
88005	SafeSeal Duet 3 1/2" x 5 1/4"	200	5
88010	SafeSeal Duet 3 1/2" x 9"	200	6
88015	SafeSeal Duet 2 1/4" x 4"	200	10
88025	SafeSeal Duet 5 1/4" x 10"	200	10
88030	SafeSeal Duet 7 1/2" x 13"	200	5
88035	SafeSeal Duet 10" x 14"	200	5
88040	SafeSeal Duet 12" x 17"	200	5





Various sized reels allows for customized sterilization of a wide variety of items.

- Individual and flexible formats save time and lower costs. Ideal for specialty items of various sizes.
- Special grade, heavy-weight paper provides superior strength.
- Blue tinted transparent material for easy viewing and possible instrument breaching.

Item No.	Description	Unit/Pkg	Pkg/Case
9510+	Peelable Reels 2" x 100'	10	1
9520+	Peelable Reels 3" x 100'	10	1
9530+	Peelable Reels 4" x 100'	10	1
9540+	Peelable Reels 6" x 100'	10	1

MEDICOM

* Color change guide available on medicom.ca







Protect Yourself this Cold Season!

Microban[®] Hand Sanitizer reduces the risk of infection and contamination. It contains 70% (v/v) ethyl alcohol so it is effective against the bacteria and viruses commonly found on hands such as E.Coli, MRSA, Rhinovirus and Noroviruses. It also contains glycerin, aloe and vitamin E so it leaves hands feeling soft, smooth and non-sticky.

Effective. Fast. Gentle.

14 MEDICOM[®] US CATALOG 2012




Microban® Hand Sanitizer



Alcohol-based hand sanitizer which offers a fast and simple means of reducing bacteria and conditioning hands in many healthcare settings.

- High alcohol concentration (70% Ethyl alcohol) to help reduce the risk of infection and contamination from bacteria and viruses commonly found on hands (such as E. coli, MRSA, rhinoviruses and noroviruses).
- Kills 99.9% of germs
- Aloe vera and vitamin E ingredients serve to moisturize and rejuvenate hands and leave them feeling soft, smooth and non-sticky.

Item No.	Description	Unit/Pkg	Pkg/Case
1260-125R	Microban Sanitizer 4.2 fl. oz Original	1	12
1260-500R	Microban Sanitizer 17 fl. oz Original	1	12
1260-100R	Microban Sanitizer 33.8 fl. oz Original	1	6
1260-130	Microban Sanitizer 44 fl. oz Original	1	12
1269-000	Microban 44 fl. oz Wall Mounted Dispenser	1	6



WHAT YOU SHOULD KNOW.

Waterless, alcohol-based hand rubs are the preferred method for hand hygiene in most situations due to the superior efficacy of these agents in rapidly reducing bacterial counts on hands and their ease of use.¹

Some viruses and bacteria can live from **20 minutes up to 2 hours** on various surfaces.²

According to the Center for Disease Control, there are more than 52 million cases of the common cold each year in the U.S. 3

In 2002, approximately **10 million adults** in the U.S. (ages 18-69) were unable to work due to health problems.⁴

Failure to wash, or insufficiently washing hands, contributes to almost **50%** of all food borne illness outbreaks.²

Alcohol-based hand rubs using ethanol as the main active ingredient demonstrate activity against many non-enveloped viruses and are good predictors of activity against enveloped viruses and vegetative bacteria, including fungi.

¹ CDC, APIC, SHEA. How-to Guide: Improving Hand Hygiene. A Guide for Improving Practices among Health Care Workers___www.IHI.org

- ² C. Condayan. Handwashing Fact Sheet. MicrobeWorld. September 2009 ___http://microbeworld.org ³ CDC. Fact Stats A to Z: Common cold 1996
- ⁴ CDC. Ounce of Prevention Campaign. April 24, 2008 ____http://www.cdc.gov/ounceofprevention/






MEDICOM. Denti-Gares

For all your preventive needs, use Denti-Care[®] regularly, the complete line of preventive products including prophy pastes, fluoride gels and foams whose proven results are recognized throughout the dental









Safe and compatible gel that is the optimal choice for treating patients with cosmetic and porcelain restorations as well as patients suffering from root caries or tooth sensitivity.

- Smooth, thixotropic consistency and excellent tasting flavors for patient pleasing formula.
- 4-minute application (contains 9,050 ppm fluoride).
- Ideal for patients with reduced salivary flow or acid intolerance.

Item No.DescriptionUnit/PkgPkg/Case10023-MUN2% Neutral Na FL Gel M1110023-RBUN2% Neutral Na FL Gel RB1110023-STUN2% Neutral Na FL Gel ST11





Optimized patient comfort with reduced gagging and risk of ingestion due to stable, airy foam and no tray overflow.

- 60-second express foam application for saving time and money. Low pH and 12,300 ppm fluoride content for effective and rapid results.
- Inhibit caries formation by reducing enamel solubility and enhancing remineralisation.
- Unique bottle designed to maximize the use of the entire content — for up to 130 applications.

Item No.	Description	Unit/Pkg	Pkg/Case
10037-BGUN	Acidic Topical Foam FL BG	1	1
10037-GRUN	Acidic Topical Foam FL GR	1	1
10037-MUN	Acidic Topical Foam FL M	1	1
10037-OUN	Acidic Topical Foam FL O	1	1
10037-RBUN	Acidic Topical Foam FL RB	1	1
10037-STUN	Acidic Topical Foam FL ST	1	1





Due to the gentle neutral pH, foam is ideal for use on patients with cosmetic or porcelain dental restorations, crowns or bridges.

- Stable, airy formula and no tray overflow for optimized patient comfort and reduced gag reflex. 4-minute application and formula containing 9,050 ppm fluoride.
- Innovative bottle design maximizes the use of the entire content for up to 130 applications.
- Delightful, great tasting flavors for better patient acceptance.

Item No.	Description	Unit/Pkg	Pkg/Case
10038-MUN	Dentifoam 2% NSF Top Foam M	1	1
10038-RBUN	Dentifoam 2% NSF Top Foam RB	1	1

Bubble gum = BG • Cherry = CH • Grape = GR • Mint = M • Orange = O • Raspberry = RB • Strawberry = ST • Unflavored = UN



1

Bubble gum = BG • Cherry = CH • Grape = GR • Mint = M • Orange = O • Raspberry = RB • Strawberry = ST • Unflavored = UN



at home.

Item No.

10022-UN



Denti-Powder Prophylaxis Powder



Ideal for quick stain removal and deep cleaning into pits, fissures and inter proximal areas for sub-gingival plaque removal and selective polishing procedures.

- Reduced powder particle size keeps clogging to a minimum and allows for optimal use of most airflow devices.
- Convenient bottle design allows for no-spill, easy pouring and reduces waste.
- Flip cap maintains powder freshness.





Patients will enjoy this great tasting and easy-to-use rinse application with no trays required.

- Esthetically compatible with composite and porcelain dental restorations.
- Provides increased protection for patients with high risk of decay.
- Large, 67.6 fl. oz./2 litre format.





Ideal for professional use and in conjunction with home care for gingivitis, periodontitis, oral irrigation and postoperative healing.

- Provides an extra boost in the reduction of tissue inflammation, bleeding and plaque accumulation.
- Enhanced light mint flavor reduces any bitter aftertaste and increases patient acceptance.
- Convenient 16 oz. size for daily home care rinsing.





Save time with these convenient and comfortable-fit foam trays, created with a typodont design to reduce overflow and provide even coverage of fluoride to all surfaces.

Available in hinged and single alternatives:

Hinged

- Easy, removable tabs for fast and simple clean up.
- Color-coded sizes: Small (Yellow), Medium (White), Large (Blue).

Single

- Provides a more comfortable experience, one arch at a time.
- Ideal for children and patients with a hyperactive gag reflex.

Item No.	Description	Unit/Pkg	g Pkg/Case	Item No.	Description	Unit/Pkg	Pkg/Case	Item No.	Description	Unit/Pkg	Pkg/Case	Item No.	Description	Unit/Pkg	g Pkg/Ca
10032-MBUN	I Denti-Care Prophy Powder M	1	1	10044-BUN	2% Neutral Na FL Rinse Berry	1	1	10025-H	0.12% Chg Oral Rinse 500ml	12	1	10029-LUN	Disposable Gel Tray - L Hinged	50	1
10032-RBUN	Denti-Care Prophy Powder RB	1	1	10044-MUN	2% Neutral Na FL Rinse M	1	1					10029-MUN	Disposable Gel Tray - M Hinger	d 50	1
												10029-SUN	Disposable Gel Tray - S Hinged	50	1
												10030-LUN	Disposable Gel Tray - L Singles	50	1
												10030-MUN	Disposable Gel Tray - M Singles	50	1
												10030-SUN	Disposable Gel Tray - S Singles	50	1

 $Bubble gum = BG \bullet Cherry = CH \bullet Grape = GR \bullet Mint = M \bullet Orange = O \bullet Raspberry = RB \bullet Strawberry = ST \bullet Unflavored = UN$







Duraflor[®], a name you've trusted since 1984, treats hypersensitivity and postoperative sensitivity. The 5% Sodium Fluoride Varnish formulas also help remineralize tooth enamel and protects teeth from acid erosion.







Duraflor 5% Sodium Fluoride Varnish in Tube



Satisfying healthcare professionals and their clients for more than 30 years, this product comes in a handy and economical tube format and plays an important role in preventive dental health. Ideal for hypersensitivity treatments.

- Proven fluoride treatment.
- Sweetened with Xylitol for increased patient acceptance.
- 0.34 fl. oz tube allows for multiple treatments.



Duraflor 5% Sodium Fluoride Varnish in Unit Dose



A classic fluoride treatment available in convenient, unit dosing. Features the highest concentration of fluoride available—and fast and easy application.

- Proven, 30-year, in-practice use for dependable, effective results.
- Quick application and long-lasting results.
- 0.5 mL unit doses.



Duraflor Halo 5% Sodium Fluoride Varnish in Unit Dose




Releases more than 3X the fluoride than the leading competitive brand, over a 4-hour period.* Recommended for the treatment of hypersensitive teeth, post-operative sensitivity and to seal dentinal tubules for cavity preparation. Prophy paste ring included in each box.

- "Peel and Apply" pre-measured, one-step package for easy, consistent application which eliminates messes and the risk of cross contamination.
- Highest allowable fluoride concentration for fast action. Quickly dries to natural white tooth color.
- Sweetened with Xylitol for increased patient acceptance.
- Free of nut, sesame, corn, shellfish, penicilin, egg, dairy products and gluten.
- 0.5 mL unit doses.

Item No.	Description	Unit/Pkg	Pkg/Case
1015-SM32	Duraflor Halo SM 32	32	6
1015-WB32	Duraflor Halo WB 32	32	6
1015-SM250	Duraflor Halo SM 250	250	1
1015-WB250	Duraflor Halo WB 250	250	1

* Study FVR 07-218-2: "Short-term Fluoride Release from a Fluoride Varnish," Thermametric Technologies, Inc. 11/09/07. Varnish is a registered trademark of 3M or 3M ESPE.



Bubble gum = BG • Spearmint = SM • Wild Berry = WB • Raspberry = RB





From bibs to head rest covers, Medicom[®] has a wide selection of high quality products for all your daily needs.

- Dental Bibs
- Dental Tray Covers
- Head Rest Covers
- Paper Cups
- Plastic Cups





Durable patient protection thanks to high absorption. Ideal for cleanings and minor procedures.

- Horizontal embossed pattern and a unique water-repellant edge for enhanced shielding from fluids.
- Strong, tear-resistant plastic backing allows for increased clothing protection.
- Many color options to coordinate with Medicom[®] plastic cups.

Item No.	Description	Unit/Pkg	Pkg/Case
8280	Dental Bibs 2 Ply - Lavender	500	1
8281	Dental Bibs 2 Ply - Grey	500	1
8282	Dental Bibs 2 Ply - Blue	500	1
8283	Dental Bibs 2 Ply - White	500	1
8284	Dental Bibs 2 Ply - Yellow	500	1
8285	Dental Bibs 2 Ply - Green	500	1
8286	Dental Bibs 2 Ply - Dusty Rose	500	1
8287	Dental Bibs 2 Ply - Beige	500	1
8288	Dental Bibs 2 Ply - Peach	500	1
8289	Dental Bibs 2 Ply - Aqua	500	1





Enhanced patient protection with additional layer of paper and plastic backing. Suitable for minor procedures, like cleanings, where extra protection is desired.

- Superior absorption and clothing protection with additional layer of paper and plastic backing.
- Excellent protection thanks to horizontal embossed pattern and a unique water-repellant edge.
- Many color options to coordinate with Medicom[®] plastic cups.

Item No.	Description	Unit/Pkg	Pkg/Case
8380	Dental Bibs 3 Ply - Lavender	500	1
8381	Dental Bibs 3 Ply - Grey	500	1
8382	Dental Bibs 3 Ply - Blue	500	1
8383	Dental Bibs 3 Ply - White	500	1
8384	Dental Bibs 3 Ply - Yellow	500	1
8385	Dental Bibs 3 Ply - Green	500	1
8386	Dental Bibs 3 Ply - Dusty Rose	500	1
8387	Dental Bibs 3 Ply - Beige	500	1
8388	Dental Bibs 3 Ply - Peach	500	1
8389	Dental Bibs 3 Ply - Aqua	500	1










Covers help prevent cross contamination and unwanted water seepage to your equipment and tray surfaces. Contribute to overall cleanliness and patient comfort.

- Available in a variety of colors, shapes and sizes to accommodate assortment of equipment and decors.
- Made with high-quality paper stock.
- Single use, latex free.

Item No.	Description	Unit/Pkg	Pkg/Case
5593	Dental Tray Cover 8 1/2" x 12 1/4"	1000	1
5593-BL	Dental Tray Cover 8 1/2" x 12 1/4" - Blue	1000	1
5593-DR	Dental Tray Cover 8 1/2" x 12 1/4" - Dusty Rose	1000	1
5593-LAV	Dental Tray Cover 8 1/2" x 12 1/4" - Lavender	1000	1
5594	Dental Tray Cover 9 1/2" x 12 1/4"	1000	1
5595	Dental Tray Cover 9" x 13 1/2"	1000	1
5595-BL	Dental Tray Cover 9" x 13 1/2" - Blue	1000	1
5595-DR	Dental Tray Cover 9" x 13 1/2" - Dusty Rose	1000	1
5595-LAV	Dental Tray Cover 9" x 13 1/2" - Lavender	1000	1
5596	Dental Tray Cover 10 1/4" x 15 3/4"	1000	1
5597	Dental Tray Cover 11" x 17 1/2"	1000	1
5598	Dental Tray Cover Round 13 1/2"	1000	1
5599	Dental Tray Cover 11" x 11"	1000	1



Increase patient comfort with fresh covers to protect them and your equipment at the same time.

- Easy-to-place and remove covers use fluid repellant materials to save time and reduce cross contamination.
- Soft poly-blend fabric is highly breathable and contributes to patient comfort.
- Protect the appearance and quality of your dental chair.

3014 Headrest Cover 10" x 10" - White	500	4
	000	1
3015 Headrest Cover 10" x 10" - Blue	500	1
3016 Headrest Cover 10" x 10" - Dusty Rose	500	1
3020 Headrest Cover 10" x 10" - Lavender	500	1
3017 Headrest Cover 10" x 13" - White	500	1
3018 Headrest Cover 10" x 13" - Blue	500	1
3019 Headrest Cover 10" x 13" - Dusty Rose	500	1
3021 Headrest Cover 10" x 13" - Lavender	500	1





Rinse-cup outside designs put patients at ease and contribute to their overall comfort. "Rolled lip" for better lip feel and control.

- Made from biodegradable material to minimize harm to the environment.
- Happy Teeth, Bubbles and Hello Kitty® designs to choose from.
- Reinforced, coated, interior avoids leakage.

Item No.	Description	Unit/Pkg	Pkg/Case
113-CH	Paper Cups - Bubbles Design 4oz	100	10
114-CH	Paper Cups - Healthy Teeth 4oz	100	10
115-CH	Paper Cups - Bubbles Design 5oz	100	10
116-CH	Paper Cups - Healthy Teeth 5oz	100	10
120-HK	Paper Cups - Hello Kitty 5oz-Mix	100	10





Pastel cup colors coordinated with Dry-Back® dental bibs. Plastic material for solid feel.

• Double coating with built-in grip control outer-wall design for increased strength and stable feel.

- Special rolled-rim design for improved "lip feel" and control.
- Attractively priced and high quality.

Item No.	Description	Unit/Pkg	Pkg/Case
104	Plastic Cups 5oz - Peach	100	10
105	Plastic Cups 5oz - Green	100	10
106	Plastic Cups 5oz - Lavender	100	10
107	Plastic Cups 5oz - Dusty Rose	100	10
108	Plastic Cups 5oz - Yellow	100	10
109	Plastic Cups 5oz - Beige	100	10
110	Plastic Cups 5oz - White	100	10
111	Plastic Cups 5oz - Grey	100	10
112	Plastic Cups 5oz - Blue	100	10





MEDICOM[®]

Cover up from head to toe with a variety of quality products.

- Isolation Gowns
- Lab Gowns
- Shoe Covers
- Bouffant Caps







For procedures requiring protection from spray and splatter, fluid-repellent design provides added protection against cross contamination.

- Elastic cuffs, waist and neckties for comfortable fit.
- Sewn seams provide for superior strength. Available in a variety of colors and sizes.
- Single-use, disposable and latex free.

Item No.	Description	Unit/Pkg	Pkg/Case
8009	Isolation Gown Reg - Yellow	10	5
8010	Isolation Gown XL - Yellow	10	5
8023	Isolation Gown Reg - Blue	10	5





Front opening snaps, knitted cuffs and collar for comfortable, easy in-and-out fit.

- 3-ply SMS fabric for increased breathability.
- Protection against cross contamination.
- Single-use, disposable and latex free.
- Available in sky blue.

ltem No.	Description	Unit/Pkg	Pkg/Case
8016	Lab Gown S - Sky Blue	10	5
8017	Lab Gown M - Sky Blue	10	5
8018	Lab Gown L - Sky Blue	10	5
8019	Lab Gown XL - Sky Blue	10	5







Footwear protection with fluid repellant, slip-on covers.

- Sewn seams provide increased tear strength.
- Elastic opening for easy on-and-off convenience.
- Non-conductive, non-skid and plain options.

Item No.	Description	Unit/Pkg	Pkg/Case
8006	Shoe cover N-Conduct N-Skid-XL - Blue	100	3
8007	Shoe cover N-Conduct N-Skid-Reg - Blue	100	3
8008	Shoe cover N-Conduct Reg - Blue	100	3





Hair and head fluid protection with breathable, low-linting design.

- Latex free, elastic headband for reduced allergic reactions.
- Fluid repellant material for optimal protection.
- Available in sky blue and white.

Item No.	Description	Unit/Pkg	Pkg/Case
8000	Bouffant Caps 21" White - In Box	100	5
8001	Bouffant Caps 21" Blue - In Box	100	5
8003	Bouffant Caps 24" White - In Box	100	5
8004	Bouffant Caps 24" Blue - In Box	100	5







MEDICOM[®]

Medicom[®] is dedicated to offering the widest selection of high quality products.

- Evacuation Products
- Air-water Syringe Tips
- Tongue Depressors
- Barriers
- Dental Dams
- Scalpels and Blades



Saliva Evacu	a Ejectors ation Prod	S ucts	Pro-Tip T Single-Us Syringe T Evacuation	e Air/W	/ater	2	High Vo Evacuato	lume ors n Proc	ucts		Tongue	Depr	essor	2	Barrier F Barriers	ilm	
	S.		all a	5		111		1	P			2	3				
 Provides opt maintaining you work effice Non-rusting a forms into dee for ease of use Soft bonded, aspirations fo safety and con World-renown 	ts shape t ciently and o lloy wire that sired configu e. smooth tip to r increased p mfort.	to help easily. ration o avoid atient	All the benefits of a water syringe tip w popular 3-in-1 m range of converter • Unique dual chamk separate air from w instant dry air and p spray when and wh • Single-use design p contamination. • Rigid design is effer retraction and pliak edges increases pa	vhich fits odels - s availab per desigr ater, guar perfect mi ere you n prevents c ctive for ti ple, round	s most wide I.e. to antees isted eed it. cross issue ed	requ acce "S" • 45° and asp • Str cor wo • Str vor	eased suction fouring it. Enjoys and visibility shaped design. Angled tip for redideal provide the form of the subscription of	y inc with u duced s ncrease ccess to ole tips i	reased inique, splatter d	patie smo cons	l for oral exa ent care. Preci oth polished f istent performant ex free.	ision cu finish. I	ut with Delivers	 Cont hass conv Un cou use for Ad ap 	e time and i tamination wit sle-free film wh venient rolls. ique packaging fo unter-top dispense e, tear-off sheets convenience. hesive coating all plication and rem ailable in blue and	h sing ich con eatures er with s ows for oval.	le-use, mes in single-
Item No. Description	Unit/Pk	g Pkg/Case	Item No. Description	Unit/Pkg	Pkg/Case	Item No.	Description	Unit/Pkg	Pkg/Case	Item No.	Description	Unit/Pkg	Pkg/Case	Item No.	Description	Unit/Pkg	Pkg/Case
0500 Saliva Ejector Clear White /		10	1007072 Pro-Tip Turbo Syringe Tip - Bag o	250	12	505	High Vol Evac Non-Vented-White	50	20	40003	5.5" Non Sterile Tongue Depressor	500	10	5050	Barrier Film - Blue	1200	8
501 Saliva Ejector Clear White /	100	10	1007073 Pro-Tip Turbo Syringe Tip - Bulk	250	6	0510	High Vol Evac Vent 1 side-White	50	20	40004	5.5" Sterile Tongue Depressor	500	10	5051	Barrier Film - Clear	1200	8
502 Saliva Ejector White White /	100 White Tip	10	1007080 Pro-Syringe	1	1					40006	6" Non Sterile Tongue Depressor	500	10				
vvriite vvnite /	white tip		1007090 Pro-Tip Turbo Adec Converter	1	2					40007	6" Sterile	500	10				
			1007091 Pro-Tip Turbo Engle Converter	1	2						Tongue Depressor						
			1007092 Pro-Tip Turbo DentalEZ Converte	1 r	2												
			1007093 Pro-Tip Turbo DCI Converter	1	2												







Designed to keep common surfaces in the operatory free from harmful bacteria caused by airborne aerosols, splash and spatter.

- Intelligent design allows for easy and fast changing of barriers between patients to save valuable time.
- Wide ranging shapes and sizes to cover virtually everything in the dental operatory— from tray covers to syringe sleeves.
- Simple, grab-and-go dispenser boxes for convenience.

Item No.	Description	Unit/Pkg	Pkg/Case
9551	SafeTBarrier TR SLV 11 5/8" x 14 1/2"	500	4
9552	SafeTBarrier HD RST CVER 10" x 14"	250	4
9553	SafeTBarrier TR SLV 10 1/2" x 14"	500	4
9554	SafeTBarrier TR SLV 11 1/2" x 16 1/2"	500	4
9555	SafeTBarrier XR SLV 15 1/2" x 26"	250	4
9556	SafeTBarrier T-LIGHT 4" x 5 3/4"	500	12
9557	SafeTBarrier SYR SLV 2 1/2" x 10"	500	12
9558	SafeTBarrier SYR SLV w ADH 2.5" x 11.5"	500	12
9559	SafeTBarrier SYR SLV 2 1/2" x 10" Blue	500	12
9560	SafeTBarrier SYR SLV No Open 2.5" x 10"	500	12
9564	SafeTBarrier SYR SLV w ADH & No Open	500	12





Extreme strength and flexibility allow for manipulation in the most difficult to reach areas.

- Maximum barrier protection and isolation of operating worksite reduces cross contamination by controlling moisture.
- Next generation, built-in strength and improved tear resistance save time and avoid potentially costly wastage.
- Made with natural rubber latex.

ltem No.	Description	Unit/Pkg	Pkg/Case
1255HB	Dentdam 5" x 5" Heavy Blue - Unscent	52	6
1255HG	Dentdam 5" x 5" Heavy Green - Mintscen	t 52	6
1255MB	Dentdam 5" x 5" Med Blue - Unscent	52	6
1255MG	Dentdam 5" x 5" Med Green - Mintscent	52	6
1255TB	Dentdam 5" x 5" Thin Blue - Unscent	52	6
1255TG	Dentdam 5" x 5" Thin Green - Mintscent	52	6
1266HB	Dentdam 6" x 6" Heavy Blue - Unscent	36	6
1266HG	Dentdam 6" x 6" Heavy Green - Mintscent	t 36	6
1266MB	Dentdam 6" x 6" Med Blue - Unscent	36	6
1266MG	Dentdam 6" x 6" Med Green - Mintscent	36	6
1266TB	Dentdam 6" x 6" Thin Blue - Unscent	36	6
1266TG	Dentdam 6" x 6" Thin Green - Mintscent	36	6





Single-use scalpels avoid cross contamination and are suitable for simple surgical procedures.

- Precision molding of blade into handle, with a replaceable blade guard, for complete user protection.
- Comfort fit handle design for good grip and easy manipulation.
- Made of high-quality stainless steel.

Item No.	Description	Unit/Pkg	Pkg/Case
90000-10	Sterile Disposable Scalpel #10	10	1
90000-11	Sterile Disposable Scalpel #11	10	1
90000-12	Sterile Disposable Scalpel #12	10	1
90000-15	Sterile Disposable Scalpel #15	10	1
90000-20	Sterile Disposable Scalpel #20	10	1
90000-21	Sterile Disposable Scalpel #21	10	1
90000-22	Sterile Disposable Scalpel #22	10	1
90000-23	Sterile Disposable Scalpel #23	10	1
90000-24	Sterile Disposable Scalpel #24	10	1





Blade-only option allows for cost savings and avoids cross contamination.

- Precision, cutting-edge design ideal for multiple incisions.
- Sterile package sterilized by gamma radiation of 2.5 MRads.
- Made of high-quality stainless or carbon steel.

Item No.	Description	Unit/Pkg	Pkg/Case
90010-10	Sterile Blades Stainless #10	100	1
90010-11	Sterile Blades Stainless #11	100	1
90010-12	Sterile Blades Stainless #12	100	1
90010-15	Sterile Blades Stainless #15	100	1
90010-20	Sterile Blades Stainless #20	100	1
90010-21	Sterile Blades Stainless #21	100	1
90010-22	Sterile Blades Stainless #22	100	1
90010-23	Sterile Blades Stainless #23	100	1
90010-24	Sterile Blades Stainless #24	100	1
90011-10	Sterile Blades Carbon Steel #10	100	1
90011-11	Sterile Blades Carbon Steel #11	100	1
90011-12	Sterile Blades Carbon Steel #12	100	1
90011-15	Sterile Blades Carbon Steel #15	100	1
90011-21	Sterile Blades Carbon Steel #21	100	1
90011-22	Sterile Blades Carbon Steel #22	100	1



MEDICOM[®] Safe+Mask[®]

MASKS ADDITIONAL INFORMATION

Latex free and fiberglass free • Dye free inner layers • All single use

Regulated Facial Protection You Can Trust

To assure adequate facial protection for any type of procedure, the American Society for Testing Materials (ASTM) recommends to **only select ASTM-rated masks**. The ASTM plays a leadership role in addressing the standardization needs of the global marketplace. ASTM-2100 is a standard which covers the testing and requirements for materials used in the design of medical face masks. Depending on the testing performance results, the standard is able to categorize each mask according to 3 barrier levels: low, moderate and high.

ASTM Requirements	High Barrier	Moderate Barrier	Low Barrier
Resistance to penetration by synthetic blood (mm Hg)	160	120	80
Bacterial Filtration Efficiency (BFE %)	≥ 98	≥ 98	≥ 95
Sub-micron particulate filtration efficiency at 0.1micron (PFE %)	≥ 98	≥ 98	≥ 95
Differential pressure ΔP (mm H ₂ O/cm ²)	< 5.0	< 5.0	< 4.0
Flame Spread	Class1	Class1	Class1

High Barrier masks are ideal for:

Moderate Barrier masks are ideal for:

- Oral exams for high-risk patients
- Procedures involving high levels of spray, spatter and airborne particles
- Ultrasonic scalers and high-speed aerosols
- Air polishers and air-water syringes
- Oral surgery (tooth extractions, maxo/facial reconstruction)

- Oral exams for known higher risk patients
- Procedures involving moderate levels of spray, spatter and airborne particles
- High or low-speed handpieces
- Dry, low aerosol procedures
- Bleaching, orthodontic and periodontic procedures
- Lab work with exposure to blood or contaminated materials

Low Barrier masks are ideal for:

- Oral exams for low-risk patients
- Procedures involving low levels of spray and spatter
- Orthodontic treatment with little or no spray
- Hand instrumentation with little moisture or spray
- Fluoride treatments
- Basic lab work





www.medicom.ca

A.R. Medicom Inc.

1200, 55th Avenue, Montreal (Quebec) H8T 3J8 Canada

Tel: 514 636-6262 Email: sales@medicom.com

A.R. Medicom Inc. (Asia) Ltd.

Rm. 1511-B, Kailey Industrial Centre 12 Fung Yip Street, Chai Wan, Hong Kong

Tel: 011 852-2147-4968 Email: sales@medicom-asia.com Customer Service Tel: 514 636-6131

Toll Free (North America) Tel: 1 800 361-2862

A.R. Medicom Inc. (Japan) Ltd.

Edomachi 96, Strong Bldg. 8F Chuo-Ku, Kobe, Hyogo, Japan 650 0033

Tel: 011 81-78-327-2333 Email: sales@medicom-japan.com Medicom USA

4049 Allen Station Rd. Augusta, GA 30906

Tel: 706 790-9115 Email: sales@medicom.com

A. R. Medicom Inc. (Taiwan) Ltd.

4F-3, No. 163, Sec. 1, Keelung Road Xinyi District, Taipei City, Taiwan

Tel: +886 2 2766 8655 Email: info@medicom-taiwan.com

Medicom Healthcare B.V.

Parallelweg 80a 3931 MT Woudenberg, The Netherlands

Tel: 011 (31-33) 286-6460 Email: nlsales@medicom.com

A.R. Medicom Inc. Healthcare (Shanghai) Ltd.

Rm. 503, 728 Xinhua Road, Changhing District Shanghai 200050, China

Tel: +86 21 5273 9363 Email: info@medicom-china.com