

Alaris™ Pump Module Set Compatibility Card

Alaris Pump Module Model 8100



WARNING:

Use only administration sets labeled for use with Alaris Pump module. Use of any set other than those stated may cause improper instrument operation, resulting in inaccurate fluid delivery or other potential hazard.

NOTE:

Not all sets are available in all countries.

Model	Administrative Set Description
10010453	GEM 20DP V/NV NTG 1.2MF 1SS DEHP FREE
10010454	GEM 20D V/NV NTG .2MF 1SS DEHP FREE
10010483	GEMINI™ SYRINGE VTD SET, DEHP FREE
10010541	GEM V/NV MC 20DP 4SS CKV
10010547	GEMINI 20D V/NV 1.2MF, NITRO TUBE, 1SS
10010761	GEM V/NV 20DP W/YELLOW STRIPE TBG
10010871	GEM V/NV MC 20DP 2SMS
10011301	GEM V/NV 20DP CKV W/YELLOW STRIPE TBG
10011462	GEM 60DP CKV 4VS
10011652	GEM V/NV 2SS
10012144	GEM 20DP V/NV CKV 3SS
10012182	GEM 10DP CKV 4WY STCK 4SS
10012283	GEM V/NV QS MC 60D 3SS 0.2MF
10012293	GEM V/NV QS MC 60D CKV 3SS
10012645	GEM IN-LINE DC 2QS CKV 1SS
10012802	GEM V/NV 20DP CKV MNF 3SS
10013034	GEM V/NV 20DP CKV 3SS 2G 4W STPCK
10013037	GEM NV BLD 1SS
10013072	GEM 20DP CKV 2SS LOW SORB DEHP FREE
10013186	GEM V/NV 60DP CKV 3SS
10013361	GEMINI INFUSION SET, 2SS Y & BAG PORT
10013854	GEM W/SS BAP 20 15MF LOW SORBING

Model	Administrative Set Description
10013889	GEM 20DP CKV 3SS
10013890	GEM 20DP LOW SORB DEHP FREE
10014035	GEM 60DP MC CKV 4SS
10015012	GEM 20DP MC .2MF 2SS DEHP-FREE
10015294	GEM V/NV 2CKV 3SMS 20DP
10015312	GEM V/NV 20DP CKV 3SS 2G 3W STPCK
10015414	GEM NV BLD 15DP W/180 FLTR
10015489	GEM 20DP 3SS CKV DEHP FREE
10015645	GEM V/NV 20DP 2SS DEHP FREE
10015862	GEM 20DP W/SS DC
10015896	GEM 10DP CKV 3SS 2G 4WAY STCK
10031773	GEM SET 20DP 2CKV 3SS W/EXT SET 1SS
10062818	GEM NV BLD SET W/180 FLT DEHP-FFP
10137405	GEM 20DP CKV 3SS
10142568	GEM SET 20DP 4SS W/EXT 4G-4W STPCK
10178573	GEM 20DP V/NV CKV 4SS
10321213	GEM W/SS DC 20DP
10382963	GEM SET W/CV 1SS 3P-MNFLD & 3W-STPCK EXT
10404198	GEM 20DP CV 3SS DEHP-FREE
10406194	GEM 20DP CV 3SS & 0.2M FILTER DEHP-FREE
10408730	GEM V/NV 20DP CV 3SS
10561554	GEM 20DP W/15MF CKV 3SS
10779313	GEM 20DP V/NV MC 4VS DEHP-FREE
10813621	GEM 20DPCKV 4SS 3G 4-WAY STPCK DEHP FREE
10817613	GEM 10DP 2CKV 3SS DEHP FREE
10817920	GEM BURETTE 60DP PS 2SS DEHP FREE
10821753	GEM V/NV MC 4SS 20DP 10PK DEHP FREE
10863358	GEM V/NV 20DP WITH 15M FILTER NO PORTS

Model	Administrative Set Description
10879047	GEM V/NV 10DP 1CV 2SS DEHP FREE
10879391	GEM V/NV 20DP 2CV 3SS
10927766	GEM 20DP V/NV CKV 3SS
10933805	GEM 20DP V/NV CKV 2SS
10934345	GEM 20DP V/NV CKV 3SS 2G-3W STPCK
10942011	GEM V/NV 20DP DEHP FREE
11045328	GEM V/NV 20DP MC 0.2 MICRON FILTER
11059729	GEM SET CKV 20DP 4W STPCK 4SS
11063735	GEM V/NV 20DP MC 1SS
11118122	GEM V/NV MC 60DP 4SS 0.2MF
11130148	GEM 20DP V/NV 0.2MF 1VS
11133698	GEM V/NV 2VS 20DP W/15MF
11147991	GEM V/NV 1VS 20DP
11147992	GEM V/NV MNF 2VS 20DP
11171447	GEM V/NV 20DP 2CV 3SS DEHP FREE
11202653	GEM DUAL V/NV 20DP 0.2MF CKV VS
11258588	GEM HALF SET VS
11287205	GEM V/NV 20DP CKV 4SS DEHP FREE
11401542	GEM MC V/NV 20DP 4 VS
11404930	GEM V/NV 10DP CKV 3SS
11419365	GEM V/NV 20DP CKV 5SS 3P MNFLD & 3W STCK
11426964	GEM V/NV 20DP LOW SORBING NO PORTS
11426965	GEM V/NV 20DP CKV 5SS DEHP FREE
11484001	GEM V/NV 20DP CKV 1SS DEHP FREE
11522558	GEM V/NV 20DP PS CKV 3SS DEHP FREE
11532269	GEM V/NV 20DP CKV 2SS 0.2MF LOW SORBING
11582773	GEM V/NV 20DP CKV 3SS 2G-3WAY STPCK
11590100	GEM V/NV 60DP CKV 2SS DEHP FREE

Model	Administrative Set Description
11607704	GEM 20DP 2SS 120" DEHP FREE
11607787	GEM V/NV 60DP MC 2SS DEHP FREE
11609481	GEM 20DP CKV 0.2MF 2SS LS-SEG DEHP FREE
11612593	GEM 20DP NO PORTS 0.2MF
11613191	GEM BURETTE 60DP 0.2 MF 2SS DEHP FREE
10010454-0006	GEM 20D V/NV NTG .2MF 1SS DEHP FREE
10010541-07	GEM V/NV MC 20DP 4SS CKV DEHP Free
10013361-0006	GEM W/SS DC 20DP 1CV 2SS
10013361T	GEM INFUSION SET, 2SS & BAG PORT TEX
10013854-0006	GEM W/SS BAP 20 15MF LOW SORBING
10014855A	GEM DC 20D 2SS CHK VL AND DEHP FREE TBG
10015861A	GEM DC 20D 3SS 1 CHK VL & DEHP FREE TBG
10061661A	GEM DC 20D 4SS CHK VL AND DHEP FREE TBG
10074281A	GEM BURETE 60D 1SS DHEP FREE TBG
10321213T	GEMINI SMARTSITE™ DC 20DP & TEXIUM™
11063735-07	GEM V/NV 20DP MC 1SS DEHP Free
2106-0500	GEM 20DP 2VS 3 STPCK MANIFOLD
2107-0500	GEM 20DP V/NV CKV 4VS 3P-MNFLD 4W-STPCK
2110-0500	GEM V/NV 20 2 INJ STS SPLT SEP
2111-0500	GEM V/NV 2YS 60DP 20PK
2120-0007	GEM 20DP V/NV 1CKV 2SPLT SPTM DEHP FREE
2120-0500	GEM V/NV CV 2VS 20DP 20PK
2121-0500	GEM V/NV CKV .2MF 1VS 20D 20P
2122-0007	107-IN GEM V/NV 20-DP 3 VS CHECK VALVE
2123-0007	102-IN GEM V/NV 20DP 2VS CHCK VALVE
2125-0500	GEM V/NV CKV 2VS 20PS 20PK
2126-0007	GEM 20DP CKV 3 SPLIT SEPTUM PORTS DEHP F
2126-0500	GEM V/NV CKV 3VS 20DP 20PK

Model	Administrative Set Description
2129-0500	GEM 20DP V/NV CKV 3VS
2130-0500	GEM V/NV CKV .2MF 2VS 20D 20P
2131-0500	GEM V/NV 1.2MF 1VS 20D 20P
2140-0600	GEM V/NV MC 3VS 20D 10P
2141-0600	GEM V/NV MC 60DP
2147-0600	GEM V/NV QS MC 3VS 60PS 10P
2171-0007	108-IN GEM NV BLD DEHP FREE 2SPKS MF
2177-0000	GEM NV BLOOD SET SPLT SEPTUM
2200-0006	GEM V/NV 20DROP 20PACK
2200-0500	GEM V/NV PRIMARY 20DP
22000-B007T	ALARIS SYSTEM W/TEXIUM &BAG ACCESS PORT
22003-0004	ALS SET DC 20D DEHP FREE
2202-0007	GEM 20DP V/NV 1.2MF DEHP FREE
2202-0500	GEM V/NV 1.2MF 20DP 20PK
2203-0006	GEMINI LIGHT RESISTANT INFUSION SET
2203-0500	GEM V/NV 20DP LGT SN STV
2204-0007	GEM V/NV 20DP DEHP-FREE NO PORTS
2205-0000	GEMINI ENTERAL N/V SPIKE IN-LINE D/CH
2206-0007	GEM 20DP YELLOW STRIPED TUBING DEHP FREE
2210-0006	GEM PRIMARY V 24HR EUROPE
2210-0008	GEM V/NV 2YS 20D 20P
2210-0500	GEM V/NV PRI W/2 INJ SITES 20
2211-0500	GEM V/NV 2YS 60D 20P
2214-0006	GEMINI INFUSION SET, 2 Y, 15ÁM
2220-0006	GEMINI V/NV CKV 2YS 20DROP 20PACK
2220-0500	GEM V/NV CKV 2YS 20DP 20PK
2230-0006	GEM V/NV CKV .2MF 20DP 20PK
2230-0008	GEM V/NV CKV .2MF 20DP 20PK

Model	Administrative Set Description
2230-0500	GEM V/NV CKV .2MF 20DP 20PK
22300-0004	ALS SET DC 20D .2 MICFILTER DEHP FREE
2231-0006	GEMINI FILTER SET, 1.2µM, 1 Y, PE
2231-0008	GEM V/NV NTG 1.2MF 1YS 20D 20P
2231-0500	GEM V/NV NTG 1.2MF 1YS 20D 20P
22401-0007	GEM BURETTE 60DP PS .2MF
2241-0006	GEM V/NV METERED CHAMBER 60DROP
2241-0008	GEM V/NV MC 60D 10P
2241-0600	GEM V/NV MC 60D 10P
2260-0006	GEMINI NITRO SET 24HR EUROPE
2260-0008	GEM V/NV NTG 20DP 20PK
2260-0500	GEM V/NV NTG 20DP 20PK
22600-0007T	47" ALARIS SYSTEM 20DP W/TEXIUM
22601-B007T	86" ALARIS 20DP W/TEXIUM BAG ACCESS PRT
22602-B007T	118"DED-SET W/TEXIUM AND BAG ACCES PORT
22603-B007T	47"DED-SET 20D W/TEXIUM & BAG ACCESS VLV
22604-B007	118" ALARIS 20DP W/BAG ACCESS PORT
2264-0006	GEM NV/V NTG 0.2MF 1YS 20D 20P
2264-0008	GEM V/NV NTG .2MF 1YS 20D 20P
2264-0500	GEM V/NV NTG .2MF 1YS 20D 20P
2270-0500	GEMINI BLOOD SET 15DP NV 180MF NO PORTS
22700-0004	ALS BLOOD SET DC 20D 200MIC FILTER DEHP
2277-0006	GEMINI BLOOD SET 24HR EUROPE
2277-0008	GEM NV BLD 1YS 15DP 10PK
2278-0500	GEMINI BLOOD SET 20DP NV 200MF
2279-0006	GEMINI BLOOD SET, 1 SPIKE, DEHP FREE
2280-0000	GEMINI SYRINGE ADSET PRI 20PACK
2280-0006	GEMINI SYRINGE ADAPTER SET

Model	Administrative Set Description
2280-0008	GEM V SYR 20PK
24001-0007	47" GEM SET 20DP CKVLV NDLFREE SS VLV
24006-0004	ALS SET DC 20D 15MICF 3SS CKV DEHP FREE
24007-0004	ALS SET DC 20D 15MICF CKV 5SS DEHP FREE
24008-0007	ALS SET 20D V/NV DC CHCK VLV 2SS 3P/MANI
24009-0007T	ALRS SET 20D CHCK VLV 3 SS & TXUM DEHP F
2401-0006	GEM V/NV 1SMS INTL 20DP 20PK
2401-0500	GEMINI SET W/ONE SS VLV PORT
24010-0007T	ALRS SET 20D CHCK VLV 3SS & TXUM DEHP F
2403-0000	GEMINI HALF SET 1SS 20PK
2403-0007	GEM HALF SET 1SS 87IN DEHP FREE
2406-0500	GEM V/NV 2SMS 20DP 20PK
2408-0500	GEM 10DP CKV 4WAY-STPCK 4SS
2409-0500	GEM SET 20DP CKV 4SS 4STPCK AND EXT SET
2410-0500	GEMINI INFUSION SET, 2SS Y
2411-0500	GEM V/NV SS 60DP 20PK
2412-0500	GEM 20DP .2MF 1SS PORT 113IN
2413-0007	GEM20DP CKV 5SS 4STPCK4-GANGED DEHPFREE
2414-0006	GEMINI INFUSION SET, 2 SS Y
2414-0500	GEM 20DP 15MF 2SS
2416-0007	GEM 20DP V/NV 2CV 3SS DEHP FREE
2419-0007	GEM 20DP V/NV CKV 3SS
2420-0006	GEM V/NV CKV 2SMS INTL 20DP 20PK
2420-0007	GEM V/NV 20D 1CV 2SS DEHP FREE
2420-0500	GEM V/NV 20DP CKV 2SS 117-IN 20PK
2421-0500	GEM V/NV 2SMS 60DP 20PK
2422-0007	GEM 20DP V/NV CV 4SS 4WY STCK DEHP FREE
2423-0007	GEM 10DP CKV 4SS 2G 4WAY STCK DEHP FREE

Model	Administrative Set Description
2424-0007	GEM 20DP V/NV CKV 4SS DEHP FREE
2425-0500	GEM V/NV CKV 2SS 20PS 20PK
2426-0007	GEM 20DP CKV 3SS DEHP FREE
2426-0500	GEM V/NV CKV 3 SS 20DP 20PK
2428-0007	GEM 20DP V/NV 2CKV 4SS DEHP FREE
2429-0007	GEM 20DP V/NV 1CKV 4SS DEHP FREE
2429-0500	GEM 20DP V/NV CKV 3SS
2430-0500	GEM V/NV CKV .2MF 2SMS 20D 20P
24301-0007T	GEM 20DP 0.2MF LOW SORB DEHP FREE TEX
24302-0004	ALS SET DC 20D 15MICF 2SS DEHP FREE
2432-0007	GEM 20DP V/NV CKVLV 3SS 0.2MF DEHP FREE
2433-0007	GEM 20DP V/NV 1CKV 2SS DEHP FREE
2434-0007	GEM SET 20DP DEHP FREE CKV .2 FLTR 1SS
2435-0007	GEM 20DP V/NV 1CKV 4SS DEHP FREE
2436-0500	GEM 20 DP V/NV 2 CKV AND 5SS
2438-0500	GEMINI 10DP CKV 4SS DEHP
2440-0006	GEM V/NV MC 4SMS 20DP 10PK
2440-0500	GEM V/NV CKV 4SS 20DP 2GNGD-4WSTPCK
2440-0600	GEMINI MC W/3SS VLV PTS
2441-0006	GEM V/NV MC 4SS PORTS 60DP 10PK
2441-0007	GEM V/NV MC 4SS 60DP DEHP FREE
2441-0600	GEM V/NV MC 3 SMS 60DP 10PK
2442-0007	GEM V/NV MC VS 60 DP 3SS DEHP FREE
2443-0600	GEM V/NV MC 0.2MF 3SMS 60DP 10PK
2444-0007	GEMINI 20DP BURETTE 1SS LOW SORBING
2447-0007	GEM BURETTE V/NV 60DP PS 3SS DEHP FREE
2447-0600	GEM V/NV QS MC 3SMS 60PS 10PK
2448-0007	GEM V/NV MC 3SS PORTS 20DP10PK DEHP Free

Model	Administrative Set Description
2448-0600	GEM V/NV MC 3SS PORTS 20DP 10PK
2449-0600	GEM V/NV MC 3SS PORTS 60DP 10PK
2450-0007	GEM60DP CV 3SS 3WSTCK 2 1WAYVLV DEHPFREE
2450-0500	GEMINI 10DP CKV 3SS 2G-3W STOPCOCK
2451-0007	GEM V/NV 20DP 5SS 3CV DEHP FREE
2452-0007	GEM 10DP 3SS VLV 2CK VLV, DEHP-FREE
2455-0500	GEM V/NV SS 20DP 20PK
24600-0007	ALS SET DC 20DP 3SS CKVLV LOW SRB DEHP F
24601-B007T	118" DED SET 20DP W/TEXIUM CHCK VLV 2SS
24602-B007T	118" ALARIS 20DP W/0.2 FILTER TEXIUM 2SS
24603-0004	ALS SET DC 20D 15MICF 2SS CHK-VLV DEHP F
24604-0004	ALS SET DC 20D 1.2MICF DEHP FREE
2461-0007	GEMINI 20DP CKV 3SS DEHP FREE
2462-0007	GEM HALF SET 1SS LOW SORBING DEHP FREE
2463-0007	GEM V/NV 20DP 1SS DEHP FREE
2465-0007	GEM 20DP 0.2MF LOW SORBING DEHP FREE
2466-0007	GEM LOW SORB W/SS 20DP DEHP-FREE
2477-0000	GEM NV BLD 1SMS 10PK
2477-0006	GEM NV BLD 1SMS 10PK
2477-0007	GEM NV BLD 1SS DEHP FREE
2478-0000	GEM NV BLD 1SMS BLUNT SPK 10PK
2929-0500	GEM V/NV 1VS 20DP 20PK
2941-0500	GEM V/NV 60DP 20PK
2942-0500	GEM V/NV 20DP 20PK
9943-0008	GEM 1YS 20PK
9943T	GEM HALF SET LOW VOL MICROBORE 84-IN
C24101E	GEM DC 20DP W/15MF .2MF 1SS AMBER TBG
C24103E	GEM 20DP W/CV 3SS PORTS & 0.2M FILTER

Model	Administrative Set Description
C24104E	GEM 20DP CKV 4SMS 1STPCK
C24105E	GEM V/NV CKV 5SS 20DP 1 ROLLER CLAMP
C24106E	GEM 20DP CKV 2SMS 2STPCK T-CONCTR
C24107E	GEM 60DP CKV 2SMS 2STPCK T-CONNECTOR
C24111E	GEM 20DP 2SS CKV W/33 IN /10 IN EXT 1SS
C24112	GEM SP HALF SET
C24112E-0007	GEM60DP 2YS 3PRT MFLD 3WAYSTCK DEHPFREE
C24116E	GEM V/NV CKV 2 SMS 20 DP 1 SLD CLP
C24117	GEM V/NV NTG 20DP 1 SLD CLP
C24119E	GEM V/NV NTG .2MF 1SS 20DP 2SLD

© 2008, 2016 CareFusion or one of its affiliates. All rights reserved. Alaris, Gemini, Texium, and SmartSite are trademarks or registered trademarks of CareFusion or one of its affiliates. All other trademarks are the property of their respective owners.