

INSTRUCTIONS FOR USE

Treatment Accessories



Cup Dispenser

Instructions for use:

1. Pull off top of the dispenser.
2. Load cups with opening pointing up—dispenser holds 50 cups.
3. Replace top.



Ref.#'s PCC



Not made with natural rubber latex.



CROSSTEX[®]
A CANTEL MEDICAL COMPANY

Manufactured for: Crosstex International, Inc.
10 Ranick Road | Hauppauge, NY 11788 | 631.582.6777 | 888.276.7783 | Crosstex@crosstex.com

Made in China

PCC-IFU
Rev. 0 – 11/16