



STERISIL<sup>®</sup> *Safe, Simple, Sterisil*


Citrisil<sup>™</sup>

Cleans and inhibits bacteria growth in dental unit waterlines


## THE ONLY ALL-IN-ONE TABLET

### PRODUCT HIGHLIGHTS

- Cleans and maintains dental unit waterlines to  $\leq 10$  CFU/ml
- Shock tablet included to kill pre-existing odor causing bacteria
- Eliminates the need to empty bottles and purge lines at the end of each day
- Safe for dental equipment, patients, and clinicians
- Does not adversely affect composite bond strength
- Maintains antimicrobial tubing
- Does not contain oxidizers like iodine or chlorine
- Available for .70 to 1L bottles & 2L bottles




White Tablet



Blue Tablet

**48 MAINTENANCE  
TABLETS**



Shock Treatment

**2 SHOCK TABLETS**  
INCLUDED WITH EACH BOX

Visit us online at [Sterisil.com](http://Sterisil.com)





**STERISIL®** *Safe, Simple, Sterisil*



## PROTECT YOUR PATIENTS AND STAFF AGAINST BACTERIA

Citrisil™ is the only all-in-one shock and daily maintenance dental waterline tablet. Use the orange shock tablet in your lines monthly; allowing it to work overnight to kill pre-existing odor causing bacteria. Use the white or blue maintenance tablet in your bottle daily to continuously clean and maintain dental waterlines. Drop in a new maintenance tablet every time a bottle is refilled. The tablet effervesces and continually stays in the lines to prevent bacteria growth while also being safe for patients and staff.


Part #	White Maintenance Tablets
<b>C50-W</b>	(48) White Maintenance Tablets & (2) Shock Tablets, for .70–1L Bottles
<b>C50-2W</b>	(48) White Maintenance Tablets & (2) Shock Tablets, for 2L Bottles
<b>CS-50-W</b>	Case of (50) C50-W Boxes, White Tablets, for .70–1L Bottles
<b>CS-50-2W</b>	Case of (50) C50-2W Boxes, White Tablets, for 2L Bottles
Part #	Blue Maintenance Tablets
<b>C50-B</b>	(48) Blue Maintenance Tablets & (2) Shock Tablets, for .70–1L Bottles
<b>C50-2B</b>	(48) Blue Maintenance Tablets & (2) Shock Tablets, for 2L Bottles
<b>CS-50-B</b>	Case of (50) C50-B Boxes, Blue Tablets, for .70–1L Bottles
<b>CS-50-2B</b>	Case of (50) C50-2B Boxes, Blue Tablets, for 2L Bottles
Part #	Shock Tablets
<b>ECS-20</b>	(20) Shock Tablets
<b>CS-50-ECS20</b>	Case of (50) ECS-20 Boxes



Proudly made  
in the USA.



0120  
Sterisil, Inc.  
835 S. Hwy 105 Ste. D  
Palmer Lake, CO USA 80133  
+1 719 622 7200



EMERGO EUROPE  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands



For more information contact your **Sterisil® sales representative or local dealer.**

719 622 7200 or online at [Sterisil.com](http://Sterisil.com)



GP1005.B