

Econo Rolls



DIAMETER:
Medium

LENGTH:
1-1/2"

Satisfaction in Savings. Richmond Econo Rolls are for practices whose first consideration is price. Made with 100% cotton with a methylcellulose binder, the absorbent roll is smooth and rigid.

- Economical
- 100% cotton
- Smooth and rigid
- Non-Sterile
- Durable infection control roll dispensers are available; they are fully autoclavable and will keep rolls sterile until use; available in blue and white.



CODE	PRODUCT	DESCRIPTION	PACKAGING	CASE COUNT	ORIGIN
□ 216206	1-1/2"	MED. DIA. (3/8") - NON-STERILE	2000/BOX	12 BOXES	USA
□ 200224	INFECTION CONTROL ROLL DISPENSER	BLUE - AUTOCLAVABLE	1 DISP + 200 ROLLS*	6 DISPENSERS	USA
□ 200424	INFECTION CONTROL ROLL DISPENSER	WHITE - AUTOCLAVABLE	1 DISP + 200 ROLLS*	6 DISPENSERS	USA

*Dispensers are shipped with Braided Rolls only.

All products must be ordered in full case quantities. Please contact Richmond Dental customer service for information about quantity discounts and product category discounts.