

EVA

Digital X-Ray System

- Exceptional Image Quality
- Ideal For Any Practice
- Full Integration
- Best Supported Sensor
- Great Value








Go Digital With EVA

Affordable Digital Value

Reliable. EVA digital x-ray system delivers big value. Backed by the proven reliability of tens of thousands of installations ranging from single chair clinics to dental chains, EVA offers a rugged sensor design coupled with the highest standard of technology.

Functional. Just a few millimeters thin, EVA delivers crisp, clear radiographs with virtually any intraoral x-ray head. The CMOS imaging chip is designed for maximum durability, and the sensor sizes are nearly identical to Pedo and Adult size films. Fast USB connectivity delivers images to the computer monitor in seconds.

Value. EVA's value makes it an ideal choice for both large dental chains and value-conscious single practices. EVA comes with EVAsoft software, unlimited licenses, a 2-year warranty, lifetime free technical support, positioning and training videos and the EVA Loyalty Program.



EVA Digital Intraoral X-Ray System

- ✓ Affordable Digital Technology
- ✓ More Than 10 Years of EVA Imaging
- ✓ Consistent Image Capture, Fewer Retakes
- ✓ Integrates With All Major Practice Management Software
- ✓ TWAIN Driver Integrates With Most Imaging Software
- ✓ Auto-Save Stores Images For You
- ✓ Attach EVA X-Rays Directly Into Electronic Claims

The Dent-X Advantage

- ✓ Free Lifetime Support and Log-On Help
- ✓ Standard Two-Year Manufacturer's Warranty
- ✓ Service Multiple Operatories With One Sensor
- ✓ Made in the USA



EVA Digital X-Ray Systems

EVA systems are available in combo packs with 1, 2, or 3 sensors. Create the system to suit your practice with any size sensor combination.


Each EVA system comes with:

- ✓ One EVA USB Station
- ✓ One Box of Protective Sheaths
- ✓ RINN Positioning Aids
- ✓ EVAsoft Software
- ✓ 2-Year Warranty
- ✓ Free Lifetime Support
- ✓ In-Office Training
- ✓ Twain Driver

EVA Universal TWAIN System

- ✓ Make EVA work with your existing imaging software
- ✓ Digital solution
- ✓ Always free Lifetime Support
- ✓ EVA's TWAIN driver allows EVA sensors to acquire in most imaging software systems

Maximize Your ROI



With the purchase of one EVA sensor and a few inexpensive EVA USB Stations, EVA can easily service multiple operatories. Take advantage of the value and capability that EVA offers your entire practice.

EVAsoft Software



EVAsoft Dental Image Management software meets the imaging demands of today's digital practices. Every EVA System ships with EVAsoft Imaging Software.

EVAsoft's many features make it the best choice for imaging software. EVAsoft easily integrates with every popular Practice Management Software, and useful features such as Keep-On-Top allows viewing x-rays while simultaneously editing a patient's chart. The ProIntegrate feature links to share patient information. A sophisticated software that is easy to learn and simple to operate, EVAsoft has all of the advanced features necessary to deliver rapid, effective diagnoses from multiple imaging modalities.


From the simplicity of navigation to the outstanding value it offers, EVAsoft Dental Image Management software is the choice for any practice!

EVAsoft:

- ✓ EVA Digital X-Ray Sensor, Digital Pan & Intraoral Camera Software Suite
- ✓ Unlimited Multi-User Software per Network
- ✓ Stand-Alone or Network Environments
- ✓ Integrates with All Major Practice Management Software using ProIntegrate
- ✓ Developed and Supported in America, by ImageWorks

The EVAsoft Experience:

- ✓ Easy-to-Use with Simple Tab Navigation
- ✓ Automatic X-Ray Sequence Selection
- ✓ Automatic Image Save Function
- ✓ Endodontic Exposure Mode
- ✓ Decay and Abscess Detection Tool
- ✓ Automatic Image Contrast Enhancement
- ✓ View X-Rays While Modifying Practice Management Charts
- ✓ Easily Send X-Rays to Insurance Companies and Referring Physicians



EVA System Specifications


Specifications	EVA Size #2	EVA Size #1
Part number	9992410200	9992410100
Active area	1" x 1.4" (25.8 x 36 mm)	0.8" x 1.2" (20 x 30 mm)
Capsule dimensions	1.2" x 1.74" x 0.2" (30.8 x 44 x 5 mm)	1" x 1.5" x 0.2" (25.6 x 38.4 x 5 mm)
Image size (unpacked)	2.1 MB	1.3 MB
Sensor type	CMOS	
Gray shades	4096 (12 bits)	
Actual imager resolution	16.7 lp/mm	
Pixel size	30 µm x 30 µm	
Sensor cable length	6 ft. (2 meters)	
Power supply	USB powered device—no external power supply needed	
Interface	USB	
Shipping weight	3.7 lbs.	
Warranty	2 Years Parts and Factory Service in the United States and Canada	
Safety certifications	FDA, CSA, (US and Canada) CE, IEC 60601-1-2	
Recommended EVA computer requirements*	Windows 7 (32 or 64 bit) (Professional or Ultimate), Core 2 Duo Processor (equivalent or better); 4+ GB RAM; USB 2.0 PORTS; DVD-RW; 250 GB Hard Drive; High Resolution LCD Monitor with DVI connection.	
Recommended EVAsoft computer requirements*	Windows 7 (32 or 64 bit); 4+ GB RAM; USB 2.0 PORT1; CD-R / DVD-RW; Hard Drive - 150 GB; Monitor with 1000:1 Contrast Ratio and 0.255 dot pitch; Minimum resolution of 1600x1200 @ 32 bits color depth.	



Copyright © 2010 ImageWorks. All rights reserved.
Reproduction of the whole or any part of the contents without written permission is prohibited.
EVA: Registered Trademark of ImageWorks. Printed in the U.S.A.

SM-D155-06

*Computer not included. While not required, Internet access will greatly improve Technical Support assistance. It is highly recommended that you consult your IT professional for a Data Back-up system. Specifications are subject to change without notice.



ImageWorks
250 Clearbrook Road, Suite 240
Elmsford, NY 10523 USA
www.ImageWorksCorporation.com

1.800.592.6666
P 914.592.6100
F 914.592.6148

