

TECHNICAL DATA SHEET

Medicom SafeGauze Non-Woven Sponge Items # 4522, 4522-5 and 4544

Intended Use

A hydrophilic non-sterile sponge intended to cover a wound, such as to be placed directly on a patient's wound to absorb exudates.

Characteristics	Specification			
Description	Soft non-woven material, folded			
Material Composition	65% Viscose / 35% Polyester (Latex free)			
Color	White			
Ply	4 ply			
Dimensions	Size	Item Code	Width & Length (cm)	
	2" x 2"	4522, 4522-5	5,1	
	4" x 4"	4544	10,2	
Packaging	Item Code	Sponges/Sleeve	Sleeves/Case	Sponges/Case
	4544	200	10	2000
	4522		20	4000
	4522-5		25	5000

Item Code	Description	Color	Specific Attribute
4522	Medicom SafeGauze Non-Woven Sponges, 4 Ply, 2" x 2" (5 x 5 cm), 200 Sponges/Sleeve, 20 Sleeves/Case, 4000 Sponges/Case	White	Non-Sterile
4522-5	Medicom SafeGauze Non-Woven Sponges, 4 Ply, 2" x 2" (5 x 5 cm), 200 Sponges/Sleeve, 25 Sleeves/Case, 5000 Sponges/Case	White	Non-Sterile
4544	Medicom SafeGauze Non-Woven Sponges, 4 Ply, 4" x 4" (10 x 10 cm), 200 Sponges/Sleeve, 10 Sleeves/Case, 2000 Sponges/Case	White	Non-Sterile

Visual

Inner Sleeve



TECHNICAL DATA SHEET**Medicom SafeGauze Non-Woven Sponge
Items # 4522, 4522-5 and 4544****Regulatory**

Health Canada Class 1
FDA Class 1, Product Code NAC

Storage Conditions

Store under cool dry place.

Country of Origin

USA