



Headrest Covers

Poly Coated – manufactured with a combination of soft facial tissue and polypropylene. Fluid resistant. Latex-free.

FEATURES

Soft facial tissue

Polypropylene

BENEFITS

Comfort, absorbency

Strong, fluid resistant

Plastic – clear soft polyethylene with a 3" gusset to lie flat on the headrest. Latex-free.

FEATURES

Clear, soft polyethylene

1 mil thick

3" gusset

BENEFITS

Comfort

Strength

Cover lies flat on headrest

CODE	DESCRIPTION	COLOR	QTY
Poly Coated Headrest Covers			
LOWHA	10" x 10"	White	500/case
L3WHA	10" x 13"	White	500/case
Plastic Headrest Covers			
L0CPA	9.5" x 11"	Clear	250/box
L3CPA	9.5" x 14"	Clear	250/box

Distributed by: