

TECHNICAL DATA SHEET

Medicom SafeBasics Tissue/Poly Head Rest Covers Items # 3014 to 3021

Intended Use

A headrest cover is intended to be mounted to the top of a dental chair backrest and helps keep the dental chair sanitary and dry. The fluid repellent material protects patients against cross-contamination.

| Characteristic | | Specification | | | |
|----------------|-------------|---|----------------------|---------------------|-------------------|
| Description | | Flat pouch closed on three sides made of embossed material composed of tissue with plastic backing. Open side borders are uneven. | | | |
| Material | Composition | Tissue Layer : 100% Wood Pulp Poly Backing : Polyethylene | | | |
| Dimensions | | Size | Interior Length (cm) | Interior Width (cm) | Extra Length (cm) |
| | | 10" x 10" (25 x 25 cm) | 25,4 | 25,4 | 0,5 |
| | | 10" x 13" (25 x 33 cm) | 25,4 | 33,0 | 0,5 |
| Packaging | | 500 Head Rest Covers/Case | | | |

| Item Code | Color | Size |
|-----------|------------|-----------|
| 3014 | White | 10" x 10" |
| 3015 | Blue | |
| 3016 | Dusty Rose | |
| 3020 | Lavender | |
| 3017 | White | 10" x 13" |
| 3018 | Blue | |
| 3019 | Dusty Rose | |
| 3021 | Lavender | |

| Item Code | Description | Specific Attributes |
|-----------|--|-------------------------|
| 3014 | Medicom SafeBasics Tissue/Poly Head Rest Covers, 10" x 10" (25 x 25 cm), 500 Head Rest Covers/Case | White, Non-Sterile |
| 3015 | Medicom SafeBasics Tissue/Poly Head Rest Covers, 10" x 10" (25 x 25 cm), 500 Head Rest Covers/Case | Blue, Non-Sterile |
| 3016 | Medicom SafeBasics Tissue/Poly Head Rest Covers, 10" x 10" (25 x 25 cm), 500 Head Rest Covers/Case | Dusty Rose, Non-Sterile |
| 3017 | Medicom SafeBasics Tissue/Poly Head Rest Covers, 10" x 13" (25 x 33 cm), 500 Head Rest Covers/Case | White, Non-Sterile |
| 3018 | Medicom SafeBasics Tissue/Poly Head Rest Covers, 10" x 13" (25 x 33 cm), 500 Head Rest Covers/Case | Bleu, Non-Sterile |
| 3019 | Medicom SafeBasics Tissue/Poly Head Rest Covers, 10" x 13" (25 x 33 cm), 500 Head Rest Covers/Case | Dusty Rose, Non-Sterile |
| 3020 | Medicom SafeBasics Tissue/Poly Head Rest Covers, 10" x 10" (25 x 25 cm), 500 Head Rest Covers/Case | Lavender, Non-Sterile |
| 3021 | Medicom SafeBasics Tissue/Poly Head Rest Covers, 10" x 13" (25 x 33 cm), 500 Head Rest Covers/Case | Lavender, Non-Sterile |

TECHNICAL DATA SHEET**Medicom SafeBasics Tissue/Poly Head Rest Covers
Items # 3014 to 3021****Visual****Regulatory**

FDA Class 1, Product Code NRU
Health Canada – Not regulated

Storage Conditions

Store in a cool dry place

Country of Origin

China