

GUARDING AGAINST INFECTION™

TIDI® focuses on you, the healthcare professional, continually evaluating your needs. So we understand that it's imperative to create an exam environment that presents a professional image and makes patients feel comfortable. It's essential that exam products be easy and convenient to use. And it's critical to minimize any chance for cross contamination.

In response, we engineer specific single-use infection prevention solutions to meet these needs. You can be confident in our reliable, consistent product quality, assured by our certification in ISO 9001 and ISO 13485.

With these assurances, healthcare professionals confidently reach for TIDI's® single-use infection control products to maintain a hygienic environment, put staff and patients at ease, and focus on care.



Clear Choices for Confident Care

Our simplified product organization includes three distinct levels that clearly describe our differentiated, value-added products, based on insight gleaned from healthcare professionals. This makes it easy to choose the best TIDI® products for your procedures, your environment, and the comfort and convenience of you and your patients.

PRODUCT LEVELS



High-performance products designed with thoughtful features for the highest-quality care, making the professional's job easier while enhancing the patient experience.



A wide variety of convenient, highquality products that offer comfort and protection on demand.



Economical single-use products that provide functional everyday protection.

TABLE OF CONTENTS

2–7	Bibs + Towels					
8-9	Headrest Covers					
10–11	Tray Covers					
12–15	Sheaths + Barriers					
16	Specialty Gloves					
17	Eye Shields + Face Shields					
18-19	Facemasks					
20-21	Cups + Saliva Ejectors					
22	Cotton + Gauze + Nonwovens					
23	Sterilization + Autoclave					
24	Disposal + Office Sundry					
25	Glossary					



Sustainable Solutions

TIDI® Products is focused on reducing resource consumption, minimizing environmental impact and acting in an environmentally responsible manner, in order to provide effective solutions without compromising the well-being of future generations.

Product Samples

Experience the softness of a headrest cover, the absorbency of a towel, the convenience of a barrier, or the performance of any TIDI® product with a no-charge sample pack. **Call 1.800.521.1314 or visit TIDIProducts.com.**

Online Information


Explore our broad range of product areas and learn more about how our solutions can enhance the care experience. Visit **TIDIProducts.com**.

Value-Added Packaging

Newly designed packaging provides clear, concise, instant communication. Quickly and simply verify product identity and description for ease in ordering and selecting the right TIDI® products for your needs.







Our most absorbent choices are excellent for prophy and minor procedures BIBS + TOWELS when you need to protect patients from fluids.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
919401	White	13" x 18"	500	3-ply tissue, poly-back, diamond-embossed
919402	Green	13" x 18"	500	3-ply tissue, poly-back, diamond-embossed
919403	Blue	13" x 18"	500	3-ply tissue, poly-back, diamond-embossed
919404	Yellow	13" x 18"	500	3-ply tissue, poly-back, diamond-embossed
919405	Gray	13" x 18"	500	3-ply tissue, poly-back, diamond-embossed
919406	Mauve	13" x 18"	500	3-ply tissue, poly-back, diamond-embossed
919407	Beige	13" x 18"	500	3-ply tissue, poly-back, diamond-embossed
919408	Peach	13" x 18"	500	3-ply tissue, poly-back, diamond-embossed
919410	Teal	13" x 18"	500	3-ply tissue, poly-back, diamond-embossed
919400	Lavender	13" x 18"	500	3-ply tissue, poly-back, diamond-embossed
917401	White	13" x 18"	500	3-ply tissue, poly-back, waffle-embossed
917402	Green	13" x 18"	500	3-ply tissue, poly-back, waffle-embossed
917403	Blue	13" x 18"	500	3-ply tissue, poly-back, waffle-embossed
917404	Yellow	13" x 18"	500	3-ply tissue, poly-back, waffle-embossed
917405	Gray	13" x 18"	500	3-ply tissue, poly-back, waffle-embossed
917406	Mauve	13" x 18"	500	3-ply tissue, poly-back, waffle-embossed
917407	Beige	13" x 18"	500	3-ply tissue, poly-back, waffle-embossed
917408	Peach	13" x 18"	500	3-ply tissue, poly-back, waffle-embossed
917410	Teal	13" x 18"	500	3-ply tissue, poly-back, waffle-embossed
917400	Lavender	13" x 18"	500	3-ply tissue, poly-back, waffle-embossed
917411	White	17" x 18"	500	3-ply tissue, poly-back, waffle-embossed
917413	Blue	17" x 18"	500	3-ply tissue, poly-back, waffle-embossed







Choice BIBS + TOWELS

Choice
The right bibs to choose for reliable patient protection, suited for exams, consultations and x-rays.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
919461	White	13" x 18"	500	2-ply tissue, poly-back, diamond-embossed
919462	Green	13" x 18"	500	2-ply tissue, poly-back, diamond-embossed
919463	Blue	13" x 18"	500	2-ply tissue, poly-back, diamond-embossed
919464	Yellow	13" x 18"	500	2-ply tissue, poly-back, diamond-embossed
919465	Gray	13" x 18"	500	2-ply tissue, poly-back, diamond-embossed
919466	Mauve	13" x 18"	500	2-ply tissue, poly-back, diamond-embossed
919467	Beige	13" x 18"	500	2-ply tissue, poly-back, diamond-embossed
919468	Peach	13" x 18"	500	2-ply tissue, poly-back, diamond-embossed
919470	Teal	13" x 18"	500	2-ply tissue, poly-back, diamond-embossed
919459	Lavender	13" x 18"	500	2-ply tissue, poly-back, diamond-embossed
917461	White	13" x 18"	500	2-ply tissue, poly-back, waffle-embossed
917462	Green	13" x 18"	500	2-ply tissue, poly-back, waffle-embossed
917463	Blue	13" x 18"	500	2-ply tissue, poly-back, waffle-embossed
917464	Yellow	13" x 18"	500	2-ply tissue, poly-back, waffle-embossed
917465	Gray	13" x 18"	500	2-ply tissue, poly-back, waffle-embossed
917466	Mauve	13" x 18"	500	2-ply tissue, poly-back, waffle-embossed
917467	Beige	13" x 18"	500	2-ply tissue, poly-back, waffle-embossed
917468	Peach	13" x 18"	500	2-ply tissue, poly-back, waffle-embossed
917470	Teal	13" x 18"	500	2-ply tissue, poly-back, waffle-embossed
917459	Lavender	13" x 18"	500	2-ply tissue, poly-back, waffle-embossed
917471	White	17" × 18"	500	2-ply tissue, poly-back, waffle-embossed
917473	Blue	17" x 18"	500	2-ply tissue, poly-back, waffle-embossed







Choice

BIBS + TOWELS TIDI® offers a variety of relaxing and soothing prints for patients of all ages.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
983813	Cat in the Hat [†]	13" x 18"	500	Specialty printed 2-ply tissue, poly-back, waffle-embossed
981410	Bugs & Things	10" x 13"	250	Specialty printed 2-ply tissue, poly-back, waffle-embossed
981413	Bugs & Things	13" x 18"	250	Specialty printed 2-ply tissue, poly-back, waffle-embossed
981913	Schooltime	13" x 18"	500	Specialty printed 2-ply tissue, poly-back, waffle-embossed
917490	TIDI® Tooth	13" x 18"	500	Specialty printed 2-ply tissue, poly-back, waffle-embossed
981610	Under the Sea	10" x 13"	250	Specialty printed 2-ply tissue, poly-back, waffle-embossed
981613	Under the Sea	13" x 18"	500	Specialty printed 2-ply tissue, poly-back, waffle-embossed
982413	Watercolors	13" x 18"	500	Specialty printed 2-ply tissue, poly-back, waffle-embossed
982433	Watercolors	13" x 18"	500	Specialty printed 3-ply tissue, poly-back, waffle-embossed
983910	Horton Hears a Who [†]	10" x 13"	250	Specialty printed 2-ply tissue, poly-back, waffle-embossed
983913	Horton Hears a Who [†]	13" x 18"	500	Specialty printed 2-ply tissue, poly-back, waffle-embossed

 $^{^{\}dagger}\text{Dr.}$ Seuss Properties TM & @ 2004 Dr. Seuss Enterprises, L.P. All Rights Reserved.



CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
7400	Clear	7" x 5" x 6"	1	Acrylic countertop dispenser fits 13" x 18" bibs and towels



Convenient acrylic countertop dispenser fits all TIDI® 13" x 18" bibs and towels.







The highest level of absorbency for procedures such as root canals, extractions, and ultrasonic scaling. The contour neck assures a perfect SPECIALTY BIBS fit and comfort for patients.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
917901	White	18" x 22"	400	Contour bib, 3-ply tissue, poly-back
917902	Green	18" x 22"	400	Contour bib, 3-ply tissue, poly-back
917903	Blue	18" x 22"	400	Contour bib, 3-ply tissue, poly-back
917906	Mauve	18" x 22"	400	Contour bib, 3-ply tissue, poly-back
917909	Lavender	18" x 22"	400	Contour bib, 3-ply tissue, poly-back
917441	White	17" x 18"	500	Contour bib, 3-ply tissue, poly-back
917443	Blue	17" x 18"	500	Contour bib, 3-ply tissue, poly-back
917677	Blue	18" x 25"	250	Oral surgery bib, 3-ply tissue, poly-back, shell pattern
917678	Mauve	18" x 25"	250	Oral surgery bib, 3-ply tissue, poly-back, shell pattern

SPECIALTY BIBS and comfort during exams.

High absorbency plus a contour neck provide patients with confidence

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
917341	White	17" x 18"	500	Contour bib, 2-ply tissue, poly-back
917343	Blue	17" x 18"	500	Contour bib, 2-ply tissue, poly-back
917344	Mauve	17" x 18"	500	Contour bib, 2-ply tissue, poly-back
917345	Lavender	17" x 18"	500	Contour bib, 2-ply tissue, poly-back





Everyday SPECIALTY BIBS A selection of bibs providing to help protect patients from and oral surgery procedures.

A selection of bibs providing substantial coverage and poly barriers to help protect patients from bleach and fluids during endodontic and oral surgery procedures.

CODE	COLOR/ PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
917673	Blue	18" x 25"	250	Chainless bib with ties, tissue/poly
917676	Mauve	18" x 25"	250	Chainless bib with ties, tissue/poly
917480	Lavender	29" x 21"	100	Endo bib, waist-length with contour neck, tissue/poly/tissue
917483	Blue	29" x 21"	100	Endo bib, waist-length with contour neck, tissue/poly/tissue
917486	Mauve	29" x 21"	100	Endo bib, waist-length with contour neck, tissue/poly/tissue
917493	Blue	29" x 42"	50	Oral surgery bib, knee-length with contour neck, tissue/poly/tissue
917670	White	36" x 40"	100	Patient apron, cutout neck, poly



BIB CHAIN STUDIES

The following are outcomes of several small studies recently conducted on bib holders. One study, conducted by Noel Kelsch, RDHAP, assessed the role bib holders play in harboring microorganisms. In comparing levels of bacterial found on used bib holders with those found in airport restrooms, similar levels of bacteria were observed - as well as increased levels of microbial contamination in the crevices of the bib holders, including underneath smooth tubing. Of note, the bib holders sampled were reported to have already been disinfected with a hospital grade disinfectant.¹

Kelsch's findings are enclosed in a study conducted by the University of North Carolina that led to the conclusion that bib holders are potential sources of contamination, with one in five bib clips sampled having significant levels of microorganisms. Importantly, bacterial growth on the bib holders was largely attributed to saliva and plaque as well as skin and flesh, and included E. coli and Staph aureus.²





917676 Tissue/poly bib with ties




range of colors to enhance any office décor.



¹ Kelsch N. Don't clip that crud on me. RDH, 2011

² Oral microbiology lab, Chapel Hill School of Dentistry, University of North Carolina. Data on file.



Facial grade tissue provides a pleasing appearance and our softest comfort, along with absorbency and strength. A poly HEADREST COVERS liner also prevents leak-through.



CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
919511	White	10" × 10"	500	FABRICEL® headrest cover, poly-lined
919512	White	13" x 10"	500	FABRICEL® headrest cover, poly-lined
919513	White	13" x 13"	500	FABRICEL® headrest cover, poly-lined
919514	White	15" x 11-1/2"	500	FABRICEL® headrest cover, poly-lined
919515	White	10" x 12", Chayes	500	FABRICEL® headrest cover, poly-lined

Soft tissue provides comfort for the patient and absorbency **HEADREST COVERS** for infection control. Poly-lined to prevent leak-through.




CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
919711	White	10" x 10"	500	Tissue/poly headrest cover
919712	Green	10" x 10"	500	Tissue/poly headrest cover
919713	Blue	10" x 10"	500	Tissue/poly headrest cover
919716	Mauve	10" x 10"	500	Tissue/poly headrest cover
919719	Lavender	10" x 10"	500	Tissue/poly headrest cover
919811	White	13" x 10"	500	Tissue/poly headrest cover
919812	Green	13" x 10"	500	Tissue/poly headrest cover
919813	Blue	13" x 10"	500	Tissue/poly headrest cover
919816	Mauve	13" x 10"	500	Tissue/poly headrest cover
919819	Lavender	13" x 10"	500	Tissue/poly headrest cover



Everyday HEADREST COVERS protection for equipment.

A reliable barrier offers peace of mind for patients and

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
919611	White	10" x 10"	500	All-poly headrest cover
919612	Clear	10" x 10"	500	All-poly headrest cover
919621	White	13" x 10"	500	All-poly headrest cover
919622	Clear	13" x 10"	500	All-poly headrest cover
919631	White	13" x 13"	500	All-poly headrest cover




Choice Use our tray also minimize on the tray.

Use our tray covers as a convenient infection prevention barrier. They also minimize noise and improve visibility of tools, burs and files resting on the tray.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
917511	White	8-1/2" x 12-1/4"	1000	Ritter (B) tray cover
917512	Green	8-1/2" x 12-1/4"	1000	Ritter (B) tray cover
917513	Blue	8-1/2" x 12-1/4"	1000	Ritter (B) tray cover
917514	Yellow	8-1/2" x 12-1/4"	1000	Ritter (B) tray cover
917516	Mauve	8-1/2" x 12-1/4"	1000	Ritter (B) tray cover
917519	Lavender	8-1/2" x 12-1/4"	1000	Ritter (B) tray cover
917520	TIDI® Tooth	8-1/2" x 12-1/4"	1000	Ritter (B) tray cover
917521	White	9-1/2" x 12-3/8"	1000	Chayes (A) tray cover
917523	Blue	9-1/2" x 12-3/8"	1000	Chayes (A) tray cover
917531	White	9-1/2" x 13-1/2"	1000	Midwest (E) tray cover
917533	Blue	9-1/2" x 13-1/2"	1000	Midwest (E) tray cover
917536	Mauve	9-1/2" x 13-1/2"	1000	Midwest (E) tray cover
917541	White	10-1/4" x 15-3/4"	1000	S.S. White (D) tray cover
917543	Blue	10-1/4" x 15-3/4"	1000	S.S. White (D) tray cover
917551	White	11" × 17-1/4"	1000	Weber (C) tray cover
917553	Blue	11" × 17-1/4"	1000	Weber (C) tray cover
917561	White	13-1/2" diameter	1000	Circular (G) tray cover
917571	White	5" x 8"	1000	Mini (F) tray cover
917573	Blue	5" x 8"	1000	Mini (F) tray cover
917581	White	9" x 13-1/2"	1000	Certified tray cover
917591	White	11" x 11"	1000	Cox tray cover

For extra protection and infection prevention, use our tray covers with our poly tray sleeves, found on page 15





SHEATHS + BARRIERS

Aid in the protection from cross contamination for equipment, instruments and other surfaces in the exam environment.

PATENT PENDING



Custom fit for ergonomic handling. Soft, rounded edges for greater patient comfort.





with *The SureClear Window*™



COMPETITOR SHEATH



A competitor's camera sheath gives the dentist a blurry image, making the view of the tooth difficult

TIDISHIELD™ BARRIER



The TIDIShield™ Custom
Barrier, with it's tight fit around
the camera, provides a perfect
clear view and superior
infection prevention

NO SHEATH



Without a protective sheath, the intra-oral camera offers an unobstructed view, but with the threat of cross-contamination





CONFIDENT CURES

The TIDIShield™ Curing Light Sleeve with The SureCure Window™ is proven to produce surfaces as well cured as those where no barrier is used.

INFECTION PREVENTION

Unique individually wrapped, easy-peel outer wrapping protects inner sleeve from contamination. Eliminates the potential for bacteria build-up on curing light surfaces.

COMFORT

Custom fit with soft, rounded edges for greater patient comfort.

COST-EFFECTIVE

Extends the life of equipment by reducing the need for harsh disinfectants on your curing light. Prevents hard to remove composite build up on the curing light tip.







- Custom fit that eliminates distortion
- SureClear Window[™] for superior image quality
- Individually wrapped for enhanced infection prevention
- Extends the life of equipment by eliminating the need for harsh disinfectants
- Unique easy peel design
- Not made from natural latex rubber
- Instructions for use printed directly on the packaging
- Reorder part number printed on the sleeve

CODE	QUANTITY	MANUFACTURER	CAMERA MODEL
21001	500 (5 boxes of 100)	Carestream/Kodak	1500
20808	500 (5 boxes of 100)	Schick	USB
20917	500 (5 boxes of 100)	Camsight	Optum II
20989	500 (5 boxes of 100)	Camsight	Opticam
20829	500 (5 boxes of 100)	Carestream	STV Pro I & II Kodak 1000
20862	500 (5 boxes of 100)	Cygnus	IDT Oral Vision 5000, 6000
20834	500 (5 boxes of 100)	Dental Vision	Ultracam
20816	500 (5 boxes of 100)	Dentsply	AcuCam 90-degree lens, AcuCam Plus, Polo
20893	500 (5 boxes of 100)	Dentsply	AcuCam Concept III
20856	500 (5 boxes of 100)	Dentsply	Insight
21013	500 (5 boxes of 100)	Digital Doc	Iris
20883	500 (5 boxes of 100)	Digital Doc	NCC Model
20911	500 (5 boxes of 100)	Digital Doc	Aurora
20968	500 (5 boxes of 100)	Digital Doc	Net Design
20990	500 (5 boxes of 100)	Digital Doc	Icon
20810	500 (5 boxes of 100)	DMD	Telicam
20800	500 (5 boxes of 100)	Elca	Elcacam
20818	500 (5 boxes of 100)	Expanded Optics	Explorer
20981	500 (5 boxes of 100)	GC Corp	Crystalcam
20812	500 (5 boxes of 100)	ICS	Vidicam, Intracam, EasyCam, Prospective, Celebrity, Accucure 300
20868	500 (5 boxes of 100)	ICS	Easy Cam, Easy Doc, True Vision
20884	500 (5 boxes of 100)	ICS	Easy Cam II (Tom Cat Supercam and Slimwand, Easy Doc)
20802	500 (5 boxes of 100)	IDT	DT Oral Vision H-3, H-7, H-8, H-9
20987	500 (5 boxes of 100)	Imagin	Hot Rod
21026	500 (5 boxes of 100)	Imagin Systems	Image Master
20991	500 (5 boxes of 100)	Miras	Mirrocam
20984	500 (5 boxes of 100)	Planmeca	AF Camera Sheath
20831	500 (5 boxes of 100)	Pro-Den Systems / Dent-X	Pro-Scope 1000 and Oral Scan, Easy Doc
20806	500 (5 boxes of 100)	Prodentec	Prism
20916	500 (5 boxes of 100)	DE Systoms	
	000 (0 boxes of 100)	RF Systems	DP 7S or DP 5/6
20964	500 (5 boxes of 100)	RF Systems	DP 7S or DP 5/6 SS21
21006	,		
	500 (5 boxes of 100)	RF Systems	SS21
21006	500 (5 boxes of 100) 500 (5 boxes of 100)	RF Systems RF Systems	SS21 Einstein
21006 20867	500 (5 boxes of 100) 500 (5 boxes of 100) 500 (5 boxes of 100)	RF Systems RF Systems Schick	SS21 Einstein Schick CDR Cam 2000
21006 20867 20855	500 (5 boxes of 100) 500 (5 boxes of 100) 500 (5 boxes of 100) 500 (5 boxes of 100)	RF Systems RF Systems Schick Siemens	SS21 Einstein Schick CDR Cam 2000 Ceracam/Minicam Ultra
21006 20867 20855 20970	500 (5 boxes of 100) 500 (5 boxes of 100) 500 (5 boxes of 100) 500 (5 boxes of 100) 500 (5 boxes of 100)	RF Systems RF Systems Schick Siemens Sirona	SS21 Einstein Schick CDR Cam 2000 Ceracam/Minicam Ultra Sirocam – 3
21006 20867 20855 20970 20899	500 (5 boxes of 100) 500 (5 boxes of 100)	RF Systems RF Systems Schick Siemens Sirona Sometech	SS21 Einstein Schick CDR Cam 2000 Ceracam/Minicam Ultra Sirocam – 3 Dr. Camscope
21006 20867 20855 20970 20899 20915	500 (5 boxes of 100) 500 (5 boxes of 100)	RF Systems RF Systems Schick Siemens Sirona Sometech Sometech	SS21 Einstein Schick CDR Cam 2000 Ceracam/Minicam Ultra Sirocam – 3 Dr. Camscope Oracam
21006 20867 20855 20970 20899 20915 20928	500 (5 boxes of 100) 500 (5 boxes of 100)	RF Systems RF Systems Schick Siemens Sirona Sometech Sometech Sometech	SS21 Einstein Schick CDR Cam 2000 Ceracam/Minicam Ultra Sirocam – 3 Dr. Camscope Oracam Oracam
21006 20867 20855 20970 20899 20915 20928 20993	500 (5 boxes of 100) 500 (5 boxes of 100)	RF Systems RF Systems Schick Siemens Sirona Sometech Sometech Sometech Sometech	SS21 Einstein Schick CDR Cam 2000 Ceracam/Minicam Ultra Sirocam – 3 Dr. Camscope Oracam Oracam ST-121
21006 20867 20855 20970 20899 20915 20928 20993 20929	500 (5 boxes of 100) 500 (5 boxes of 100)	RF Systems RF Systems Schick Siemens Sirona Sometech Sometech Sometech Sometech Sometech Sometech Sometech Sometech	SS21 Einstein Schick CDR Cam 2000 Ceracam/Minicam Ultra Sirocam – 3 Dr. Camscope Oracam Oracam ST-121 Suni-Cam II
21006 20867 20855 20970 20899 20915 20928 20993 20929	500 (5 boxes of 100)	RF Systems RF Systems Schick Siemens Sirona Sometech Sometech Sometech Sometech Sometech Sometech Suni Suni	Einstein Schick CDR Cam 2000 Ceracam/Minicam Ultra Sirocam – 3 Dr. Camscope Oracam Oracam ST-121 Suni-Cam II Suni-Cam





Continued

CODE	QUANTITY	MANUFACTURER	CAMERA MODEL
20805	500 (5 boxes of 100)	Vision Tools	Spectra Vu
20852	500 (5 boxes of 100)	Welch Allyn	Reveal
20932	500 (5 boxes of 100)	Win US Tech	WIN 100D
20833	500 (5 boxes of 100)	Yoshida Dental	Pict Cam
20876	500 (5 boxes of 100)	Yoshida Dental	Pict Cam
20905	500 (5 boxes of 100)	Yoshida Dental	Pict Cam



- Custom fit for ergonomic handling
- Soft rounded edges for greater patient comfort
- Individually wrapped for enhanced infection
- Extends the life of equipment by eliminating the need for harsh disinfectants
- Unique easy peel design
- Not made from natural latex rubber
- Instructions for use printed directly on the packaging
- Reorder part number printed on the sleeve

CODE	QUANTITY	MANUFACTURER	SIZE
20977	500 (5 boxes of 100)	Carestream/Kodak 6100	Size 0
20978	500 (5 boxes of 100)	Carestream/Kodak 6100	Size 1
20979	500 (5 boxes of 100)	Carestream/Kodak 6100	Size 2
20841	500 (5 boxes of 100)	Cygnus Ray MPS	Size 1
20999	500 (5 boxes of 100)	Dexis	Universal
20957	500 (5 boxes of 100)	Palodex	Size 1
20861	500 (5 boxes of 100)	Regam	Size 2
20998	500 (5 boxes of 100)	Schick	Size 0
20824	500 (5 boxes of 100)	Schick	Size 1
20825	500 (5 boxes of 100)	Schick	Size 2
20819	500 (5 boxes of 100)	Suni	Size 1
20997	Bulk (1 box of 500)	Suni	XL – X-Ray Over Bag
№ 21040	500 (5 boxes of 100)	Progeny	Size 2
21041	500 (5 boxes of 100)	Progeny	Size 1



Soft X-ray sensor material provides comfort for patients.

CUSTOM CAMERA COVERS FROM **EQUIPMENT MANUFACTURERS**

TIDI® custom-fit barriers, designed for a perfect fit, are available direct from the equipment manufacturer. For barriers designed to fit the following models, contact the manufacturers at the numbers listed below.

CODE	MANUFACTURER	CAMERA	CONTACT INFORMATION
20927	Air Techniques	Acclaim	Air Techniques Customer Service: 516-433-7676 (Option 1)
20836	Air Techniques	Vistacam	Air Techniques Customer Service: 516-433-7676 (Option 1)
20908	Gendex	AcuCam IV	Dental Imaging Technologies: 800-205-3570
21003	Dexis	DexCam 3	Dental Imaging Technologies: 800-205-3570
21004	Gendex	GXC-300	Dental Imaging Technologies: 800-205-3570
20967	Sirona	Sirocam	Sirona USA: 800-659-5977
20892	Sirona	Sirocam 2	Sirona USA: 800-659-5977
20941	Sopro	S717	Sopro Customer Service: 800-289-6367
20972	Sota	Claris	Sota Customer Service: 888-266-7682
20809	Schick Technologies	USB	Schick Customer Service: 1-888-4-SCHICK
21007	Carestream	Kodak 1500	Carestream: 800-944-6365

CUSTOM X-RAY SENSOR SHEATHS FROM EQUIPMENT MANUFACTURERS



CODE	MODEL	SIZE	CONTACT INFORMATION
20890	Sirona	Size 2	Sirona USA: 800-659-5977
20904	Sirona	Size 1	Sirona USA: 800-659-5977



- Custom fit SureCure Window™reduces light degradation* for confident cures
- Keeps LED tip clean
- Individually wrapped for enhanced infection prevention
- Extends the life of equipment by eliminating the need for harsh disinfectants
- Unique easy peel design
- Not made from natural latex rubber
- Instructions for use printed directly on the packaging
- Reorder part number printed on the sleeve

*As compared to other barriers

CODE	QUANTITY	MANUFACTURER	CURING LIGHT MODEL
21100	500 (5 boxes of 100)	Kerr	Demi Ultra
21101	500 (5 boxes of 100)	3M	Elipar S10
21102	500 (5 boxes of 100)	Ivoclar Vivadent	Bluephase
21103	500 (5 boxes of 100)	Dentsply	SmartLite Max
21105	500 (5 boxes of 100)	Ultradent	Valo
21106	500 (5 boxes of 100)	SDI	Radii Plus
21109	500 (5 boxes of 100)	Coltene-Whaledent	SPEC3



Visit www.TIDIProducts.com for a complete list of curing light sleeves available.

STANDARD FIT SURFACE BARRIERS

Conveniently protect surfaces against cross contamination.

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
915001	Clear	27-1/2" x 24"	225/Roll, 4 Rolls/Case	Poly chair cover sleeve, on a roll
915002	Clear	2-1/2" x 10"	500/Box, 24 Boxes/Case	Poly air/water syringe sleeve, closed end
915003	Clear	4" × 6"	500/Box, 24 Boxes/Case	T-style light handle cover
915004	Clear	15" x 26"	250/Box, 4 Boxes/Case	Poly X-ray equipment sleeve
915005	Clear	10-1/2" x 14"	500/Box, 12 Boxes/Case	Poly tray sleeve
915006	Clear	2-1/2" x 10"	500/Box, 24 Boxes/Case	Poly air/water syringe sleeve, open end
915007	Clear	11-5/8" x 14-1/2"	500/Box, 4 Boxes/Case	Poly tray sleeve
915008	Clear	11-1/2" x 16"	500/Box, 4 Boxes/Case	Poly tray sleeve



915003 T-style light handle cover



Simplified surface protection. Additional headrest covers on page 9.





NITRILE 12" GLOVES

Nitrile Gloves provide a non-latex alternative that provide extra protection while handling instruments during the cleaning, rinsing, drying, packaging and sorting procedures that take place during Instrument Reprocessing.

- Latex Free
- Powder-Free
- 8 mils at fingertips
- Ambidextrous for quick and easy donning
- Textured for added grip
- Beaded cuff
- 12 inches long for additional protection

PRODUCT NUMBER	SIZE	UNIT QUANTITY
940001-1	Small	10 dispenser boxes of 50 gloves; 500 gloves/case
940002-1	Medium	10 dispenser boxes of 50 gloves; 500 gloves/case
940003-1	Large	10 dispenser boxes of 50 gloves; 500 gloves/case
940004-1	X-Large	10 dispenser boxes of 50 gloves; 500 gloves/case



POWDERFREE LATEX 18 ML GLOVES

- 18 mils thick at fingertips, eliminates need for double gloving
- Fitted for right and left hands, reduces hand fatigue
- Textured for added grip
- Beaded cuff
- 12 inches long for additional protection

PRODUCT NUMBER	SIZE	UNIT QUANTITY
932477-1	Small	10 dispenser boxes of 50 pair; 1000 gloves/case
932480-1	Medium	10 dispenser boxes of 50 pair; 1000 gloves/case
932483-1	Large	10 dispenser boxes of 50 pair; 1000 gloves/case
932486-1	X-Large	10 dispenser boxes of 50 pair; 1000 gloves/case

POWDERFREE LATEX 14 ML GLOVES

- 14 mils thick at fingertips, eliminates need for double gloving
- Ambidextrous for quick and easy donning
- Textured for added grip
- Beaded cuff
- 12 inches long for additional protection
- UL Classified as meeting NFPA requirements

PRODUCT NUMBER	SIZE	UNIT QUANTITY
BS0450-1	Small	20 dispenser boxes of 50 gloves; 1000 gloves/case
BS0460-1	Medium	20 dispenser boxes of 50 gloves; 1000 gloves/case
BS0470-1	Large	20 dispenser boxes of 50 gloves; 1000 gloves/case
BS0480-1	X-Large	20 dispenser boxes of 50 gloves; 1000 gloves/case





TIDISHIELD™ ASSEMBLE 'N GO™ EYE SHIELD

- Weighs only 3/10 of an ounce!
- "Wrap-around" design provides both front and side eye protection
- Replacement optical-grade replacement lenses available for maximized cost-effect replenishment
- Frames are autoclavable and reusable
- Perfect for both patient and provider

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
9210-100-6	Clear	One Size	600 (6 boxes of 100)	Disposable lens for use with 9452-10 frames
9452-10	Combo Pack	One Size	5 frames and 25 lenses (10 boxes)	Reusable frames with disposable eye shields
9452P-10	Combo Pack	One Size	5 pink frames and 25 lenses (10 boxes)	Reusable frames with disposable eye shields

TIDISHIELD™ ASSEMBLE 'N GO™ FACE SHIELD

- Weighs less than 1.5 ounces!
- Designed to maximize airflow and prevent fogging
- Provides over-the-top, side and front face protection
- Disposable, optical grade, anti-static replacement shields
- Includes a clear, snap-in nosepiece for support
- Replacement shields available for maximized cost-effective replenishment
- Visors are autoclavable and reusable



CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
2210-40-6	Clear	13" x 7"	240 (6 boxes of 40)	Disposable optical grade shields for use with 2452-6 visors
2210-100-5	Clear	13" x 7"	500 (5 boxes of 100)	Disposable optical grade shields for use with 2452-6 visors
2452-6	Mini Pack	One Size	3 visors and 25 shields (6 boxes)	Reusable visors with disposable face shields

TIDISHIELD™ DISPOSABLE FACE SHIELD

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
FS150-100	Clear	Standard 11.5" wide x 9" long	100 (4 boxes of 25)	One piece face shield with foam headband, elastic strap, and optical grade lens
FS150-50	Clear	Standard 11.5" wide x 9" long	50 (bulk)	One piece face shield with foam headband, elastic strap, and optical grade lens
FS160F-50	Clear	Full size 14" wide x 11.5" long	50 (bulk)	One piece face shield with foam headband, elastic strap, and optical grade lens



ACCESSORIES

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
9810-50	Side Shield Clear	One size	50 Side Shields (25 pairs)	Disposable Side Shields, use with prescription glasses

TIDISHIELD FLIP 'N GO™ EYE SHIELD

- · Perfect for patient or provider
- FLAT for compact storage, FLIP to wear
- Extremely lightweight and comfortable

- Wrap-around design for front and side eye protection
- · Optical-grade, anti-static, lens

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
9300FE-50	5 assorted colored frames/clear lens	One size fits most	50 (2 boxes of 25)	One piece disposable eyewear; stores flat and assembles with one simple step













TIDISHIELD™ PROCEDURE FACEMASKS

Experience the difference with TIDIShield™ Fluid Resistant Procedure facemasks and compare the quality to traditional tissue lined facemasks typically used in dental offices. Proudly made in the USA to maximize comfort and protection.



9040DB Ocean Blue



9040DSA Pink Punch



9040DY Lemon Drop



9040DL Lavender Swirl

MAXIMUM-COMFORT - Extremely soft and lightweight material; perfect for longer wear.

IMPROVED EARLOOP DESIGN - Longer, rounded earloops that will not pull or slice into ears.

SUPERIOR BREATHABILITY - Lightweight material that is latex-free, fiberglass-free, odorless and non-irritating allowing you to breath easier.

FLUID RESISTANT - Water repellent SMS material prevents liquids from passing through the mask.

UNIQUE DESIGN - 3-ply downward pleat design provides an exact fit and eliminates the chance for fluid pooling.

CODE	COLOR/PRINT	UNIT QUANTITY	LEVEL	DESCRIPTION
9040DB	Ocean Blue	10 boxes of 50/500 facemasks/case	Level 2	TIDIShield™ Fluid Resistant Procedure Facemask
9040DSA	Pink Punch	10 boxes of 50/500 facemasks/case	Level 2	TIDIShield™ Fluid Resistant Procedure Facemask
9040DY	Lemon Drop	10 boxes of 50/500 facemasks/case	Level 1	TIDIShield™ Fluid Resistant Procedure Facemask
9040DL	Lavender Swirl	10 boxes of 50/500 facemasks/case	Level 1	TIDIShield™ Fluid Resistant Procedure Facemask
9010	Blue	10 boxes of 50/500 facemasks case	Level 3	Fluid Resistant Facemasks, High Barrier Protection
9020	Blue	10 boxes of 40/400 facemasks case	Level 3	Fog Free Fluid Resistant Facemasks, High Barrier Protection

ASTM F2100-07 Performance Classes for Facemasks

TIDIShield™ Facemasks are tested with the ASTM F2100-07 "Standard for Performance of Materials Used in Medical Facemasks." This is the only specific standard available that provides performance classes for facemasks to fall into.

REQUIREMENTS	9040DB 9040DSA	9040DY 9040DL	9010 9020
ASTM F1862 Resistance to penetration by synthetic blood (mm HG)	120	80	160
ASTM F2101 Bacterial Filtration Efficiency (BFE)	≥ 95	≥ 99	≥ 99
ASTM F2299 Particulate Filtration Efficiency (PFE)	≥ 99	≥ 99	≥ 99
Delta P Breathability	< 2.0	< 2.0	< 5.0
Flame Spread	Class 1	Class 1	Class 1

The performance classes rate from low, moderate and high barrier.



TIDISHIELD™ PROCEDURE FACEMASKS

Reliable Protection with TIDIShield™-

- OSHA recommends changing masks during patient treatments if the mask becomes wet to maintain high filterabiltiy.
- Made from a water repellent SMS material that prevents liquid from passing through.
- Provides lasting superior protection during extended procedures compared to tissue-lined facemasks that can absorb fluid and lose their protective barrier properties once wet.

Superior Protection with TIDIShield™

- Exceeds the required ASTM Level 1 BFE & PFE filtration testing standards; offering you superior protection compared to traditional Level 1 procedure facemasks.
- Ocean Blue and Pink Punch masks are more resistant than Level 1 procedure facemasks to fluid penetration and meet ASTM Level 2 requirements.





















Not Required for Level 1 Facemasks










Match colors and prints to our tissue and paper products for a friendly PATIENT CUPS and coordinated exam environment.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
9211	White	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9212	Green	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9213	Blue	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9214	Yellow	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9215	Gray	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9216	Mauve	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9217	Beige	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9218	Peach	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9210	Lavender	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9241	White	3.5 oz.	1000, 100/bag	Plastic cup, ribbed design
9243	Blue	3.5 oz.	1000, 100/bag	Plastic cup, ribbed design
9220	TIDI® Tooth	5 oz.	1000, 100/bag	Waxed paper cup
9221	TIDI® Tooth	5 oz.	500, 100/bag	Waxed paper cup
9226	Jazz	5 oz.	1000, 100/bag	Waxed infused paper cup
9240	Jazz	3 oz.	5000, 100/bag	Waxed infused paper cup
9231	White	6 oz.	2500, 250/bag	Evacuation cup liner, ADEC





CUP DISPENSER

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
9246	Chrome	3 or 5 oz., 16-1/4" x 3-3/4" x 3-3/4"	1	Adjustable cup dispenser

SALIVA EJECTORS securely fastened.

Saliva ejectors contain copper-plated wire, molded in, allowing them to bend easily and then hold their shape. Smooth tips are

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
9306	White	6"	1000, 100/bag	White ejector with white tip
9307	Clear	6"	1000, 100/bag	Clear ejector with blue tip
9304	White	Standard	1000, 100/bag	High-volume evacuator tip
9333	White	Standard	1000, 100/bag	High-volume evacuator tip, vented



COTTON, GAUZE + NONWOVENS

Choose the products that meet your needs, whether pure, natural cotton or nonwovens for maximum absorbency.

COTTON-FILLED GAUZE SPONGES

Cotton-filled sponges feature a woven gauze facing with an absorbent cotton fill, in both sterile and non-sterile packaging.

BRAND	CODE	SIZE	UNIT QUANTITY	DESCRIPTION
VENTURE™	919000	2" x 2"	5000 (25 sleeves of 200)	Gauze sponges with cotton filler, non-sterile
VENTURE™	908223	3" x 3"	4000 (40 sleeves of 100)	Gauze sponges with cotton filler, non-sterile
VENTURE™	908224	4" x 4"	2000 (20 sleeves of 100)	Gauze sponges with cotton filler, non-sterile
TIDI®	919010	2" x 2"	5000 (25 sleeves of 200)	Gauze sponges with cotton filler, sterile



made

ALL-GAUZE SPONGES Gauze sponges are ideal for general cleaning and packing.

VENTURE™	908282	2" x 2"	5000 (25 sleeves of 200)	8-ply gauze sponges, non-sterile
VENTURE™	908294	2" x 2"	8000 (40 sleeves of 200)	12-ply gauze sponges, non-sterile
VENTURE™	908206	3" x 3"	4000 (20 sleeves of 200)	8-ply gauze sponges, non-sterile
VENTURE™	908293	3" x 3"	4000 (20 sleeves of 200)	12-ply gauze sponges, non-sterile
VENTURE™	908292	4" × 4"	4000 (20 sleeves of 200)	8-ply gauze sponges, non-sterile
VENTURE™	908295	4" × 4"	2000 (10 sleeves of 200)	12-ply gauze sponges, non-sterile
VENTURE™	908291	4" × 4"	2000 (10 sleeves of 200)	16-ply gauze sponges, non-sterile
VENTURE™	908200	4" x 8"	2000 (10 sleeves of 200)	12-ply gauze sponges, non-sterile

COTTON ROLLS

TIDI non-linting non x-ray detectable cotton rolls are machine-rolled into soft, pliable cylinders. Venture rolls offer quality at an economical price.

TIDI®	969120	1-1/2" x 5/16"	2000 (40 bundles of 50)	Non-sterile, #1 small diameter
TIDI®	919121	1-1/2" x 3/8"	2000 (40 bundles of 50)	Non-sterile, #2 medium diameter
TIDI®	969122	3/8" x 6"	500 (10 bundles of 50)	Non-sterile, #2 medium diameter
TIDI®	969123	1-1/2" x 1/2"	2000 (40 bundles of 50)	Non-sterile, #3 large diameter
TIDI®	969126	1-1/2" x 3/8"	2000 (40 bundles of 50)	Sterile, #2 medium diameter

NONWOVEN SPONGES

Nonwoven sponges are a highly absorbent, low-linting option for prepping and packing. Also convenient for wiping down instruments and workspace.

OMNISORB®	942022	2" x 2"	5000 (25 sleeves of 200)	4-ply (closed weave), non-sterile
OMNISORB®	942033	3" x 3"	4000 (20 sleeves of 200)	4-ply (closed weave), non-sterile
OMNISORB®	942044	4" × 4"	2000 (10 sleeves of 200)	4-ply (closed weave), non-sterile
VENTURE™	912022	2" x 2"	5000 (25 sleeves of 200)	4-ply (closed weave), non-sterile
VENTURE™	912033	3" x 3"	4000 (20 sleeves of 200)	4-ply (closed weave), non-sterile
VENTURE™	912044	4" × 4"	2000 (10 sleeves of 200)	4-ply (closed weave), non-sterile
VENTURE™	260002	2" x 2"	8000 (40 sleeves of 200)	4-ply (open weave), non-sterile
VENTURE™	260003	3" x 3"	4000 (20 sleeves of 200)	4-ply (open weave), non-sterile
VENTURE™	260004	4" × 4"	2000 (10 sleeves of 200)	4-ply (open weave), non-sterile



DURAWICK™ COUNTER TOWEL

DuraWick™ counter towels were specially created to be used to last through several rounds of tools being placed on them from your ultrasonic cleaning. The top layer of material wicks away moisture, which is then absorbed into the middle layer of the product, keeping your sterilization area clean and sanitary.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
983914	Blue	13" x 18"	100	Wicking towel for ultrasonlic tool cleaning, 3-ply Rayon/Tissue/Poly

CSR WRAP

Sterilization wrap conforms easily around various shapes and maintains its wrapped integrity to provide barrier protection against saline, blood and water while allowing effective penetration for a variety of sterilization methods.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
960172	Blue	12" x 12"	1000	Sterilization wrap
960173	Blue	15" x 15"	1000	Sterilization wrap
960174	Blue	18" x 18"	500	Sterilization wrap
960175	Blue	20" x 20"	500	Sterilization wrap
960176	Blue	24" x 24"	500	Sterilization wrap
960177	Blue	30" x 30"	400	Sterilization wrap
960178	Blue	36" x 36"	200	Sterilization wrap

STERILIZATION POUCHES

Self-sealing and have external color indicators for peace of mind.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
5549	Clear one side	2-3/4" x 10"	300	Sterilization pouches, dual indicators, latex free
5550	Clear one side	3-1/4" x 6-1/2"	200	Sterilization pouches, dual indicators, latex free
5551	Clear one side	3-1/2" x 10"	200	Sterilization pouches, dual indicators, latex free
5552	Clear one side	5-1/4" x 11"	200	Sterilization pouches, dual indicators, latex free
5553	Clear one side	7-1/2" x 12"	200	Sterilization pouches, dual indicators, latex free
5543	Blue transparent	2-1/2" x 8-1/4"	1000	Autoclave envelopes, contain latex in adhesive strip
5544	Blue transparent	3-1/2" x 9-3/4"	1000	Autoclave envelopes, contain latex in adhesive strip
5545	Blue transparent	6" x 9-3/4"	1000	Autoclave envelopes, contain latex in adhesive strip
30203	White	2-1/2" x 1-1/2" x 10-1/2"	1000	Autoclave envelopes, contain latex in adhesive strip

AUTOCLAVABLE TOWELS

A poly layer, sandwiched between layers of tissue, allows for use in kits and in autoclaves, while maintaining the integrity of the towels.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
8253	White	13" x 18"	500	Tissue/poly/tissue autoclavable towel
8252	White	19" x 22"	300	Tissue/poly/tissue autoclavable towel
8251	White	19" x 30"	300	Tissue/poly/tissue autoclavable towel



INFECTIOUS WASTE DISPOSAL BAGS

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
8653	Red	24" x 24" (10 gal)	250	High-density (12 microns) infectious waste (Biohazard) bags
8654	Red	24" x 33" (16 gal)	250	High-density (12 microns) infectious waste (Biohazard) bags
8656	Red	33" x 40" (33 gal)	250	High-density (17 microns) infectious waste (Biohazard) bags
8657	Red	40" x 48" (40 gal)	250	High-density (17 microns) infectious waste (Biohazard) bags

CHAIRSIDE DISPOSAL BAGS

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
950232	White Plain	6-1/2" x 3-1/8" x 11-3/8"	2000	Flame-retardant paper bag, chairside
950231	Blue Floral	6-1/2" x 3-1/8" x 11-3/8"	2000	Flame-retardant paper bag, chairside
950251	Blue Soft Shell	6-1/2" x 3-1/8" x 11-3/8"	2000	Non-flame-retardant paper bag, chairside
950252	White Plain	6-1/2" x 3-1/8" x 11-3/8"	2000	Non-flame-retardant paper bag, chairside
950234	Blue Floral	6-1/2" x 3-1/8" x 11-3/8"	2000 (20 bags of 100)	Non-flame-retardant poly bag, chairside

STORAGE ENVELOPES

CODE	BRAND	SIZE	UNIT QUANTITY	DESCRIPTION
950222	Kraft	8-1/2" x 10-1/2"	1000/case, 20 bundles/case, rotated	X-ray storage envelopes
950221	Kraft	10-1/2" x 12-1/2"	1000/case, 20 bundles/case, rotated	X-ray storage envelopes
950218	Kraft	14-1/2" x 17-1/2"	500/case, 10 bundles/case, rotated	X-ray storage envelopes, heavyweight
950220	Kraft	14-1/2" x 17-1/2"	500/case, 10 bundles/case, rotated	X-ray storage envelopes, lightweight
950219	Kraft	15" x 18"	50/case, bulk packed	X-ray mailer envelopes with insert



KAY-PEES® PROFESSIONAL TOWELS For sensitive skin, lint-free, low friction. made in USA

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
911898	White	13-1/2" x 17-1/2"	500 (10 boxes of 50)	4-ply tissue, super soft, highly absorbent, edge-embossed
911900	Blue	13-1/2" x 17-1/2"	500 (10 boxes of 50)	4-ply tissue, super soft, highly absorbent, edge-embossed

T/P/T: tissue/poly/tissue

A layer of polyethylene film, sandwiched between two sheets of tissue, adds strength and resistance against fluid strike-through to the soft, absorbent tissue.

Tissue (2-ply, 3-ply and 4-ply)

TIDI® offers tissue products in configurations to suit a range of needs. Multiple tissue plies add absorbency and strength, while two-ply tissue provides economical softness and absorbency that is well suited for shorter procedures.

Poly

A layer of premium linear polyethylene provides high tensile strength and reliable barrier properties.

FABRICEL®

A proprietary material that fuses a sheet of soft, supple facial-grade tissue and a layer of polyethylene film. This heat lamination process creates a very soft yet super-strong, tear-resistant material that prevents soak-through.

TIDI® Cust om Barriers

TIDI® Products is proud to partner with manufacturers of dental and medical devices that require a custom-fit infection control barrier. Our design team works with manufacturers to create custom sheaths that aid in the protection from cross contamination while allowing healthcare professionals to use the instrument without the barrier interfering with the procedure. If you are a manufacturer and would like more information on how we can develop a custom barrier for your device, please **contact TIDI® Customer Service at 800-521-1314.**















TIDI® supports the Forest Stewardship Council and responsible forestry.

Certified in ISO 9001 and EN/ISO 13485













TIDI® Products, LLC 570 Enterprise Drive Neenah, WI 54956 USA

Customer Service 800.521.1314 +1 920.751.4138 800.837.7770 FAX

www.TIDIProducts.com