

## TECHNICAL DATA SHEET

### Medicom SafeBasics Bib

**Items # 8280, 8281, 8282, 8283, 8284, 8285, 8286, 8287, 8288, 8289**

#### Intended Use

Dental bibs are made of paper intended to absorb saliva or moisture from the oral cavity during dental procedures

Characteristic	Specification		
Description	Folded piece of paper with a plastic backing. Embossed horizontal pattern on the body of the bib and vertical embossed pattern on the vertical borders.		
Material	<b>Composition</b>	Tissue Paper (10.5 – 11.5 lbs/3000 sq.ft for 1 ply) Polyethylene Film (0.40 mil +/- 10%)	
Color	<b>Item code</b>		<b>Color</b>
	8280		Lavender
	8281		Silver Grey
	8282		Blue
	8283		White
	8284		Yellow
	8285		Green
	8286		Dusty Rose
	8287		Beige
	8288		Peach
	8289		Aqua
Dimensions	<b>Unfolded Size</b>	<b>Unfolded Width (cm)</b>	<b>Unfolded Length (cm)</b>
	13" x 18"	33	45,7
Ply	2 tissue ply		
Packaging	500 Bibs/Case		

#### Visual



## TECHNICAL DATA SHEET

### Medicom SafeBasics Bib

**Items # 8280, 8281, 8282, 8283, 8284, 8285, 8286, 8287, 8288, 8289**

Item Code	Description	Color	Specific Attributes
8280	Medicom SafeBasics Dry-Back Bibs, 13" x 18" (33 x 45,7 cm), 500 Bibs/Case	Lavender	2 Ply of Tissue + 1 Ply Poly, Non-Sterile
8281	Medicom SafeBasics Dry-Back Bibs, 13" x 18" (33 x 45,7 cm), 500 Bibs/Case	Silver Grey	2 Ply of Tissue + 1 Ply Poly, Non-Sterile
8282	Medicom SafeBasics Dry-Back Bibs, 13" x 18" (33 x 45,7 cm), 500 Bibs/Case	Blue	2 Ply of Tissue + 1 Ply Poly, Non-Sterile
8283	Medicom SafeBasics Dry-Back Bibs, 13" x 18" (33 x 45,7 cm), 500 Bibs/Case	White	2 Ply of Tissue + 1 Ply Poly, Non-Sterile
8284	Medicom SafeBasics Dry-Back Bibs, 13" x 18" (33 x 45,7 cm), 500 Bibs/Case	Yellow	2 Ply of Tissue + 1 Ply Poly, Non-Sterile
8285	Medicom SafeBasics Dry-Back Bibs, 13" x 18" (33 x 45,7 cm), 500 Bibs/Case	Green	2 Ply of Tissue + 1 Ply Poly, Non-Sterile
8286	Medicom SafeBasics Dry-Back Bibs, 13" x 18" (33 x 45,7 cm), 500 Bibs/Case	Dusty Rose	2 Ply of Tissue + 1 Ply Poly, Non-Sterile
8287	Medicom SafeBasics Dry-Back Bibs, 13" x 18" (33 x 45,7 cm), 500 Bibs/Case	Beige	2 Ply of Tissue + 1 Ply Poly, Non-Sterile
8288	Medicom SafeBasics Dry-Back Bibs, 13" x 18" (33 x 45,7 cm), 500 Bibs/Case	Peach	2 Ply of Tissue + 1 Ply Poly, Non-Sterile
8289	Medicom SafeBasics Dry-Back Bibs, 13" x 18" (33 x 45,7 cm), 500 Bibs/Case	Aqua	2 Ply of Tissue + 1 Ply Poly, Non-Sterile

#### Regulatory

FDA, Not a Medical Device

Health Canada, Not a Medical Device

#### Storage Conditions

Store in a cool dry place

#### Country of Origin

USA