

The lifetime of this product is 240 liters or 365 calendar days (if water usage records are not monitored) or if iodine levels fall to less than 0.5 ppm. To maximize the life of the filter, shut off water to the dental unit when not in use.

[illegible][illegible]

To reorder, call 1-800-HU-FRIEDY or order through your dental distributor.

Product Name: Hu-Friedy Waterline Filter

Part Code: IMS-1440

Call 1-800-Hu-Friedy or contact your authorized Hu-Friedy representative for more information.

©2018 Hu-Friedy Mfg. Co., LLC. All rights reserved. [980] IMS-3219/0818

Hu-Friedy