



RESIN 8 COLORS

A SUPERIOR RESIN INSTRUMENT

How the best perform



YOUR PERFORMANCE IS OUR PASSION.

At Hu-Friedy, every instrument, from its proprietary metal and textured design, to the color for organization, is designed with one purpose—to help you perform at your best.

When we created Resin 8 Colors scalars and curettes, we combined resin handles in 8 eye-catching colors with features that maximize your comfort and efficiency to achieve the best clinical outcomes—every patient, every day.



RESIN 8 COLORS POINTS OF PERFORMANCE

EVEREDGE® TECHNOLOGY

Stays sharper longer, which means less time sharpening

STEEL INNER CORE DESIGN

For added strength and tactile sensitivity

LARGE DIAMETER RESIN HANDLE

Optimum weight for a more comfortable grip and reduced pinch force

UNIQUE TEXTURED HANDLE DESIGN

Wave grooves and unique knurling texture increase rotational control and provide a light secure grasp

SIGNATURE SERIES® GRIPS

Silicone grip provides more comfort and color identification

PERFECT ANGLES

Consistently hand crafted blades for the most efficient calculus removal

COLOR CODED SYSTEM

To easily identify scalars and curettes that correspond to anatomical areas of the mouth



RESIN 8 COLORS GRACEY COLOR CODING

Resin 8 Colors Gracey curettes are organized by areas of the mouth, providing easy identification.



GRACEYS: **1/2, 3/4**

1/2 and 3/4 Graceys available in Standard, Rigid, After Five®, After Five Rigid, Mini Five® and Mini Five Rigid varieties



GRACEYS: **5/6**

5/6 Graceys available in Standard, Rigid, After Five, After Five Rigid, Mini Five and Mini Five Rigid varieties



GRACEYS: **7/8, 9/10**

7/8 Graceys available in Standard, Rigid, After Five, After Five Rigid, Mini Five and Mini Five Rigid varieties
9/10 Graceys available in Standard



GRACEYS: **11/12, 15/16**

11/12 and 15/16 Graceys available in Standard, Rigid, After Five, After Five Rigid, Mini Five and Mini Five Rigid



GRACEYS: **13/14, 17/18**


13/14 Graceys available in Standard, Rigid, After Five, After Five Rigid, Mini Five and Mini Five Rigid
17/18 Graceys available in Standard and Rigid




11/14 and 12/13 Graceys are available in the Standard variety and come in a Resin 8 Colors red handle

USE OF GRACEY CURETTES

GRACEY COLOR CHART



SYSTEMATIC INSTRUMENTATION



Lingual Maxillary Left
Second Molar
SG11/12C8



Lingual Mandibular Right
First Premolar
SG5/6C8

RESIN 8 COLORS
SCALERS & CURETTES



ANTERIOR SICKLE SCALERS

Purple Resin 8 Colors available with the following working ends:

Nevi 1 & Nevi 1/H5	O/00 Morse
H6/7 Hygienist	137
H5/33 Hygienist	U15/30 and U15/33 Towner - Jacquette
30/33 Jacquette	



POSTERIOR SICKLE SCALERS

Red Resin 8 Colors available with the following working ends:

204 IUFW	Nevi 2	31/32 Standard Jacquette
204S	Nevi 3	34/35 Miniature Jacquette
204SD	Nevi 4	
135		




UNIVERSAL CURETTES

Lavender Resin 8 Colors available with the following working ends:

Barnhart 1/2, 1/2R, 5/6, 5/6R	Younger-Good 7/8	Langer 1/2, 3/4, 5/6, 17/18
Columbia 2R/2L, 4R/4L, 13/14	McCall 13/14, 13/14S, 17/18	

CARING FOR YOUR
RESIN INSTRUMENTS



CLEANING AND CARE OF YOUR RESIN INSTRUMENTS

In order to maintain the integrity of the resin instruments, detergents or disinfectants containing phenols or iodophors cannot be used. In addition, dry heat is not compatible with instruments with resin handles (Hu-Friedy handle #8), with resin or silicone components, inserts on any instruments, or with resin cassettes. The sterilizer equipment manufacturer’s compatibility with specific materials must be observed.



MAINTAINING THE SHARP EDGE

Although Hu-Friedy’s Resin Color Instruments are made with EverEdge® technology which keeps the working ends sharper longer, at some point every scaler needs sharpening. Hu-Friedy offers over thirty years experience in instrument finishing and sharpening to the clinician who prefers a “hands-off” style of sharpening.

Through this service, our instrument artisans recreate factory cutting blades on a daily basis, sharpening your instrument to the performance level of a brand new scaler! With Hu-Friedy’s Professional Sharpening Service, you’ll save valuable time and effort sharpening and have more time to provide better patient care!



REPLACING YOUR RESIN INSTRUMENT

Don’t forget about Hu-Friedy’s Envirodent® Program! The longest running program in the dental industry, Envirodent provides an environmentally-friendly way for dental professionals to dispose of their old, broken, or re-tipped instruments, and earn free instruments.

No landfills, no refurbishing—instruments recycled through the Envirodent program are sent to metal scrap recyclers so they can find new life in other items such as bridges, car bumpers and street lights.

Visit [Hu-Friedy.com/Envirodent](https://www.hu-friedy.com/Envirodent) for more information.

RESIN 8 COLORS PART CODES

GRACEY CURETTES

GRACEY 1/2

■ Standard.....	SG1/2C8
■ Rigid.....	SG1/2RC8
■ After Five®.....	SRPG1/2C8
■ After Five Rigid.....	SRP1/2RC8
■ Mini Five®.....	SAS1/2C8
■ Mini Five Rigid.....	SS1/2RC8

GRACEY 3/4

■ Standard.....	SG3/4C8
■ Rigid.....	SG3/4RC8
■ After Five.....	SRPG3/4C8
■ After Five Rigid.....	SRP3/4RC8
■ Mini Five.....	SAS3/4C8
■ Mini Five Rigid.....	SS3/4RC8

GRACEY 5/6

■ Standard.....	SG5/6C8
■ Rigid.....	SG5/6RC8
■ After Five.....	SRPG5/6C8
■ After Five Rigid.....	SRP5/6RC8
■ Mini Five.....	SAS5/6C8
■ Mini Five Rigid.....	SS5/6RC8

GRACEY 7/8

■ Standard.....	SG7/8C8
■ Rigid.....	SG7/8RC8
■ After Five.....	SRPG7/8C8
■ After Five Rigid.....	SRP7/8RC8
■ Mini Five.....	SAS7/8C8
■ Mini Five Rigid.....	SS7/8RC8

GRACEY 9/10

■ Standard.....	SG9/10C8
-----------------	----------

GRACEY 11/12

■ Standard.....	SG11/12C8
■ Rigid.....	SG11/12RC8
■ After Five.....	SRP11/12C8
■ After Five Rigid.....	SP11/12RC8
■ Mini Five.....	SAS11/12C8
■ Mini Five Rigid.....	SS11/12RC8

GRACEY 13/14

■ Standard.....	SG13/14C8
■ Rigid.....	SG13/14RC8
■ After Five.....	SRP13/14C8
■ After Five Rigid.....	SP13/14RC8
■ Mini Five.....	SAS13/14C8
■ Mini Five Rigid.....	SS13/14RC8

GRACEY 15/16

■ Standard.....	SG15/16C8
■ Rigid.....	SG15/16RC8
■ After Five.....	SRP15/16C8
■ After Five Rigid.....	SP15/16RC8
■ Mini Five.....	SAS15/16C8
■ Mini Five Rigid.....	SS15/16RC8

GRACEY 17/18

■ Standard.....	SG17/18C8
■ Rigid.....	SG17/18RC8

GRACEY 11/14

■ Standard.....	SG11/14C8
-----------------	-----------

GRACEY 12/13

■ Standard.....	SG12/13C8
-----------------	-----------

SICKLE SCALERS

ANTERIOR SICKLE SCALERS

■ Nevi 1.....	SCNEVI1C8
■ Nevi 1/H5.....	SNEV1/H5C8
■ H5/33.....	SH5/33C8
■ H6/7.....	SH6/7C8
■ 30/33 Jacquette.....	SJ30/33C8
■ 0/00 Morse.....	SM0/00C8
■ 137.....	SN137C8
■ U15/30 Jacquette ...	SU15/30C8
■ U15/33 Jacquette.....	SU15/33C8

POSTERIOR SICKLE SCALERS

■ 204 IUFW.....	SIUFW204C8
■ 204S.....	S204SC8
■ 204SD.....	S204SDC8
■ Nevi 2.....	SCNEVI2C8
■ Nevi 3.....	SCNEVI3C8
■ Nevi 4.....	SCNEVI4C8
■ 31/32 Jacquette.....	SJ31/32C8
■ 34/35 Jacquette	SJ34/35C8
■ 135.....	SN135C8

UNIVERSAL CURETTES

BARNHART

■ 1/2.....	SBH1/2C8
■ 1/2 Rigid.....	SBH1/2RC8
■ 5/6.....	SBH5/6C8
■ 5/6 Rigid.....	SBH5/6RC8

COLUMBIA

■ 2R/2L.....	SC2R/2LC8
■ 4R/4L.....	SC4R/4LC8
■ 13/14.....	SC13/14C8

YOUNGER-GOOD

■ 7/8.....	SYG7/8C8
------------	----------

McCALL

■ 13/14.....	SM13/14C8
■ 13S/14S.....	SM13/14SC8
■ 17/18.....	SM17/18C8

LANGER

■ 1/2.....	SL1/2C8
■ 3/4.....	SL3/4C8
■ 5/6.....	SL5/6C8
■ 17/18.....	SL17/18C8



Denotes handle color.