

## Saliva Ejectors Non-Sterile

### Intended Use

Saliva ejectors are used in dentistry in order to remove saliva, blood or other debris from the mouth during dental procedures. This disposable device is connected to an active device, normally a pneumatic system suction circuit. The tip is attached via a melting process that ensures safety and protection from detachment, to avoid any vacuum effect during aspiration and also prevents the back flow.

Item Code	Description	Type
700	Medicom Safebasics SE clear with blue tip	Regular
701	Medicom Safebasics SE clear with white tip	Regular
702	Medicom Safebasics SE Blue with Blue tip	Regular
703	Medicom Safebasics SE white with white tip	Regular
710	Medicom Safebasics SE clear blue with blue tip	Regular
711	Medicom Safebasics SE clear pink with pink tip	Regular
712	Medicom Safebasics SE clear Lavender with Lavender tip	Regular
713	Medicom Safebasics SE mint scented clear green with green tip	Scented
714	Medicom Safebasics SE strawberry scented clear red with red tip	Scented



### Packaging

105 units per bag  
10 bags per case  
1050 units per case

### Size

Length of Saliva Ejector 150 mm  
Cap Length 14 mm  
Diameter of Cannula 6.2 mm

### Material

PVC Cannula  
Latex Free  
Reinforced with zinc or copper covered wire or anodized steel wire core.

### Usage

For single-use only

### Sterility

Non-Sterile

### Country of Origin

Italy