



Seizure Pad Selection Guide

A SIMPLE SOLUTION FOR
HARD SURFACES



5756 Progressa Bed
Seizure Pad



5757 Hill-Rom® Centrella
Bed Seizure Pad



5767 Stryker® Spirit Bed
Seizure Pad

Posey Seizure Pads are form fitted to fit snugly over the existing bed side rails. The nylon cover is backed with a urethane water resistant coating that can be easily wiped clean with liquid disinfectant. Strategically located openings allow for easy access to controls on most bed models.



HELP MINIMIZE ACCIDENTAL PATIENT INJURIES RESULTING FROM STRIKING BED SIDE RAILS

- Contour fitted design features 1" medium density foam.
- Seizure pad easily attaches with hook-and-loop fasteners.
- Allows unobstructed lowering of side rails without removal.
- Seizure pads wipe clean, water resistant covers.

BED DESCRIPTION	RIGHT/LEFT UPPER DIMENSIONS	RIGHT/LEFT LOWER DIMENSIONS	NO. OF DS	ITEM #
Hill-Rom® Beds				
Centrella™ Smart+ Bed	29 1/2"L x 10"H x 2 2/3" D	29 1/2"L x 10"H x 2 2/3" D	4	5757
Advanta™ and Century Plus® Beds	39"L x 11 1/4"H x 3"D	28 1/4"L x 8 3/4"H x 3"D	4	5725
Advanta™ (Pre-02/2004), Advance® Century 8500, and Century CC Beds	Large Left: 36"L x 11"H x 3"D Large Right: 36"L x 11"H x 3"D	Right/Left: 27 1/2"L x 11 1/2"H x 3"D Small Right/Left: 18 1/4"L x 10"H x 3"D	4	5729
CareAssist®/Advanta™ 2 Beds	42"L x 15 1/2"H x 2 1/2"D	27"L x 15"H x 2 1/2"D	4	5742
VersaCare® Med/Surg Beds	31 1/2"L x 11"H x 2 1/2"D	34 1/2"L x 10 1/2"H x 2 1/2"D	4	5743
Affinity III® and IV® Maternity Beds	31"L x 11"H x 2 1/2"D	N/A	2	5744
TotalCare® Med/Surg Beds	24 1/2"L x 11 3/4"H x 3"D	39 1/2"L x 11 3/4"H x 3"D	4	5746
Resident® LTC Beds	22"L x 7 1/2"H x 3"D	N/A	2	5747
Progressa® ICU Beds	25 1/2"L x 9 1/2"H x 2 1/2"D	40 1/2"L x 9 1/2"H x 2 1/2"D	4	5756
Low® Beds	28 1/2"L x 12"H x 2 1/2"D	N/A	2	5760
Affinity II® Maternity Beds	31"L x 12 1/2"H x 3"D	N/A	2	5761
Stryker® Beds				
Spirit Select® and Spirit Plus® Beds	23 1/2"L x 8"H x 2 1/3"D	31 1/2"L x 7 1/2" H x 2 1/3"D	4	5767
Secure I and II® Beds with rounded side rails and MPS Beds	35 1/2"L x 12"H x 3"D	27 1/2"L x 10 3/4"H x 3"D	4	5755
Secure I and II® Beds with rectangular side rails	29 1/2"L x 10 3/4"H x 3"D	19"L x 10"H x 3"D	4	5762
GO Bed II® with rectangular side rails	Center: 22"L x 6"H x 1"D Right/Left: 49"L x 12"H x 2 1/2"D	Center: 20 1/2"L x 4"H x 1"D Right/Left: 39"L x 12 1/2"H x 2 1/2"D	4	5763
In-Touch™ Critical Care Beds	29"L x 11"H x 3"D	Center: 22"L x 6"H x 1"D Right/Left: 48"L x 13"H x 3"D	4	5764
S3™ Beds	36"L x 11"H x 1"D	31"L x 11"H x 1"D	4	5766
Other Manufacturers				
Carroll® Beds with Care Assist Arms	28 1/2"L x 12"H x 2 1/2"D	N/A	2	5760
Med-Mizer® Long Term Care Beds	N/A	22"L x 6"H x 1"D	2	5765