



# Patient Safety Portfolio



# Table of Contents

---

**FALLS MANAGEMENT**

<b>4-5</b>	Alarms
<b>6-7</b>	Sensors and Sensor Pads
<b>8-10</b>	Adapters, Cables, and Brackets
<b>11</b>	Floor Cushions
<b>12-14</b>	Hipsters
<b>15-16</b>	Socks and Slippers
<b>17</b>	Fall Risk Indicators
<b>18-19</b>	Gait Belts

**PATIENT SAFETY AND PROTECTION**

<b>20</b>	Activity Products
<b>21-22</b>	Skin Sleeves
<b>23-24</b>	Mitts
<b>25-26</b>	SecureSleeves
<b>27-29</b>	Soft Limb Holders
<b>30-33</b>	Twice-As-Tough Cuffs
<b>34-38</b>	Belts
<b>39-45</b>	Jackets and Vests
<b>46</b>	Accessories

**BED SAFETY**

<b>47-49</b>	Seizure Pads
--------------	--------------

**POSITIONING**

- 50-51**      Bed Positioning
- 52**         Chair Positioning

**WOUND MANAGEMENT**

- 53**         Heel Offloading Devices
- 54**         Heel Padding Devices
- 55**         Footcare
- 56**         Heel/Elbow Protection

**SPECIALTY AND PEDIATRICS**

- 57**         Respiratory
- 58**         Pediatrics and NICU
- 59**         Hand Therapy

**POSEY BED**

- 60**         Posey Bed
- 61**         Posey Bed Accessories



8345, 8645, 8374 | Alarms

## ALARMS

Posey® Alarm Systems help reduce in-hospital patient falls and are an essential part of falls prevention.

CODE	DESCRIPTION	RECOMMENDED USE	NURSE CALL CONNECTION	CUSTOM RECORDING	AC POWER	SENSOR PORTS	MAGNET PLATE
8645	Sitter On Cue	Premium Alarm with Improved Caregiver Workflow	Yes	Yes	Yes	2	No
8345	Sitter Elite	Acute Care - Single Port Posey Fall Alarm	Yes	Yes	Yes	1	Yes
8374	KeepSafe Deluxe	Acute Care - Compact and Cost-Effective	Yes	No	Yes	1	No
8323	KeepSafe Cadet	Post Acute - Economy Solution with Magnet Cords	No	No	No	0	Yes
8324	KeepSafe Scout	Post Acute - Economy Solution	No	No	No	1	No
8373	KeepSafe Essential	Post Acute - Magnet Cord Alarm	No	No	Yes	0	Yes



8621 | Wireless Transmitter

## WIRELESS TRANSMITTERS

Compatible with the Sitter Alarm series, the battery powered Posey Wireless Transmitter and Receiver are add-on reusable components that enable Bluetooth transmission from sensor to alarm.

CODE	DESCRIPTION
8621	Wireless Transmitter
8630	Wireless Receiver



SENSORS AND SENSOR PADS

Posey Sensors and Sensor Pads, in conjunction with Posey Alarms, are important parts of your falls management protocols.

8309 | Chair Sensor Pad, Single Patient Use

CODE	DESCRIPTION	SENSOR PAD SIZE	LIFESPAN	CORD LENGTH	CASE QTYS. (EACHES)
8309	Chair Sensor Pad, Single Patient Use	13" W x 13" L (33 cm x 33 cm)	30 days	3ft	30
8309EL	Chair Sensor Pad, Single Patient Use	13" W x 13" L (33 cm x 33 cm)	30 days	8ft	30
8283	Bed Sensor Pad, Single Patient Use	12" W x 120" L (30 cm x 3 m)	30 days	8ft	30
8316	Stretcher Sensor Pad	4.5" W x 24" L (11 cm x 61 cm)	15 days	3ft	6
8316EL	Stretcher Sensor Pad	4.5" W x 24" L (11 cm x 61 cm)	15 days	8ft	6
8317	Low Air-loss Sensor Pad, Single Patient Use	6½" W x 31½" L (16½ cm x 80 cm)	30 days	8ft	12
8308	Chair Sensor, Single Resident	13" W x 13" L (33 cm x 33 cm)	6 months	3ft	30
8307	Bed Sensor, Single Resident	12" W x 120" L (30 cm x 3 m)	6 months	8ft	30





8399 | HeadStart® Notification Chair Sensor

**HEADSTART® NOTIFICATION CHAIR SENSOR**

CODE	DESCRIPTION	STRAP LENGTH	CORD LENGTH
8399	HeadStart® Notification Chair Sensor	Adjusts from 31" to 77" (79 cm to 196 cm)	8 ft
8399L	HeadStart® Notification Chair Sensor (Large)	Adjusts up to 89" (226 cm)	8 ft
8399SL	HeadStart® Notification Chair Sensor, Short Cord	Adjusts from 31" to 77" (79 cm to 196 cm)	3 ft
8399LSL	HeadStart® Notification Chair Sensor, Short Cord (Large)	Adjusts up to 89" (226 cm)	3 ft



8332 | Toilet Seat Sensor

**BATHROOM SENSORS**

CODE	DESCRIPTION	SENSOR PAD SIZE
8332	Toilet Seat Sensor	3"W x 5 $\frac{1}{2}$ "L x $\frac{3}{4}$ "D (8 cm x 14 cm x 2 cm)
8333	Bedside Patient Commode	2"W x 4"L x $\frac{1}{2}$ "D (5 cm x 10 cm x 2 cm)



**ADAPTERS, CABLES, AND BRACKETS**

Posey Alarm Adapters and Brackets make it easy to integrate alarms within your facility.

**BRACKET**

CODE	DESCRIPTION
8288	Alarm Bracket, IV Pole or Wheelchair
8289	Bracket Alarm, Wheelchair
8276	Alarm Mounting Bracket with Wire Hanger
8208	Alarm Mounting Bracket, Standard Wall Bracket
8208CS	Alarm Mounting Bracket, Standard Wall Bracket, Command Strips
8208G	Alarm Mounting Bracket, GXH Track System
8208H	Alarm Mounting Bracket, Modular Horizontal Track System
8208L	Alarm Mounting Bracket, Locking
8269	Alarm Mounting Bracket, Biothane Strap
8278N	Alarm Mounting Bracket, Saddle Narrow
8208HR	Alarm Mounting Bracket, Hillrom Horizontal Track System





8282 | Posey Nurse Call Cable



8208H | Posey Alarm Wall Bracket

## CABLES AND ADAPTERS

CODE	DESCRIPTION
8282	Nurse Call Cable, 12 ft Cord
8282SL	Nurse Call Cable, 3 ft Cord
8235NCS	1/4" Phone Plug Y-Adapter, 12 ft Cord
8235NCA	Executone Adapter (C/P Not Wired), 12 ft Cord
8235NCASL	Executone Adapter (C/P Not Wired), 3 ft Cord
8235NCC	Rauland 11-Pin M-F to S12f Adapter, 12 ft Cord
8235NCSL	Rauland 11-Pin M-F to S12f Adapter, 3 ft Cord
8235NCD	Philips/Sylvania 20-Pin M-F to S12f Adapter, 12 ft Cord
8235NCDSL	Philips/Sylvania 20-Pin M-F to S12f Adapter, 3 ft Cord
8235NCF	Executone 7-Pin Futura M-F Adapter to 1/4" Phono Jack, 12 ft Cord
8235NCFSL	Executone 7-Pin Futura M-F Adapter to 1/4" Phono Jack, 3 ft Cord
8235NCG	Zettler 9-Pin, Sentinel 400, Y-Adapter, Zettler 9-Pin M-F & 1/4" Jack, 12 ft Cord
8235NCGSL	Zettler 9-Pin, Sentinel 400, Y-Adapter, Zettler 9-Pin M-F & 1/4" Jack, 3 ft Cord
8235NCHH	Dukane 18-Pin Adapter, 12 ft Cord
8235NCHHSL	Dukane 18-Pin Adapter, 3 ft Cord
8235NCJ	Dukane 12-Pin K97 M-F to S12f Adapter, 12 ft Cord
8235NCJSL	Dukane 12-Pin K97 M-F to S12f Adapter, 3 ft Cord
8235NCJJ	Hillrom 37-Pin Adapter, 12 ft Cord
8235NCJJSL	Hillrom 37-Pin Adapter, 3 ft Cord
8235NCK	Fisher Berkeley 8-Pin Adapter C54m to S54f & S12f, 12 ft Cord
8235NCKSL	Fisher Berkeley 8-Pin Adapter C54m to S54f & S12f, 3 ft Cord
8235NCL	Executone 21-pin M-F to S12F Adapter Assy, 12 ft cord

**CABLES AND ADAPTERS**

CODE	DESCRIPTION
8235NCLSL	Executone 21-pin M-F to S12F Adapter Assy, 3 ft cord
8235NCM	Jeron 8-pin Adapter (Din Plug), 12 ft cord
8235NCMSL	Jeron 8-pin Adapter (Din Plug), 3 ft cord
8235NCMM	Rauland Responder 8-pin Adapter, 12 ft cord
8235NCMMSL	Rauland Responder 8-pin Adapter, 3 ft cord
8235NCNN	Hillrom Adapter Box D35M to T35F & 1/4" Jack w/ Dummy, 12 ft cord
8235NCNNSL	Hillrom Adapter Box D35M to T35F & 1/4" Jack w/ Dummy, 3 ft cord
8235NCP	Fisher Berkeley Adapter Box D51M to T51F & 1/4" Jack w/ Dummy, 12 ft cord
8235NCPSL	Fisher Berkeley Adapter Box D51M to T51F & 1/4" Jack w/ Dummy, 3 ft cord
8235NCPP	Dukane Box Adapter DB37P Male to DB37P Female & 1/4" Bed Exit Jack, 12 ft cord
8235NCPPSL	Dukane Box Adapter DB37P Male to DB37P Female & 1/4" Bed Exit Jack, 3 ft cord
8235NCQQ	Intercall 8-pin, 12 ft cord
8235NCQQSL	Intercall 8-pin, 3 ft cord
8235NCR	Wescom Adapter C59M to S59F & 1/4" Jack w/ Dummy, 12 ft cord
8235NCRSL	Wescom Adapter C59M to S59F & 1/4" Jack w/ Dummy, 3 ft cord
8235NCRR	Dukane B27 17-pin w/ Phone Plug Adapter, 12 ft cord
8235NCRRSL	Dukane B27 17-pin w/ Phone Plug Adapter, 3 ft cord
8235NCSSL	1/4" Phone Plug Y-Adapter, 3 ft cord
8235NCVV	Executone 16-pin Lifesave Box Adapter D16M to T16F & S12F, 12 ft cord
8235NCVVSL	Executone 16-pin Lifesave Box Adapter D16M to T16F & S12F, 3 ft cord
8235NCWW	Wesom 10-pin Modular Adapter M-F w/ 1/4" Jack & Dummy, 12 ft cord
8235NCWWSL	Wesom 10-pin Modular Adapter M-F w/ 1/4" Jack & Dummy, 3 ft cord
8235NCXX	Intercall Adapter C17M to S17F and S12F with 1/4" Dummy, 12 ft cord
8235NCXXSL	Intercall Adapter C17M to S17F and S12F with 1/4" Dummy, 3 ft cord
8235NCZZ	Zettler 9-pin Sentinel 500, 12 ft cord
8235NCZZSL	Zettler 9-pin Sentinel 500, 3 ft cord

## FLOOR CUSHIONS

Posey Floor Cushions can help reduce the risk of fall-related injuries.

6027R | Reflective Bevel Cushion, Tri-fold

CODE	DESCRIPTION	DIMENSIONS	UOM
6020	Floor Cushion, Tri-fold	67" L x 26" W x 1" H	Each
6021	Floor Cushion, Two-tone, Tri-fold	67" L x 26" W x 1" H	Each
6023	Beveled Floor Cushion, Tri-fold	70" L x 38" W x 1" H	Each
6024	Economy Floor Cushion, Bi-fold	72" L x 24" W x 1 1/2" H	Each
6025	Bevel Floor Cushion Wide, Tri-fold	70" L x 44" W x 1" H	Each
6026	Beveled Floor Cushion, Tri-fold	70" L x 29" W x 1" H	Each
6027	Beveled Floor Cushion, Tri-fold	70" L x 38" W x 1" H	Each
6027R	Reflective Bevel Cushion, Tri-fold	70" L x 38" W x 1" H	Each
6028R	Reflective Narrow Beveled, Tri-fold	70" L x 29" W x 1" H	Each



## HIPSTERS

Posey Hipsters can help reduce injury for those at high risk of hip fractures.

## HIPSTERS

CODE	DESCRIPTION	SIZE	UOM
6016S	Standard Brief	Small	Each
6016M	Standard Brief	Medium	Each
6016L	Standard Brief	Large	Each
6016XL	Standard Brief	XLarge	Each
6016XXL	Standard Brief	XXLarge	Each
6017S	Incontinent Brief Removable Pads	Small	Each
6017M	Incontinent Brief Removable Pads	Medium	Each
6017L	Incontinent Brief Removable Pads	Large	Each
6017XL	Incontinent Brief Removable Pads	XLarge	Each
6017XXL	Incontinent Brief Removable Pads	XXLarge	Each
6018M	Male Fly Brief Removable Pads	Medium	Each
6018L	Male Fly Brief Removable Pads	Large	Each
6019S	Ez-On Brief Removable Pads	Small	Each
6019M	Ez-On Brief Removable Pads	Medium	Each
6019L	Ez-On Brief Removable Pads	Large	Each
6019XL	Ez-On Brief Removable Pads	XLarge	Each
6019XXL	Ez-On Brief Removable Pads	XXLarge	Each
6030XL	Ladies' Community Hipster Removable Pads	XLarge	Each
6031S	Men's Community Hipster Removable Pads	Small	Each





6018M | Male Fly Brief Removable Pads



6017L | Incontinent Brief Removable Pads



6019S | EZ-On Brief Removable Pads



6030XL | Ladies' Community Hipster Removable Pads



## SOCKS AND SLIPPERS

Double-sided non-skid socks and slippers help prevent patient falls.



## SOCKS AND SLIPPERS

CODE	DESCRIPTION	COLOR/PATTERN	SIZE	UOM
6239G	Fall Management Socks	Green	Standard	Pair
6239LG	Fall Management Socks	Green	Large	Pair
6239LO	Fall Management Socks	Orange	Large	Pair
6239LP	Fall Management Socks	Purple	Large	Pair
6239LR	Fall Management Socks	Red	Large	Pair
6239LY	Fall Management Socks	Yellow	Large	Pair
6239O	Fall Management Socks	Orange	Standard	Pair
6239P	Fall Management Socks	Purple	Standard	Pair
6239R	Fall Management Socks	Red	Standard	Pair
6239Y	Fall Management Socks	Yellow	Standard	Pair
6240S	Non-Skid Slippers	White	Small	Pair
6240M	Non-Skid Slippers	Blue	Medium	Pair
6240M-L	Non-Skid Slippers	Brown	Med-Large	Pair
6240L	Non-Skid Slippers	Beige	Large	Pair
6240XL	Non-Skid Slippers	Navy	XLarge	Pair
6242S	Quick Dry Slipper	White	Small	Pair
6242M	Quick Dry Slipper	White	Medium	Pair
6242L	Quick Dry Slipper	White	Large	Pair
6242XL	Quick Dry Slipper	White	XLarge	Pair
6242XXL	Quick Dry Slipper	White	XXLarge	Pair
6243M	Fall Management Slippers	Red	Medium	Pair
6243L	Fall Management Slippers	Red	Large	Pair
6243XL	Fall Management Slippers	Red	XLarge	Pair
6245S	Fall Management Slippers	Purple	Small	Pair
6245M	Fall Management Slippers	Purple	Medium	Pair
6245L	Fall Management Slippers	Purple	Large	Pair
6250M	Fall Management Slippers	Yellow	Medium	Pair
6250M-L	Fall Management Slippers	Yellow	Med-Large	Pair
6250L	Fall Management Slippers	Yellow	Large	Pair
6250XL	Fall Management Slippers	Yellow	XLarge	Pair
6252Y	Fall Management Pediatric Socks	Yellow	Pediatric	Pair



6236Y | Falls Risk Indicators

## FALLS RISK INDICATORS

Fall risk indicators provide a discreet way for identifying if a patient is at risk of a fall.

CODE	DESCRIPTION	COLOR/PATTERN	UOM
6248Y	Blanket	Yellow	Each
6248BY	Blanket	Yellow	36/case
6247Y	Bracelet	Yellow Star	50/box
6251Y	Bracelet	Yellow "Fall Risk"	50/box
6233Y	Magnet	Yellow	12/box
6235Y	Essential Kit with Regular Size Socks	Yellow	1 Each
6237Y	Essential Kit with Large Size Socks	Yellow	1 Each
6236Y	Deluxe Kit with Regular Size Socks	Yellow	1 Each
6238Y	Deluxe Kit with Large Size Socks	Yellow	1 Each

Essential Kit contains: Socks, Magnet, and Bracelet

Deluxe Kit contains: Socks, Magnet, Bracelet, and Blanket



GAIT BELTS

Gait belts help reduce the risk of unwanted expenses and injury, and can serve as a low-cost solution for facilities.

6556 | Gait Belt



6531 | Gait Belt, Pastel



6528 | Gait Belt, Navy

CODE	COLOR/PATTERN	DIMENSIONS	UOM
6524	White	54 "	Each
6524L	White	74 "	Each
6556	White	66 "	50/Case
6556S	White	54 "	50/Case
6556L	White	78 "	50/Case
6525	Coral Reef	54 "	Each
6525L	Coral Reef	74 "	Each
6527	Fiesta	54 "	Each
6527L	Fiesta	74 "	Each
6528	Navy	54 "	Each
6528L	Navy	74 "	Each
6529	Southwest	54 "	Each
6529L	Southwest	74 "	Each
6531	Pastel	54 "	Each
6531L	Pastel	74 "	Each
6549	USA Flag	54 "	Each
6549L	USA Flag	74 "	Each
65252	Lavender	54 "	Each
65252L	Lavender	74 "	Each
6547Y	Yellow	54 "	Each
6547YL	Yellow Long Gait Belt	74 "	Each





6556Q | Quick-Release Belts



6546B | EZ Clean Gait Belt, Blue

## QUICK-RELEASE BELTS

CODE	COLOR/PATTERN	DIMENSIONS	UOM
6524Q	White	52"	Each
6524QL	White	71 $\frac{1}{2}$ "	Each
6556Q	White	66"	50/Case
6556QS	White	54"	50/Case
6556QL	White	78"	50/Case
6528Q	Navy	52"	Each
6528QL	Navy	71 $\frac{1}{2}$ "	Each
6531Q	Pastel	52"	Each
6531QL	Pastel	71 $\frac{1}{2}$ "	Each

## EZ CLEAN GAIT BELTS

CODE	COLOR/PATTERN	DIMENSIONS	UOM
6546	Black	63"	Each
6546B	Blue	63"	Each
6546BL	Blue	84"	Each
6546L	Black	84"	Each
6546O	Orange	63"	Each
6546OL	Orange	84"	Each
6546P	Pink	63"	Each
6546PL	Pink	84"	Each
6546Y	Yellow	63"	Each
6546YL	Yellow	84"	Each



7405 | Posey Activity Apron II

ACTIVITY PRODUCTS

Posey Activity Products are designed for patients with diminished cognitive function. The bright colors and texture variation provide visual and tactile stimulation, and the buttons and shoestring treading may aid improvement of fine motor skills.

CODE	DESCRIPTION	UOM
7400	Posey Activity Apron	Each
7405	Posey Activity Apron II	Each
7410	Happy Hands	Each



7400 | Posey Activity Apron



6000 | Skin Sleeves Arm-Medium

## SKIN SLEEVES

Posey Skin Sleeves are versatile, non-compression skin coverings that offer a protective barrier to shield your patient's fragile skin from tears, bruises, and abrasions. They help prevent patient disruption of medical treatments by camouflaging IV and wound sites.

SKIN SLEEVES

CODE	DESCRIPTION	DIMENSIONS	SKIN COLOR	COLOR/SIZE/BAND	UOM
6000S	Skin Sleeves Arm-Small	Small-Length: 15½", Circumference: 8½"	Light	Red	Pair
6000	Skin Sleeves Arm-Medium	Medium-Length: 16½", Circumference: 11"	Light	Green	Pair
6000L	Skin Sleeves Arm-Large	Large-Length: 18", Circumference: 11"	Light	Yellow	Pair
6000XL	Skin Sleeves Arm-XLarge	X-Large-Length: 20", Circumference: 13"	Light	Blue	Pair
6002	Skin Sleeves Padded Arm	Medium-Length: 16½", Circumference: 11"	Light	Green	Pair
6003S	Skin Sleeves Arm Dark-Small	Small-Length: 16½", Circumference: 11"	Dark	Green	Pair
6003	Skin Sleeves Arm Dark-Medium	Medium-Length: 16½", Circumference: 11"	Dark	Green	Pair
6001	Skin Sleeves Leg	Medium-Length: 13½", Circumference: 11"	Light	Green	Pair
6005	Skin Sleeves Padded Leg	Medium-Length: 13½", Circumference: 11"	Light	Green	Pair
6000B	Skin Sleeves Bulk Packed	Medium-Length: 16½", Circumference: 11"	Light	Green	25 Pairs



6000 | Posey Skin Sleeves Arm



6001 | Posey Skin Sleeves



6003S | Posey Skin Sleeves Arm Dark - Small





2811 | Peek-A-Boo® Mitt

## MITTS

Posey Mitts are used as an option to help protect patients who are prone to disrupting medical treatment or self-harm. A mitt hinders picking and tube pulling and is considered to be a less restrictive device that should be considered and tried before more restrictive methods are utilized.

CODE	DESCRIPTION	STYLE	TYPE	UOM
2811	Peek-A-Boo® Mitt	Padded with Inspection Flap	Fiber	Pair
2816	Posey Finger Control	Mesh Back	Fiber	Pair
2819	Posey Double Padded Mitt	Double Padded	Fiber	Pair
2810	Diversinary Mitt	Single Padded	Polysyrene Bead	Pair
2814	Posey Double-Security	Mesh Back	Fiber	Pair
2826	Posey Economy Mitt	Mesh Back	Foam	Pair
2813	Posey Secure Mitts	Mesh Back	Fiber	Pair
2815	Posey Finger Control	Open End	Fiber	Pair
2809	Mitt Connecting Strap	Connecting Strap	N/A	Pair





2810 | Diversionary Mitt

PEDIATRIC MITTS

CODE	ITEM DESCRIPTION	STYLE	FILL	UOM
2911S	Peek-A-Boo Mitt, Small	Padded with Inspection Flap	Fiber	Pair
2911M	Peek-A-Boo Mitt, Medium	Padded with Inspection Flap	Fiber	Pair
2911L	Peek-A-Boo Mitt, Large	Padded with Inspection Flap	Fiber	Pair



2911M | Peek-A-Boo Mitt, Medium



8168S, 8168M, 8168L | SecureSleeves

## SECURESLEEVES

The Posey SecureSleeve® is a multi-purpose external arm splint that helps protect tubes, IV sites, and wounds by reducing patient arm movement.

ADULT

CODE	DESCRIPTION	DIMENSIONS	LIMB CIRCUMFERENCE	UOM
8168S	Small	11"L x 14½"W (28 cm x 37 cm)	3"-14" (8-36 cm)	Each
8168M	Medium	13"L x 18"W (33 cm x 46 cm)	6"-17" (15-43 cm)	Each
8168L	Large	14"L x 25"W (36 cm x 64 cm)	8"-21" (20-53 cm)	Each
8169	Leg	18"L x 27"W (46 cm x 69 cm)	14"-27" (36-69 cm)	Each
8162	X-Ray Translucent Clip	N/A	N/A	Each

PEDIATRIC

CODE	DESCRIPTION	DIMENSIONS	LIMB CIRCUMFERENCE	UOM
8168I	Infant	3"L x 4"W (8 cm x 10 cm)	1½"-4" (4-10 cm)	Each
8168IL	Infant, Large	4½"L x 4"W (11 cm x 10 cm)	1½"-4" (4-10 cm)	Each
8168P	Pediatric	7"L x 7"W (18 cm x 18 cm)	3"-10" (8-25 cm)	Each
8168XS	Pediatric, Large	10"L x 12"W (25 cm x 30 cm)	3"-12" (8-30 cm)	Each



8168P | Pediatric



## SOFT LIMB HOLDERS

Posey Limb Holders are designed to limit the movement of limbs in patients who are at risk of disrupting life-saving treatment, such as pulling lines or tubes, and causing injury to themselves or others.

2532 | Soft Limb Holders, Quick-Release Buckle

## FOAM LIMB HOLDERS

CODE	CUFF CLOSURE	CUFF DIMENSION	CONNECTING STRAP	STRAP ATTACHMENT	STRAP LENGTH	UOM
2510	D-Ring	13" L x 3" W (33 cm x 8 cm)	Single Strap Tie	Tie	55" (140 cm)	Pair
25281	Slide Buckle	11" L x 2½" W (28 cm x 6 cm)	Double Strap Tie	Tie	56" (142 cm)	Pair
2530	Slide Buckle	13" L x 3" W (33 cm x 8 cm)	Double Strap Tie	Tie	56" (142 cm)	Pair
2531	Quick-Release Buckle	13" L x 3" W (33 cm x 8 cm)	Double Strap Tie	Tie	33" (84 cm)	Pair
2532	Quick-Release Buckle	13" L x 3" W (33 cm x 8 cm)	Quick-Release: Two-Piece	Quick-Release Buckle	31" (79 cm)	Pair
2532A	Quick-Release Buckle	13" L x 3" W (33 cm x 8 cm)	Quick-Release: Two-Piece w/ Plastic Slider	Quick-Release Buckle	31" (79 cm)	Pair
2533	Quick-Release Buckle	13" L x 3" W (33 cm x 8 cm)	Quick-Release: One-Piece	Quick-Release Buckle	38" (97 cm)	Pair
2534	Quick-Release Buckle	13" L x 3" W (33 cm x 8 cm)	Quick-Release: Two-Piece Double Strap	Quick-Release Buckle	38" (97 cm)	Pair



2531 | Double Strap Tie

2552 | Quick-Release Buckle

QUILTED LIMB HOLDERS

CODE	CUFF CLOSURE	CUFF DIMENSION	CONNECTING STRAP	STRAP ATTACHMENT	STRAP LENGTH	UOM
2540	Slide Buckle	12¼" L x 3¾" W (31 cm x 10 cm)	Double Strap Tie	Tie	55" (140 cm)	Pair
2541	Quick-Release Buckle	12¼" L x 3¾" W (31 cm x 10 cm)	Double Strap Tie	Tie	55" (140 cm)	Pair
2550	Quick-Release Buckle	11½" L x 3" W (29 cm x 8 cm)	Double Strap Tie	Tie	55" (140 cm)	Pair
2551	Quick-Release Buckle	13" L x 3" W (33 cm x 8 cm)	Quick-Release: Two-Piece	Quick-Release Buckle	31" (79 cm)	Pair
2551A	Quick-Release Buckle	13" L x 3" W (33 cm x 8 cm)	Quick-Release: Two-Piece w/ Plastic Slider	Quick-Release Buckle	31" (79 cm)	Pair
2552	Quick-Release Buckle	13" L x 3" W (33 Cm X 8 Cm)	Quick-Release: One-Piece	Quick-Release Buckle	38" (97 Cm)	Pair
2554	Quick-Release Buckle	13" L x 3" W (33 cm x 8 cm)	Quick-Release: Two-Piece Double Strap	Quick-Release Buckle	38" (97 cm)	Pair

DELUXE LIMB HOLDER

CODE	CUFF CLOSURE	CUFF DIMENSION	CONNECTING STRAP	STRAP ATTACHMENT	STRAP LENGTH	UOM
2631	Wrap	13½" L x 4¼" W (34 cm x 11 cm)	Double Strap Tie	Tie	54 1/2" (138 cm)	Pair



## PEDIATRIC LIMB HOLDERS

CODE	DESCRIPTION	STYLE	STRAP LENGTH	UOM
4733	Infant-2 years	5½"L x 1"W (14 cm x 3 cm)	26" (66 cm)	Dozen
4734	2-5 years	7"L x 1½"W (18 cm x 4 cm)	32" (81 cm)	Dozen



4733 | Infant-2 years



2700Q | Twice-As-Tough Cuffs, Three Point Buckle

## **TWICE-AS-TOUGH® CUFFS**

Posey Twice-as-Tough® Cuffs are used to protect patients assessed as being in extreme danger of harming themselves or others. They are designed to limit the movement of limbs in patients assessed to be at risk of disrupting life-saving treatment, line or tube pulling, exacerbating skin conditions, self-injury, or compromising wound site integrity.

## TWICE-AS-TOUGH CUFFS

CODE	CUFF CLOSURE	CUFF DIMENSION	CONNECTING STRAP	STRAP ATTACHMENT	STRAP LENGTH	UOM
2789Q	Quick-Release Buckle	12"L x 2½"W (30 cm x 6 cm)	Single Two-Piece	Quick-Release	45" (114 cm)	Pair
2790Q	Quick-Release Buckle	12"L x 2½"W (30 cm x 6 cm)	Double Strap	Quick-Release	45" (114 cm)	Pair
2791Q	Quick-Release Buckle	14½"L x 2½"W (37 cm x 6 cm)	Single Strap	Quick-Release	46" (117 cm)	Pair
2790	D-Ring	12"L x 2½"W (30 cm x 6 cm)	Single Strap	D-Ring	48½" (123 cm)	Pair
2791	D-Ring	14½"L x 2½"W (37 cm x 6 cm)	Double Strap	D-Ring	50" (127 cm)	Pair
2700Q	Quick-Release Buckle	12"L x 2½"W (30 cm x 6 cm)	Single Strap	Three-Point Buckle	100" (254 cm)	Pair
2700QL	Quick-Release Buckle	14½"L x 2½"W (37 cm x 6 cm)	Single Strap	Three-Point Buckle	100" (254 cm)	Pair
2792	D-Ring	12"L x 2½"W (30 cm x 6 cm)	Single Strap	Lock	46" (117 cm)	Pair
2792Q	Quick-Release Buckle	12"L x 2½"W (30 cm x 6 cm)	Single Strap	Lock	46" (117 cm)	Pair
2792QXL	Quick-Release Buckle	12"L x 2½"W (30 cm x 6 cm)	Double Strap	Lock	90" (229 cm)	Pair
2793	D-Ring	14½"L x 2½"W (37 cm x 6 cm)	Single Strap	Lock	44½" (113 cm)	Pair
2798	Lock	12"L x 2½"W (30 cm x 6 cm)	Single Strap	Lock	72" (183 cm)	Pair
2799	Double Locking	14½"L x 2½"W (37 cm x 6 cm)	Single Strap	Lock	72" (183 cm)	Pair
2900	Biothane Cuff	14½"L x 2½"W (37 cm x 6 cm)	Single Strap	Lock	Varies	Pair



2789Q | Quick-Release Buckle



2799 | Double Locking



2900 | Biothane Cuff

## 2900 CONNECTING STRAPS

CODE	ITEM DESCRIPTION	UOM
2346	Biothane Connecting Strap	Each
2350	Polypropylene Connecting Strap	Each
2351	Biothane Connecting Strap	Each

Available in multiple lengths: 18", 24", 36", 48", 60", or 72" (46 cm, 61 cm, 91 cm, 122 cm, 152 cm, or 183 cm)

CONNECTED TWICE-AS-TOUGH CUFFS

CODE	DESCRIPTION	DIMENSIONS	STRAP LENGTH	UOM
2794	Connected Wrist, Bed	12" L x 2½" W (30 cm x 6 cm)	32" L (81 cm)	Each
2795	Connected Ankle, Bed	14½" L x 2½" W (37 cm x 6 cm)	32" L (81 cm)	Each
2795L	Connected Ankle, Bed	14½" L x 2½" W (37 cm x 6 cm)	42" L (107 cm)	Each
2796	Connected Wrist, Stretcher	12" L x 2½" W (30 cm x 6 cm)	Adj. 46"-86" (117-218 cm)	Each
2797	Connected Ankle, Stretcher	14½" L x 2½" W (37 cm x 6 cm)	Adj. 46"-86" (117-218 cm)	Each



2796 | Connected Wrist, Stretcher



2794 & 2795 | Connected Twice-as-Tough Cuffs, Bed

TWICE-AS-TOUGH CUFFS STRETCHER

CODE	DESCRIPTION	DIMENSIONS	STRAP LENGTH	UOM
2750	Stretcher, Regular (Blue)	12" L x 2½" W (30 cm x 6 cm)	9" L x 2" W (23 cm x 5 cm)	Pair
2755	Stretcher, Large (Red)	14½" L x 2½" W (37 cm x 6 cm)	9" L x 2" W (23 cm x 5 cm)	Pair



2750 | Stretcher, Regular (Blue)

SPECIALTY TWICE-AS-TOUGH CUFFS

CODE	DESCRIPTION	CUFF DIMENSIONS	WAIST BELT DIMENSIONS	PAD DIMENSIONS	UOM
2800	Twice-as-Tough Ambulatory Belt	12" L x 2½" W (30 cm x 6 cm)	48" L x 2½" W (122 cm x 6 cm)	N/A	Each
2801	Twice-as-Tough Leg Hobble	14½" L x 2½" W (37 cm x 6 cm)	N/A	13½" L x 2" W (34 cm x 5 cm)	Each





1135 | Posey Self-Releasing Roll Belt

## SELF-RELEASING BELTS

Posey Self-Releasing Belts help reduce the risk of forward sliding and act as a reminder for patients to ask for assistance with ambulation.



1135 | Posey Self-Releasing Roll Belt



3656 | The Posey Torso Support

## SELF-RELEASING BELTS

Posey Self-Releasing Roll Belts allow a patient to sit up or roll from side to side in bed.

CODE	DESCRIPTION	DIMENSIONS	UOM
1135	Posey Self-Releasing Roll Belt	4 Straps: 33" L x 2" W, Waist Circumference: 36"-54"	Each
1135QR	Posey Self-Releasing Roll Belt - Quick-Release	2 Straps: 34" L x 2" W, Waist Circumference: 36" x 54"	Each
1135QXL	Posey Self-Releasing Roll Belt, XL	4 Straps: 41" L x 2" W, Waist Circumference: 54" x 68"	Each

## TORSO SUPPORT

The Posey Torso Support provides extra support around the chest and over the shoulders to help prevent leaning and forward sliding.

CODE	DESCRIPTION	DIMENSIONS	UOM
3656S	Posey Torso Support	Small, Waist: 54" L x 5"W, Shoulder: 34" x 3"	Each
3656M	Posey Torso Support	Medium, Waist: 54" L x 5"W, Shoulder: 37" x 3"	Each
3656L	Posey Torso Support	Large, Waist: 60" L x 5" W, Shoulder: 37" x 3"	Each
3656XL	Posey Torso Support	XLarge, Waist: 72" L x 5" W, Shoulder: 41" x 3"	Each
3656XXL	Posey Torso Support	XXLarge, Waist: 84" L x 5" W, Shoulder: 46" x 3"	Each

## SELF-RELEASING WRAP-AROUND

The Self-Releasing Wrap-Around is ideal for residents with minimal upper trunk control who require positioning support.

CODE	DESCRIPTION	DIMENSIONS	UOM
3658S-M	Posey Self-Releasing Wrap-Around	Medium, 66" L x 5" W	Each
3658M-L	Posey Self-Releasing Wrap-Around	Large, 84" L x 5" W	Each
3658XL	Posey Self-Releasing Wrap-Around, XL	XLarge, 100" L x 5" W	Each
3658XXL	Posey Self-Releasing Wrap-Around, XXL	XXLarge, 121" L x 5" W	Each

## UNIVERSAL WRAP-AROUND STRAP

Padded belts help add extra comfort and are positioned low over thighs to help prevent sliding.

CODE	DESCRIPTION	DIMENSIONS	UOM
3659L	Universal Wrap-Around Strap	84" L x 5" W	Each
3659XXL	Universal Wrap- Around Strap	121" L x 5" W	Each

## SELF-RELEASING OMNI BELT

The Posey Self-Releasing Omni Belt helps remind patients to call for assistance before exiting a bed or wheelchair.

CODE	DESCRIPTION	DIMENSIONS	UOM
4116PQ	Posey Self-Releasing Omni Belt	92" L x 2" W	Each
4126Q	Posey Self-Releasing Padded Belt	2 Pads: 13" L x 4" W	Each
4126V	Posey Self-Releasing Padded Belt	W Pads: 12.5" L x 4" W	Each

DISPOSABLE KNEE AND BODY STRAP

Posey Disposable Knee and Body Straps are designed to assist in patient positioning for medical treatment on an examination or operating table.

CODE	DESCRIPTION	DIMENSIONS	UOM
5551	Disposable Knee and Body Strap	Overall: 33" L x 3" W	Each



3658L | Posey Self-Releasing Wrap-Around



4126Q | Padded belts



5551 | Disposable Knee and Body Strap



## RESTRAINT BELTS

Posey offers several types of restraint belts, including the Roll Belt and Pelvic Soft Belt, which help prevent sliding, and the General Purpose Belt, which serves as a fifth point restraint.

1201 | Restraint Belt

### RESTRAINT BELTS

CODE	DESCRIPTION	DIMENSIONS	UOM
1231	Posey Roll Belt Cotton	17½" with small ½"	Each
1231Q	Quick-Release Roll Belt	Pad: 17½" L x 6" W; Strap: Long 59" L, Short 33" W	Each
1231QXL	Posey Roll Belt, Size XL	Pad: 28½" L x 6" W; Strap: Long 82" L, Short 50" W	Each
4125	Posey Soft Belt w/ Ties	Pad: 21½" L x 4½" W; Straps: 73"	Each
4125L	Posey Long Soft Belt	Pad: 32" L x 4½" W; Straps: 73"	Each
4125Q	Quick-Release Soft Belt	Pad: 21½" L x 4½" W; Straps: 68"	Each
4125C	Pelvic Belt	Pad: 21½" L x 5" W; Strap: Pelvic Belt 72", 2 Venter Straps, 72" & 45"	Each
5550	Posey General Purpose	66" L x 3" W	Each
5550B	Posey General Purpose	90" L x 3" W	Each



4125 | Posey Soft Belt w/Ties



5550 | Posey General Purpose



1334 | Posey Key Lock Belt



1337 | Posey Swedish Belt

LOCKING BELTS

The Posey Key Lock Belt allows the patient to sit up or move from side to side in the bed. The interlocking harness helps prevent the patient from sliding out of the belt.

The Swedish Belt is a waist restraint belt with a key lock that is independent of the side straps.

CODE	DESCRIPTION	DIMENSIONS	STRAP ATTACHMENT	UOM
1334	Posey Key Lock Belt	12" x 4", Pad Waist Sizes 24-56", Straps 42"	2 Bed Attachment Straps: 42"	Each
1337	Posey Swedish Belt	53½"	2 Side Straps: 30"	Each






## SLEEVED JACKET

The Posey Sleeved Jacket is for patients at risk for falls or those requiring a restraint to assist with medical treatment. The sleeved jacket is used for positioning in bed, a wheelchair, or geri-chair.

## SLEEVED JACKET

CODE	DESCRIPTION	STRAP LENGTH	MATERIAL	STRAP TYPE	DIMENSIONS	UOM
3350XS	Posey Sleeved Jacket	65" (165 cm)	Poly/Cotton Plaid with Zipper	A-Single Strap Tie	XSmall	Each
3350S	Posey Sleeved Jacket	65" (165 cm)	Poly/Cotton Plaid with Zipper	A-Single Strap Tie	Small	Each
3350M	Posey Sleeved Jacket	65" (165 cm)	Poly/Cotton Plaid with Zipper	A-Single Strap Tie	Medium	Each
3350L	Posey Sleeved Jacket	65" (165 cm)	Poly/Cotton Plaid with Zipper	A-Single Strap Tie	Large	Each
3350XL	Posey Sleeved Jacket	75" (191 cm)	Poly/Cotton Plaid with Zipper	A-Single Strap Tie	XLarge	Each
3350XXL	Posey Sleeved Jacket	87" (221 cm)	Poly/Cotton Plaid with Zipper	A-Single Strap Tie	XXLarge	Each
3350XXXL	Posey Sleeved Jacket	92" (234 cm)	Poly/Cotton Plaid with Zipper	A-Single Strap Tie	XXXLarge	Each
3353S	Posey Sleeved Jacket	65" (165 cm)	Breezeline Mesh with Zipper	A-Single Strap Tie	Small	Each
3353M	Posey Sleeved Jacket	65" (165 cm)	Breezeline Mesh with Zipper	A-Single Strap Tie	Medium	Each
3353L	Posey Sleeved Jacket	65" (165 cm)	Breezeline Mesh with Zipper	A-Single Strap Tie	Large	Each
3353XL	Posey Sleeved Jacket	75" (191 cm)	Breezeline Mesh with Zipper	A-Single Strap Tie	XLarge	Each
3360S	Sleeved Jacket w/ QR	65" (165 cm)	Poly/Cotton Plaid with QR Buckle	C-Single Strap QR	Small	Each
3360M	Sleeved Jacket w/ QR	65" (165 cm)	Poly/Cotton Plaid with QR Buckle	C-Single Strap QR	Medium	Each
3360L	Sleeved Jacket w/ QR	65" (165 cm)	Poly/Cotton Plaid with QR Buckle	C-Single Strap QR	Large	Each
3360XL	Sleeved Jacket w/ QR	75" (191 cm)	Poly/Cotton Plaid with QR Buckle	C-Single Strap QR	XLarge	Each
3360XXL	Sleeved Jacket w/ QR	87" (221 cm)	Poly/Cotton Plaid with QR Buckle	C-Single Strap QR	XXLarge	Each
3360XXXL	Sleeved Jacket w/ QR	92" (234 cm)	Poly/Cotton Plaid with QR Buckle	C-Single Strap QR	XXXLarge	Each



## ROLL JACKET

The Posey Roll Jacket is for patients at risk for falls or those requiring a restraint to assist with medical treatment. The roll jacket allows the patient to roll from side to side and sit up in bed. The V-neckline and off-centered rear zipper maximize the patient's comfort.

CODE	DESCRIPTION	STRAP LENGTH	MATERIAL	STRAP TYPE	UOM
3320S	Posey Roll Jacket, Plaid	76" (193 cm)	Poly/Cotton Plaid	A-Single Strap Tie	Each
3320M	Posey Roll Jacket, Plaid	76" (193 cm)	Poly/Cotton Plaid	A-Single Strap Tie	Each
3320L	Posey Roll Jacket, Plaid	76" (193 cm)	Poly/Cotton Plaid	A-Single Strap Tie	Each
3320XL	Posey Roll Jacket, Plaid	87" (221 cm)	Poly/Cotton Plaid	A-Single Strap Tie	Each
3320XXL	Posey Roll Jacket, Plaid	87" (221 cm)	Poly/Cotton Plaid	A-Single Strap Tie	Each



IV ACCESS JACKET

Posey IV Access Jackets are for patients at risk for falls and/or who require a restraint to assist with medical treatment. The two zippers located on the back shoulders of the jacket provide access to IV sites and monitoring lines.

3100M | IV Access Jacket

CODE	DESCRIPTION	STRAP LENGTH	MATERIAL	STRAP TYPE	DIMENSIONS	UOM
3100S	IV Access Jacket	64" (163 cm)	Poly/Cotton	A-Single Strap Tie	Small	Each
3100M	IV Access Jacket	64" (163 cm)	Poly/Cotton	A-Single Strap Tie	Medium	Each
3100L	IV Access Jacket	64" (163 cm)	Poly/Cotton	A-Single Strap Tie	Large	Each
3100XL	IV Access Jacket	75" (191 cm)	Poly/Cotton	A-Single Strap Tie	XLarge	Each
3110S	IV Access Jacket	64" (163 cm)	Poly/Cotton	A-Single Strap Tie	Small	Each
3110M	IV Access Jacket	64" (163 cm)	Poly/Cotton	A-Single Strap Tie	Medium	Each
3110L	IV Access Jacket	64" (163 cm)	Poly/Cotton	A-Single Strap Tie	Large	Each
3110XL	IV Access Jacket	75" (191 cm)	Poly/Cotton	A-Single Strap Tie	XLarge	Each

SIZING TABLE

Always use the proper product. Products that are too small or large may compromise patient comfort and could result in severe injury or death.

BINDING COLOR	SIZE	WEIGHT lb. (kg.)	CHEST in. (cm)
White	X-Small	60-115 (27-52)	25-32 (64-81)
Red	Small	112-160 (51-73)	31-37 (79-94)
Green	Medium	135-203 (61-92)	35-40 (89-102)
Yellow	Large	160-225 (73-102)	38-44 (97-112)
Blue	X-Large	180-247 (82-112)	42-48 (107-122)
Black	XX-Large	220-275 (100-125)	46-55 (117-140)
Yellow/Black	XXX-Large	265-305 (120-138)	54-60 (137-152)
Blue/Black	XXXX-Large	295-340 (133-154)	58-64 (147-163)

• Posey Belts are not color-coded, but are sized according to this table. • Flame-retardant fabric is available on request. • Patient weight and size are a general indicator only. Consider individual physical characteristics to choose the right product for each patient. Refer to product label for specific sizing information.



3053 | Safety Vest

## SAFETY VEST

The Posey Safety Vest is for patients at risk for falls or those requiring a restraint to assist with medical treatment. The safety vests are sleeveless restraints that help prevent unassisted bed or chair exits.



3060 | Posey Safety Vest



3050 | Posey Safety Vest



## SAFETY VEST

CODE	DESCRIPTION	STRAP LENGTH	MATERIAL	STRAP TYPE	DIMENSIONS	UOM
3050S	Posey Safety Vest	64" (163 cm)	Poly/Cotton Plaid with Zipper	A-Single Strap Tie	Small	Each
3050M	Posey Safety Vest	64" (163 cm)	Poly/Cotton Plaid with Zipper	A-Single Strap Tie	Medium	Each
3050L	Posey Safety Vest	64" (163 cm)	Poly/Cotton Plaid with Zipper	A-Single Strap Tie	Large	Each
3050XL	Posey Safety Vest	75" (191 cm)	Poly/Cotton Plaid with Zipper	A-Single Strap Tie	XLarge	Each
3050XXL	Posey Safety Vest	83" (211 cm)	Poly/Cotton Plaid with Zipper	A-Single Strap Tie	XXLarge	Each
3050XXXL	Posey Safety Vest	88" (224 cm)	Poly/Cotton Plaid with Zipper	A-Single Strap Tie	XXXLarge	Each
3050XXXXL	Posey Safety Vest	93" (236 cm)	Poly/Cotton Plaid with Zipper	A-Single Strap Tie	XXXXLarge	Each
3060S	Posey Safety Vest	64" (163 cm)	Poly/Cotton Plaid with Zipper	C-Single Strap QR	Small	Each
3060M	Posey Safety Vest	64" (163 cm)	Poly/Cotton Plaid with Zipper	C-Single Strap QR	Medium	Each
3060L	Posey Safety Vest	64" (163 cm)	Poly/Cotton Plaid with Zipper	C-Single Strap QR	Large	Each
3060XL	Posey Safety Vest	75" (191 cm)	Poly/Cotton Plaid with Zipper	C-Single Strap QR	XLarge	Each
3060XXL	Posey Safety Vest	83" (211 cm)	Poly/Cotton Plaid with Zipper	C-Single Strap QR	XXLarge	Each
6063S	Posey Safety Vest	64" (163 cm)	Breezeline Mesh with Zipper	C-Single Strap QR	Small	Each
3063M	Posey Safety Vest	64" (163 cm)	Breezeline Mesh with Zipper	C-Single Strap QR	Medium	Each
3063L	Posey Safety Vest	64" (163 cm)	Breezeline Mesh with Zipper	C-Single Strap QR	Large	Each
3063XL	Posey Safety Vest	75" (191 cm)	Breezeline Mesh with Zipper	C-Single Strap QR	XLarge	Each
3063XXL	Posey Safety Vest	83" (211 cm)	Breezeline Mesh with Zipper	C-Single Strap QR	XXLarge	Each
3070S	Posey Safety Vest	64" (163 cm)	Poly/Cotton Plaid with Hook & Loop	A-Single Strap Tie	Small	Each
3070M	Posey Safety Vest	64" (163 cm)	Poly/Cotton Plaid with Hook & Loop	A-Single Strap Tie	Medium	Each
3070L	Posey Safety Vest	64" (163 cm)	Poly/Cotton Plaid with Hook & Loop	A-Single Strap Tie	Large	Each
3070XL	Posey Safety Vest	75" (191 cm)	Poly/Cotton Plaid with Hook & Loop	A-Single Strap Tie	XLarge	Each
3070XXL	Posey Safety Vest	83" (211 cm)	Poly/Cotton Plaid with Hook & Loop	A-Single Strap Tie	XXLarge	Each

## CRISS-CROSS VEST

The Posey Criss-Cross Vest is for patients at risk for falls or those requiring a restraint to assist with medical treatment. The criss-cross vest is a sleeveless restraint that can be used in either a bed or a chair.

3311M | Criss-Cross Vest

CODE	DESCRIPTION	DIMENSIONS	STRAP LENGTH	MATERIAL	STRAP TYPE	UOM
3311S	Posey Criss-Cross Vest	58" (147 cm)	Breezeline Mesh	A-Single Strap Tie	Small	Each
3311M	Posey Criss-Cross Vest	61" (155 cm)	Breezeline Mesh	A-Single Strap Tie	Medium	Each
3311L	Posey Criss-Cross Vest	64" (163 cm)	Breezeline Mesh	A-Single Strap Tie	Large	Each
3311XL	Posey Criss-Cross Vest	75" (191 cm)	Breezeline Mesh	A-Single Strap Tie	XLarge	Each

ACCESSORIES

CODE	DESCRIPTION	UOM
2785	Posey Crisis Bag	Each
1074	Posey Key	Each
2370	Universal Anchor Strap	Each



2370 | Universal Anchor Strap





## SEIZURE PADS

Seizure pads are form-fitted to fit snugly over existing bed side rails to help minimize accidental patient injuries resulting from striking bedside rails.

5757 | Centrella Smart + Bed

### SEIZURE PADS - HILL-ROM BEDS

CODE	DESCRIPTION	RIGHT/LEFT UPPER DIMENSIONS	RIGHT/LEFT LOWER DIMENSIONS	NO. OF PADS
5757	Centrella Smart + Bed	29 $\frac{1}{2}$ " L x 10" W x 2 $\frac{2}{3}$ " H	29 $\frac{1}{2}$ " L x 10" W x 2 $\frac{2}{3}$ " H	4
5725	Advanta and Century Plus Beds	39" L x 11 $\frac{1}{4}$ " W x 3" H	28 $\frac{1}{4}$ " L x 8 $\frac{3}{4}$ " W x 3" H	4
5729	Advanta (Pre-02/2004), Advance Century 8500, and Century CC Beds	Large Left: 36" L x 11" W x 3" H Large Right: 36" L x 11" W x 3" H	Right/Left: 27 $\frac{1}{2}$ " L x 11 $\frac{1}{2}$ " W x 3" H Small Right/Left: 18 $\frac{1}{4}$ " L x 10" W x 3" H	4
5742	Care Assist/Advanta 2 Beds	42" L x 15 $\frac{1}{2}$ " W	27" L x 15" W x 2 $\frac{1}{2}$ " H	4
5743	VersaCare Med/Surg Beds	31 $\frac{1}{2}$ " L x 11" W x 2 $\frac{1}{2}$ " H	34 $\frac{1}{2}$ " L x 10 $\frac{1}{2}$ " W x 2 $\frac{1}{2}$ " H	4
5744	Affinity III and IV Maternity Beds	31" L x 11" W x 2 $\frac{1}{2}$ " H	N/A	2
5746	TotalCare Med/Surg Beds	24 $\frac{1}{2}$ " L x 11 $\frac{3}{4}$ " W x 3" H	39 $\frac{1}{2}$ " L x 11 $\frac{3}{4}$ " W x 3" H	4
5747	Resident LTC Beds	22" L x 7 $\frac{1}{2}$ " W x 3" H	N/A	2
5756	Progressa ICU Beds	25 $\frac{1}{2}$ " L x 9 $\frac{1}{2}$ " W x 2 $\frac{1}{2}$ " H	40 $\frac{1}{2}$ " L x 9 $\frac{1}{2}$ " W x 2 $\frac{1}{2}$ " H	4
5760	Low Beds	28 $\frac{1}{2}$ " L x 12" W x 2 $\frac{1}{2}$ " H	N/A	2
5761	Affinity II Maternity Beds	31" L x 12 $\frac{1}{2}$ " W x 3" H	N/A	2

## SEIZURE PADS - STRYKER BEDS

CODE	DESCRIPTION	RIGHT/LEFT UPPER DIMENSIONS	RIGHT/LEFT LOWER DIMENSIONS	NO. OF PADS
5767	Spirit Select and Spirit Plus Beds	23 $\frac{1}{2}$ " L x 8" W x 2 $\frac{1}{3}$ " H	31 $\frac{1}{2}$ " L x 7 $\frac{1}{2}$ " W x 2 $\frac{1}{3}$ " H	4
5755	Secure I and II Beds with rounded side rails and MPS Beds	35 $\frac{1}{2}$ " L x 12" W x 3" H	27 $\frac{1}{2}$ " L x 10 $\frac{3}{4}$ " W x 3" H	4
5762	Secure I and II Beds with rectangular side rails	29 $\frac{1}{2}$ " L x 10 $\frac{3}{4}$ " W x 3" H	19" L x 10" W x 3" H	4
5763	GO Bed II with rectangular side rails	Center: 22" L x 6" W x 1" H Right/Left: 49" L x 12" W x 2 $\frac{1}{2}$ " H	Center: 20 $\frac{1}{2}$ " L x 4" W x 1" H Right/Left: 39" L x 12 $\frac{1}{2}$ " W x 2 $\frac{1}{2}$ " H	4
5764	In-Touch Critical Care Beds	29" L x 11" W x 3" H	Center: 22" L x 6" W x 1" H Right/Left: 48" L x 13" W x 3" H	4
5766	S3 Beds	36" L x 11" W x 1" H	31" L x 11" W x 1" H	4



5767 | Spirit Select and Spirit Plus Beds



5756 | Progressa ICU Beds


## SEIZURE PADS - OTHER MANUFACTURERS

CODE	DESCRIPTION	RIGHT/LEFT UPPER DIMENSIONS	RIGHT/LEFT LOWER DIMENSIONS	NO. OF PADS
5760	Carroll Beds with Care Assist Arms	28 $\frac{1}{2}$ " L x 12" W x 2 $\frac{1}{2}$ " H	N/A	2
5765	Med-mizer Long Term Care Beds	N/A	22" L x 6" W x 1" H	2



## SIDE RAIL PADS

CODE	DESCRIPTION	RIGHT/LEFT UPPER DIMENSIONS	UOM
5705	Soft Split Side Rail Protectors	36"L x 16"H x 1/2" Thick (91 cm x 41 cm x 1 cm)	Pair
5705LS	Soft Split Side Rail Protectors	36"L x 16"H x 1/2" Thick 27 1/2"L x 10 3/4"W x 3"H	Pair
5708	Side Rail Wedges	35"L x 17"H x 2"D (89 cm x 43 cm x 5 cm)	Pair
5605	Deluxe Guard Rail Pads - 72"	72"L x 14"H (183 cm x 36 cm)	Pair
5606	Deluxe Guard Rail Pads - 60"	60"L x 14"H (152 cm x 36 cm)	Pair
5607	Deluxe Guard Rail Pads - 43"	43"L x 14"H (109 cm x 36 cm)	Pair
5709	Deluxe Half Guard Rail Pads	40"L x 16"H x 1"D (102 cm x 41 cm x 3 cm)	Pair
5730	Horseshoe Wedge	43"L x 14"H (109 cm x 36 cm)	Each
5715	Gap Filler	15" L x 38" W x 28" H (38cm x 97cm x 71cm)	Pair
5707	Elastic Mesh Side Rail Protectors	35"L x 6"H x 2"D (89 cm x 15 cm x 5 cm)	Each
5707H	Elastic Mesh Side Rail Protectors - Half	Fit rails up to 72"L (183 cm)	Each



5705 | Soft Split Side Rail Protectors



5607 | Deluxe Guard Rail Pads - 43"



5700 | Roll Guard

## BED POSITIONING

Posey positioning products provide upper torso support and are ideal for patients who lean laterally.





6309M | Lateral Wedge, Medium



5700 | Roll Guard

CODE	DESCRIPTION	DIMENSIONS	UOM
5700	Roll Guard	Overall: 25"L x 43"W (64 cm x 109 cm)	1 Each
5700SC	Roll Guard w/Vinyl Cover	Overall: 25"L x 43"W (64 cm x 109 cm)	1 Each
5716	Soft Rails Double Bolster w/ Brushed Polyester Cover	Foundation: 33"L x 33"W (84 cm x 84 cm)	1 Each
5716SC	Soft Rails Double Bolster w/ Vinyl Cover	Foundation: 33"L x 33"W (84 cm x 84 cm)	1 Each
5750	Foam Mattress Cover	6"H to 8"H (15 cm x 20 cm), 34"W to 39"W (86 cm x 99 cm), 80"L (203 cm)	1 Each
5750NC	Perimeter Mattress Cover	6"H to 8"H (15 cm x 20 cm), 34"W to 39"W (86 cm x 99 cm), 80"L (203 cm)	1 Each
5751	Mattress Cover, Bariatric	6"H to 10"H (15 cm x 25 cm), 48"W (122 cm), 80"L to 86"L (203 cm x 218 cm)	1 Each
6302S	Abduction Wedge, Small	8" (20 cm) Deep by 5" (13 cm) Thick with a 7" (18 cm) Front Width Tapering to 5" (13 cm) at the Rear	1 Each
6302L	Abduction Wedge, Large	22"L x 5"H (56 cm x 13 cm) and Tapers from 15" (38 cm) at the Feet to 7" (18 cm) at the Thighs	1 Each
6306	Posey Bed Fellow	64"L x 14"W (163 cm x 36 cm)	1 Each
6307	Posey Separator	7"H x 4"W x 5"D (18 cm x 10 cm x 13 cm) with Two 20" (51 cm) Straps	1 Each
6309M	Lateral Wedge, Medium	16"L x 6"W x 7"H (41 cm x 15 cm x 18 cm)	1 Each
6309L	Lateral Wedge, Large	19"L x 7"W x 9"H (48 cm x 18 cm x 23 cm)	1 Each
6309XL	Lateral Wedge, XLarge	23"L x 7"W x 9"H (58 cm x 18 cm x 23 cm)	1 Each



6542 | E-Z Release Hugger



8215 | Urine Drainage Bag Holder/Cover Canvas with Straps

## CHAIR POSITIONING

Posey offers seating and accessories to make sitting in a chair or wheelchair easier for patients and residents.

CODE	DESCRIPTION	DIMENSIONS	UOM
8215	Posey Urine Drainage Bag Holder/Cover Canvas with Straps	13 1/2" L x 10 1/2" W (34 cm x 27 cm)	Each
8215V	Posey Urine Drainage Bag Holder/Cover Vinyl with Straps	14 1/2" L x 10 1/2" W (34 cm x 27 cm)	Each
8270	Posey Urine Bottle Holder	11 1/2" L x 7" W (29 cm x 18 cm)	Each
8275	Urine Bag Cover w/Window, No Straps	11" L x 10" W (28 cm x 25 cm)	Each
6351	Foot Hugger 1"	13"W x 9"D x 1"H (33 cm x 23 cm x 3 cm)	Each
6542	E-Z Release Hugger	Fits 18-22" (46 cm x 56 cm) Standard Full or Desk Arm Wheelchair	Each
6526G	Posey Comfy Seat For Geri Chair	Seat 22"W x 20"D (56 cm x 51 cm), Back 22"W x 31"H (56 cm x 79 cm), Side 20"D x 20"H (51 cm x 51 cm)	Each
6355	Leg Hugger	15 1/2"W x 9"H (39 cm x 23 cm)	Each
6526	Posey Comfy Seat For Wheelchair	Seat 20"W x 18"D (51 cm x 46 cm), Back 20"W x 17"H (51 cm x 43 cm), Side 15"D x 14"H (38 cm x 36 cm)	Each
6320	Posey Deluxe Lateral Arm Support	16" L x 7" W (41 cm x 18 cm)	Each



6128 | PRO-heeLx Heel Protector

## HEEL OFFLOADING DEVICES

PRO-heeLx® is designed to float the heels above the bed surface and helps protect against heel pressure injuries. The Posey Premium Heel Guard is a high resilience foam heel offloading device that floats the heel and helps reduce the risk of foot drop.

CODE	DESCRIPTION	DIMENSIONS	WOMEN'S SHOE SIZE	MEN'S SHOE SIZE	UOM
6145	Posey Premium Heel Guard	14" L x 11" H	N/A	N/A	Each
6218S	PRO-heeLx Small Heel Protector	N/A	Up to 6	N/A	Each
6128	PRO-heeLx Heel Protector	N/A	7–11	Up to 11	Each
6128L	PRO-heeLx Large Heel Protector	N/A	N/A	12+	Each
6129	PRO-heeLx Stabilizer Wedge	Fits all PRO-heeLx sizes	N/A	N/A	Each





6148SM | Podus w/ Ambulation Sole



6145 | Posey Premium Heel Guard

HEEL PADDING DEVICES

The multi-purpose foot boot helps in the prevention of heel and toe ulcers and safeguards against foot drop. Heel guards help isolate the heel to help protect against pressure injuries.

CODE	DESCRIPTION	DIMENSIONS	UOM
6147SM	Deluxe Podus Boot, S-M	Calf Circ: 13–16", Foot Circ: 7–10"	Each
6147ML	Deluxe Podus Boot, M-L	Calf Circ: 15–18", Foot Circ: 9–14"	Each
6148SM	Podus w/ Ambulation Sole	Calf Circ: 13–16", Foot Circ: 7–10"	Each
6148ML	Podus w/ Ambulation Sole	Calf Circ: 15–18", Foot Circ: 9–14"	Each
6127	Posey Heel Guard	12" L x 7½" H	Pair

## FOOTCARE

Posey footcare products help reduce patient discomfort by minimizing friction and shearing, offering calf support, keeping sheets and blankets off of the patient's feet, and raising the heel to reduce the risk of pressure injuries.

6113 | Quilted Long Boot, Plaid

CODE	DESCRIPTION	DIMENSIONS	UOM
6113	Quilted Long Boot, Plaid	14" L x 9" H	Pair
6115	Posey Quilted Heel	10" L x 9" H	Pair
6115L	Quilted Heel Pillow, Large	11½" L x 11" H	Pair
6432	Posey Foot Elevator 4"	Accommodates limbs up to 12" circumference	Each
6433	Posey Foot Elevator 6"	Accommodates limbs up to 12" circumference	Each
6430	Bed Cradle & Foot Support	36" L x 21½" H	Each
6436	Foot Support, Cover Only for 6430	36" L x 21½" H	Each
6435	Bed Cradle, High Top Frame Only	36" L x 24" H	Each
6437	Foot Support, Cover Only for 6435	36" L x 24" H	Each



6113 | Quilted Long Boot, Plaid



6115 | Posey Quilted Heel





## HEEL-ELBOW PROTECTION

This sock-like covering helps provide comfort and helps protect against limb friction and shear forces.

6224M | Heel-Elbow Protector, Medium

CODE	DESCRIPTION	DIMENSIONS	UOM
6224S	Heel-Elbow Protector, Small	Small, up to 8½" Circumference	Pair
6224M	Heel-Elbow Protector, Medium	Medium, up to 11" Circumference	Pair
6224L	Heel-Elbow Protector, Large	Large, up to 13½" Circumference	Pair
6224XL	Heel-Elbow Protector, XLarge	XLarge, up to 16" Circumference	Pair



8197M | Posey Foam Trach Ties



8199 | Posey Cufflator™

## RESPIRATORY THERAPY

Posey offers easy-to-use methods for securing patient tracheal tubes and pulse oximeters. The Posey Cufflator helps to eliminate the need for syringes, electricity, or complicated connections required to check cuff pressure.

CODE	DESCRIPTION	DIMENSIONS	PATIENT TYPE	UOM
8199	Posey Cufflator	N/A	N/A	Each
8197S	Posey Foam Trach Ties	9" L x 1" W (23 cm x 3 cm)	Neonatal and Infants 7-9" (18-23 cm)	12
8197M	Posey Foam Trach Ties	18½" L x 1" W (47 cm x 3 cm)	Adolescent and Adult 9-14" (23-35 cm)	12
8197L	Posey Foam Trach Ties	23" L x 1" W (60 cm x 3 cm)	Adult 13-19" (33-48 cm)	12
8197XL	Posey Foam Trach Ties	30" L x 1" W (70 cm x 3 cm)	Bariatric and Adult 16-24" (41-61 cm)	12
8196S	Posey Trach Tube Secure	9" L x 1" W (23 cm x 3 cm)	Neonatal and Infants 7-9" (18-23 cm)	12
8196M	Posey Trach Tube Secure	18½" L x 1" W (47 cm x 3 cm)	Adolescent and Adult 9-14" (23-35 cm)	12
8196L	Posey Trach Tube Secure	23" L x 1" W (60 cm x 3 cm)	Adult 13-19" (33-48 cm)	12





4648 | Posey Movable ID Bracelets, Newborn



4644 | Posey Eye Protector, Small



8141 | Posey IV Shield

PEDIATRICS & NICU

Posey offers several products to help monitor and protect pediatric patients.

CODE	DESCRIPTION	UOM
4644	Posey Eye Protector, Small	Dozen
4645	Posey Eye Protector, Premie	Dozen
4646	Posey Eye Protector, Newborn	Dozen
4648	Posey Movable ID Bracelets, Newborn	Dozen
4649	Posey Movable ID Bracelets, Premie	Dozen
4749	Posey Circumstraint Board Straps	Dozen
8141	Posey IV Shield	Dozen



6507 | Posey Palm Grip-Small

## HAND THERAPY

Posey offers products that can help provide anti-contracture hand therapy without restricting finger movement, and can act as a diversionary device for individuals prone to pulling at IV or feeding tubes.

CODE	DESCRIPTION	DIMENSIONS	COLOR/SIZE/ BAND	UOM
5608	Posey Palm Cone	5"L (13 cm)	Yellow	1 each
6507	Posey Palm Grip-Small	5"L x 2½"W (13 cm x 6 cm)	Red	1 each
6508	Posey Palm Grip-Large	5"L x 3"W (13 cm x 8 cm)	Green	1 each
6560	Posey Finger Contracture	5"L x 3"W (13 cm x 8 cm)	Blue	1 each



8070 | Posey Bed All-In-One

**POSEY BED**

The Posey Bed is a patented bed enclosure system and is a less restrictive restraint alternative to physical restraints such as belts, vest, or jackets. The Posey Bed helps to reduce the risk of unassisted bed exits and patient falls and provides a safe patient environment.

**POSEY BED**

CODE	DESCRIPTION	DIMENSIONS (SIZE)	UOM
8060	Posey Bed A-Frame Canopy + Frame	51" L x 36" W x 80" H	Each
8060T	8060 Replacement Tent	Fits 8060	Each
8070	Posey Bed All-In-One	48" L x 36" W x 60" H	Each
8070T	8070 Replacement Canopy	Fits 8070	Each

As of December 2019, TIDI Products no longer sells Posey Beds into the European Union.

POSEY BED ACCESSORIES

Posey offers several products to be used with the Posey Bed, a patented bed enclosure system that is a less restrictive restraint alternative to physical restraints such as belts, vests, or jackets.

POSEY BED ACCESSORIES

CODE	DESCRIPTION	DIMENSIONS	UOM
8002	Head & Foot Board Padding	35" H x 17" L x 2" Thick	Each
8007	BCS Mattress Pad	82" L x 36" W x 1" Thick	Each
8018	Posey Attachment Straps	N/A	Each
8021	Posey Bed Filler Cushions	Fits 8060 & 8070	Each
8025	Posey Torso Cushion	Fits 8060 & 8070	Each
8071	Transport Cover For 8070	48" L x 36" W x 60" H	Each



8007 | BCS Mattress Pad



8021 | Posey Bed Filler Cushions



8021 | Posey Bed Filler Cushions



## ONLINE INFORMATION

Explore our broad range of product areas and learn more about how our solutions support caregivers and protect patients at **tidiproducs.com**.

---

Call **1.800.521.1314**, visit **tidiproducs.com**, or email us at **excellence@tidiproducs.com**.



TIDI Products, LLC • 570 Enterprise Drive • Neenah, WI 54956 USA  
Customer Service 800.521.1314 • International +1.920.751.4036 • Fax 800.837.7770  
[www.tidiproducs.com](http://www.tidiproducs.com)

Copyright (C)2019 TIDI Products, LLC. All Rights Reserved. Posey is a registered trademark property of TIDI Products, LLC.  
For U.S. and foreign patent information, see [tidiproducs.com/virtual-patent-marking](http://tidiproducs.com/virtual-patent-marking). Not all products are available in all markets.