

# Keeping Your Patients Safe and Comfortable



Posey has an array of innovative, easy-to-use solutions that keep your patients safe and comfortable

Whenever possible, non-restrictive interventions should be attempted before applying restraints, according to both federal and state regulations, along with The Joint Commission (formerly JCAHO) standards. Avoid line pulling, wound picking, or any other disruption to treatment by disorientated or discomforted patients through the use of effective alternatives, without the documentation required when using restraint devices.



## **Peek-A-Boo® Mitts**

- Prevents patients from grasping or dislodging tubes or other devices
- Enables full mobility
- Allows for circulation checks without mitt removal



## **SecureSleeves®**

- Prevents patients from grabbing and dislodging tubes
- Provides patient comfort with unique beaded technology
- Enables full hand mobility



## **SkinSleeves® Protector**

- Serves as a least-restrictive option
- Hides delicate sites
- Provides an extra layer of padding to protect fragile skin from tears and bruises
- Reduces skin tears by 53%\*

\*Bank, D, Nix, D. Preventing skin tears in a nursing and rehabilitation center: an interdisciplinary effort. *Ostomy Wound Management*. 2006; 52(9):38–46.

# KEEPING YOUR PATIENTS SAFE AND COMFORTABLE



## SkinSleeves® Protector

- Camouflage IV and wound sites to help prevent patient disruption of medical treatments
- Help protect thin, fragile skin tissue
- Provide warmth for patients' arms and legs
- Color-coded sizing bands to speed visual searches
- Can be machine washed

Item	Size	Color Band	Skin Color	Style
<b>6000S</b>	Small	Red	Light	Arm
<b>6000</b>	Medium	Green	Light	Arm
<b>6000L</b>	Large	Yellow	Light	Arm
<b>6000XL</b>	X-Large	Blue	Light	Arm
<b>6000B</b>	Medium	Green (Bulk)	Light	Arm
<b>6003S</b>	Small	Red	Dark	Arm
<b>6003</b>	Medium	Green	Dark	Arm
<b>6002</b>	Medium	Green (Padded)	Light	Arm
<b>6001</b>	Medium	Green	Light	Leg
<b>6005</b>	Medium	Green (Padded)	Light	Leg



## SecureSleeves®

- Apply easily using a hook-and-loop closure
- Help prevent sliding by pinning or clipping to the patient's sleeve
- Splint can be left in place during X-Ray because of its translucent fabric
- Facilitate comfort with tiny, flame retardant polystyrene bead filling
- Can be machine washed

ADULT SIZES			
Item	Description	Dimensions	Limb Circumference
<b>8168S</b>	Small	11"L x 14½"W (28 cm x 37 cm)	3"-14" (8-36 cm)
<b>8168M</b>	Medium	13"L x 18"W (33 cm x 46 cm)	6"-17" (15-43 cm)
<b>8168L</b>	Large	14"L x 25"W (36 cm x 64 cm)	8"-21" (20-53 cm)
<b>8169</b>	Leg	18"L x 27"W (46 cm x 69 cm)	14"- 27" (36-69 cm)

  

PEDIATRIC SIZES			
Item	Description	Dimensions	Limb Circumference
<b>8168I</b>	Infant	3"L x 4"W (8 cm x 10 cm)	1½"-4" (4-10 cm)
<b>8168IL</b>	Infant, Large	4½"L x 4"W (11 cm x 10 cm)	½"-4" (4-10 cm)
<b>8168P</b>	Pediatric	7"L x 7"W (18 cm x 18 cm)	3"-10" (8-25 cm)
<b>8168XS</b>	Pediatric, Large	10"L x 12"W (25 cm x 30 cm)	3"-12" (8-30 cm)



## Hand Mitts

- Fit on either hand via an ambidextrous design
- Give patients padding for extra comfort
- Secure in place with hook-and-loop straps or quick-release buckles
- Help facilitate skin checks and allow skin to breathe with available mesh backing material
- Can use optional connecting straps
- Can be machine washed

ADULT MITTS				
Item	Description	Style	Fill	Finger Separators
<b>2810</b>	Soft Hand Mitt	Breathable spun PET	Polystyrene bead	NA
<b>2811</b>	Peek-A-Boo Mitt	Double Padded	Fiber Fill	NA
<b>2813</b>	Secure Mitt	Mesh Back	Fiber Fill	5
<b>2814</b>	Double Security Mitt	Mesh Back	Fiber Fill	NA
<b>2816</b>	Finger Control Mitt	Mesh Back	Fiber Fill	5
<b>2819</b>	Double Security Mitt	Double Padded	Fiber Fill	NA
<b>2826</b>	Economy Mitt	Mesh Back	Foam	3

  

PEDIATRIC PEEK-A-BOO® MITTS	
Item	Description
<b>2911S</b>	Infant Mitts (0-1 year old)
<b>2911M</b>	Toddler Mitts (1-3 years old)
<b>2911L</b>	Youth Mitts (3-7 years old)