



# Rest Easy

Knowing That They're  
Resting Easy

THE POSEY SITTER ON CUE™ SYSTEM IS:

- › Intuitive
- › Reliable
- › Easy to use

*Take comfort knowing that the industry leader in falls  
management is looking after your patients.*





# Easy to Use Easy to Integrate

## MAKING FALLS MANAGEMENT EASIER ON CAREGIVERS

The Posey Sitter On Cue Alarm System is a restraint-free solution engineered to alert caregivers of a potential falls risk, giving them more time to react and intervene when every second counts.

The new features help you integrate the alarm system seamlessly into your current process. Training your team is easy with the intuitive interface, and the alarms are compatible with all existing Posey sensors, covering your needs from the ED to ICU.



## POSEY OFFERS A WIDE VARIETY OF COMPATIBLE SENSORS.



CHAIR SENSORS



BED SENSORS



ALARM BELTS

## SPECIALITY SENSORS



8317  
LOW AIR LOSS  
SENSORS



8316  
STRETCHER  
SENSORS



8332  
TOILET  
SENSORS

# Designed with Caregivers in Mind

## + STATUS LIGHT

Status is highly visible and can be viewed without disturbing the patient.

**GREEN** flash **GO**  
**RED** flash **NO**



## + AUDIO CUES

Gentle audio cues help guide a caregiver through setup.




## + HOLD BUTTON

Pressing hold disables the alarm for 30 seconds, allowing you to freely move your patient from their bed to a chair without setting off the alarm.



## + DUAL PORTS

Transition patients between two locations without having to reconfigure the alarm system each time.



"I like that it doesn't rely only on a visual but tells you that it's working."

"The sound cues letting you know a patient is being monitored is extremely helpful"

"Large light allows for better monitoring"

— ANONYMOUS NURSES\*

\*DATA ON FILE



# A Blueprint For Reducing Falls

PRACTICE	QUESTIONS TO CONSIDER
MOUNT ALARMS AT THE POINT OF CARE	<ul style="list-style-type: none"> <li>• Are your alarms there when you need them?</li> <li>• Do your alarms get lost in the shuffle?</li> </ul>
CONNECT ALARMS TO NURSE CALL	<ul style="list-style-type: none"> <li>• Would it be easier for your staff if alerts were integrated into the nurse call system?</li> <li>• Does your staff waste precious time trying to figure out which room an alarm is coming from?</li> </ul>
STOCK SENSOR PADS IN EVERY DEPARTMENT THAT HAS ALARMS	<ul style="list-style-type: none"> <li>• Have you needed a sensor pad but couldn't find one?</li> <li>• Is there a dedicated space on your floor that you know will have sensor pads stocked?</li> <li>• Would it be easier if a fresh sensor pad was stocked in every patient room?</li> </ul>
REVIEW FALLS PROTOCOL REGULARLY	<ul style="list-style-type: none"> <li>• Do you know when your falls protocol was last updated?</li> <li>• How often is a sensor pad intended to be used on a patient that is a fall risk?</li> <li>• Do you know how many sensors you used last year compared to how many patients were deemed a high fall risk?</li> </ul>

DESCRIPTION	TYPE	PART #
Sitter On Cue	Alarm	8645
Chair Sensor Pad, Single Patient Use, 30-day, 3-ft cord	Sensor Pad	8309
Chair Sensor Pad, Single Patient Use, 30-day, 8-ft cord	Sensor Pad	8309EL
Bed Sensor Pad, Single Patient Use, 30-day, 8-ft Cord	Sensor Pad	8283
Stretcher Sensor Pad, Disposable, 15-day, 3-ft Cord	Sensor Pad	8316
Low Pressure Mattress Sensor Pad, Single Patient Use, 30-day, 8-ft Cord	Sensor Pad	8317
Disposable self-releasing alarm belt, adjusts from 31" to 77" (79cm to 196cm)	Sensor Pad	8399
Disposable self-releasing alarm belt (Large), adjusts up to 89" (226cm)	Sensor Pad	8399L
Toilet Seat Sensor	Sensor Pad	8332
Bedside Commode Sensor	Sensor Pad	8333
Posey Nurse Call Cable, 12-ft	Accessory	8282
Posey Nurse Call Cable, 3-ft	Accessory	8282SL



TIDI Products, LLC • 570 Enterprise Drive • Neenah, WI 54956 USA  
 Customer Service 800.521.1314 • International +1.920.751.4036 • Fax 800.837.7770  
[www.tidiproductions.com](http://www.tidiproductions.com)

© TIDI Products, LLC. All rights reserved. Posey is a registered trademark of TIDI Products, LLC. On Cue is a trademark of TIDI Products, LLC.

MK024-19