



Nitrile

These hazmat gloves offer dual-layer protection against opioids and hazardous substances that makes blood and bodily fluids more visible

- Tested against both fentanyl and gastric acid (vomit) to simulate hazardous, “real world” overdose situations
- Feature a dual layer, dual color design (white exterior, black interior)
- Deliver advanced barrier protection (AQL 0.65)
- Triple washed, and have a true pH balanced interior coating for improved hand health



Industries

- Correctional Facilities
- Emergency Medical Services
- Healthcare
- Law Enforcement

Recommended For

- Administering drugs
- Inspection, selecting, checking parts
- Sample taking and processing
- Transferring liquids and solids
- Emergency services
- Overdose situations
- Standard patient examinations



TECHNICAL DATA SHEET

PRODUCT INFORMATION

	ApexPro™ LC100 AP12
Material	Nitrile
Color	White Exterior, Black Interior
Glove Design	Powder-Free, Chlorinated, Fully Textured
Cuff	Beaded
Regulatory/Standards Compliance	ASTM D6319, FDA21 CFR 177-2600, ISO 13485, ISO 9001, NFPA 1999:2018
Packaging	100 gloves per box; 10 boxes per case; 2XL and 3XL-90 gloves per box; 10 boxes per case
Storage	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Country of Origin	China
User Needs Segment	High Risk
Available sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10), XXL (10.5 - 11), XXXL (11.5 - 12)
Anti-static	Not Tested

PHYSICAL PROPERTIES

	Typical Values		Testing Method
Length (mm/inches)	300 / 11.8		ASTM D6319
Freedom from Holes (Inspection level I)	0.65 AQL		ISO 2859-1
Palm Thickness (mm/mils)	0.15 / 5.9		ASTM D6319
Finger Thickness (mm/mils)	0.21 / 8.3		ASTM D6319
	BEFORE AGING	AFTER AGING	
Ultimate Tensile Strength (MPa)	20-25	20-25	ASTM D412 & D573
Elongation at Break (%)	≥550	≥500	ASTM D412

ORDERING INFORMATION

Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)	XXL (10.5 - 11)	XXXL (11.5 - 12)
Product Code	AP12-1X	AP12-1	AP12-2	AP12-3	AP12-4	AP12-5	AP12-6

Performance Standards and Regulatory Compliance



For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV
T: +32 (0) 2 528 74 00
F: +32 (0) 2 528 74 01

Asia Pacific Region

Ansell Global Trading Center
(Malaysia) Sdn Bhd
T: +603 8310 6688
F: +603 8310 6699

North America Region

Ansell Healthcare Products LLC
T: +1 800 800 0444
F: +1 800 800 0445

Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.
T: +52 442 248 1544 / 248 3133

Australia

Ansell Limited
T: +61 1800 337 041
F: +61 1800 803 578

Russia

Ansell PYC
Ten. +7 495 258 13 16

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking © 2021 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.