



FrontLine™ 9 FL9

Nitrile Exam Gloves

NBR

Nitrile

An effective barrier against blood-borne pathogens

- Advanced non-latex nitrile polymer offers superior stretch and comfort
- Meets or exceeds the emergency medical glove requirements of NFPA 1999
- Fentanyl Permeation Testing: zero breakthrough with Fentanyl Citrate 100mcg/2ml for 240 minutes
- Non-stick exterior surface prevents tapes from sticking to the glove
- AQL of 0.65 exceeding industry standards
- Human skin pH 5.5 balanced coating for reduced risk of irritation and perspiration



Industries

- Emergency Medical Services
- Healthcare
- Law Enforcement

Recommended For

- Administering drugs
- Standard patient examinations
- Emergency services





FrontLine™ 9 FL9

Nitrile Exam Gloves

TECHNICAL DATA SHEET

PRODUCT INFORMATION

	FrontLine™ 9 FL9
Material	Nitrile
Color	Blue
Glove Design	Powder-Free, Chlorinated, Fully Textured
Cuff	Beaded
Manufacturing/QMS Audit Standards	ISO 13485, ISO 9001
Regulatory/Standards Compliance	ASTM D6319, FDA21 CFR 177-2600, ISO 13485, NFPA 1999:2013
Storage	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Country of Origin	China
User Needs Segment	Robust
Available sizes	S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10), XXL (10.5 - 11), XXXL (11.5 - 12)
Anti-static	Not Tested

PHYSICAL PROPERTIES

	Typical Values		Testing Method
Length (mm/inches)	240 / 9.5		ASTM D3578, ASTM D3577
Freedom from Holes (Inspection level I)	0.65 AQL		ASTM D5151
Palm Thickness (mm/mils)	0.14 / 5.5		ASTM D3578, ASTM D3577
Finger Thickness (mm/mils)	0.16 / 6.3		ASTM D3578, ASTM D3577
	BEFORE AGING	AFTER AGING	
Ultimate Tensile Strength (MPa)	20	18	ASTM D412 & D573
Elongation at Break (%)	550	500	ASTM D412 & D573

ORDERING INFORMATION

Size	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)	XXL (10.5 - 11)	XXXL (11.5 - 12)
Product Code	FL9-1	FL9-2	FL9-3	FL9-4	FL9-5	FL9-6

Performance Standards and Regulatory Compliance



ISO 13485

For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV
T: +32 (0) 2 528 74 00
F: +32 (0) 2 528 74 01

Asia Pacific Region

Ansell Global Trading Center
(Malaysia) Sdn Bhd
T: +603 8310 6688
F: +603 8310 6699

North America Region

Ansell Healthcare Products LLC
T: +1 800 800 0444
F: +1 800 800 0445

Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.
T: +52 442 248 1544 / 248 3133

Australia

Ansell Limited
T: +61 1800 337 041
F: +61 1800 803 578

Russia

Ansell PYC
Ten. +7 495 258 13 16

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking © 2021 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

