

GAMMEX® PI Radiation Attenuation

M E D I C A L S P E C I A L T Y G L O V E S

Polyisoprene, personal protective shield protects hands from scattered radiation exposure during fluoroscopic procedures¹

KEY FEATURES AND BENEFITS:

- This polyisoprene, radiation attenuation personal protective shield reduces skin-dose exposure by up to 53.4% (beam energy level 60 kVp) versus double gloving with conventional surgical gloves that attenuate under 1% of scattered radiation²
- Eliminates the risk of developing Type I latex allergy and sensitization
- This glove's thickness is optimized for the ideal balance of protection and comfort
- Black in color, this glove features a micro-textured finish for good grip control


IMPORTANT: Classified as a personal protective shield and are not a medical device. The gloves must be used in combination with medical gloves if they come in contact with blood or other bodily fluids.

1. CAUTION: Radiation attenuation gloves are designed to protect hands from scattered radiation exposure originating from the X-ray beam during fluoroscopic procedures. They are not intended for use in direct or primary X-ray beams.

2. Wagner LK, Mulhern OR. Radiation-attenuating surgical gloves: effects of scatter and secondary electron production. Radiology. 1996;200(1):45-48.

Recommended For

- Protecting hands from image guided interventional procedures (fluoroscopy)
- Protection from Type I latex allergy in HCW's or patients



GAMMEX® Innovation that Protects

Driven by innovation, GAMMEX® surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.



GAMMEX® PI Radiation Attenuation

MEDICAL SPECIALTY GLOVES

PRODUCT INFORMATION

Material	Polyisoprene with Tungsten
Powder Content	Powder-Free
Color	Black
Thickness	Standard
Cuff Style	Straight
External Glove Surface	Textured Fingers
Internal Glove Surface	Chlorinated
Allergy Prevention	Latex (Type I)
Grip Level	Moderate
Tested For Use With Chemotherapy Drugs	No, not tested in accordance with ASTM D6978
Freedom from Holes (Inspection level I)	1.5 AQL
Sterile	Yes
Product Standards	ASTM F2547, EN 61331-1, EN ISO 374-5:2016, EN 388:2016, EN 420:2003 + A1:2009, EN 421:2010
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
Packaging	5 pairs per box
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

ORDERING INFORMATION

Size	6.5	7	7.5	8	8.5	9
Product Code	20873565	20873570	20873575	20873580	20873585	20873590

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 844 494 7854
Email: insidesalesus@ansell.com
Sandel US Tel: 866 764 3327
Sandel email: sandel-cs@ansell.com

Central & South America
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil
Tel: +55-11-3356-3100
Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

India
Tel: +91 22 38124500
Email: info.india@ansell.com
Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com
China
Tel: +86(21) 3827 5005
Email: infochina@ansell.com

ansell.com

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending:
www.ansell.com/patentmarking.
© 2021 Ansell Limited. All Rights Reserved.