



TECHNICAL DATA SHEET

PRODUCT NAME SaniTex Plus™ Germicidal Disposable Wipe

PRODUCT OVERVIEW An EPA registered, low-level disinfectant wipe ideal for use on a wide range of hard, non-porous surfaces and equipment.

COMPOSITION

- n-Alkyl (68% C₁₂, 32% C₁₄) dimethyl ethylbenzyl ammonium chlorides, 0.125%
- n-Alkyl (60% C₁₄, 30% C₁₆, 5% C₁₂, 5% C₁₈) dimethyl benzyl ammonium chlorides, 0.125%
- Isopropanol, 14.85%

DETAILS

- A broad spectrum, low-level disinfectant (bactericidal and virucidal) effective against 16 microorganisms including bloodborne pathogens, viruses and bacteria.
- Total disinfection time is 3 minutes.
- Ready-to-use, pre-moistened wipes saturated with the precise amount of germicidal solution to effectively clean and disinfect, which saves time.
- Ideal for use in dental and alternate care settings on a wide range of hard, non-porous surfaces and equipment.
- The unique deep-well lid provides plenty of room to store the next wipe and a secure seal to prevent moisture loss.
- White wipe with mild alcohol scent.

PACKAGING Product is available in the following configurations:

Product Code	Description
SANLC	SaniTex Plus™ Germicidal Disposable Wipe (Large) 160 wipes/canister (6 in x 6.75 in / 15.2 cm x 17.1 cm)
SANXLC2	SaniTex Plus™ Germicidal Disposable Wipe (Extra Large) 65 wipes/canister (7.5 in x 15 in / 19.0 cm x 38.1 cm)

STORAGE REQUIREMENTS Do not store near heat or open flame. When not in use, keep center cap of lid closed to prevent moisture loss.

SHELF LIFE 2 years from date of manufacture

REGULATORY INFORMATION

US EPA Registered Pesticide
 EPA Reg. No.: 9480-6-64285
 EPA Est. No.: 9480-NY-1
 Medical Device Classification: not a medical device
 FDA Product code: n/a
 FDA Regulation: n/a
 Medical Device Listing Number: n/a