

2019 Dental Catalog



Dedicated to delivering innovative infection prevention products & services for patients, caregivers, and other healthcare providers, which improve outcomes and help save lives.

Table of Contents

| | |
|-------------------------------|----|
| Chemistries | 4 |
| Dental Water Management | 12 |
| Equipment | 18 |
| PPE | 26 |
| Preventives | 40 |
| Single-Use Disposables | 48 |
| Sterilization | 74 |
| Index | 86 |

Crosstex offers a broad range of powerful solutions for cleaning, surface disinfection and instrument reprocessing.

Proper reprocessing of reusable devices is critical to ensuring patient safety and avoiding cross-contamination. Cleaning solutions, including high-, intermediate- and low-level instrument and surface disinfectants as well as enzymatic detergents for cleaning sensitive instruments and evacuation lines are among the many products available. Many of these solutions are non-corrosive and pH-neutral, for maximum compatibility with instruments, equipment and surfaces.



Surface Disinfection Microorganism Chart



SaniTex Plus™ 2
Cleaner,
Sanitizer &
Disinfectant

SaniTex Plus™
Cleaner,
Disinfectant &
Deodorizer

**Sani-Cloth®
AF3**
Germicidal
Disposable
Wipe

**Super
Sani-Cloth®**
Germicidal
Disposable
Wipe

**Sani-Cloth®
Bleach**
Germicidal
Disposable
Wipe

SaniTex Plus™
Germicidal
Disposable
Wipe

**Sani-Cloth®
Plus**
Germicidal
Disposable
Cloth

**Sani-Cloth®
HB**
Germicidal
Disposable
Wipe

| Efficacy level | Intermediate | Intermediate | Intermediate | Intermediate | Intermediate | Low | Low | Low |
|---|--------------|--------------|--------------|--------------|--------------|-----------|-----------|------------|
| Total disinfection contact time | 2 minutes | 6 minutes | 3 minutes | 2 minutes | 4 minutes | 3 minutes | 3 minutes | 10 minutes |
| <i>Acinetobacter baumannii</i> | ● | ● | ● | ● | ● | | | |
| <i>Adenovirus, Type 5</i> | | | ● | ● | | | | |
| <i>Aspergillus niger</i> | | ● | | | | | | ● |
| Avian Influenza A Virus | ● | ● | ● | | ● | | | |
| <i>Campylobacter jejuni</i> | | | ● | ● | ● | ● | ● | ● |
| <i>Candida albicans</i> | | ● | ● | ● | ● | | | ● |
| <i>Clostridium difficile</i> | | | | | ● | | | |
| Cytomegalovirus (CMV) | | | ● | | ● | | | |
| <i>Enterobacter aerogenes</i> | ● | | ● | | | | | ● |
| <i>Enterococcus faecalis</i> , Vancomycin resistant (VRE) | ● | ● | ● | ● | ● | ● | ● | ● |
| <i>Escherichia coli</i> | | | ● | ● | ● | ● | ● | ● |
| <i>Escherichia coli</i> O157:H7 | | ● | ● | ● | ● | ● | ● | ● |
| <i>Escherichia coli</i> , Drug Resistant | ● | | ● | ● | ● | | | ● |
| Hepatitis B Virus (HBV) | ● | ● | ● | ● | ● | ● | ● | ● |
| Hepatitis C Virus (HCV) | ● | ● | ● | ● | ● | ● | ● | |
| Herpes Simplex Virus, Type 2 | ● | ● | ● | ● | ● | ● | ● | ● |
| Human Coronavirus | | | ● | ● | ● | | | |
| Human Immunodeficiency Virus (HIV-1) | ● | ● | ● | ● | ● | ● | ● | ● |
| Influenza A Virus | ● | ● | ● | ● | ● | ● | ● | ● |
| Influenza B Virus | | | ● | | ● | | | ● |
| <i>Klebsiella pneumoniae</i> | | | ● | ● | ● | | | ● |
| <i>Klebsiella pneumoniae</i> , Drug Resistant | ● | | ● | ● | ● | | | ● |
| <i>Listeria monocytogenes</i> | ● | | ● | | ● | | | ● |
| <i>Mycobacterium bovis</i> , BCG (TB) | ● | ● | ● | ● | ● | | | |
| Norovirus | | | | | ● | | | |
| <i>Pseudomonas aeruginosa</i> | ● | ● | ● | ● | ● | ● | ● | ● |
| <i>Pseudomonas aeruginosa</i> , Drug Resistant | | | | | | | | ● |
| Respiratory Syncytial Virus (RSV) | | | ● | ● | ● | ● | ● | |
| Rhinovirus | ● | ● | | ● | ● | | | |
| Rotavirus (Strain WA) | ● | ● | ● | ● | ● | | | |
| <i>Salmonella enterica</i> | ● | ● | ● | ● | ● | ● | ● | ● |
| <i>Serratia marcescens</i> | ● | | ● | | ● | | | ● |
| <i>Staphylococcus aureus</i> | ● | ● | ● | ● | ● | ● | ● | ● |
| <i>Staphylococcus aureus</i> , Methicillin Resistant (MRSA) | ● | ● | ● | ● | ● | ● | ● | ● |
| <i>Staphylococcus aureus</i> , Vancomycin Resistant (VRSA) | | | ● | | ● | | | |
| <i>Streptococcus pneumoniae</i> , Penicillin resistant | | | ● | | ● | | | |
| <i>Streptococcus pyogenes</i> | | | ● | | ● | | | ● |
| <i>Trichophyton interdigitale</i> | ● | ● | | | ● | | | ● |
| Vaccinia Virus | | | | ● | | | | ● |

● = Effective against

For a comprehensive list of all microorganisms each product is effective against, please refer to the technical documentation.



Surface Disinfectants

● SaniTex Plus™ 2 Cleaner, Sanitizer & Disinfectant

- **EPA registered** – Intermediate-level disinfectant (tuberculocidal, bactericidal, virucidal and fungicidal).
- **Low alcohol/quaternary ammonium formulation** - Ready-to-use solution provides exceptional cleaning and disinfection on a wide array of hard, non-porous surfaces and equipment in dental settings.
- **Broad spectrum with rapid action** – effective against 22 microorganisms in 2 minutes including TB, bloodborne pathogens, bacteria, viruses and fungi.
- **3-in-1 product** – Clean, sanitize, and disinfect with one product.

EPA REG NO.: 70144-1-64285



SaniTex Plus™ 2 Cleaner, Sanitizer & Disinfectant

| Ref. # | Description | Quantity |
|--------|--------------|-------------------|
| JSSDP2 | Spray Bottle | 24 fl oz (710 mL) |
| JSSGP2 | Bottle | 1 gal (3.8 L) |

● SaniTex Plus™ Cleaner, Disinfectant & Deodorizer

- **EPA registered** – Intermediate-level disinfectant (tuberculocidal, bactericidal, virucidal and fungicidal).
- **Low alcohol/quaternary ammonium formulation** – Ready-to-use solution provides excellent cleaning and disinfecting capabilities on a broad range of hard, non-porous surfaces and equipment in dental settings.
- **Broad spectrum** – Effective against 21 microorganisms in 6 minutes including TB, bloodborne pathogens, bacteria, viruses and fungi.

EPA REG NO.: 1130-15-64285



SaniTex Plus™ Cleaner, Disinfectant & Deodorizer

| Ref. # | Description | Quantity |
|--------|--------------|----------------|
| JSSDP | Spray Bottle | 1 qt (0.95 L) |
| JSSGP | Bottle | 1 gal (3.78 L) |

DISINFECTANT LEVEL KEY

- High-level
- Intermediate-level
- Low-level

● Sani-Cloth® AF3 Germicidal Disposable Wipe

- **EPA registered** – Intermediate-level disinfectant (tuberculocidal, bactericidal and virucidal).
- **Unique, alcohol-free/quaternary ammonium formulation** – designed for disinfecting a broad range of surfaces and equipment in dental settings.
- **Broad spectrum** – Effective against 45 microorganisms in 3 minutes, including TB, bloodborne pathogens, bacteria and viruses.
- **Reliable** – Ready-to-use, thick, strong pre-moistened wipe has high fluid capacity and is pliable to get into crevices and corners of equipment.



EPA REG NO.: 9480-9

● Super Sani-Cloth® Germicidal Disposable Wipe

- **EPA registered** – Intermediate-level disinfectant (tuberculocidal, bactericidal and virucidal).
- **Fast-acting alcohol/quaternary ammonium formulation** – Ideal for use in fast-paced environments requiring short contact time and broad coverage of microorganisms, allowing for quick room/operator turnover on a broad range of surfaces and equipment in dental settings.
- **Broad spectrum** – Effective against 30 microorganisms in 2 minutes including TB, bloodborne pathogens, bacteria and viruses.
- **Reliable** – Ready-to-use, thick, strong pre-moistened wipe has high fluid capacity and is pliable to get into crevices and corners of equipment.



EPA REG NO.: 9480-4

● Sani-Cloth® Bleach Germicidal Disposable Wipe

- **EPA registered** – Intermediate-level disinfectant (sporicidal, tuberculocidal, bactericidal, virucidal and fungicidal).
- **High risk prevention** – Stabilized dilution of sodium hypochlorite (bleach) ideal for disinfecting high risk areas endemic with *Clostridium difficile* (*C. diff*) spores and Norovirus.
- **Broad spectrum** – Effective against 50 microorganisms in 4 minutes including *C. diff* spores, TB, bloodborne pathogens, bacteria, viruses and fungi.
- **Reliable** – Ready-to-use, thick, strong pre-moistened wipe has high fluid capacity and is pliable to get into crevices and corners of equipment.



EPA REG NO.: 9480-8

All wipes not made with natural rubber latex.

Sani-Cloth® AF3 Germicidal Disposable Wipe

| Ref. # | Description | Quantity |
|-------------------------------|--|--------------|
| SANAF3 (P13872) | Large Wipe 6 in x 6.75 in (15.2 cm x 17.1 cm) | 160/Canister |
| SANAF3L23 (P63884) | Extra Large Wipe 7.5 in x 15 in (19.1 cm x 38.1 cm) | 65/Canister |
| SANAF3I (H59200) | Large Individual Packet 5 in x 8 in (12.7 cm x 20.3 cm) | 50/Box |
| SANAF3L23I (U27500) | Extra-Large Individual Packet 11.5 in x 11.75 in (29.2 cm x 29.8 cm) | 50/Box |

Super Sani-Cloth® Germicidal Disposable Wipe

| Ref. # | Description | Quantity |
|----------------------------|--|--------------|
| SANLS (Q55172) | Large Wipe 6 in x 6.75 in (15.2 cm x 17.1 cm) | 160/Canister |
| SANLSL2 (Q86984) | Extra Large Wipe 7.5 in x 15 in (19.1 cm x 38.1 cm) | 65/Canister |
| SANLSI (H04082) | Large Individual Packet 5 in x 8 in (12.7 cm x 20.3 cm) | 50/Box |
| SANLSLI (U87295) | Extra-Large Individual Packet 11.5 in x 11.75 in (29.2 cm x 29.8 cm) | 50/Box |

Sani-Cloth® Bleach Germicidal Disposable Wipe

| Ref. # | Description | Quantity |
|------------------------------|--|-------------|
| SANLB2 (P54072) | Large Wipe 6 in x 10.5 in (15.2 cm x 26.6 cm) | 75/Canister |
| SANLB2L2 (P25784) | Extra Large Wipe 7.5 in x 15 in (19.1 cm x 38.1 cm) | 65/Canister |
| SANLB (H58195) | Large Individual Packet 5 in x 7 in (12.7 cm x 17.7 cm) | 40/Box |
| SANLB2L2I (U26595) | Extra-Large Individual Packet 11.5 in x 11.75 in (29.2 cm x 29.8 cm) | 40/Box |



Surface Disinfectants

● SaniTex Plus™ Germicidal Disposable Wipe

- **EPA registered** – Low-level disinfectant (bactericidal and virucidal).
- **Low alcohol/quaternary ammonium formulation** – Ready-to-use, pre-moistened wipes saturated with the precise amount of germicidal solution to effectively clean and disinfect.
- **Broad spectrum** – Effective against 16 microorganisms in 3 minutes including bloodborne pathogens, bacteria and viruses.
- **Multiple uses** – Ideal for use in dental settings on a broad range of hard, non-porous surfaces and equipment.



EPA REG NO.: 9480-6-64285

SaniTex Plus™ Germicidal Disposable Wipe

| Ref. # | Description | Quantity |
|----------------|---|--------------|
| SANLC | Large Wipe 6 in x 6.75 in (15.2 cm x 17.1 cm) | 160/Canister |
| SANXLC2 | Extra Large Wipe 7.5 in x 15 in (19.1 cm x 38.1 cm) | 65/Canister |

● Sani-Cloth® Plus Germicidal Disposable Cloth

- **EPA registered** – Low-level disinfectant (bactericidal and virucidal).
- **Low alcohol/quaternary ammonium formulation** – Ready-to-use, pre-moistened wipes ideal for disinfecting a broad range of surfaces and equipment in dental settings.
- **Broad spectrum** – Effective against 16 microorganisms in 3 minutes including bloodborne pathogens, bacteria and viruses.



EPA REG NO.: 9480-6

Sani-Cloth® Plus Germicidal Disposable Cloth

| Ref. # | Description | Quantity |
|------------------------|---|--------------|
| SANL (Q89072) | Large Wipe 6 in x 6.75 in (15.2 cm x 17.1 cm) | 160/Canister |
| SANXL2 (Q85084) | Extra Large wipe 7.5 in x 15 in (19.1 cm x 38.1 cm) | 65/Canister |

● Sani-Cloth® HB Germicidal Disposable Wipe

- **EPA registered** – Low-level disinfectant (bactericidal, virucidal and fungicidal).
- **Alcohol-free/quaternary ammonium formulation** – Ready-to-use, pre-moistened wipes ideal for disinfecting sensitive hard, non-porous surfaces and equipment in dental settings.
- **Broad spectrum** – Effective against 100+ microorganisms in 10 minutes including bloodborne pathogens, bacteria, viruses and fungi.



EPA REG NO.: 61178-4-9480

Sani-Cloth® HB Germicidal Disposable Wipe

| Ref. # | Description | Quantity |
|--------------------------|---|--------------|
| SANHBL (Q08472) | Large Wipe 6 in x 6.75 in (15.2 cm x 17.1 cm) | 160/Canister |
| SANHBXL2 (Q85484) | Extra Large Wipe 7.5 in x 15 in (19.1 cm x 38.1 cm) | 65/Canister |

DISINFECTANT LEVEL KEY

- High-level ● Intermediate-level ● Low-level

High-Level Disinfectants/Sterilants & Accessories

● Rapicide™ OPA/28 High-Level Disinfectant

- **Innovative/safe formulation** – Ready-to-use, high-level disinfectant used to reprocess heat-sensitive semi-critical instruments with an active ingredient of 0.575% ortho-phthalaldehyde (OPA).
- **Fast & economical** – High-level disinfection in just 10 minutes at room temperature (68 °F/20 °C) with a reuse period of up to 28 days.[†]
- **Excellent materials compatibility**, non-corrosive to metal.

[†]Solution needs to be tested with Rapicide™ OPA/28 Test Strips prior to each use to verify MRC (minimum recommended concentration).



Rapicide™ OPA/28 High-Level Disinfectant

| Ref. # | Description | Quantity |
|-----------|-------------|---------------|
| ML02-0127 | Bottle | 1 gal (3.8 L) |

Rapicide™ OPA/28 Test Strips

- **Be compliant** – Test strips are designed specifically for use with Rapicide OPA/28 High-Level Disinfectant and provide MRC (minimum recommended concentration) verification.
- **User friendly** – Large, easy-to-read indicator pad with an intuitive pass/fail reference color chart on the bottle.
- **Fast** – 3-second dip time with a 90-second read time.



Rapicide™ OPA/28 Test Strips

| Ref. # | Description | Quantity |
|-----------|-------------|-----------|
| ML02-0137 | Test Strips | 50/Bottle |

● Sani Glut™ Glutaraldehyde 3%

- **Multi-purpose formula** – Biodegradable and formaldehyde-free reusable manual liquid chemical sterilant and high-level disinfectant, for use with immersible and heat sensitive instruments and/or devices.
 - 10 hours for sterilization (at 25 °C)
 - 25 minutes for high-level disinfection (at 25 °C)
- **Activation** – Requires addition of activator (included with bottle).
- **Economical** – Provides up to 28 days of reuse.

Use appropriate test strip prior to each use to verify MRC: Micro-Scientific Micro-MEC™ 1.8% Glutaraldehyde Monitor Test Strips (M60054)



Sani Glut™ Glutaraldehyde 3%

| Ref. # | Description | Quantity |
|--------|-------------|---------------|
| JGLUT | Bottle | 1 gal (3.8 L) |

Sani HLD Holder

- **Convenient** – Use with cleaners, disinfectants or cold sterilants
- **Easy-to-use** – Inner basket is easily removed and allows for convenient draining, while hinged lid contains calendar to notate when solution needs to be changed
- **High-capacity** – Holds up to 1 gal (3.8 L) of solution



Sani HLD Holder

| Ref. # | Description | Quantity |
|--------|--|----------|
| JHOLD | 13 ⁷ / ₈ in x 9 ⁵ / ₈ in x 5 ¹ / ₈ in (35.2 cm x 24.4 cm x 13.0 cm) | Each |



Cleaning Solutions & Detergents

●●● Sani ProZyme™ Enzymatic Detergent

- **Reliable, powerful formulation** – Protease enzymes break down organic soils such as blood, tissue and protein rich fluids.
- **Cost effective and safe** – Highly concentrated, low-foaming, biodegradable, phosphate-free, rust-inhibiting properties with a neutral pH; perfect for manual or automated washing.
- **Usage** – Dilute 0.5 oz to 2 oz per gallon of water; each bottle makes up to 256 gallons of solution (depending upon debris level).
- **Manual soak time:** 2 to 3 min (15 to 20 min for heavy debris).



Sani ProZyme™ Enzymatic Detergent

| Ref. # | Description | Quantity |
|--------|-------------|---------------|
| JED | Bottle | 1 gal (3.8 L) |

●●● CrossZyme™ Concentrated Ultrasonic Cleaner

- **Low-foaming, enzymatic formulation** – Provides effective cleaning to remove blood and organic debris immediately from instruments and easy rinsing over a wide range of concentrations and temperatures; powerfully chelated to improve performance in hard water conditions.
- **Biodegradable** - free from phosphate and chloride components to enhance the removal of blood and proteinaceous fluids.
- **Neutral pH** – safe for use on instruments and endoscopes, plastic, stainless steel, glass and aluminum; contains rust inhibitors and is non-corrosive.
- **Usage** – 1:20 dilution (general cleaning); 1:10 dilution (heavy cleaning).



CrossZyme™ Concentrated Ultrasonic Cleaner

| Ref. # | Description | Quantity |
|--------|-------------|----------------|
| JEZ | Bottle | 1 gal (3.78 L) |

●●●● EfferZyme™ Enzymatic Cleaning Tablets

- **Enzymatic formula** – Biodegradable, non-chlorine formulation contains fast-acting protease enzymes that enhance the removal of blood and other soils that may be resistant to other chemical detergents.
- **Neutral pH formula** – Contains no harsh alkaline components and is non-corrosive.
- **Cost effective** – Tablet formulation reduces shipping, handling and storage costs.
- **Usage** – 2 tablets per 1 gallon of water; each box makes 26 gallons of solution.



EfferZyme™ Enzymatic Cleaning Tablets

| Ref. # | Description | Quantity |
|--------|----------------------|----------|
| JET | Effervescent Tablets | 52/Box |

USAGE KEY

● Ultrasonic Cleaner

● Instrument Cleaning

● Instrument Presoak

● Evacuation System Cleaning

●● General Purpose Ultrasonic Cleaner, Non-ammoniated Concentrate

- **Multi-purpose formula** – Used for the removal of buffing compounds, pumice, rouge, tripoli, food particles, foreign matter and soils of all types.
- **Safe** – Use on dentures, burs, bridges, crowns, orthodontic appliances, instruments, diamonds and glass slides.
- **Usage** – Dilute 1 part concentrate with 10 parts water; use as an ultrasonic carrier bath when an auxiliary beaker is used for other cleaning tasks.
- **Cleaning time** – 2 to 10 minutes.



General Purpose Ultrasonic Cleaner, Non-ammoniated Concentrate

| Ref. # | Description | Quantity |
|--------|-------------|-----------------|
| JEZNA | Bottle | 1 gal (3.785 L) |

●● Tartar & Stain Remover

- **Powerful** – Proven biodegradable and phosphate-free ultrasonic acidic cleaner* for the removal of tartar, tobacco and food stains, permanent cement, certain temporary cements, solder, rust and hard water scale from bridges, dentures and orthodontic appliances.
- **Ready-to-use** – No dilution required.
- **Cleaning time** – 5 to 10 minutes.

*Not to be used directly in an ultrasonic bath, rather place solution in a beaker first before adding to ultrasonic bath.



Tartar & Stain Remover

| Ref. # | Description | Quantity |
|--------|-------------|-----------------|
| JEZTS | Bottle | 1 gal (3.785 L) |

● Sani Vac™ Evacuation System Cleaner

- **Superior cleaning** - Unique, low-foaming, multi-proteolytic enzymatic formula efficiently removes dental debris by breaking down saliva, mucous and blood as well as neutralizes odor in evacuation lines.
- **Compatible** - Biodegradable, pH-neutral and non-oxidizing solution meets the EPA's BMP's for Dental Amalgam Waste[^] to promote continuous, high-performance suction in wet/dry vacuum systems and amalgam separators.
- **Reduced risk** – Routine cleaning will help prevent downtime and costly line repairs by ensuring evacuation lines remain unclogged for optimal suction and increased patient safety.[†]
- **Economical** - The highly concentrated formula offers a low cost-per-use when diluted and treats up to 3 lines per ounce of concentrate and 384 lines per gallon.



Sani Vac™ Evacuation System Cleaner

| Ref. # | Description | Quantity |
|--------|-------------|----------------|
| JVAC | Bottle | 1 gal (3.78 L) |

[^] EPA Best Management Practices for Dental Amalgam Waste (www.epa.gov/eg/dental-effluent-guidelines)

[†] "Possibility of Cross-Contamination Between Dental Patients by Means of the Saliva Ejector", The Journal of the American Dental Association, April 1993, pp 77-80.



Dental Water Management

Water management and compliance in the dental office

Mechanical Room




Syclone™ Amalgam Separator

Rapid-velocity separation combined with high capacity canister to deliver better operational efficiency



VistaClear™ Centralized Water Filtration System

Filters city water for direct operator connections and provides a way to introduce periodic antimicrobial treatments.



Stericenter



VistaPure™ Water Purification System

Produces high-purity process water for filling sterilizers, dental bottles and more—no more buying bottled water!



VistaCool™ Direct-to-Drain Cooling System for Autoclave Wastewater

Cools and drains autoclave wastewater automatically—saving staff time while protecting cabinets and plumbing.

Operator



DentaPure™ Municipal Cartridge for units plumbed to city water

DentaPure™ Independent Water Bottle Cartridge

Provides 365 days (or 240L of water) of safe, compliant dental unit water



Liquid Ultra™ Solution

Multiple DUWL treatment options for bottle fed systems

Your complete dental water partner for safety, compliance and efficiency.

From waterline treatment to wastewater management, water quality affects nearly every function of modern dental offices. Poor water quality can harm patients, staff and equipment—not to mention the reputation of your practice.

Our comprehensive line of dental water products can help to improve patient care, staff productivity and practice profitability, while helping to maintain compliance with regulatory guidelines.





Municipally Plumbed Systems

DentaPure™ 365-Day Municipal Cartridge*

Easily adapted for installation onto incoming waterlines plumbed directly to city water via included poly-flo fittings.

- EPA Registered / FDA 510K Cleared.
- Cartridge reduces bacterial exposure by eluting non-allergenic, elemental iodine into dental water.
- Change every 365 days, or 240L of water.
- Crosstex reminds you when it is time to replace your cartridge.
- Not restricted by the EPA rule BMP for dental amalgam waste.

**One installation provides 365 days (or 240L of water if usage records are kept) of safe, compliant dental unit water*



For J-box installation



DentaPure™ 365-Day Municipal Cartridge

| Ref. # | Description | Quantity |
|--------|---------------------|----------|
| DP365M | Municipal Cartridge | 1 Each |

Bottle Systems

DentaPure™ 365-Day Independent Water Bottle Cartridge*

Installs onto your unit's existing water bottle pickup tube via included luer fitting.

- EPA Registered / FDA 510K Cleared.
- Can use either tap or distilled water.
- Cartridge reduces bacterial exposure by eluting non-allergenic, elemental iodine into dental water.
- Change every 365 days, or 240L of water.
- Crosstex reminds you when it is time to replace your cartridge.
- Not restricted by the EPA rule BMP for dental amalgam waste.

**One installation provides 365 days (or 240L of water if usage records are kept) of safe, compliant dental unit water*



For bottle installation



DentaPure™ 365-Day Independent Water Bottle Cartridge

| Ref. # | Description | Quantity |
|--------|------------------------------------|----------|
| DP365B | Independent Water Bottle Cartridge | 1 Each |

DentaPure™ 60-Day Independent Water Bottle Cartridge*

Installs onto your unit's existing water bottle pickup tube via included luer fitting.

- EPA Registered / FDA 510K Cleared.
- Perfect for short-term use.
- Cartridge reduces bacterial exposure by eluting non-allergenic, elemental iodine into dental water.
- Change every 60 days, or 40L of water.
- Crosstex reminds you when it is time to replace your cartridge.
- Not restricted by the EPA rule BMP for dental amalgam waste.

**One installation provides 60 days (or 40L of water if usage records are kept) of safe, compliant dental unit water*



For bottle installation



DentaPure™ 60-Day Independent Water Bottle Cartridge

| Ref. # | Description | Quantity |
|--------|------------------------------------|----------|
| DP40B | Independent Water Bottle Cartridge | 1 Each |

Testing

Iodine Test Strips

- Test strips may be used to verify adequate levels (higher than 0.5 ppm) of iodine release.
- These tests strips are only to be used to determine the level of IODINE present in dental unit waterlines (DUWLs), and cannot be used to measure the presence of Colony Forming Units (CFUs).



| Iodine Test Strips | | |
|--------------------|-----------------------|----------|
| Ref. # | Description | Quantity |
| DPTEST | 50 Iodine test strips | 1 Bottle |

Shock/In-Line Treatment

Liquid Ultra™ Solution

The powerful and versatile dental unit waterline shock treatment for external bottle systems.

- EPA Registered
- **Powerful** – Kills biofilm bacteria; removes existing biofilm and prevents/suppresses formation of biofilm in dental unit waterlines.
- **Versatile** – May be used as a weekly treatment; in conjunction with a tablet protocol and to rapidly treat waterlines prior to initiating a new treatment protocol, or to quickly treat DUWLs with extremely high CFU conditions.



| Liquid Ultra™ Solution | | |
|------------------------|---|---------------------|
| Ref. # | Description | Quantity |
| DS503L | Liquid Ultra™ Solution Each set contains: 1 bottle of solutions 1 & 2 (3 fl oz/bottle) | 1 box of 10 sets |



Dental Water Systems

VistaPure™ Water Purification System

- Provides two grades of high-purity water for steam sterilizers, dental bottles, ultrasonic cleaners and the rinse cycle of instrument washers.
 - Steam sterilizer-ready, distilled-quality (deionized) water eliminates the need to buy distilled water or own a distiller.
 - Bottle-ready, hyperfiltered water eliminates the need to buy bottled water and ensures high-quality water for mixing with bottle and residual cleaners, and for the rinse cycle of instrument washers.
- System connects directly to municipal water supply and sends purified water anywhere in the steri-center for easy, on-demand use.
- Included with the VistaPure™ System:
 - Handheld, point-and-shoot wand makes it easy for anyone to fill steam sterilizers, no matter where they're located.
 - Mounted faucet makes it easy to fill bottles over a sink.
 - In-line and handheld TDS meters monitor water quality — and make it easy to know when it's time to replace your filters.
 - 4.25 gal water storage tank (only one storage tank required). A larger tank is available if needed.
- Low ongoing cost of system ownership.
- No electricity needed.



VistaPure™ Water Purification System

| Ref. # | Description | Quantity |
|---------|--|----------|
| CVV3000 | VistaPure™ Water Purification System | 1 Each |
| CVS3010 | VistaPure™ Instrument Washer Adapter Kit | 1 Each |
| CVS3020 | VistaPure™ AutoFill Kit for SciCan STATIM® Autoclave | 1 Each |
| CVS3030 | VistaPure™ AutoFill Kit for SciCan BRAVO™ Autoclave | 1 Each |
| CVR3251 | VistaPure™ Replacement Filter – Prefilter (#1) | 1 Each |
| CVR3252 | VistaPure™ Replacement Filter – Deionization Filter (#2) | 1 Each |
| CVR3253 | VistaPure™ Replacement Filter – Deionization Filter (#3) | 1 Each |
| CVR3254 | VistaPure™ Replacement Filter – Hyperfilter (#4) | 1 Each |
| CVR3255 | VistaPure™ Replacement Filter – Polishing Filter (#5) | 1 Each |
| CVR2551 | VistaPure™ Legacy Filter – Polishing Filter (#5) | 1 Each |
| CVS9111 | VistaPure™ Water Storage Tank (4.25 gal) | 1 Each |
| CVS9113 | VistaPure™ Water Storage Tank (9.50 gal) | 1 Each |

VistaClear™ Centralized Water Filtration System

- Provides dental practices with centralized submicron water filtration, backflow prevention and a convenient method for delivering periodic antimicrobial treatments.
- Eliminates the need for dental bottles.
- Eliminates wait-time and the need for storage tanks and bottle systems.
- Centralized, patented cleaning method makes regular waterline maintenance easier and more efficient.
- Built-in centralized, CSA-certified backflow prevention.
 - No need for costly additional backflow preventers.
 - No ongoing cost of annual inspections of reduced pressure zone (RPZ) devices.
- Eliminates the need for purchased water and the daily waterline tablets necessary for bottle systems.
- Enables dedicated waterlines for up to eight operatories.
 - Multiple systems can be networked into zones.
 - Dedicated waterlines can easily be turned on and off individually.




VistaClear™ Centralized Water Filtration System

| Ref. # | Description | Quantity |
|---------------|---|----------|
| CVV1000-28-HP | VistaClear™ HP Dental Waterline Filtration System | 1 Each |
| CVR5450 | VistaClear™ Replacement Filter for Legacy Systems | 1 Each |
| CVR9280 | VistaClear™ Retrofit Kit for Upgrade to HP Design | 1 Each |
| CVR9720 | VistaClear™ HP Replacement Filter Annual Kit (Both Filters) | 1 Each |
| CVR9721 | VistaClear™ HP Replacement Sediment Filter | 1 Each |
| CVR9722 | VistaClear™ HP Replacement Carbon Filter | 1 Each |
| CVDM125-10C | VistaClear™ HP 10-Op Distribution Manifold | 1 Each |

VistaCool™ Direct-to-Drain Cooling System for Autoclave Wastewater

- Direct-to-drain thermal reduction system that dramatically reduces the temperature of autoclave wastewater — eliminating condenser bottles and tanks.
- Eliminates the need for office staff to empty steaming-hot water bottles manually – no user maintenance required.
- Prevents autoclave wastewater from melting plumbing.
- Built-in backflow prevention meets or exceeds local plumbing code.
- Helps prevent mold & mildew caused by spills and humid environments created by exhaust from pitchers and bottles.
- Helps protect cabinets from rot and rust by sending cooled wastewater directly and safely to the drain.
- Beyond periodic visual inspections of the tubing, there is nothing more to do to ensure proper operation.



VistaCool™ Direct-to-Drain Cooling System for Autoclave Wastewater

| Ref. # | Description | Quantity |
|---------|--|----------|
| CVV7501 | VistaCool™ Single-Autoclave Thermal Reduction System | 1 Each |
| CVV7502 | VistaCool™ Double-Autoclave Thermal Reduction System | 1 Each |
| CVS7545 | VistaCool™ Adjustable Tank Stand for Single-Autoclave System | 1 Each |
| CVS7547 | VistaCool™ Adjustable Tank Stand for Double-Autoclave System | 1 Each |
| CVS7560 | VistaCool™ Adapter Kit for Midmark M9 & M11 Autoclaves | 1 Each |

Dental Water System Accessories

DrainHub™ Multi-Port Drain Adapter

- CSA-certified, patent-pending device accepts up to eight drain connections in one small hub—saving space, cost and time.
- Easy to install: top and bottom slip hubs accept 1-1/2 in pipe; bottom hub glues into 2 in coupling.
- 1/2" NPT threaded ports accept a variety of push-to-connect, compression, and hose barb adapter fittings.



DrainHub™ Multi-Port Drain Adapter

| Ref. # | Description | Quantity |
|-----------|--|--------------------|
| CVDH158 | DrainHub™ Multi-Port Drain Adapter (8-Port) | 1 Each |
| CVDH8200 | DrainHub™ Fittings Accessory Kit | 20 Asstd. Fittings |
| CVDH15882 | DrainHub™ (8-Port) with Fittings Accessory Kit | 1 Each |

VistaCheck™ Dual Check Valve Backflow Preventer

- CSA-certified dual check valve backflow prevention device can eliminate the need for reduced pressure zone (RPZ) valves, vacuum breakers, and air gaps.
- Made from medical-grade components suitable for food contact.
- Generally accepted by state and local code authorities as suitable for backflow and cross-connection protection.



VistaCheck™ Dual Check Valve Backflow Preventer

| Ref. # | Description | Quantity |
|-----------|---|----------|
| CVVC250-A | VistaCheck™ Dual Check Backflow Preventer 1/4 in (0.635 cm) | 1 Each |



Equipment

Flexible options for your practice to deliver conscious sedation



Portable

- Can move between operatories for increased flexibility.
- Occupies less space and no requirement for central storage unit.
- Utilizes E cylinders and optimal solution when frequency of administration is low.
- Available in enclosed option to keep cylinders out of sight and make disinfection faster.



Centralized

- Optimal solution when nitrous oxide/oxygen sedation is used frequently.
- Requires central storage area for cylinders and initial investment in gas storage system.
- Utilizes G and H cylinders so more cost-effective over time.



Chairside

- Eliminate long mixed gas and vacuum hoses draping from adjacent cabinets.
- Create uncluttered, safer work zone for doctor and staff.
- Minimize displacement of patient mask caused by drag/pull and weight of long hoses.
- Kits available for a wide variety of dental chair brands.

Please visit the website at accutron-inc.com to view our full product catalog or contact our Customer Care team with questions on how to best set up nitrous oxide/oxygen sedation in your practice.



For more information, please visit:



Accutron-inc.com

Complete solutions for nitrous oxide/oxygen sedation

With over 45 years of experience, Accutron by Crosstex provides industry-leading solutions for nitrous oxide/oxygen sedation. The wide portfolio of products are designed to meet the needs of the clinician and patient for conscious sedation.

Nitrous oxide/oxygen sedation can provide a valuable additional revenue stream for your practice while offering patients relaxed, comfortable experiences. The innovative portfolio of Accutron Nitrous Oxide/Oxygen Sedation equipment and accessories, includes single-use masks, scavenging systems, manifolds, flowmeter systems and associated components.





Flowmeters

Digital Ultra™ Flushmount Flowmeter

- Automatic compensation to easily adjust total flow or gas ratio while the other stays constant.
- Both digital and electronic flowtube displays for ease of use.
- Maximum nitrous oxide ratio at 70% for patient safety.
- Audible and visual alarms to notify of low oxygen.
- Oxygen fail-safe and flush for additional safety.
- Emergency oxygen resuscitation connector.
- Sealed flush surface makes it easy to disinfect or barrier-protect.
- Reliable and accurate with minimal maintenance.
- Available with remote control and printer options, as well as multiple configurations including cabinet mount and portable.



Digital Ultra™ Flushmount Flowmeter

| Ref. # | Description |
|-------------|---|
| 56600 | Digital Ultra™ Flushmount Flowmeter Pkg. G (includes mounting hardware) |
| 56800 | Digital Ultra™ Flushmount Flowmeter Pkg. H (includes mounting hardware, Remote Flow System & RFS Flush Mount Faceplate) |
| 56800-CHAIR | Digital Ultra™ Flushmount Flowmeter Pkg. H CHAIR (includes Remote Flow System A, 25 ft of Vacuum Hose and 25 ft of Mixed Gas Hoses - requires purchase of appropriate chair brand RFS chairmount kit) |
| 56950 | Digital Ultra™ Flushmount Flowmeter Pkg. I (includes mounting hardware, Accu-Vac & Accu-Vac Flush Panel Mount Faceplate) |
| 51000 | Digital Newport™ Flowmeter System (enclosed 4-Cylinder) |

For full list, visit accutron-inc.com

Manifolds

Digi-Flo™ Automatic Switching Manifold

- Automatic cylinder changeover.
- Alarm with audio and visual warnings.
- Single or multiple zone valve and/or alarm panel capacity.
- Clear and precise LED numeric line pressure readings.
- Easily accessible components.
- Optional conduit adapter box available if local codes require wiring inside conduit.



Digi-Flo™ Automatic Switching Manifold

| Ref. # | Description |
|--------|---|
| 49150 | Digi-Flo™ Automatic Switching Manifold/ Alarm Pkg. A, Desk Alarm Panel (includes pre-install kit) |
| 49250 | Digi-Flo™ Automatic Switching Manifold/ Alarm Pkg. B, Wall Alarm Panel (includes pre-install kit) |
| 49350 | Digi-Flo™ Automatic Switching Manifold/ Alarm Pkg. C, Desk Alarm Panel & Security System (includes pre-install kit) |
| 49450 | Digi-Flo™ Automatic Switching Manifold/ Alarm Pkg. D, Wall Alarm Panel & Security System (includes pre-install kit) |

For full list, visit accutron-inc.com

Nasal Masks

Axess™ Low Profile Nasal Mask

- Low profile design provides unencumbered access to the oral cavity.
- Fits snugly without tape or uncomfortable tubes protruding into the patient's nostril.
- Lightweight, autoclavable scavenging circuit minimizes drag and pull, enhancing patient comfort.
- Available in scented and unscented options.
- Translucent design allows visual monitoring of patient's breathing.
- Individually wrapped for greater infection protection.
- Lower cost than competing single-use disposable nasal masks.
- Single-use minimizes cross-contamination; saves staff time.

Not made with natural latex rubber.



Axess™ Low Profile Nasal Mask

| Ref. # | Description | Quantity |
|----------|---------------------------|----------|
| 53037-9 | Small Clear Unscented | 24/Box |
| 53037-12 | Small Birthday Bubblegum | 24/Box |
| 53037-16 | Small Fresh Mint | 24/Box |
| 53035-9 | Medium Clear Unscented | 24/Box |
| 53035-12 | Medium Birthday Bubblegum | 24/Box |
| 53035-16 | Medium Fresh Mint | 24/Box |
| 53034-9 | Large Clear Unscented | 24/Box |
| 53034-12 | Large Birthday Bubblegum | 24/Box |
| 53034-16 | Large Fresh Mint | 24/Box |

ClearView™ Classic Nasal Mask

- Mask-in-mask design minimizes ambient nitrous oxide in operatory.
- Clear outer mask allows visual monitoring of patient's breathing.
- Soft inner mask creates improved facial seal.
- Low profile expands clinical field of view.
- Appealing colors and scents engage and relax patients.
- Single-use minimizes cross-contamination; saves staff time.
- Compatible with other scavenging circuits, including PIP+™ Scavenging Circuits.

Not made with natural latex rubber.



ClearView™ Classic Nasal Mask

| Ref. # | Description | Quantity |
|----------|------------------------------|----------|
| 33037-9 | Pediatric Unscented Grey | 12/Box |
| 33037-10 | Pediatric Outlaw Orange | 12/Box |
| 33037-11 | Pediatric Sassy Strawberry | 12/Box |
| 33037-12 | Pediatric Birthday Bubblegum | 12/Box |
| 33037-17 | Pediatric Groovy Grape | 12/Box |
| 33035-9 | Adult Unscented Grey | 12/Box |
| 33035-10 | Adult Outlaw Orange | 12/Box |
| 33035-12 | Adult Birthday Bubblegum | 12/Box |
| 33035-14 | Adult French Vanilla | 12/Box |
| 33035-16 | Adult Fresh Mint | 12/Box |
| 33034-9 | Large Adult Unscented Grey | 12/Box |
| 33034-10 | Large Adult Outlaw Orange | 12/Box |
| 33034-14 | Large Adult French Vanilla | 12/Box |
| 33034-16 | Large Adult Fresh Mint | 12/Box |



Nasal Masks

Personal Inhaler Plus+™ Nasal Masks

- Sets the standard when it comes to conscious sedation therapy.
- Provides patients with reliable and comfortable nitrous oxide/oxygen delivery.
- Single-use minimizes cross-contamination; saves staff time.
- Lightly scented and brightly colored nasal masks.
- Unscented option also available.
- Available in small, medium and large sizes to provide options for optimal fit.

Not made with natural latex rubber.



Personal Inhaler Plus+™ Nasal Masks

| Ref. # | Description | Quantity |
|----------|--|----------|
| 33017-1 | Small Variety Pack 1 (Orange, Strawberry & Bubblegum) | 24/Box |
| 33017-9 | Small Unscented Grey | 24/Box |
| 33017-10 | Small Outlaw Orange | 24/Box |
| 33017-11 | Small Sassy Strawberry | 24/Box |
| 33017-12 | Small Birthday Bubblegum | 24/Box |
| 33017-17 | Small Groovy Grape | 24/Box |
| 33016-1 | Medium Variety Pack 1 (Orange, Strawberry & Bubblegum) | 24/Box |
| 33016-2 | Medium Variety Pack 2 (Pina Colada, Vanilla & Orange) | 24/Box |
| 33016-3 | Medium Variety Pack 3 (Mint, Peach & Vanilla) | 24/Box |
| 33016-9 | Medium Unscented Grey | 24/Box |
| 33016-10 | Medium Outlaw Orange | 24/Box |
| 33016-11 | Medium Sassy Strawberry | 24/Box |
| 33016-12 | Medium Birthday Bubblegum | 24/Box |
| 33016-13 | Medium Pina Colada | 24/Box |
| 33016-14 | Medium French Vanilla | 24/Box |
| 33016-15 | Medium Perfect Peach | 24/Box |
| 33016-16 | Medium Fresh Mint | 24/Box |
| 33015-3 | Large Variety Pack 3 (Mint, Peach & Vanilla) | 24/Box |
| 33015-9 | Large Unscented Grey | 24/Box |
| 33015-14 | Large French Vanilla | 24/Box |
| 33015-15 | Large Perfect Peach | 24/Box |
| 33015-16 | Large Fresh Mint | 24/Box |

Scavenging Circuits

Axess™ Scavenging Circuits

- Steam-sterilizable circuit is environmentally friendly and cost-effective.
- Slender, lightweight circuit is more comfortable for patients and allows unobstructed access for clinician.
- Minimizes drag and pull, so nasal mask stays in place.

Not made with natural latex rubber.



Axess™ Scavenging Circuits

| Ref. # | Description | Quantity |
|--------|---|-----------|
| 53005 | Axess™ Scavenging Circuit for In-line vacuum controller | 1 Package |
| 53006 | Axess™ Scavenging Circuit for Remote Flow System (RFS™) and Accu-Vac™ | 1 Package |

ClearView™ Scavenging Circuits

- Autoclavable.
- X-Ray compatible.
- Meets NIOSH and ADA flow rate recommendations.
- Fits most brands of flowmeters.

Not made with natural latex rubber.



ClearView™ Scavenging Circuits for Standard Bag Tee

| Ref. # | Description | Quantity |
|--------|--|-----------|
| 43003 | ClearView™ Scavenging Circuit II (includes Scavenging Circuit with Vacuum Control Valve, Corrugated Tube, Reservoir Bag and complimentary ClearView™ Single-Use Nasal Masks) | 1 Package |
| 43005 | ClearView™ Scavenging Circuit I (includes Scavenging Circuit with Vacuum Control Valve and complimentary ClearView™ Single-Use Nasal Masks) | 1 Package |

ClearView™ Scavenging Circuits for Remote Flow Systems & Accu-Vac™

| Ref. # | Description | Quantity |
|--------|--|-----------|
| 43004 | ClearView™ Scavenging Circuit II (includes Scavenging Circuit, Corrugated Tube, Reservoir Bag and complimentary ClearView™ Single-Use Nasal Masks) | 1 Package |
| 43006 | ClearView™ Scavenging Circuit I (includes Scavenging Circuit and complimentary ClearView™ Single-Use Nasal Masks) | 1 Package |



Scavenging Circuits

PIP+™ Scavenging Circuits

- Autoclavable.
- X-ray compatible.
- Meets NIOSH and ADA flow rate recommendations.
- Flexible system accommodates most brands of flowmeters.

Not made with natural latex rubber.



PIP+™ Scavenging Circuits for Standard Bag Tee

| Ref. # | Description | Quantity |
|--------|---|-----------|
| 32007 | PIP+™ Scavenging Circuit II (Includes Scavenging Circuit with Vacuum Control Valve, Corrugated Tube, Reservoir Bag and complimentary PIP+™ Single-Use Nasal Masks) | 1 Package |
| 33010 | PIP+™ Scavenging Circuit I (Includes Scavenging Circuit with Vacuum Control Valve and complimentary PIP+™ Single-Use Nasal Masks) | 1 Package |

PIP+™ Scavenging Circuits for Remote Flow Systems & Accu-Vac™

| Ref. # | Description | Quantity |
|--------|--|-----------|
| 32009 | PIP+™ Scavenging Circuit II (Includes Scavenging Circuit, Corrugated Tube, Reservoir Bag and complimentary PIP+ Single-Use Nasal Masks) | 1 Package |
| 33008 | PIP+™ Scavenging Circuit I (Includes Scavenging Circuit and complimentary PIP+™ Single-Use Nasal Masks) | 1 Package |

Accessories

Axess™ Capnography Adapters

- Attaches between Axess™ Nasal Mask and scavenging circuit.
- Capnography sampling takes place within the nasal mask and directly measures patient's exhalations at the source.
- Female Luer-Lok® port provides a secure connection to capnograph monitor sample line.
- Connects to either right or left side of nasal mask for additional flexibility.
- Single-use and individually wrapped to minimize cross-contamination and save staff time.

Not made with natural latex rubber.



| Axess™ Capnography Adapters | | |
|-----------------------------|-----------------------------|----------|
| Ref. # | Description | Quantity |
| 53055 | Axess™ Capnography Adapters | 24/Box |

ClearView™ Capnography Adapters

- Attaches between ClearView™ Nasal Mask and scavenging circuit.
- Capnography sampling takes place within the nasal mask and directly measures patient's exhalations at the source.
- Female Luer-Lok® port provides a secure connection to capnograph monitor sample line.
- Connects to either right or left side of nasal mask for additional flexibility.
- Single-use and individually wrapped to minimize cross-contamination and save staff time.
- Bundled option with ClearView™ Nasal Masks available for ease of ordering.

Not made with natural latex rubber.



| ClearView™ Capnography Adapters | | |
|---------------------------------|---------------------------------|----------|
| Ref. # | Description | Quantity |
| 43055 | ClearView™ Capnography Adapters | 12/Box |

| ClearView™ Nasal Mask with Capnography Adapter | | |
|--|------------------------------|----------|
| Ref. # | Description | Quantity |
| 33037-9-CAP | Pediatric Unscented Grey | 12/Box |
| 33037-10-CAP | Pediatric Outlaw Orange | 12/Box |
| 33037-11-CAP | Pediatric Sassy Strawberry | 12/Box |
| 33037-12-CAP | Pediatric Birthday Bubblegum | 12/Box |
| 33037-17-CAP | Pediatric Groovy Grape | 12/Box |
| 33035-9-CAP | Adult Unscented Grey | 12/Box |
| 33035-10-CAP | Adult Outlaw Orange | 12/Box |
| 33035-12-CAP | Adult Birthday Bubblegum | 12/Box |
| 33035-14-CAP | Adult French Vanilla | 12/Box |
| 33035-16-CAP | Adult Fresh Mint | 12/Box |
| 33034-9-CAP | Large Adult Unscented Grey | 12/Box |
| 33034-10-CAP | Large Adult Outlaw Orange | 12/Box |
| 33034-14-CAP | Large Adult French Vanilla | 12/Box |
| 33034-16-CAP | Large Adult Fresh Mint | 12/Box |



Masks & PPE

| <div> <div>MASKENOMICS</div> <div> <div>N95</div> <div>LEVEL 3</div> <div>LEVEL 2</div> <div>LEVEL 1</div> </div> <div>FILTRATION SCALE</div> </div> | ASTM F2100-11 Standards | High Fluid Barrier ASTM Level 3 | Moderate Fluid Barrier ASTM Level 2 | Light Fluid Barrier ASTM Level 1 |
|--|--|------------------------------------|--|-------------------------------------|
| | Fluid resistance, mm Hg | 160 | 120 | 80 |
| | BFE, @ 3 microns | ≥98% | ≥98% | ≥95% |
| | PFE, @ 0.1 micron | ≥98% | ≥98% | ≥95% |
| | Delta P, mm H ₂ O/cm ² | <5.0 | <5.0 | <4.0 |
| | Flame spread | Class 1 | Class 1 | Class 1 |
| | | | | |

ASTM specification definitions

- **Fluid resistance** measures the mask's ability to prevent fluids from traveling through the material. The higher the fluid resistance (filtration), the better the protection. (Measured in mm Hg.)
- **BFE (bacterial filtration efficiency)** measures the mask's ability to prevent bacteria from traveling through the material. (Measured at 3 microns.)
- **Delta P (differential pressure)** determines breathing resistance—the higher the Delta P, the less the breathability, but the better the filtration. (Measured in mm H₂O/cm².)
- **PFE (particulate filtration efficiency)** measures the mask's ability to prevent sub-micron particulates from traveling through the material. The size of the particles filtered is critical. (Measured at 0.1 micron.)
- **Flame spread** is a ranking derived by laboratory standard test methodology of a material's propensity to burn rapidly and spread flames.

ASTM Level 3

- Ideal for procedures where moderate to high amounts of fluid, spray and/or aerosols are produced.
- Complex oral surgery.
- Crown preparation.
- Implant placement.
- Periodontal surgery.
- Use of ultrasonic scalers (Magnetostriptive and Piezo).
- Laser based applications.*
- Root planing & Scaling.

ASTM Level 2

- Ideal for procedures where light to moderate amounts of fluid, spray and/or aerosols are produced.
- Limited oral surgery.
- Endodontics.
- Prophylaxis.
- Restoratives/ composites.
- Sealants.

ASTM Level 1

- Ideal for procedures where light amounts of fluid, spray and/or aerosols are produced.
- Patient exams.
- Operatory cleaning/maintenance.
- Impressions.
- Lab trimming, finishing & polishing.
- Orthodontics.

*Masks are considered a secondary control and are not meant to replace recommended primary engineering controls for laser plume exposure.

CAUTION: Like patients, every procedure is unique. Use your clinical judgment to determine the appropriate level of barrier protection based on the length of the procedure, the amount of fluid or aerosol, and standard precautions.

Personal Protection Equipment (PPE) is at the heart of any infection prevention program.

Crosstex manufactures an extensive portfolio of PPE to meet all your practice needs, including our Made-in-USA face masks that meet all three ASTM levels. Our award-winning Secure Fit™ Technology face masks were the first clinical face masks to feature adjustable aluminum nose and chin closures that minimize the potential of inhaling airborne contaminants by reducing gapping around the borders of the mask and provide a custom fit regardless of face shape or size.

Selecting the appropriate mask for a particular procedure is a critical component of your Personal Protective Equipment (PPE) protocol. Although masks may look similar, each mask has notable differences affecting its level of protection. Understanding ASTM performance levels can help ensure your mask will provide appropriate protection to minimize the spread of potentially infectious diseases.





Masks with Secure Fit™ Technology

● Ultra™ Sensitive Earloop Mask with Secure Fit™ Mask Technology

- Fluid resistant outer layer.
- Fluid resistant, extra soft white hypoallergenic inner cellulose layer – will not lint, tear or shred.
- Free of chemicals, inks and dyes.
- Extra-long aluminum nose piece.

Ultra™ Sensitive FogFree™ Earloop Mask

- Provides all the features of Ultra™ Sensitive Earloop Mask with Secure Fit™ Mask Technology with the addition of a FogFree™ Strip on mask inner layer.

Ultra™ Sensitive FogFree™ Earloop Mask with Shield

- Provides all the features of Ultra™ Sensitive FogFree™ Earloop Mask with Secure Fit™ Mask Technology with addition of the attached optically clear, distortion-free shield.



FogFree™ Mask



FogFree™ Mask with shield

Ultra™ Sensitive Earloop Mask with Secure Fit™ Mask Technology

| Ref. # | Description | Quantity |
|----------|-------------|----------|
| GCFCXSSF | White | 50/Box |

Ultra™ Sensitive FogFree™ Earloop Mask with Secure Fit™ Mask Technology

| Ref. # | Description | Quantity |
|-----------|-------------|----------|
| GCFCXSFSF | White | 40/Box |

Ultra™ Sensitive FogFree™ Earloop Mask with Shield with Secure Fit™ Mask Technology

| Ref. # | Description | Quantity |
|---------|-------------|----------|
| GCPWSSF | White | 25/Box |

● Ultra™ Earloop Mask with Secure Fit™ Mask Technology

- Fluid resistant outer layer.
- Fluid resistant white inner layer.
- Extra-long aluminum nose piece.



Ultra™ Earloop Mask with Secure Fit™ Mask Technology

| Ref. # | Description | Quantity |
|----------|-------------|----------|
| GCFCXUSF | Blue | 50/Box |

● Procedural Earloop Mask with Secure Fit™ Mask Technology

- Fluid resistant outer layer.
- Fluid resistant white inner layer.
- Extra-long aluminum nose piece.



Procedural Earloop Mask with Secure Fit™ Mask Technology

| Ref. # | Description | Quantity |
|---------|-------------|----------|
| GCPBLSF | Blue | 50/Box |
| GCPLVSF | Lavender | 50/Box |
| GCPPKSF | Pink | 50/Box |

FILTRATION LEVEL KEY

● N95

● ASTM
Level 3

● ASTM
Level 2

● ASTM
Level 1

● Low
Performance

● Minimum
Performance

● Isofluid™ Plus Pure Earloop Mask with Secure Fit™ Mask Technology

- Fluid resistant outer layer.
- Fluid resistant, extra soft white hypoallergenic inner cellulose layer – will not lint, tear or shred.
- Free of chemicals, inks and dyes.
- Extra-long aluminum nose piece.



Isofluid™ Plus Pure Earloop Mask with Secure Fit™ Mask Technology

| Ref. # | Description | Quantity |
|------------------|-------------|----------|
| GPLUSWHSF | White | 50/Box |

● Isofluid™ Earloop Mask with Secure Fit™ Mask Technology

- Fluid resistant outer layer.
- White medical grade tissue inner layer.
- Extra-long aluminum nose piece.



Isofluid™ Earloop Mask with Secure Fit™ Mask Technology

| Ref. # | Description | Quantity |
|----------------|-------------|----------|
| GCIBLSF | Blue | 50/Box |
| GCILVSF | Lavender | 50/Box |
| GCIPKSF | Pink | 50/Box |

Isofluid™ FogFree™ Earloop Mask

- Provides all the features of Isofluid™ Earloop Mask with Secure Fit™ Mask Technology with the addition of a FogFree™ Strip on mask inner layer.

Isofluid™ FogFree™ Earloop Mask with Secure Fit™ Mask Technology

| Ref. # | Description | Quantity |
|-----------------|-------------|----------|
| GCICXBSF | Blue | 40/Box |

Isofluid™ FogFree™ Earloop Mask with Shield

- Provides all features of Isofluid™ FogFree™ Earloop Mask with Secure Fit™ Mask Technology with the addition of attached optically clear, distortion-free face shield.

Isofluid™ FogFree™ Earloop Mask with Shield with Secure Fit™ Mask Technology

| Ref. # | Description | Quantity |
|-----------------|-------------|----------|
| GCIPWBSF | Blue | 25/Box |



Standard Masks

● Ultra™ Sensitive Earloop Mask

- Fluid resistant outer layer.
- Fluid resistant, extra soft white hypoallergenic inner cellulose layer – will not lint, tear or shred.
- Free of chemicals, inks and dyes.
- Extra-long aluminum nose piece.



Ultra™ Sensitive FogFree™ Earloop Mask

- Provides all the features of Ultra™ Sensitive Earloop Mask with the addition of a FogFree™ Strip on mask inner layer.

Ultra™ Sensitive FogFree™ Earloop Mask with Shield

- Provides all the features of Ultra™ Sensitive FogFree™ Earloop Mask with the addition of attached optically clear, distortion-free face shield.

Ultra™ Sensitive Earloop Mask

| Ref. # | Description | Quantity |
|--------|-------------|----------|
| GCFXCS | White | 50/Box |

Ultra™ Sensitive FogFree™ Earloop Mask

| Ref. # | Description | Quantity |
|---------|-------------|----------|
| GCFXCSF | White | 40/Box |

Ultra™ Sensitive FogFree™ Earloop Mask with Shield

| Ref. # | Description | Quantity |
|--------|-------------|----------|
| GCPWS | White | 25/Box |

● Ultra™ Earloop Mask

- Fluid resistant outer layer.
- Fluid resistant white inner layer.
- Extra-long aluminum nose piece.



Ultra™ FogFree™ Earloop Mask

- Provides all the features of Ultra™ Earloop Mask with the addition of a FogFree™ Strip on mask inner layer.

Ultra™ FogFree™ Earloop Mask with Shield

- Provides all the features of Ultra™ FogFree™ Earloop Mask with the addition of attached optically clear, distortion-free face shield.

Ultra™ Earloop Mask

| Ref. # | Description | Quantity |
|--------|-------------|----------|
| GCFXCU | Blue | 40/Box |

Ultra™ FogFree™ Earloop Mask

| Ref. # | Description | Quantity |
|--------|-------------|----------|
| GCFXCU | Blue | 40/Box |

Ultra™ FogFree™ Earloop Mask with Shield

| Ref. # | Description | Quantity |
|--------|-------------|----------|
| GCPW | Blue | 25/Box |

● Procedural Earloop Mask

- Fluid resistant outer layer.
- Fluid resistant white inner layer.
- Extra-long aluminum nose piece.



Procedural Earloop Mask

| Ref. # | Description | Quantity |
|----------|-----------------|----------|
| GCPBL | Blue | 50/Box |
| GCPLV | Lavender | 50/Box |
| GCPPK | Pink | 50/Box |
| GCPYE | Yellow | 50/Box |
| GCPBL100 | Blue, bulk pack | 100/Box |

FILTRATION LEVEL KEY

● N95

● ASTM
Level 3

● ASTM
Level 2

● ASTM
Level 1

● Low
Performance

● Minimum
Performance

● Surgical Tie-on Mask

- Fluid resistant outer layer.
- Fluid resistant white inner layer.
- Extra-large panel for surgical procedures.
- Extra-long aluminum nose piece.



Surgical Tie-on Mask

| Ref. # | Description | Quantity |
|--------|-------------|----------|
| GCS | Blue | 50/Box |

● Isofluid™ Plus Earloop Mask

- Fluid resistant outer layer.
- Fluid resistant, extra soft white hypoallergenic inner cellulose layer – will not lint, tear or shred.
- Extra-long aluminum nose piece.



Isofluid™ Plus Earloop Mask

| Ref. # | Description | Quantity |
|---------|--------------|----------|
| GPLUSBL | Blue | 50/Box |
| GPLUSKA | Kaleidoscope | 50/Box |
| GPLUSPK | Pink | 50/Box |

● Isofluid™ Earloop Mask

- Fluid resistant outer layer.
- White medical grade tissue inner layer.
- Extra-long aluminum nose piece.



Isofluid™ Earloop Mask

| Ref. # | Description | Quantity |
|----------|-----------------|----------|
| GCIBL | Blue | 50/Box |
| GCIGR | Green | 50/Box |
| GCILV | Lavender | 50/Box |
| GCIPK | Pink | 50/Box |
| GCISA | Sapphire | 50/Box |
| GCITE | Teal | 50/Box |
| GCITQ | Turquoise | 50/Box |
| GCIBL100 | Blue, bulk pack | 100/Box |

Isofluid™ FogFree™ Earloop Mask

- Provides all the features of Isofluid™ Earloop Mask with the addition of a FogFree™ Strip on mask inner layer.

Isofluid™ FogFree™ Earloop Mask with Shield

- Provides all the features of Isofluid™ FogFree™ Earloop Mask with the addition of attached optically clear, distortion-free face shield.

Isofluid™ FogFree™ Earloop Mask

| Ref. # | Description | Quantity |
|--------|-------------|----------|
| GCICXB | Blue | 40/Box |
| GCICXS | Sapphire | 40/Box |
| GCICXT | Turquoise | 40/Box |

Isofluid™ FogFree™ Earloop Mask with Shield

| Ref. # | Description | Quantity |
|--------|-------------|----------|
| GCIPWB | Blue | 25/Box |
| GCIPWS | Sapphire | 25/Box |
| GCIPWT | Turquoise | 25/Box |



Additional Masks

● Isolator™ Plus N95 Particulate Respirator

- NIOSH approved N95 Respirator
- Meets guidelines for TB exposure control
- Resistance to penetration by synthetic blood: 160 mmHg
- Inner/outer fluid resistant layers
- Chin pocket construction

Although the N95 closely resembles a surgical mask, it is a respirator and must be used in accordance with all OSHA regulations regarding respiratory protection.



Isolator™ Plus N95 Particulate Respirator

| Ref. # | Description | Quantity |
|--------|-------------------|----------|
| GPRN95 | Blue/White Stripe | 28/Box |

● Le Petit Masque

- Fluid resistant pink outer layer
- Soft white spunbound inner layer
- ASTM Level 1
- Smaller than traditional earloop masks, for a better fit



Le Petit Masque

| Ref. # | Description | Quantity |
|---------|-------------|----------|
| GCAPTPK | Pink | 50/Box |

● Surgical Molded Mask

- Splash resistant barrier
- Flared edge prevents irritation
- Soft, flexible nose piece

Does not meet ASTM classifications.



Surgical Molded Mask

| Ref. # | Description | Quantity |
|--------|-------------|----------|
| GCLBL | Blue | 50/Box |
| GCPK | Pink | 50/Box |

● Isolite™ Earloop Mask

- Outer medical grade tissue layer
- White medical grade tissue inner layer
- Extra-long aluminum nose piece

Does not meet ASTM classifications.



Isolite™ Earloop Mask

| Ref. # | Description | Quantity |
|--------|-------------|----------|
| GCLBL | Blue | 50/Box |
| GCLPK | Pink | 50/Box |
| GCLYE | Yellow | 50/Box |

FILTRATION LEVEL KEY

● N95

● ASTM
Level 3

● ASTM
Level 2

● ASTM
Level 1

● Low
Performance


● Minimum
Performance

Additional Face Protection

Patient Safety Mask with Shield

- Fluid resistant.
- Covers patient's eyes and nose.
- Protects fragile eye and mucosal tissue from flying particles, spray and splatter.
- Makes patients feel safe, but not "in the dark," with attached clear-vision shield.
- Form-fitting profile, flexible materials provide maximum access to mouth.
- White spunbond inner and outer layers.

Not made with natural rubber latex.



Patient Safety Mask with Shield

| Ref. # | Description | Quantity |
|--------|-------------|----------|
| GCPAT | White | 25/Box |

Face Shield

- Optically clear, distortion-free wrap-around face shield.
- 1.5 in (3.81 cm) foam headband holds shield away from face; "floats" lightly on forehead, with no pressure on temples; vented for increased air flow.
- Protects mask and face from direct splatter; may prolong mask life.
- Sonically welded elastic headband for added strength.
- Anti-fog treatment on inside and outside of shield.



Face Shield

| Ref. # | Description | Quantity |
|--------|--|----------|
| GCSS9 | Full length 9 in Clear Face Shield | 18/Box |
| GCSS9B | Full length 9 in Clear Face Shield, bulk | 72/Box |
| GCSS | ¾ length 7 in Clear Face Shield | 18/Box |
| GCSSB | ¾ length 7 in Clear Face Shield, bulk | 72/Box |

See Clear® Lens Cleaning Wipes

- Effective for cleaning eye protection gear, eyewear, face shields and loupes.
- Ideal for glass and polycarbonate lenses.
- Safe for use on anti-reflective, UV and glare-free coated lenses.
- Anti-fog, anti-static and anti-streaking.
- Will not scratch lenses and dries quickly.
- Light, fresh fragrance.




See Clear® Lens Cleaning Wipes

| Ref. # | Description | Quantity |
|--------|---|----------|
| SEEC | Individual Wipe 6 in x 5 in (15.2 cm x 12.7 cm) | 120/Box |



Mask Feature Guide



- A** Extra-long, enclosed aluminum nose closure forms a strong seal around nose and cheekbones for maximum protection.
- B** Optional attached wrap-around face shield is optically clear and distortion-free.
- C** The elastic, flat ear loops minimize pressure to the ears and are attached to the outside of the mask to eliminate skin irritation and provide the most comfortable fit.
- D** Enclosed aluminum chin closure significantly reduces gapping at the sides and bottom of the mask, minimizing the potential of inhaling airborne contaminants.
- E** Unique vapor barrier on the outside layer blocks moisture and prevents fogging.
- F** Fluid resistant, spunbond outer layer for maximum protection.
- G** Ultrasoft, odorless, medical grade facial tissue outer layer.
- H** Highest quality meltblown filter media for best filtration efficiency and breathability.
- I** Extra soft, white hypoallergenic inner cellulose layer for sensitive skin is free of chemicals, inks and dyes; will not lint tear or shred for ultimate comfort.
- J** Fluid resistant, white, spunbond inner layer.
- K** Soft, white medical-grade facial tissue inner layer.
- L** FogFree™ Strip blocks and absorbs moisture; forms a strong seal preventing fogging of eyewear while it cushions the nose ridge.
- M** Sof-Bond™ Ultrasonic Mask assembly ensures strongest construction and highest quality, eliminating holes and defects.

| | A | B | C | D | E | F | G | H | I | J | K | L | M |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| Ultra™ Sensitive Earloop Mask with Secure Fit™ Mask Technology | • | | • | • | | • | | • | • | | | | • |
| Ultra™ Sensitive FogFree™ Earloop Mask with Secure Fit™ Mask Technology | • | | • | • | • | • | | • | • | | | • | • |
| Ultra™ Sensitive FogFree™ Earloop Mask with shield with Secure Fit™ Mask Technology | • | • | • | • | • | • | | • | • | | | • | • |
| Ultra™ Earloop Mask with Secure Fit™ Mask Technology | • | | • | • | | • | | • | | • | | | • |
| Ultra™ Sensitive Earloop Mask | • | | • | | | • | | • | • | | | | • |
| Ultra™ Sensitive FogFree™ Earloop Mask | • | | • | | • | • | | • | • | | | • | • |
| Ultra™ Sensitive FogFree™ Earloop Mask with shield | • | • | • | | • | • | | • | • | | | • | • |
| Ultra™ Earloop Mask | • | | • | | | • | | • | | • | | | • |
| Ultra™ FogFree™ Earloop Mask | • | | • | | • | • | | • | | • | | • | • |
| Ultra™ FogFree™ Earloop Mask with shield | • | • | • | | • | • | | • | | • | | • | • |
| Procedural Earloop Mask with Secure Fit™ Mask Technology | • | | • | • | | • | | • | | • | | | • |
| Procedural Earloop Mask | • | | • | | | • | | • | | | | | • |
| Surgical Tie-on Mask | • | | • | | | • | | • | | | | | • |
| Isofluid™ Plus Pure Earloop Mask with Secure Fit™ Mask Technology | • | | • | • | | • | | • | • | | | | • |
| Isofluid™ Earloop Mask with Secure Fit™ Mask Technology | • | | • | • | | • | | • | | | • | | • |
| Isofluid™ FogFree™ Earloop Mask with Secure Fit™ Mask Technology | • | | • | • | | • | | • | | | • | • | • |
| Isofluid™ FogFree™ Earloop Mask with Shield with Secure Fit™ Mask Technology | • | • | • | • | | • | | • | | | • | • | • |
| Isofluid™ Plus Earloop Mask | • | | • | | | | | • | • | | | | • |
| Isofluid™ Earloop Mask | • | | • | | | | | • | | • | | | • |
| Isofluid™ FogFree™ Earloop Mask | • | | • | | | | | • | | • | | • | • |
| Isofluid™ FogFree™ Earloop Mask with Shield | • | • | • | | | | | • | | • | | • | • |
| Isolite™ Earloop Mask | • | | • | | | | • | • | | | • | | • |

Mask Ordering Information

| Product Name, Quantity | | Ref. #, Color | | |
|------------------------------------|---|--|---|-------------------|
| | N95 Isolator™ Plus N95 Particulate Respirator 28/Box | GPRN95, Blue/White Stripe | | |
| | ASTM LEVEL 3 Ultra™ Sensitive Earloop Mask with Secure Fit™ Mask Technology 50/Box | GCFCXSSF, White | | |
| | Ultra™ Sensitive FogFree™ Earloop Mask with Secure Fit™ Mask Technology 40/Box | GCFCXSFSF, White | | |
| | Ultra™ Sensitive FogFree™ Earloop Mask with shield with Secure Fit™ Mask Technology 25/Box | GCPWSSF, White | | |
| | Ultra™ Earloop Mask with Secure Fit™ Mask Technology 50/Box | GCFCXUSF, Blue | | |
| | Ultra™ Sensitive Earloop Mask 50/Box | GCFCXS, White | | |
| | Ultra™ Sensitive FogFree™ Earloop Mask 40/Box | GCFCXS, White | | |
| | Ultra™ Sensitive FogFree™ Earloop Mask with Shield 25/Box | GCPWS, White | | |
| | Ultra™ Earloop Mask 40/Box | GCFCXU, White | | |
| | Ultra™ FogFree™ Earloop Mask 40/Box | GCFCX, White | | |
| | Ultra™ FogFree™ Earloop Mask with Shield 25/Box | GCPW, White | | |
| | | ASTM LEVEL 2 Procedural Earloop Mask with Secure Fit™ Mask Technology 50/Box | GCPBLSF, Blue | GCPLVSF, Lavender |
| Procedural Earloop Mask 50/Box | | GCPBL, Blue GCPLV, Lavender | GCPPK, Pink GCPYE, Yellow | GCPPK, Pink |
| Procedural Earloop Mask 100/Box | | GCPBL100, Blue | | |
| Surgical Tie-On Mask 50/Box | | GCS, Blue | | |
| | ASTM LEVEL 1 Isofluid™ Plus Pure Mask with Secure Fit™ Mask Technology 50/Box | GPLUSWHSF, White | | |
| | Isofluid™ Earloop Mask with Secure Fit™ Mask Technology 50/Box | GCIBLSF, Blue | GCILVSF, Lavender | GCIPKSF, Pink |
| | Isofluid™ FogFree™ Earloop Mask with Secure Fit™ Mask Technology 40/Box | GCICXBSF, Blue | | |
| | Isofluid™ FogFree™ Earloop Mask with Shield with Secure Fit™ Mask Technology 25/Box | GCIPWBSF, Blue | | |
| | Isofluid™ Plus Earloop Mask 50/Box | GPLUSBL, Blue | GPLUSKA, Kaleidoscope | GPLUSPK, Pink |
| | Isofluid™ Earloop Mask 50/Box | GCIBL, Blue GCIGR, Green GCILV, Lavender | GCIPK, Pink GCISA, Sapphire GCITE, Teal | GCITQ, Turquoise |
| | Isofluid™ Earloop Mask 100/Box | GCIBL100, Blue | | |
| | Isofluid™ FogFree™ Earloop Mask 40/Box | GCICXB, Blue | GCICXS, Sapphire | GCICXT, Turquoise |
| | Isofluid™ FogFree™ Earloop Mask with Shield 25/Box | GCIPWB, Blue | GCIPWS, Sapphire | GCIPWT, Turquoise |
| | Le Petit Masque 50/Box | GCAPTPK, Pink | | |
| | LOW PERFORMANCE Surgical Molded Mask 50/Box | GCBL, Blue | GCPK, Pink | |
| | MINIMUM PERFORMANCE Isolite™ Earloop Mask 50/Box | GCLBL, Blue | GCLPK, Pink | GCLYE, Yellow |



Patient Towels, Bibs, and Liners

The Crosstex success story began with the development of our first product, the disposable patient towel. For over 60 years, we have been the benchmark of quality in patient towels, bibs, and liners. Our unique manufacturing process:

- Provides the often imitated, but never duplicated, Crosstex embossing pattern, which allows moisture to flow horizontally across the towel, protecting your patient's clothing.
- Offers our exclusive reinforced edge embossing design for added strength and durability where towel clips are placed.
- Bonds the layers of tissue (and poly, when applicable) securely together, eliminating separation of the layers.
- Creates the strength of the towel horizontally through the combination of our proprietary manufacturing equipment and tissue formation. This prevents tearing when a towel clip is applied.



Sani-Tab™ Chain-Free Patient Towels

Don't let infection spread from dangerous clip!

Dental bib clips and chains have been shown to pick up dangerous microorganisms from skin, saliva and plaque, presenting a potential source of cross-contamination in the dental operator. Reduce the risk with Sani-Tab™ Chain-Free Patient Towels — the only chain-free patient towel with adhesive tabs eliminating the need for bib chains and clips. Sani-Tab maximizes infection control and patient comfort while saving time and money!

- Easy to use – simply position towel anywhere on the patient and apply pressure to the tab.
- Bib stays where placed to maximize coverage, unlike bib clips and other disposable towel holders.
- Can be placed high on patient to protect clothing.



Towel Dispenser

- Made from polished acrylic.
- Sits on counter/shelf or easily mounts to wall with supplied mounting kit.
- Accommodates folded towels up to 19 in x 13 in (48.3 cm x 33.0 cm).
- Holds approximately 50 towels.
- Front opening allows for easier grasping of towels.



Towel Dispenser

| Ref. # | Description | Quantity |
|--------|---|----------|
| PTC | 5 in x 6 in x 7 in (12.7 cm x 15.2 cm x 17.7 cm) | Each |

3-Ply Tissue with 1-Ply Poly

| | Size | Fold | Embossing Pattern | Ref. #, Color | | |
|---|--------------------------------------|----------|-------------------|--|--|---|
| Polyback™ Towels 500/Case | 19 in x 13 in (48.3 cm x 33.0 cm) | Crosstex | Horizontal | WPXAAQ, Aqua WPXBL, Blue WPXGR, Green WPXPE, Peach | WPXBG, Beige WPXWH, White WPXDR, Dusty Rose WPXLV, Lavender | WPXSI, Silver Grey WPXYE, Yellow WPXFP, Happy Feet WPXTB, Teddy Bear |
| Poly-Gard™ Towels 500/Case | 19 in x 16 in (48.3 cm x 40.6 cm) | Crosstex | Horizontal | WGXBL, Blue | WGXWH, White | |
| Poly-Gard™ Deluxe Towels 300/Case | 18 in x 30 in (45.7 cm x 76.2 cm) | Vantage | Pebble | W3D1930BL, Blue | | |
| Polyback™ Contour Bib Towels 400/Case | 18 in x 22 in (45.7 cm x 55.9 cm) | Vantage | Pebble | WERUECBL, Blue WERUECGR, Green | WERUEC, White | WERUECDR, Dusty Rose |
| Advantage Plus™ Towels 500/Case | 18 in x 13 in (45.7 cm x 33.0 cm) | Vantage | Vertical | WPXAAQ, Aqua WPXABL, Blue WPXAGR, Green WPXAPE, Peach | WPXABG, Beige WPXAWH, White WPXADR, Dusty Rose WPXALV, Lavender | WPXASI, Silver Grey WPXAYE, Yellow |
| Sani-Tab™ Chain-Free Patient Towel 400/Case | 19 in x 13 in (48.3 cm x 33.0 cm) | Crosstex | Horizontal | WPXBLT, Blue | WPXWHT, White | WPXDRT, Dusty Rose |

2-Ply Tissue with 1-Ply Poly

| | Size | Fold | Embossing Pattern | Ref. #, Color | | |
|--|--------------------------------------|----------|-------------------|---|---|--|
| Econoback™ Towels 500/Case | 19 in x 13 in (48.3 cm x 33.0 cm) | Crosstex | Horizontal | WEXAAQ, Aqua WEXBL, Blue WEXGR, Green WEXPE, Peach WEXBG, Beige | WEXWH, White WEXDR, Dusty Rose WEXLV, Lavender WEXSI, Silver Grey WEXYE, Yellow | WEXFP, Happy Feet WEXTB, Teddy Bear WEXPP, Pink Ribbon |
| Ultra-Gard™ Towels 500/Case | 19 in x 16 in (48.3 cm x 40.6 cm) | Crosstex | Horizontal | WUXBL, Blue | WUXWH, White | WUXYE, Yellow |
| Advantage™ Towels 500/Case | 18 in x 13 in (45.7 cm x 33.0 cm) | Vantage | Vertical | WEXAAQ, Aqua WEXABL, Blue WEXAGR, Green WEXAPE, Peach | WEXABG, Beige WEXAWH, White WEXADR, Dusty Rose WEXALV, Lavender | WEXASI, Silver Grey WEXAYE, Yellow |
| Sani-Tab™ Chain-Free Patient Towels 400/Case | 19 in x 13 in (48.3 cm x 33.0 cm) | Crosstex | Horizontal | WEXBLT, Blue WEXWHT, White | WEXDRT, Dusty Rose | WEXLVT, Lavender |

Extra Heavy Tissue with 1-Ply Poly

| | Size | Fold | Embossing Pattern | Ref. #, Color | | |
|---|--------------------------------------|----------|-------------------|--|--|--------------------------------------|
| Proback™ Disposable Paper Towels 500/Case | 19 in x 13 in (48.3 cm x 33.0 cm) | Crosstex | Horizontal | WEPBAQ, Aqua WEPBBL, Blue WEPBGR, Green WEPBPE, Peach | WEPBWH, White WEPBDR, Dusty Rose WEPBLV, Lavender WEPBSI, Silver Grey | WEPBYE, Yellow WEPBTB, Teddy Bear |

3-Ply Tissue

| | Size | Fold | Embossing Pattern | Ref. #, Color | | |
|--|--------------------------------------|----------|-------------------|---|---|------------------------------------|
| Professional™ Regular Towels 500/Case | 19 in x 13 in (48.3 cm x 33.0 cm) | Crosstex | Horizontal | WTXAAQ, Aqua WTXBL, Blue WTXGR, Green | WTXBG, Beige WTXWH, White WTXDR, Dusty Rose | WTXYE, Yellow WTXFP, Happy Feet |

2-Ply Tissue

| | Size | Fold | Embossing Pattern | Ref. #, Color | |
|---|--------------------------------------|----------|-------------------|---------------|--------------|
| Professional™ Thrift Towels 500/Case | 19 in x 13 in (48.3 cm x 33.0 cm) | Crosstex | Horizontal | WT5BL, Blue | WT5WH, White |

All towels, bibs, and liners not made with natural rubber latex.



Aprons & Gowns

Barrier Gowns

- Touch-tab Velcro® Closure at neckline.
- Hand-sewn ties at waist.
- Elastic cuffs to fit over gloves.
- Open back keeps wearer cooler.
- Made of 100% non-woven spun bonded polypropylene with a polyethylene coating.

Disposable barrier gowns are not intended for use as a general chemical barrier or as a surgical gown.



Barrier Gowns

| Ref. # | Description | Quantity |
|------------|-------------------|----------|
| NBG | One Size Fits All | 10/Bag |

Patient Barrier Aprons

Plastic & Nylon

- Protects the patient from liquids.
- Features two alligator clips attached to the front of the neck band to hold a disposable patient towel/bib and a touch-tab Velcro® Closure for fast and easy fastening.
- Available in three convenient sizes, barrier aprons may be wiped clean quickly and easily with any antibacterial product or warm soap and water.



Nylon

- Lighter in weight.
- Greater patient comfort.

Patient Barrier Aprons, Waist Length 27 in x 30 in (68.6 cm x 193.0 cm)

| Ref. # | Description | Quantity |
|--------------|----------------|----------|
| NN3BL | Blue Nylon | 10/Bag |
| NP3BL | Blue Plastic | 10/Bag |
| NP3WH | White Plastic | 10/Bag |
| NP3YE | Yellow Plastic | 10/Bag |

Patient Barrier Aprons, Knee Length 27 in x 44 in (68.6 cm x 111.8 cm)

| Ref. # | Description | Quantity |
|--------------|----------------|----------|
| NN4BL | Blue Nylon | 10/Bag |
| NP4BL | Blue Plastic | 10/Bag |
| NP4WH | White Plastic | 10/Bag |
| NP4YE | Yellow Plastic | 10/Bag |

Patient Barrier Aprons, Ankle Length 36 in x 54 in (91.4 cm x 137.2 cm)

| Ref. # | Description | Quantity |
|--------------|---------------|----------|
| NP5BL | Blue Plastic | 10/Bag |
| NP5WH | White Plastic | 10/Bag |

Gloves

Ambi-ance™ Ambidextrous Nitrile Powder-Free Gloves

- Soft stretch nitrile.
- Puncture and tear-resistant.
- Textured fingertips and beaded cuff.
- Space-saver packaging – 200 gloves/box.
- Color – Blue.

Not made with natural rubber latex.



Ambi-ance™ Ambidextrous Nitrile Powder-Free Gloves

| Ref. # | Description | Quantity |
|--------------|-------------|----------|
| 19010 | Extra Small | 200/Box |
| 19011 | Small | 200/Box |
| 19012 | Medium | 200/Box |
| 19013 | Large | 200/Box |
| 19014 | Extra Large | 200/Box |

Hand Sanitizers

Sani-Hands® Instant Hand Sanitizing Wipes

- Premoistened wipe delivery system cleans and sanitizes hands.
- Formulated with 70% Ethyl Alcohol.
- Wipe format provides friction to remove soil from hands.
- Contains soothing Aloe and Vitamin E to moisturize the skin.



Sani-Hands® Instant Hand Sanitizing Wipes

| Ref. # | Description | Quantity |
|--------------|--|--------------|
| SAND | Medium Canister 6 in x 7.5 in (15.2 cm x 19.1 cm) | 135/Canister |
| SANDK | Large Canister 6 in x 7.5 in (15.2 cm x 19.1 cm) | 220/Canister |
| SANDI | Individual Packets 5 in x 8 in (12.7 cm x 20.3 cm) | 100/Box |

SaniTyze™ Hand Sanitizer

- Contains 70% Ethyl Alcohol meeting CDC recommended guidelines.*
- Latex compatible—contains no lanolin or petroleum based emollients.
- Contains Aloe Vera and Vitamin E to moisturize, condition and soften skin.
- Quick drying, leaves no sticky residue.
- After application, product does not require rinsing.
- Nut and gluten-free.



SaniTyze™ Hand Sanitizer

| Ref. # | Description | Quantity |
|-------------|----------------------------------|----------|
| JWG | Bottle 4 fl oz (118 mL) | 1 Each |
| JWG8 | Pump Bottle 8 fl oz (237 mL) | 1 Each |
| JWGP | Pump Bottle 18 fl oz (532 mL) | 1 Each |

*CDC 2002 Guideline for Hand Hygiene in Health-Care Settings, Vol. 51, No. RR-16

SaniTyze™ Wall-mounted Dispenser

- Convenient, wall-mounted one-touch dispenser.
- Dispenses a pre-measured amount to minimize waste.
- Easily accommodates SaniTyze™ Hand Sanitizer 1000 mL pouch.



SaniTyze™ Wall-mounted Dispenser

| Ref. # | Description | Quantity |
|----------------|--|----------|
| PGEL1 | Wall-mounted Dispenser for 1000 mL Pouch | 1 Each |
| JWGEL1 | Pouch 1000 mL | 1 Each |
| JWGPRO1 | Promo Pack: 2 Dispensers, plus 6 x 1000 mL Pouches | 1 Each |



Preventives

Notes

Preventives are the bedrock of dental care, and clinicians require efficient, ergonomic tools that improve the patient experience.

Designed with both professionals and patients in mind, more hygienists choose Patient's Choice™ Dental Products for preventives that offer superior patient comfort and safety while improving workflow. Efficient, ergonomic and smooth-running prophylaxis angles team up with our patient-pleasing assortment of fluorides, varnishes and prophylaxis pastes.





Fluoride Varnish & Trays

ProGuard™ 5% Sodium Fluoride Varnish

- Convenient, easy-to-apply, fully enclosed single-dose packaging.
- Minimal set-up and clean-up time.
- Smooth consistency with excellent flow characteristics.
- Efficient fluoride release with high substantivity.
- Sets fast, even in the presence of saliva.
- Dries to a clear finish.
- Available in mint and melon flavors.
- Sweetened with Xylitol.
- Free of: soy, milk, eggs, gluten, aspartame and saccharin.

Not made with natural rubber latex.



ProGuard™ 5% Sodium Fluoride Varnish

| Ref. # | Description | Quantity |
|---------|-------------|----------|
| UFVPM | Mint | 50/Box |
| UFVPMML | Melon | 50/Box |

Sparkle V™ 5% Sodium Fluoride Varnish

- Single-dose packaging.
- Minimal set-up and clean-up time.
- Sets rapidly in the presence of saliva.
- Dries to a natural tooth shade on contact.
- Available in bubblegum and mint flavors.
- Sweetened with Xylitol.
- Free of: dairy, eggs, gluten, soy, peanuts, tree nuts, tree nut derivatives, seed products, aspartame and saccharin.

Not made with natural rubber latex.



Sparkle V™ 5% Sodium Fluoride Varnish

| Ref. # | Description | Quantity |
|--------|-------------|----------|
| UFVBG | Bubble Gum | 120/Box |
| UFVM | Mint | 120/Box |

Zap™ Fluoride Gel

- Smooth, creamy gel with 1.23% APF fluoride ion.
- Advanced thixotropic formula flows under pressure for greater coverage.
- Unique formula for longer-lasting flavor retention.
- Soft dispensing bottle with convenient side window.
- Available in 9 patient-pleasing flavors.
- Sweetened with Xylitol.
- Free of: soy, milk, eggs, gluten, aspartame, and parabens.



Zap™ Fluoride Gel

| Ref. # | Description | Quantity |
|---------|----------------|----------------|
| UFGABG | Bubble Gum | 16 oz (480 mL) |
| UFGAC | Cherry | 16 oz (480 mL) |
| UFGACC | Cotton Candy | 16 oz (480 mL) |
| UFGAG | Grape | 16 oz (480 mL) |
| UFGAMAR | Marshmallow | 16 oz (480 mL) |
| UFGAOV | Orange-Vanilla | 16 oz (480 mL) |
| UFGAS | Strawberry | 16 oz (480 mL) |
| UFGAW | Watermelon | 16 oz (480 mL) |
| UFGADFM | Dye-Free Mint | 16 oz (480 mL) |

Zap™ Neutral pH Fluoride Gel

- 0.90% fluoride ion.
- Special antibacterial formulation.
- Unique formula for longer-lasting flavor retention.
- Sweetened with Xylitol.
- Free of: soy, milk, eggs, gluten, aspartame, and parabens.




Zap™ Neutral pH Fluoride Gel

| Ref. # | Description | Quantity |
|--------|-------------|----------------|
| UFGNM | Mint | 16 oz (480 mL) |

Zap™ Dual Arch Fluoride Gel Trays

- Soft, anatomical design for comfortable, natural fit.
- Buckle ridges improve fluoride access to teeth.
- 3 color-coded sizes for easy identification.
- Long, interlocking handles for ease of application.
- Unique "space-saver" package for easy access and storage.
- Single-use.



Zap™ Dual Arch Fluoride Gel Trays

| Ref. # | Description | Quantity |
|--------|----------------|----------|
| YFTHS | Small - Yellow | 50/Box |
| YFTHM | Medium - White | 50/Box |
| YFTHL | Large - Blue | 50/Box |



Prophy Angles

Twist™ Plus Oscillating Disposable Propy Angles

- Ergonomic design improves handling.
- 120° reciprocating motion allows for continuous-contact polishing and reduces splatter.
- Robust gear mechanism provides a smooth, continuous motion.
- Recessed grip enhances the ergonomics of the angle, providing a secure and comfortable grasp.
- Smaller head improves visibility and increases maneuverability.
- Enhanced cup fluting improves interproximal application.
- Single-use.




Twist™ Plus Oscillating Disposable Propy Angles

| Ref. # | Description | Quantity |
|-----------|-----------------|----------|
| TPLUSPASC | Soft Cup, Gray | 100/Bag |
| TPLUSPAFC | Firm Cup, White | 100/Bag |
| TPLUSPATB | Tapered Brush | 100/Bag |

Sparkle™ Disposable Propy Angles

- Beveled gear design runs vibration-free/noise-free.
- Flared cup design improves subgingival access.
- Light colored angle with dark colored cup provides greater visibility.
- Single-use.



Sparkle™ Disposable Propy Angles

| Ref. # | Description | Quantity |
|---------|------------------|----------|
| TPASSC1 | Soft Cup, Purple | 100/Box |
| TPASSC5 | Soft Cup, Purple | 500/Box |
| TPASFC1 | Firm Cup, White | 100/Box |
| TPASFC5 | Firm Cup, White | 500/Box |
| TPATBC1 | Tapered Brush | 100/Box |

Sparkle™ EZ Contra Disposable Angles

- Advanced ergonomic design reduces hand fatigue and micro-traumas; allows for greater access.
- Provides all benefits of the Sparkle Propy Angle (see above).



Sparkle™ EZ Contra Disposable Angles

| Ref. # | Description | Quantity |
|--------|------------------|----------|
| TPASEZ | Soft Cup, Purple | 100/Box |
| TPAFEZ | Firm Cup, White | 100/Box |

Twist™ Disposable Propy Angles

- Ergonomically designed angle that provides a 90° reciprocating “back and forth” rotating motion for continuous contact.
- Smooth back and forth motion reduces frictional heat build-up, keeping the tooth surface cooler.
- Unique oscillating movement eliminates splatter.
- Soft web cup and ribbed flair edges, provides easy application and spreading of prophy paste.
- Single-use.



Twist™ Disposable Propy Angles

| Ref. # | Description | Quantity |
|--------|-----------------|----------|
| TPASC | Soft Cup, Gray | 100/Box |
| TPAFC | Firm Cup, White | 100/Box |
| TPATB | Tapered Brush | 100/Box |
| TPAFB | Flat Brush | 100/Box |

Prophy Paste

Sparkle™ Prophy Paste

- Contains 1.23% fluoride ion.
- Formulated to maximize stain removal with minimal enamel loss.
- Advanced formulation provides consistent, splatter-free results.
- Completely rinses out, leaving no residual grit.
- Available in multiple grits and flavors.
- Sweetened with Xylitol.
- Free of: nuts, soy, milk, eggs, gluten, and aspartame.
- Single-use.
- Berrylicious assorted flavors include: Raspberry, Blueberry, and Strawberry.
- Assorted flavors include: Cherry, Mint, and Orange Vanilla.



Sparkle™ Prophy Paste Coarse Grit

| Ref. # | Description | Quantity |
|--------|-------------------------------|-------------------|
| UPC12 | Mint | 12 oz (340 g) Jar |
| UPCM | Mint | 200/Box |
| UPCB | Berrylicious Assorted Flavors | 200/Box |
| UPCBG | Bubblegum | 200/Box |
| UPCCH | Cherry | 200/Box |
| UPCOV | Orange Vanilla | 200/Box |
| UPCA | Assorted Flavors | 200/Box |

Sparkle™ Prophy Paste Medium Grit

| Ref. # | Description | Quantity |
|--------|-------------------------------|----------|
| UPMM | Mint | 200/Box |
| UPMB | Berrylicious Assorted Flavors | 200/Box |
| UPMBG | Bubblegum | 200/Box |
| UPMCH | Cherry | 200/Box |
| UPMOV | Orange Vanilla | 200/Box |
| UPMA | Assorted Flavors | 200/Box |

Sparkle™ Prophy Paste Fine Grit

| Ref. # | Description | Quantity |
|--------|-------------|----------|
| UPFBG | Bubblegum | 200/Box |

Sparkle™ Free Prophy Paste

- All the benefits of Sparkle™ Prophy Paste in a unique fluoride/dye free formula.
- Available in multiple grits and flavors.
- Sweetened with Xylitol.
- Free of: nuts, soy, milk, eggs, gluten, and aspartame.
- Single-use.



Sparkle™ Free Prophy Paste Coarse Grit

| Ref. # | Description | Quantity |
|---------|-----------------|----------|
| UPSFCC | Cinnamon | 200/Box |
| UPSFCE | Fruity | 200/Box |
| UPSFCS | Spearmint | 200/Box |
| UPSFCWC | White Chocolate | 200/Box |

Sparkle™ Free Prophy Paste Medium Grit

| Ref. # | Description | Quantity |
|----------|-----------------|----------|
| UPSFMC | Cinnamon | 200/Box |
| UPSFME | Fruity | 200/Box |
| UPSFMS | Spearmint | 200/Box |
| UPSFMCWC | White Chocolate | 200/Box |



Topical Gel

GumNumb™ Topical Anesthetic Gel

- Contains 20% benzocaine.
- Rapid onset within 30 seconds and a short duration with virtually no systemic absorption.
- Special formulation minimizes salivation.
- Provides temporary relief of pain from local anesthetic injections, periodontal curettage, impression taking and intra-oral radiographs.
- Available in 6 flavors.
- Free of: soy, milk, eggs, gluten, and aspartame.



| GumNumb™ Topical Anesthetic Gel | | |
|---------------------------------|-------------|--------------------|
| Ref. # | Description | Quantity |
| UTBG | Bubblegum | 1 oz (29.6 mL) Jar |
| UTCH | Cherry | 1 oz (29.6 mL) Jar |
| UTG | Grape | 1 oz (29.6 mL) Jar |
| UTM | Mint | 1 oz (29.6 mL) Jar |
| UTPC | Pina Colada | 1 oz (29.6 mL) Jar |
| UTS | Strawberry | 1 oz (29.6 mL) Jar |

Notes



Single-Use Disposables

Notes

The items that run your practice, such as air/water syringe tips, articulating and occlusion materials and essentials such as HVE tips, evacuator tips and saliva ejectors.

Surgical and waste accessories such as scalpels, blades and sharps containers, operator and imaging barriers, and cottons and non-wovens.





Saliva Ejectors


Safe-Flo™ Saliva Ejectors and Valves

- The unique, built-in one-way valve provides a physical barrier to prevent backflow as described in the CDC's warning on saliva ejectors.
- The valve closes when the suction is stopped or restricted to prevent residual fluids from retracting into a patient's mouth.
- The integrity of suction and performance of the saliva ejector or HVE tip is not affected, so there is no need to change how you work in order to protect your patients.
- Available as a complete saliva ejector valve system or as an adapter valve for use with your favorite brand of saliva ejectors or HVE tips.



Safe-Flo™ Saliva Ejectors

- Fully assembled saliva ejector with Safe-Flo™ Valve replaces your current saliva ejector.
- Slightly shortens the exposed length of the saliva ejector, making it less flexible, allowing for easier tissue retraction and better access.



Safe-Flo™ SE Valve

- Use with your favorite saliva ejector brand, allowing you to continue to use your product of choice while protecting your patients against backflow.



Safe-Flo™ HV Valve

- Use with your favorite HVE tip brand, to prevent backflow in HVE lines.



Safe-Flo™ Adapter

- Allows the use of the HV suction line with a standard saliva ejector, providing a higher level of aspiration.



Safe-Flo™ Saliva Ejectors and Valves

| Ref. # | Description | Quantity |
|-----------|---------------------------------|----------|
| M1000WHBU | Safe-Flo™ Saliva Ejector, White | 100/Bag |
| M1001CLBU | Safe-Flo™ Saliva Ejector, Clear | 100/Bag |
| M1010SUSE | Safe-Flo™ SE Valve, Blue | 100/Bag |
| M1030HVSE | Safe-Flo™ HV Valve, Gray | 100/Bag |
| M1020ADAP | Safe-Flo™ Adapter, Blue/Gray | 50/Bag |

Saliva Ejectors

Comfort Plus™ Premium Saliva Ejectors

- The kinder, gentler saliva ejector...takes in fluid, not tissue!
- Tip design enhances patient comfort – will not grab mucosal tissue.
- Facilitates uninterrupted fluid and light debris removal.
- Optimizes fluid removal as tip design prevents tissue from blocking suction.
- Holds shape when formed as result of proprietary manufacturing process.
- Available in multiple color combinations, with and without scent.



Comfort Plus™ Premium Saliva Ejectors

| Ref. # | Description | Quantity |
|---------|---------------------------------|----------|
| ZWWCP | White/White | 100/Bag |
| ZCWCP | Clear/White | 100/Bag |
| ZWPCPBG | White/Pink Bubblegum scented | 100/Bag |
| ZWGCPM | White/Green Mint scented | 100/Bag |

Hygoformic™ Saliva Ejectors

- Combination saliva ejector and tongue retractor.
- Unique coil design is adjustable for maximum patient comfort.
- Smooth edges eliminate irritation – ideal for patients with sensitive tissue.
- Environmentally friendly – biodegradable.



Hygoformic™ Saliva Ejectors

| Ref. # | Description | Quantity |
|--------|-------------|----------|
| ZFF | White | 100/Bag |
| ZFFAD | Adapters | 10/Bag |

Premium™ Saliva Ejectors & Advantage™ Saliva Ejectors

- Holds shape when formed as result of proprietary manufacturing process.
- Soft, non-removable tips are bonded to the tube for maximum patient safety.
- Provides optimal suction.
- Matte finish eliminates glare.
- Available in multiple color combinations.
- Advantage Saliva Ejectors have identical manufacturing process to Premium Saliva Ejectors except with a glossy finish.



Premium™ Saliva Ejectors

| Ref. # | Description | Quantity |
|--------|---------------|----------|
| ZCB | Frosted/Blue | 100/Bag |
| ZCW | Frosted/White | 100/Bag |
| ZWW | White/White | 100/Bag |
| ZPP | Pink/Pink | 100/Bag |
| ZCMC | Multi-color | 100/Bag |

Advantage™ Saliva Ejectors

| Ref. # | Description | Quantity |
|--------|-------------|----------|
| ZBBI | Blue/Blue | 100/Bag |
| ZCBI | Clear/Blue | 100/Bag |
| ZCCI | Clear/Clear | 100/Bag |
| ZWWI | White/White | 100/Bag |



Evacuator Tips

HVE Combo-Tip™ Evacuators

- Manufactured from high density polyethylene.
- 2 tips in one – vented on one end/non-vented on the other.
- Unique “S” tip design eliminates sharp edges enhancing patient safety and comfort.
- Stiff, firm construction allows for cheek retraction.



HVE Combo-Tip™ Evacuators

| Ref. # | Description | Quantity |
|--------|-------------------|----------|
| ZET | Combo Tip™, White | 100/Bag |

MaxVac Plus™ Combo Tip™ Evacuators

- Identical to the Combo Tip in design but not in color.
- Available in 5 décor-pleasing colors.



MaxVac Plus™ Combo Tip™ Evacuators

| Ref. # | Description | Quantity |
|---------|-------------------------------------|------------------|
| ZETMVBL | Blue | 100/Bag |
| ZETMVGR | Green | 100/Bag |
| ZETMVPK | Pink | 100/Bag |
| ZETMVLV | Purple | 100/Bag |
| ZETMVYE | Yellow | 100/Bag |
| ZETMVP | Assorted (2 bags/color per case) | 10 Bags/ Case |

Non-Vented Evacuator Tips

- Identical to Combo Tip™ Evacuators except both ends are non-vented.



Non-Vented Evacuator Tips

| Ref. # | Description | Quantity |
|--------|-------------------|----------|
| ZETNV | Non-vented, White | 50/Bag |

Slotted Tip

- Manufactured from high density polyethylene.
- Slotted end for increased suction.
- Stiff, firm construction allows for cheek retraction.



Slotted Tip

| Ref. # | Description | Quantity |
|--------|----------------|----------|
| ZE550 | Slotted, White | 50/Bag |

Sparkle™ Air Water Syringe Tips

- Requires no adapter – locking groove ensures tip is securely installed.
- Quad Core air channel ensures complete air/water separation.
- Translucent colors allow for visualization of water flow/blockage.
- Beveled edge provides increased patient comfort and safety.
- Rounded non-scratch end prevents O-ring damage.
- Single-use eliminates risk of cross-contamination and cleaning of syringe between patients.
- Quick, easy, economical.
- Available in 5 translucent colors and opaque white.

Made outside of the USA.



| Sparkle™ Air Water Syringe Tips | | |
|---------------------------------|--|----------|
| Ref. # | Description | Quantity |
| BCSAWSCL | Clear | 250/Bag |
| BCSAWSWH | White | 250/Bag |
| BCSAWS | Assorted (Blue, Purple, Yellow, Amber) | 250/Bag |
| BCSAWSB | Assorted (Blue, Purple, Yellow, Amber) | 1500/Bag |

Surgical Aspirator Tips

- Manufactured from durable, high quality polyethylene with a 30° angle.
- Tips are designed to fit all high-volume evacuation lines.
- Strong construction enables tip to be used for cheek retraction as it will not bend.
- Available in 3 tip widths.

Made outside of the USA.



| Surgical Aspirator Tips | | |
|-------------------------|---|----------|
| Ref. # | Description | Quantity |
| ZSABL | Small Tip, Blue 0.0625 in (0.16 cm) | 25/Bag |
| ZSAWH | Standard Tip, White 0.125 in (0.32 cm) | 25/Bag |
| ZSAGR | Large Tip, Green 0.25 in (0.64 cm) | 25/Bag |



Cotton and Non-Wovens

Cotton Rolls

- Made from 100% cotton.
- Contains no cellulose or rayon fibers.
- Soft, pliable, non-linting.
- Will not adhere to mucous membrane.
- Maintains shape when wet.
- Absorbs up to 9 times their weight.
- Available in Premium and Advantage Plus.
- Non-sterile.



Made outside of the USA.

Premium™ Cotton Rolls

| Ref. # | Description | Quantity |
|--------|---|-----------|
| DNC | #2 Medium 1.5 in x 0.375 in (3.8 cm x 1 cm) | 2,000/Box |

Advantage Plus™ Cotton Rolls

| Ref. # | Description | Quantity |
|--------|---|-----------|
| DNAC | #2 Medium 1.5 in x 0.375 in (3.8 cm x 1 cm) | 2,000/Box |

Cotton Rolls, Braided Non-Sterile

- Made from 100%, short fiber cotton, non-sterile.
- Bound by thin strands of polyester filament.
- Starch-free for maximum softness.
- Tightly braided to maintain shape when wet.
- Highly absorbent.
- Will not adhere to mucous membrane.



Made outside of the USA.

Cotton Rolls, Braided Non-Sterile

| Ref. # | Description | Quantity |
|--------|--------------------------------------|-----------|
| DNCB | 1.5 in x 0.375 in (3.8 cm x 1 cm) | 2,000/Box |

Cotton Roll Dispenser Deluxe

- Flip-up top for easy loading.
- Holds up to 50 cotton rolls.
- Slide out drawer dispenses up to 4 rolls.
- Non-skid base for easy dispensing.



Cotton Roll Dispenser Deluxe

| Ref. # | Description | Quantity |
|--------|--|----------|
| PDCBL | Blue 4 in x 3.3 in x 2 in (10.2 cm x 8.4 cm x 5.1 cm) | 1 Each |
| PDCWH | White 4 in x 3.3 in x 2 in (10.2 cm x 8.4 cm x 5.1 cm) | 1 Each |

Cotton Tipped Applicators

- 100% cotton tip.
- All wood dowels.
- Extremely tight bud minimizes linting.
- Applicators are sealed in autoclavable bags of 100.
- Available in Premium™ Cotton Tipped Applicators (domestic) and Advantage Plus™ Cotton Tipped Applicators (imported).



Made outside of the USA.

Premium™ Cotton Tipped Applicators

| Ref. # | Description | Quantity |
|--------|----------------|-----------|
| H3C | 3 in (7.6 cm) | 1,000/Box |
| H6C | 6 in (15.2 cm) | 1,000/Box |

Advantage Plus™ Cotton Tipped Applicators

| Ref. # | Description | Quantity |
|--------|----------------|-----------|
| H3AC | 3 in (7.6 cm) | 1,000/Box |
| H6AC | 6 in (15.2 cm) | 1,000/Box |

Amalgam Squeeze Cloths

- Made from fine, closely woven white cotton.
- Minimal linting.



Amalgam Squeeze Cloths

| Ref. # | Description | Quantity |
|--------|---------------|----------|
| S5C | 3 in (7.6 cm) | 500/Box |

Premium™ All-Gauze Sponges

- Made from USP 100% gauze.
- Available in 8-ply, 12-ply or 16-ply.
- Highly absorbent.
- Non-sterile.

Made outside of the USA.



Premium™ All-Gauze Sponges

| Ref. # | Description | Quantity |
|--------|--|------------|
| ENC212 | 12-Ply 2 in x 2 in (5.1 cm x 5.1 cm) | 8,000/Case |
| ENC2 | 8-Ply 2 in x 2 in (5.1 cm x 5.1 cm) | 5,000/Case |
| ENC312 | 12-Ply 3 in x 3 in (7.6 cm x 7.6 cm) | 4,000/Case |
| ENC38 | 8-Ply 3 in x 3 in (7.6 cm x 7.6 cm) | 4,000/Case |
| ENC416 | 16-Ply 4 in x 4 in (10.2 cm x 10.2 cm) | 2,000/Case |
| ENC412 | 12-Ply 4 in x 4 in (10.2 x 10.2 cm) | 2,000/Case |
| ENC48 | 8-Ply 4 in x 4 in (10.2 x 10.2 cm) | 4,000/Case |

Ultra-Gauze™ Non-Woven Sponges

- 4-Ply, non-sterile.
- 40 g weight (weight per m² of rayon/poly fabric).
- More absorbent than traditional sponges, exclusive 70/30 polyester formed non-woven material replicates the look of woven gauze.
- Ideal for placement on incision/extraction sites as it easily releases, for increased patient comfort.
- Extremely soft, smooth, strong—virtually lint-free.
- Ideal for sustained patient use, general cleaning.



Ultra-Gauze™ Non-Woven Sponges

| Ref. # | Description | Quantity |
|---------|------------------------------------|------------|
| ENCNWU | 2 in x 2 in (5.1 cm x 5.1 cm) | 5,000/Case |
| ENCNWUB | 2 in x 2 in (5.1 cm x 5.1 cm) | 8,000/Case |
| ENC3NWU | 3 in x 3 in (7.6 cm x 7.6 cm) | 4,000/Case |
| ENC4NWU | 4 in x 4 in (10.2 cm x 10.2 cm) | 2,000/Case |



Cotton and Non-Wovens

Cotton-Filled Sponges

- Made from USP cotton gauze.
- Pure, virgin white cotton with 100% cotton-fiber filling.
- Highly absorbent.
- Available in Premium™ Cotton-Filled Sponges, Advantage Plus™ Cotton-Filled Sponges in various gram weights and sizes. Sterile and non-sterile.



Made outside of the USA.

Premium™ Cotton-Filled Sponges, Exodontia, Non-Sterile

| Ref. # | Description | Quantity |
|--------|---|------------|
| ENCC | 2 in x 2 in (5.1 cm x 5.1 cm) 125 g weight | 200/Sleeve |
| ENCC3 | 3 in x 3 in (7.6 cm x 7.6 cm) 125 g weight | 200/Sleeve |
| ENCC4 | 4 in x 4 in (10.2 cm x 10.2 cm) 125 g weight | 200/Sleeve |

Premium™ Cotton-Filled Sponges, Exodontia, Sterile

| Ref. # | Description | Quantity |
|--------|---|------------|
| ESCC | 2 in x 2 in (5.1 cm x 5.1 cm) 125 g weight | 200/Sleeve |

Advantage Plus™ Cotton-Filled Sponges, Exodontia, Non-Sterile

| Ref. # | Description | Quantity |
|--------|--|------------|
| ENACC | 2 in x 2 in (5.1 cm x 5.1 cm) 90 g weight | 200/Sleeve |

Non-Woven Sponges

- 4 Ply, Non-Sterile.
- Unique rayon/polyester blend.
- Provides maximum absorbency.
- Less adhesion to wounds for faster healing.
- Soft, smooth, strong, virtually lint-free.
- Ideal for sustained patient use, disinfecting, general cleaning.
- Available in Premium™ Non-Wovens, Advantage Plus™ Non-Wovens, and Advantage™ Non-Wovens in multiple gram weights and sizes.



Premium™ Non-Woven Sponges, Non-Sterile

| Ref. # | Description | Quantity |
|---------|--|------------|
| ENCNW4K | 2 in x 2 in (5.1 cm x 5.1 cm) 50 g weight | 200/Sleeve |
| ENCNW | 2 in x 2 in (5.1 cm x 5.1 cm) 50 g weight | 200/Sleeve |
| ENC3NW | 3 in x 3 in (7.6 cm x 7.6 cm) 40 g weight | 200/Sleeve |
| ENC4NW | 4 in x 4 in (10.2 cm x 10.2 cm) 40 g weight | 200/Sleeve |

Advantage Plus™ Non-Woven Sponges, Non-Sterile

| Ref. # | Description | Quantity |
|--------|--|------------|
| ENCNWL | 2 in x 2 in (5.1 cm x 5.1 cm) 35 g weight | 200/Sleeve |

Made outside of the USA.

Advantage™ Non-Woven Sponges, Non-Sterile

| Ref. # | Description | Quantity |
|----------|---|------------|
| ENCNWLA | 2 in x 2 in (5.1 cm x 5.1 cm) 30 g weight | 200/Sleeve |
| ENC3NWLA | 3 in x 3 in (7.6 cm x 7.6 cm) 30 g weight | 200/Sleeve |
| ENC4NWLA | 4 in x 4 in (10.2 cm x 10.2 cm) 30 gram weight | 200/Sleeve |

Made outside of the USA.

Ultra-Pure™ Non-Woven Sponges

- 4-ply, non-sterile.
- Biodegradable – a great way to reduce your carbon footprint!
- Made from 100% USA-grown cotton, providing increased comfort and faster absorbency when compared to rayon/polyester and cotton blends.
- Hypoallergenic, no toxic chemicals used in manufacturing process; all natural whitening process (H₂O₂ whitening) without use of bleach and/or chlorine.
- Exclusive automated manufacturing process reduces cross-contamination risks.
- Ideal for placement on incision/extraction sites as it easily releases, for increased patient comfort.
- Lint-free – extremely soft, smooth, strong and will not adhere to instruments.

Made outside of the USA.



Ultra-Pure™ Non-Woven Sponges

| Ref. # | Description | Quantity |
|---------|--|------------|
| ENC100 | 2 in x 2 in (5.1 cm x 5.1 cm) 45 g weight | 200/Sleeve |
| ENC4100 | 4 in x 4 in (10.2 cm x 10.2 cm) 40 g weight | 200/Sleeve |

Cotton & Non-Woven Sponge Chart

| Product | Style | Ply | Sterility | Size | Gram Weight | Quantity | Ref. # |
|------------------------------|---------------|--------|-------------|---------------------------------|-------------|------------|----------|
| Premium™ Gauze Sponges | 100% Gauze | 12-Ply | Non-Sterile | 2 in x 2 in (5.1 cm x 5.1 cm) | - | 8,000/Case | ENC212 |
| | 100% Gauze | 8-Ply | Non-Sterile | 2 in x 2 in (5.1 cm x 5.1 cm) | - | 5,000/Case | EN212 |
| | 100% Gauze | 12-Ply | Non-Sterile | 3 in x 3 in (7.6 cm x 7.6 cm) | - | 4,000/Case | ENC312 |
| | 100% Gauze | 8-Ply | Non-Sterile | 3 in x 3 in (7.6 cm x 7.6 cm) | - | 4,000/Case | ENC38 |
| | 100% Gauze | 16-Ply | Non-Sterile | 4 in x 4 in (10.2 cm x 10.2 cm) | - | 2,000/Case | ENC416 |
| | 100% Gauze | 12-Ply | Non-Sterile | 4 in x 4 in (10.2 cm x 10.2 cm) | - | 2,000/Case | ENC412 |
| | 100% Gauze | 8-Ply | Non-Sterile | 4 in x 4 in (10.2 cm x 10.2 cm) | - | 4,000/Case | ENC48 |
| Premium™ Non-Sterile Sponges | Cotton-filled | - | Non-Sterile | 2 in x 2 in (5.1 cm x 5.1 cm) | 125 | 5,000/Case | ENCC |
| | Cotton-filled | - | Non-Sterile | 3 in x 3 in (7.6 cm x 7.6 cm) | 125 | 4,000/Case | ENCC3 |
| | Cotton-filled | - | Non-Sterile | 4 in x 4 in (10.2 cm x 10.2 cm) | 125 | 2,000/Case | ENCC4 |
| Sterile Sponges | Cotton-filled | - | Sterile | 2 in x 2 in (5.1 cm x 5.1 cm) | 125 | 5,000/Case | ESCC |
| Advantage Plus™ Sponges | Cotton-filled | - | Non-Sterile | 2 in x 2 in (5.1 cm x 5.1 cm) | 90 | 5,000/Case | ENACC |
| Premium™ Sponges | Non-Woven | 4-Ply | Non-Sterile | 2 in x 2 in (5.1 cm x 5.1 cm) | 50 | 5,000/Case | ENCNW |
| | Non-Woven | 4-Ply | Non-Sterile | 3 in x 3 in (7.6 cm x 7.6 cm) | 40 | 4,000/Case | ENC3NW |
| | Non-Woven | 4-Ply | Non-Sterile | 2 in x 2 in (5.1 cm x 5.1 cm) | 50 | 4,000/Case | ENCNW4K |
| | Non-Woven | 4-Ply | Non-Sterile | 4 in x 4 in (10.2 cm x 10.2 cm) | 40 | 2,000/Case | ENC4NW |
| Advantage Plus™ Sponges | Non-Woven | 4-Ply | Non-Sterile | 2 in x 2 in (5.1 cm x 5.1 cm) | 35 | 5,000/Case | ENCNWL |
| Advantage™ Sponges | Non-Woven | 4-Ply | Non-Sterile | 2 in x 2 in (5.1 cm x 5.1 cm) | 30 | 5,000/Case | ENCNWLA |
| | Non-Woven | 4-Ply | Non-Sterile | 3 in x 3 in (7.6 cm x 7.6 cm) | 30 | 4,000/Case | ENC3NWLA |
| | Non-Woven | 4-Ply | Non-Sterile | 4 in x 4 in (10.2 cm x 10.2 cm) | 30 | 2,000/Case | ENC4NWLA |
| Ultra-Gauze™ Sponges | Non-Woven | 4-Ply | Non-Sterile | 2 in x 2 in (5.1 cm x 5.1 cm) | 40 | 5,000/Case | ENCNWU |
| | Non-Woven | 4-Ply | Non-Sterile | 2 in x 2 in (5.1 cm x 5.1 cm) | 40 | 8,000/Case | ENCNWUB |
| | Non-Woven | 4-Ply | Non-Sterile | 3 in x 3 in (7.6 cm x 7.6 cm) | 40 | 4,000/Case | ENC3NWU |
| | Non-Woven | 4-Ply | Non-Sterile | 4 in x 4 in (10.2 cm x 10.2 cm) | 40 | 2,000/Case | ENC4NWU |
| Ultra-Pure™ Sponges | Non-Woven | 4-Ply | Non-Sterile | 2 in x 2 in (5.1 cm x 5.1 cm) | 45 | 5,000/Case | ENC100 |
| | Non-Woven | 4-Ply | Non-Sterile | 4 in x 4 in (10.2 cm x 10.2 cm) | 40 | 2,000/Case | ENC4100 |



Imaging Comfort

BiteWing-Ease™ Barrier Sleeves

- All-in-one seamless barrier sleeve, comfort cushion and bite tab that adjusts snugly to fit all types of sensors.
- Enables easy sensor positioning for a more accurate image.
- Bite tab thickness protects sensor cord.
- Used with digital sensors and Rinn® XCP® system components.
- Protects valuable sensors – quick, easy and economical.



BiteWing-Ease™ Barrier Sleeves

| Ref. # | Description | Quantity |
|-------------|---|----------|
| BWEL | Large, White Adjustable to fit all size 2 sensors | 50/Box |
| BWES | Small, Blue Adjustable to fit all size 0-1 sensors | 50/Box |

Edge-Ease™ Cushions

- Provides a soft, comfortable edge to protect patient's soft tissue from digital sensor and x-ray film abrasions.
- Allows the patient to bite down harder on positioning devices for better positioning, resulting in more accurate images.
- Easily adheres to film packets, phosphor plates and digital sensor barrier sleeves.
- Economical – decreases necessity for retakes, saving both time and money.



Edge-Ease™ Cushions

| Ref. # | Description | Quantity |
|--------------|--|----------|
| BWEEB | Large, Blue Fits long edge of size 2 film and sensors | 300/Box |
| BWEEW | Small, White Fits short edge of size 1-2 film and all 0-1 sensors | 300/Box |

It's-A-Wrap™ Barriers

- Seamless, cushioned barrier sleeve to maximize patient comfort, designed for use with Schick CDR disposable bite holder tabs and Rinn-style positioners.
- Fully adjustable to all sensor types and brands.
- Protects valuable sensors – quick, easy and economical.



It's-A-Wrap™ Barriers

| Ref. # | Description | Quantity |
|---------------|---|----------|
| BWRAPL | Large, Yellow Adjustable to fit all size 2 sensors | 50/Box |
| BWRAPS | Small, White Adjustable to fit and all 0-1 sensors | 50/Box |

Sensor Slippers™ Cushions

- Cushioned, seamless sensor barrier sleeve for use with gripping style holders such as Rinn® Uni-Grip® Holder and Eezee-Grip® Holder.
- Stabilizes sensor in holder and prevents slippage enhancing the quality of the image.
- Creates softer edge that won't abrade corners of the mouth.
- Protects valuable sensors – quick, easy and economical.



Sensor Slippers™ Cushions

| Ref. # | Description | Quantity |
|--------------|-------------------------------------|----------|
| BWSSL | Large, Yellow Fits size 2 sensor | 50/Box |

Slip-Ease™ Cushions

- Unique protective cover slips over all brands of digital sensors, x-ray film, plates and positioning devices.
- Film adheres to cover with Safety Strip to prevent accidental aspiration.
- Smooths potentially sharp corners of film/sensors for increased patient comfort.
- Works best with gripping-style sensor/film holders.
- Requires just one for a full-mouth series using phosphor plates or wireless sensors.
- Protects phosphor plates from scratching.




Slip-Ease™ Cushions

| Ref. # | Description | Quantity |
|--------------|---|----------|
| BWSEL | Large, Pink Fits size 2 film and sensors | 100/Box |

Wrap-Ease™ Cushions

- All-in-one seamless barrier sleeve, comfort cushion and positioning aid that fits all brands of sensors.
- Fully adjustable, all-in-one seamless barrier sleeve and positioning aid that fits all brands of digital sensors while maximizing patient comfort.
- Ideal for vertical bite wing images.
- Eliminates the need for disposable sensor holders.
- Protects valuable sensors – quick, easy and economical.



Wrap-Ease™ Cushions

| Ref. # | Description | Quantity |
|--------------|---|----------|
| BWWEL | Large, White Adjustable to fit all size 2 sensors | 50/Box |
| BWWES | Small, Blue Adjustable to fit all size 0-1 sensors | 50/Box |

Bite Wing Tabs

- Self-adhesive, ensuring accurate, distortion-free bite wing radiographs.
- Holds securely when patients bite down ensuring precise positioning with enhanced patient comfort.
- Pop-out dispenser packaging is convenient, sanitary and economical.
- Attaches easily to both adult and pediatric size film.




Bite Wing Tabs


| Ref. # | Description | Quantity |
|------------|----------------|----------|
| YBT | Bite Wing Tabs | 500/Box |

Digital Sensor Sleeves

- Manufactured from soft polyethylene with smooth seams to eliminate patient discomfort caused by similar sleeves with rough edges.
- Convenient sizing to accommodate all digital imaging sensors.
- Provides reliable protection against cross-contamination.
- Protects valuable sensors.



Smooth seams of the Digital Sensor Sleeve by Crosstex



Rough seams of a traditional sensor sleeve

Digital Sensor Sleeves

| Ref. # | Description | Quantity |
|--------------|---|----------|
| BCXSL | Size 2 Large, Clear 1.625 in x 8.5 in (4.1 cm x 21.6 cm) | 500/Box |
| BCXSS | Size 1 Small, Clear 1.375 in x 8.25 in (3.5 cm x 21 cm) | 500/Box |



Operatory Protection

Bracket Tray Covers

- Manufactured from a unique, heavy weight paper.
- Impedes flow of moisture, minimizing seep-through.
- Available in multiple sizes and colors.



Bracket Tray Covers, Liquid-Proof

- Polyethylene-backed paper prevents seep-through.
- Available in two sizes, white only.

Dual use product:

- Poly-side down – provides positive moisture barrier.
- Poly-side up – provides ideal surface for quick mixing.



Bracket Tray Covers

| | Size | Ref. #, Color |
|--|--|--|
| Weber/Chayes Bracket Tray Covers 1000/Cartron | A 9.5 in x 12.25 in (24.1 cm x 31.1 cm) | FABL, Blue FAWH, White |
| Ritter Bracket Tray Covers 1000/Cartron | B 8.5 in x 12.25 in (21.6 cm x 31.1 cm) | FBBL, Blue FBWH, White FBDR, Dusty Rose FBGR, Green |
| Weber Hi-Lo Bracket Tray Covers 1000/Cartron | C 11 in x 17.5 in (27.9 cm x 44.5 cm) | FCBL, Blue |
| SS White #20 Bracket Tray Covers 1000/Cartron | D 10.25 in x 15.75 in (26 cm x 40 cm) | FDBL, Blue |
| Midwest Bracket Tray Covers 1000/Cartron | E 9 in x 13.5 in (22.9 cm x 34.3 cm) | FEBL, Blue |
| Mini Bracket Tray Covers 1000/Cartron | M 5 in x 8 in (12.7 cm x 20.2 cm) | FMINI, White |
| Standard Round Bracket Tray Covers 1000/Cartron | O 13.5 in diameter (34.3 cm diameter) | FOBL, Blue |
| Viking Bracket Tray Covers 1000/Cartron | V 8.5 in x 11 in (21.6 cm x 27.9 cm) | FVWH, White |
| Cox Bracket Tray Covers 1000/Cartron | X 11 in x 11 in (27.9 cm x 27.9 cm) | FCOX, White |


Bracket Tray Covers, Liquid Proof

| | Size | Ref. #, Color |
|--|--|---------------|
| Ritter Liquid Proof Bracket Tray Covers 500/Cartron | B 8.5 in x 12.25 in (21.6 cm x 31 cm) | FBC, White |
| Adec Liquid Proof Bracket Tray Covers 500/Cartron | F 8.5 in x 11 in (21.6 cm x 27.9 cm) | FF, White |

Headrest Covers, Non-Woven

- Soft, breathable non-woven fabric provides reliable protection against cross-contamination.
- Provides maximum patient comfort as fabric is quiet – does not rustle.
- Available in 2 sizes.

Made outside of the USA.



Standard Headrest Covers, Non-Woven 10 in x 10 in (25.4 cm x 25.4 cm)

| Ref. # | Description | Quantity |
|--------|-------------|----------|
| L0CB | Blue | 500/Case |
| L0C | White | 500/Case |
| L0CL | Lavender | 500/Case |

Extra Large Headrest Covers, Non-Woven 10 in x 13 in (25.4 cm x 33 cm)

| Ref. # | Description | Quantity |
|--------|-------------|----------|
| L3CB | Blue | 500/Case |
| L3C | White | 500/Case |

Headrest Covers, Plastic

- Clear, soft 1 mil thick polyethylene provides impervious barrier between the patient and the headrest for maximum protection from cross-contamination.
- 3/8 inch (0.95 cm) gusset design allows cover to lay flat on headrest for added patient comfort.
- Available in two sizes



Regular Headrest Covers, Plastic 9.5 in x 11 in (24.1 cm x 27.9 cm)

| Ref. # | Description | Quantity |
|--------|-------------|-------------------|
| L0CP | Clear | 250/Dispenser Box |
| L0CPW | White | 250/Dispenser Box |

Jumbo Headrest Covers, Plastic 9.5 in x 14 in (24.1 cm x 36.6 cm)

| Ref. # | Description | Quantity |
|--------|-------------|-------------------|
| L3CP | Clear | 250/Dispenser Box |
| L3CPW | White | 250/Dispenser Box |

Headrest Covers, Polycoated

- Soft facial tissue combined with poly to provide positive fluid barrier for maximum protection from cross-contamination.
- Strong and absorbent.
- Available in multiple sizes and colors.



Headrest Covers, Polycoated

| | Size | Ref. #, Color | | |
|---|--|---------------------------|---------------------------------|----------------|
| Standard Headrest Covers, Polycoated 500/Case | 10 in x 10 in (25.4 cm x 25.4 cm) | L0BL, Blue L0WH, White | L0DR, Dusty Rose L0GR, Green | L0LV, Lavender |
| Extra Large Headrest Covers, Polycoated 500/Case | 10 in x 13 in (25.4 cm x 33.0 cm) | L3BL, Blue L3WH, White | L3DR, Dusty Rose L3GR, Green | L3LV, Lavender |
| Jumbo Headrest Covers, Polycoated 500/Case | 13 in x 13 in (33.0 cm x 33.0 cm) | L13BL, Blue L1WH White | | |
| Extra Jumbo Headrest Covers, Polycoated 500/Case | 11.5 in x 15 in (29.2 cm x 38.1 cm) | L15WH, White | | |



Operatory Protection

Barrier Film with Finger-Lift Edge

- Unique $\frac{3}{16}$ in (0.476 cm) non-adhesive finger-lift edge allows quick, easy placement/removal.
- Light tack adhesive easily attaches to surfaces without leaving residue.
- Ideal for areas that are hard-to-reach, clean or disinfect.
- Available in blue and clear.



Barrier Film with Finger-Lift Edge
4 in x 6 in (10.2 cm x 15.2 cm)

| Ref. # | Description | Quantity |
|--------|-------------|------------------|
| BFBL | Blue | 1200 Sheets/Roll |
| BFCL | Clear | 1200 Sheets/Roll |

iBarrier® Sleeve for Tablets

- Touch sensitive sleeve reduces the risk of cross-contamination.
- Compatible with other tablets.
- Discard after each use.
- Quick, easy, economical.



iBarrier® Sleeve for Tablets

| Ref. # | Description | Quantity |
|--------|--|-----------------|
| BCIB | 8.25 in x 10 in + 1.5 in lip (21 cm x 25.4 cm + 3.8 cm lip) | 100 Sleeves/Box |

Pillow Cases

- Soft facial tissue combined with poly to provide positive fluid barrier for maximum protection from cross-contamination.
- Strong and absorbent.



Made outside of the USA.

Pillow Cases

| Ref. # | Description | Quantity |
|--------|--------------------------------------|-----------------|
| PC2130 | 21 in x 30 in (53.3 cm x 76.2 cm) | 100 Sleeves/Box |

Chair Sleeves

- Universal, clear plastic sleeves for patient chairs.
- Half sleeve design covers chair back and side chair controls.
- Reduces risk of cross-contamination for airborne aerosols and body fluids.
- Quick, easy, economical.



Made outside of the USA.

Chair Sleeves

| Ref. # | Description | Quantity |
|--------|---|------------------|
| BCCS | Clear 27.5 in x 24 in (69.9 cm x 61 cm) | 225 Sleeves/Roll |

Light Sleeves, T-Handle

- Clear plastic sleeve for T-style light handles.
- White indicator arrow for quick identification of die-cut location.
- Will not come off during use.
- Quick, easy, economical.



Made outside of the USA.

Light Sleeves, T-Handle 4 x 6" (10 x 15 cm)

| Ref. # | Description | Quantity |
|--------------|--|----------|
| BCLST | Clear with white arrow 4 in x 5.75 in (10.2 cm x 14.6 cm) | 500/Box |

X-Ray Sleeve

- Universal clear plastic sleeve fits most x-ray heads.
- Reduces risk of cross-contamination for airborne aerosols and body fluids.
- Will not adversely affect quality of the x-ray image.
- Quick, easy, economical.



Made outside of the USA.

X-Ray Sleeve

| Ref. # | Description | Quantity |
|-------------|--|----------|
| BCXR | Clear 15 in x 26 in (38.1 cm x 66 cm) | 250/Box |

Syringe Sleeves

- Clear plastic sleeve easily fits over 3-way syringe, high/low-speed hand pieces and/or air abrasion hand pieces.
- Reduces risks of cross-contamination for airborne aerosols and body fluids.



Available in 3 configurations:

- Without opening – simply place over air/water syringe or handpiece and push through.
- With opening – convenient pre-cut opening for easy insertion.
- With opening and adhesive – Convenient, pre-cut opening for easy insertion with clear adhesive strip that attaches to tubing to prevent sliding of sleeve, creating a closed barrier.



Made outside of the USA.

Syringe Sleeves without Opening

| Ref. # | Description | Quantity |
|-------------|--|----------|
| BCSS | Clear 2.5 in x 10 in (6.4 cm x 25.4 cm) | 500/Box |

Syringe Sleeves with Opening

| Ref. # | Description | Quantity |
|---------------|--|----------|
| BCSSO | Clear 2.5 in x 10 in (6.4 cm x 25.4 cm) | 500/Box |
| BCSSBO | Blue 2.5 in x 10 in (6.4 cm x 25.4 cm) | 500/Box |

Syringe Sleeves with Opening and Clear Adhesive Strip

| Ref. # | Description | Quantity |
|--------------|--|----------|
| BCSSA | Clear 2.5 in x 11.4 in (6.4 cm x 29 cm) | 500/Box |



Operatory Protection

Tray Sleeves

- Clear plastic sleeve for patient trays.
- Secure fit end-fold provides maximum protection and reduces the risk of cross-contamination from airborne aerosols and body fluids.
- Quick, easy, economical.
- Available in 4 convenient sizes.



Tray Sleeves, ADEC/Size B-Ritter

| Ref. # | Description | Quantity |
|-----------------|---|----------|
| BCTS1014 | Clear 10.5 in x 14 in (26.7 cm x 35.6 cm) | 500/Box |

Tray Sleeves, Mini Trays

| Ref. # | Description | Quantity |
|----------------|--|----------|
| BCTS710 | Clear 7.5 in x 10.5 in (19.1 cm x 26.7 cm) | 500/Box |

Made outside of the USA.

Tray Sleeves, Chayes/Size B Ritter

| Ref. # | Description | Quantity |
|-----------------|---|----------|
| BCTS1114 | Clear 11.625 in x 14.5 in (29.5 cm x 36.8 cm) | 500/Box |

Made outside of the USA.

Tray Sleeves, Midwest/Weber

| Ref. # | Description | Quantity |
|-----------------|---|----------|
| BCTS1116 | Clear 11.5 in x 16 in (29.2 cm x 40.6 cm) | 500/Box |

Made outside of the USA.

Lab Transport Bags

- Meets OSHA requirements and recommendation of NCCLS (National Committee for Clinical Laboratory Standards) for transport of medical specimens and/or dental lab cases.
- Convenient two-pouch system keeps Rx/lab slip separate from specimens/lab case as required by OSHA.
- Biohazard symbol for easy identification and OSHA compliance.
- Permanent adhesive closure does not contact gloved hand – tamper-evident.
- White front can be written on, allowing for easy recording of information.
- Quick, easy, economical.

















Lab Transport Bags

| Ref. # | Description | Quantity |
|-------------|-------------------------------------|----------|
| BAGS | 6 in x 10 in (15.2 cm x 25.4 cm) | 250/Box |






Blades and Scalpels



Bard Parker® Blades & Scalpels

| | Size | Blades (Individually packaged, 50 Blades/Box) | | | | | Disposable Scalpels (10/Box) | |
|--|------|---|---|---------------------------------|------------------------------|---------------------------|---------------------------------|------------------------------|
| | | Sterile | | | | Non-Sterile | Sterile (Disposable) | |
| | | CARBON STEEL Rib-Back® | CARBON STEEL Safetylock Rib-Back® | STAINLESS STEEL Conventional | STAINLESS STEEL Protected | CARBON STEEL Rib-Back® | STAINLESS STEEL Conventional | STAINLESS STEEL Protected |
|  | 10 | 371110 | - | 371210 | | 371310 | 37160 | - |
|  | 11 | 371111 | - | 371211 | | 371311 | 371611 | - |
|  | 12 | 371112 | - | 371212 | | 371312 | 371612 | 372612 |
|  | 15 | 371115 | 371153 | 371215 | | 371315 | 371615 | 372615 |
|  | 20 | 371120 | | 371220 | | 371320 | 371620 | |
|  | 21 | 371121 | - | 371221 | 373921 | 37321 | - | - |
|  | 22 | 371122 | - | 371222 | 373922 | 371322 | | - |
|  | 23 | | - | - | 373923 | 371323 | - | - |
|  | 24 | - | - | | | 371324 | - | - |
|  | 25 | - | - | | | 371325 | - | - |
|  | 60 | - | - | - | | 371340 | - | - |
|  | 10A | | | 371700 | | | | |
|  | 12B | | | 371712 | | | 371613 | |
|  | 15C | | | 371716 | | | | |

Handles for Scalpels & Blades

| Conventional Handles, Non-Sterile | Size | Blade Sizes | Ref. # |
|---|------|-------------|--------|
|  | 3 | 10-15 | 371030 |
|  | 4 | 20-25 | 371040 |
|  | 5 | 10-15 | 371050 |
|  | 6 | 20-25 | 371060 |
|  | 7 | 10-15 | 371070 |
|  | 9 | 10-15 | 371090 |

| Protected Handles, Non-Sterile | Size | Blade Sizes | Ref. # |
|--|------|-------------|--------|
|  | 3 | 10-15 | 374030 |



IV Catheters, Needles, and Syringes

Safety IV Catheters

- Insyte™ Autoguard™ Shielded IV Catheter made of BD Vialon™ Catheter Material with BD Instaflash™ Needle Technology.



Safety IV Catheters

| Ref. # | Gauge | Description | Quantity |
|--------|-------|------------------------|----------|
| 381423 | 22 G | 1 in, Blue | 100/Box |
| 381523 | 22 G | 1 in, Blue, with wings | 100/Box |

General Use Luer-Lok® Syringes

- Available in a variety of sizes ranging from 1 mL to 60 mL.
- Sterile.
- Single-Use Polycarbonate.



General Use Luer-Lok® Syringes

| Ref. # | Description | Quantity |
|--------|-------------|----------|
| 309628 | 1 mL | 100/Box |
| 309657 | 3 mL | 100/Box |
| 309646 | 5 mL | 100/Box |
| 302995 | 10 mL | 100/Box |
| 302830 | 20 mL | 100/Box |
| 309653 | 60 mL | 100/Box |

Hypodermic Luer-Lok® Needles

- Available in a variety of gauges and lengths.
- All gauges and lengths come in boxes of 100.
- Sterile.



Hypodermic Luer-Lok® Needles, Sterile

| Ref. # | Gauge | Length | Quantity |
|--------|-------|------------------------|----------|
| 305106 | 30 G | 0.625 in | 100/Box |
| 305128 | 30 G | 1.5 in | 100/Box |
| 305109 | 27 G | 0.625 in | 100/Box |
| 305136 | 27 G | 1.5 in | 100/Box |
| 305111 | 26 G | 0.625 in | 100/Box |
| 305122 | 25 G | 0.625 in | 100/Box |
| 305127 | 25 G | 1.5 in | 100/Box |
| 305145 | 23 G | 1 in | 100/Box |
| 305155 | 22 G | 1 in | 100/Box |
| 305156 | 22 G | 1.5 in | 100/Box |
| 305165 | 21 G | 1 in | 100/Box |
| 305167 | 21 G | 1.5 in | 100/Box |
| 305175 | 20 G | 1 in | 100/Box |
| 305176 | 20 G | 1.5 in | 100/Box |
| 305195 | 18 G | 1 in | 100/Box |
| 305196 | 18 G | 1.5 in | 100/Box |
| 305185 | 18 G | 1.5 in, Thin wall fill | 100/Box |
| 305198 | 16 G | 1.5 in | 100/Box |

Insulin Syringes with Needles

- Insulin Syringe with BD Ultra-Fine™ Needle.
- Detachable needles.
- Intradermal bevel.



Insulin Syringes with Needles

| Ref. # | Gauge | Length | Size | Quantity |
|--------|-------|--------|--------|----------|
| 328468 | 31 G | 8 mm | 0.5 mL | 100/Box |
| 328438 | 31 G | 8 mm | 0.5 mL | 100/Box |

Syringe-needle Combinations

- Available in a variety of gauges, sizes, and lengths.
- All gauges and lengths come in boxes of 100.
- Attached needle included.



1 mL Syringe-needle Combination, Slip tip

| Ref. # | Gauge | Length | Quantity |
|--------|-------|-------------|----------|
| 309623 | 27 G | 0.5 inch in | 100/Box |
| 309626 | 25 G | 0.625 in | 100/Box |

3 mL Syringe-needle Combination, Luer-Lok® Tip

| Ref. # | Gauge | Length | Quantity |
|--------|-------|----------|----------|
| 309570 | 25 G | 0.625 in | 100/Box |
| 309581 | 25 G | 1 in | 100/Box |
| 309582 | 25 G | 1.5 in | 100/Box |
| 309571 | 23 G | 1 in | 100/Box |
| 309569 | 22 G | 0.75 in | 100/Box |
| 309572 | 22 G | 1 in | 100/Box |
| 309574 | 22 G | 1.5 in | 100/Box |
| 309575 | 21 G | 1 in | 100/Box |
| 309577 | 21 G | 1.5 in | 100/Box |
| 309578 | 20 G | 1 in | 100/Box |
| 309579 | 20 G | 1.5 in | 100/Box |

3 mL Syringe-needle Combination, Luer-Lok® Tip, Thin wall

| Ref. # | Gauge | Length | Quantity |
|--------|-------|--------|----------|
| 309588 | 23 G | 1 in | 100/Box |
| 309589 | 23 G | 1.5 in | 100/Box |

5 mL Syringe-needle Combination, Luer-Lok® Tip

| Ref. # | Gauge | Length | Quantity |
|--------|-------|--------|----------|
| 309630 | 22 G | 1 in | 100/Box |
| 309631 | 22 G | 1.5 in | 100/Box |
| 309632 | 21 G | 1 in | 100/Box |
| 309633 | 21 G | 1.5 in | 100/Box |
| 309634 | 20 G | 1 in | 100/Box |
| 309635 | 20 G | 1.5 in | 100/Box |

5 mL Syringe-needle Combination, Luer-Lok® Tip

| Ref. # | Gauge | Length | Quantity |
|--------|-------|--------|----------|
| 309640 | 22 G | 1.5 in | 100/Box |
| 309644 | 20 G | 1 in | 100/Box |
| 309645 | 20 G | 1.5 in | 100/Box |



Surgical & Waste Accessories

Sharps Containers

- Convenient, vertical “point-first” drop.
- One-way funnel valve to prevent needlesticks.



Multi-Use, One Piece

| Ref. # | Description | Quantity |
|--------|---|----------|
| 305557 | 1.4 qt 7.5 in x 3 in x 3 in | 1 Each |
| 305488 | 3.3 qt, 5 in x 9.1 in x 5.3 in | 1 Each |
| 305489 | 6.9 qt, 11.5 in x 9.4 in x 5.3 in | 1 Each |
| 305490 | 8.2 qt, 13.4 in x 9.4 in x 5.3 in | 1 Each |
| 305477 | 5.0 gal, 18 in x 10.5 in x 7.5 in | 1 Each |
| 305491 | 5.0 gal, Open 18 in x 10.5 in x 7.5 in | 1 Each |

Counterbalanced Door

| Ref. # | Description | Quantity |
|--------|--|----------|
| 305426 | 3.3 qt, Red 12 in x 12 in x 4.8 in | 1 Each |
| 305427 | 5.4 qt, Clear 10.75 in X 10.75 in X 4 in | 1 Each |
| 305422 | 2.0 gal, Pearl 12.45 in x 10.72 in x 6 in | 1 Each |
| 305435 | 2.0 gal, Red 12.45 in x 10.72 in x 6 in | 1 Each |
| 305436 | 3.0 gal, Red 16.6 in x 10.7 in x 6 in | 1 Each |

Counterbalanced Side Entry 12 in x 12 in x 4.8 in

| Ref. # | Description | Quantity |
|--------|---------------|----------|
| 305443 | 5.4 qt, Red | 1 Each |
| 305425 | 5.4 qt, Pearl | 1 Each |
| 305444 | 5.4 qt, Pearl | 1 Each |

Counterbalanced Door 12 in x 12 in x 4.8 in

| Ref. # | Description | Quantity |
|--------|-------------|----------|
| 305517 | 5.4 qt, Red | 1 Each |

Phlebotomy (Dual Access) Sharps Collector

- Clear top provides visualization of collector fill-levels.



Phlebotomy (Dual Access) Sharps Collector

| Ref. # | Description | Quantity |
|--------|---|----------|
| 305635 | 1.5 qt, Red 7 in x 4.9 in x 4 in | 1 Each |
| 305487 | 1.5 qt, Red, Nestable 9.4 in x 4.9 in x 3.9 in | 1 Each |

Wallmate™ Sharps Collector

- 5.4 quart side entry pearl sharps container.
- Attached glove dispenser cabinet.
- Locking wall bracket for added security.
- Additional sizes available.



Wallmate™ Sharps Collector

| Ref. # | Description | Quantity |
|--------|---|----------|
| 305550 | Wallmate Next Generation Starter Kit, 5.4 qt | 1 Each |
| 305447 | Locking Wall Bracket for Patient/Exam Room, 5.4 qt | 1 Each |
| 305448 | Glove Box Cabinet Bracket | 1 Each |
| 305485 | Locking Wall Bracket for Multi-Use One-Piece, Fits 8 qt, 14 qt, 6 gal | 1 Each |

E-Z Scrub™ Brushes

- Single-packaged, preoperative surgical scrub brushes.
- Proprietary polyethylene soft bristles enhance mechanical action while reducing skin irritation.
- Soft sponge provides gentle cleaning on irritated areas or when film mechanical action is not required.
- Ergonomic brush design allows greater user comfort and effectiveness, ultimately improving compliance.
- Impregnated with antimicrobial formulations.



E-Z Scrub™ Brushes

| Ref. # | Description | Quantity |
|--------|--|----------|
| 371073 | Red Surgical Scrub Brush Impregnated with 4% CHG | 1 Each |
| 372053 | Brown Surgical Scrub Brush Impregnated with Povidone Iodine, Minimum Available Iodine 1% | 1 Each |



Articulating and Occlusion Testing Materials

Articulating Paper, Foil & Film

- Combined use of articulating film, foil and paper to mark high and low occlusion points provides the most accurate way to mark every point of articulation.
- Proprietary manufacturing process combined with stringent quality control guarantees a clear, positive marking without smudging of wet or dry fields.
- Color-coded packaging and product descriptions make identification easy.
- Articulating paper available in straight, curved, curved with tabs and horseshoe.
- Articulating film and foil available in straight.
- Quick, easy, economical.



Articulating Paper



Articulating Foil



Articulating Film

Articulating Paper

| Ref. # | Style | Size | Color | Quantity |
|--------|--------------------|-----------|--------------|---------------------------------|
| TPXXT | XX Thin | 32 μ | Blue | 12 Sheets/Book, 12 Books/Box |
| TPXT | Extra Thin | 38 μ | Blue | 12 Sheets/Book, 12 Books/Box |
| TPT | Thin | 71 μ | Blue | 12 Sheets/Book, 12 Books/Box |
| TPTH | Thick | 127 μ | Blue | 12 Sheets/Book, 12 Books/Box |
| TPBR | Combo | 71 μ | Blue, Red | 12 Sheets/Book, 12 Books/Box |
| TPCBRT | Thin, Curved | 71 μ | Blue, Red | 12 Sheets/Book, 12 Books/Box |
| TPCBBR | Curved with tab | 89 μ | Blue, Red | 12 Sheets/Book, 12 Books/Box |
| TPH | Horseshoe | 89 μ | Blue, Red | 12 Sheets/Book, 6 Books/Box |

Articulating Foil

| Ref. # | Style | Size | Color | Quantity |
|--------|---------|----------|--------------|---------------------------------|
| TFOB | 1 Side | 12 μ | Blue | 25 Sheets/Book, 12 Books/Box |
| TFOB2 | 2 Sides | 12 μ | Blue | 25 Sheets/Book, 12 Books/Box |
| TFOR2 | 2 Sides | 12 μ | Red | 25 Sheets/Book, 12 Books/Box |
| TFOBR | Combo | 12 μ | Blue, Red | 25 Sheets/Book, 12 Books/Box |

Articulating Film

| Ref. # | Style | Size | Color | Quantity |
|--------|-------|----------|---------------|---------------------------------|
| TFMBKR | Combo | 13 μ | Black, Red | 25 Sheets/Book, 12 Books/Box |
| TPP | Thin | 13 μ | Blue | 25 Sheets/Book, 12 Books/Box |
| TFMR | Thin | 13 μ | Red | 25 Sheets/Book, 12 Books/Box |
| TFMBR | Combo | 13 μ | Blue, Red | 25 Sheets/Book, 12 Books/Box |

Dental Dam

Dental Dam, Latex, Powder Free/Low Protein

- Manufactured from medium modulus natural latex.
- Powder-free – less than 5mg/sheet (ASTM D6124).
- Low protein – extractable protein content is below FDA allowable claim limit of 50µg/g per sheet (ASTM D5712).
- High tear-strength.
- 100% non-chlorinated.
- Available in two sizes & three gauges, both mint flavored and unflavored.

Gauge measurements:

- Thin 0.125-0.160 / Avg 0.14 mm.
- Medium 0.165-0.200 / Avg 0.18 mm.
- Heavy 0.205-0.240 / Avg 0.22 mm.

Made outside of the USA.



Dental Dam, Latex, Powder Free/Low Protein
5 in x 5 in (12.7 cm x 12.7 cm)

| Ref. # | Description | Quantity |
|--------|--------------------------|---------------|
| 19202 | Thin, Green, Mint | 52 Sheets/Box |
| 19200 | Medium, Green, Mint | 52 Sheets/Box |
| 19201 | Heavy, Green, Mint | 52 Sheets/Box |
| 19102 | Thin, Blue, Unflavored | 52 Sheets/Box |
| 19100 | Medium, Blue, Unflavored | 52 Sheets/Box |
| 19101 | Heavy, Blue, Unflavored | 52 Sheets/Box |

Dental Dam, Latex, Powder Free/Low Protein
6 in x 6 in (15.2 cm x 15.2 cm)

| Ref. # | Description | Quantity |
|--------|--------------------------|---------------|
| 19402 | Thin, Green, Mint | 36 Sheets/Box |
| 19400 | Medium, Green, Mint | 36 Sheets/Box |
| 19401 | Heavy, Green, Mint | 36 Sheets/Box |
| 19302 | Thin, Blue, Unflavored | 36 Sheets/Box |
| 19300 | Medium, Blue, Unflavored | 36 Sheets/Box |
| 19301 | Heavy, Blue, Unflavored | 36 Sheets/Box |

Dental Dam, Non-Latex

- Not made with natural rubber latex.
- Use of non-latex products drastically reduces risk of allergic reactions for both patients and staff.
- Manufactured from polyisoprene.
- Does not contain carbamate esters.
- High tear strength.
- Gauge measurements:

Medium 0.165-0.200 / Avg 0.18 mm.

Made outside of the USA.



Dental Dam, Non-Latex
6 in x 6 in (15.2 cm x 15.2 cm)

| Ref. # | Description | Quantity |
|--------|----------------------------|---------------|
| 19500 | Medium, Purple, Peppermint | 15 Sheets/Box |



Treatment Accessories

Matrix Strips

- Manufactured of mylar polyester film.
- Resists solvents – will not stretch, curl or break.
- Strong, easy, economical.



Matrix Strips

| Ref. # | Description | Quantity |
|-------------------|--|----------------------|
| KSTRIP500 | 4 in (10 cm) Long 0.375 in (0.95 cm) Wide | 500 Strips/ Tube |
| KSTRIP1000 | 0.002 in (60 µm) Thick | 1000 Strips/ Tube |

Poly Coated Mixing Pads

- Extra heavy weight paper with a high poly finish.
- Use for mixing cements, composites, rubber base materials and silicates.
- Non-skid base.
- Available in 3 sizes.



Non-skid base.

Poly Coated Mixing Pads

| Ref. # | Description | Quantity |
|---------------|------------------------------------|--------------------|
| KPAD33 | 3 in x 3 in (7.6 cm x 7.6 cm) | 100 Sheets/ Pad |
| KPAD36 | 3 in x 6 in (7.6 cm x 15.2 cm) | 100 Sheets/ Pad |
| KPAD66 | 6 in x 6 in (15.2 cm x 15.2 cm) | 100 Sheets/ Pad |

Tongue Depressors

- Manufactured from medical grade, high quality, smooth and splinter-free wood.
- Standard size.
- Self-dispensing box – quick, easy and economical.
- Available in sterile and non-sterile.




Tongue Depressors

| Ref. # | Description | Quantity |
|------------|-------------------------------|----------|
| IC | Non-Sterile 6 in (15.2 cm) | 500/Box |
| ICS | Sterile 6 in (15.2 cm) | 100/Box |

Cups

Plastic Cups

- Unique embossing for greater strength and easy dispensing.
- Double-rolled lip for easy separation and enhanced patient comfort.
- Thick-walled for added strength and patient safety.
- Double-coated for lustrous, high-gloss finish.
- Available in 9 colors.



Plastic Cups
5 oz (148 mL)

| Ref. # | Description | Quantity |
|--------|-------------------|-----------|
| CXBG | Beige | 1000/Case |
| CXBL | Blue | 1000/Case |
| CXDR | Dusty Rose | 1000/Case |
| CXGR | Green | 1000/Case |
| CXLV | Lavender | 1000/Case |
| CXPE | Peach | 1000/Case |
| CXSI | Silver Grey | 1000/Case |
| CXWH | White | 1000/Case |
| CXYE | Yellow | 1000/Case |
| CXCL | Clear/Translucent | 2500/Case |

Cup Dispenser

- Contemporary design fits décor of every office.
- Dispenses fifty 5-ounce/148 mL plastic and wax-coated paper cups—one cup at a time.
- Includes self-adhesive strips for easy wall mounting.



Made outside of the USA.

Cup Dispenser

| Ref. # | Description | Quantity |
|--------|---------------------------|----------|
| PCC | White base with clear top | 1 Each |

Medicine/Mixing Cups

- Features fluid ounce, dram, cc, mL, teaspoon and tablespoon units with clearly marked graduations.
- Convenient for dispensing or mixing.
- Unbreakable polypropylene material can be squeezed to form a pour spout.



Medicine/Mixing Cups

| Ref. # | Description | Quantity |
|--------|-----------------|------------|
| CX1 | 1 fl oz (30 mL) | 100/Sleeve |



Sterilization

Crosstex Biological Monitoring Service — the largest sterilizer testing lab in North America

Millions of spore tests...1 Mission: To protect our diverse clientele. Medical, dental, veterinary, beauty, tattoo & piercing. Offering 24/7 online biological test results & clinical consultations for all clients. Committed to protecting your patients, your staff and your practice. Learn why our clients choose crosstex for sterility assurance:

- Mail-in spore test results in as little as 24 hours for steam, 72 hours for chemical, 7 days for dry heat and EtO upon receipt of spore strips in our laboratory.
- Compatible with dry heat, chemical EtO and steam sterilizers.
- Your choice of immediate test failure notification by fax or email with a follow-up phone call the next day from one of our customer care consultants.
- Sterilization results are faxed or emailed quarterly and annually.
- Check results online 24 hours a day on our secure confidential website.
- Third party documentation of test results provides a tangible record when proof of testing must be provided.
- Highly resistant dual-species bacterial spores provide the most reliable method of testing for an effective sterilization cycle.
- One of the most cost-effective methods of biological monitoring.
- OSAP Corporate Super Sponsor.
- Corporate member of AAMI, serving on many sterilization working groups.
- Participating corporate member with DIAC, DTA, HIDA, ADAA, ADHA.

Control Strip
DO NOT OPEN
DO NOT STERILIZE


Test Strip

Customer Sterilization Information
Not responsible for envelopes not processed due to incomplete information.

| Required Information Needed on Every Test | Optional Information |
|---|--|
| Phone: () | Model #: |
| Customer ID: - (e.g. AA-123) | Operator: |
| Sterilizer #: <small>The sterilizer # identifies the sterilizer on your test report</small> | Test Time: Temp: |
| Sterilizer Brand: | Pressure: (PSI/kPa) |
| Test Date: / / | New Customer? Change of Address? <small>Please complete.</small> |
| Type: <input type="checkbox"/> Steam <input type="checkbox"/> Chemiclave <input type="checkbox"/> EtO <input type="checkbox"/> Dry Heat | Office Name: |
| <small>Please follow instructions found on box. Caution: The packaging of this product contains Dry Natural Rubber Latex. May cause allergic reactions.</small> | Address: |
| 1.800.819.3336 www.crosstexbms.com | City: |
| | State: Zip: |
| | Fax: () |
| | Email: |
| | Contact Name: |

MANUFACTURED FOR
CROSSTEX
STERILIZATION MONITORING SERVICE

For more information, please visit:

 **CrosstexBMS.com**

Our sterilization products are the ideal way to ensure sterilizer performance by providing the highest level of sterility assurance available, combining mechanical, chemical and biological process parameters.

Every load, every day, every week. Because You Can't See Sterile™, Crosstex offers world-class mail-in sterilizer monitoring at CrosstexBMS.com. For in-office biological monitoring, get results in just hours with our ConForm™ 10 Monitoring System for steam sterilizers. SteamPlus™ Type 5 Integrators can be used as an adjunct to the gold standard of biological monitoring.





Biological Monitoring

ConFirm™ 10 In-Office Biological Monitoring System

- In-office biological monitoring system provides final results in just 10 hours.
- Simply process test indicator vial in a normal steam sterilization cycle, remove, crush ampoule and place in incubator along with an unprocessed positive control vial for required 10 hour incubation period.
- Media in test vial will stay purple signifying a passed test and will turn yellow if any spores remain alive, thus signifying a failed test.
- Two incubator types are available: ConFirm™ 10 Digital Incubator and a standard dry block incubator.
- ConFirm™ 10 Incubators are equipped with a universal power supply for for 110 and 220 volts (specify when ordering).
- ISO 11138-1 and ISO 11138-3 Compliant.

Standard Dry Block Incubator made outside the USA.



ConFirm™ 10 In-Office Biological Monitoring System

| Ref. # | Description | Quantity |
|---------|--|----------|
| C10BI25 | Steam Biological Indicators | 25/Box |
| C10SK | Confirm™ 10 Digital Incubator, Record Keeper - Manual (Specify 110 V or 220 V) | 1 Each |
| CRK02 | Confirm™ Record Keeper Booklet | 1 Each |
| CCRUSH | Ampoule Crusher/Activator for In-Office System | 1 Each |
| NDB-601 | Standard Dry Block Incubator, 110 V, US Plug Includes Record Keeper | 1 Each |

ConFirm™ 24 In-Office Biological Monitoring System

- In-office biological monitoring system provides final results in just 24 hours.
- Simply process test indicator vial in a normal steam sterilization cycle, remove, crush ampoule and place in incubator along with an unprocessed control vial for required 24 hour incubation period.
- Media in test vial will stay purple signifying a passed test and will turn yellow if any spores remain alive, thus signifying a failed test.
- Two incubator types are available: ConFirm™ 24 Digital Incubator and a standard dry block incubator.
- ConFirm™ 24 Digital Incubators are equipped with a universal power supply for for 110 and 220 volts (specify when ordering).
- ISO 11138-1 and ISO 11138-3 Compliant.

Standard Incubator made outside the USA.



ConFirm™ 24 In-Office Biological Monitoring System

| Ref. # | Description | Quantity |
|---------|---|----------|
| CSBI25 | Steam Biological Indicators | 25/Box |
| CBMS10 | Starter Kit (25 Indicators, Incubator, Record Keeper, Manual) | 1 Each |
| CRK02 | Confirm™ Record Keeper Booklet | 1 Each |
| CCRUSH | Ampoule Crusher/Activator for In-Office System | 1 Each |
| NDB-601 | Standard Dry Block Incubator, 110 V, US Plug Includes Record Keeper | 1 Each |

ConFirm™ Mail-in Sterilizer Biological Monitoring Service

- The most comprehensive mail-in service in healthcare.
- Third party documentation of testing and test results provides a tangible online record when proof of testing and subsequent results must be provided.
- Results in 24 hours for steam, 72 hours for chemical vapor and 7 days for dry heat and EtO from start of incubation period at processing laboratory.
- Choice of immediate test failure notification by fax or email with a courtesy call the next day from a clinical advisor.
- Complimentary reminder service notifies practice of missed test date by either fax or email.
- Comprehensive test results are sent quarterly and annually by fax or email.
- MSM (Multi-Site Monitoring) Compliance Bridge complimentary service available for practices with multiple locations.
- Spore strips contain highly resistant dual-species bacterial spores providing the most reliable method of testing for an effective sterilization cycle.



ConFirm™ Premium Test Service

Includes pre-paid postage.

Test strips can not be split between two sterilizers.

| Ref. # | Description | Quantity |
|--------|---|----------|
| CST060 | 6 Tests (3 Strip Test, 2 Tests & 1 Control Strip) | 1 Each |
| CST120 | 12 Tests (3 Strip Test, 2 Tests & 1 Control Strip) | 1 Each |
| CST480 | 48 Tests (3 Strip Test, 2 Tests & 1 Control Strip) | 1 Each |

ConFirm™ Value Test Service

Does not include postage.


| Ref. # | Description | Quantity |
|--------|--|----------|
| CVT120 | 12 Tests (2 Strip Test, 1 Test & 1 Control Strip) | 1 Each |
| CVT520 | 52 Tests (2 Strip Test, 1 Test & 1 Control Strip) | 1 Each |



Chemical Monitoring

STEAMPlus™ Type 5 Integrators

- For use in steam sterilizers operating at 121°C to 134°C.
- Enhance sterility assurance protocol by monitoring every sterilization cycle.
- Provides a definitive pass or fail result allowing for release of all non-implant loads before results of spore strip test are received.
- Equivalent in performance to a biological indicator, and offers the highest level of sterility assurance from a chemical indicator.
- ISO 11140-1 Type 5 Compliant.



Does not replace weekly use of a biological indicator because it does not contain spores.

Note: As SteamPlus™ does not contain live organisms to prove lethality it does NOT replace the routine use/recommendations of biological monitoring

STEAMPlus™ Type 5 Integrators

| Ref. # | Description | Quantity |
|---------|---|----------|
| SSI-100 | 4 in x 0.75 in Strip (10.2 cm x 1.9 cm) | 100/Pack |
| SSI-250 | 4 in x 0.75 in Strip (10.2 cm x 1.9 cm) | 250/Pack |

Type 4 Steam Indicator Strips

- For use in steam sterilizers or chemical vapor sterilizers operating at 132°C (270°F).
- During steam sterilization the check mark turns from white to black and during chemical vapor sterilization the check mark turns from white to brown or dark brown.
- Multi-parameter – indicates if all conditions for steam sterilization have been met: time, temperature and presence of steam or chemical vapor.
- Can be used with sterilization trays, pouches, tubing, autoclave bags and CSR wraps.
- ISO 11140 -1 Type 4 Compliant.




Type 4 Steam Indicator Strips

| Ref. # | Description | Quantity |
|--------|--|----------|
| SCK | 3.75 in x 0.625 in strip (9.5 cm x 1.6 cm) | 100/Bag |

Dry Heat Chemical Indicator Strips

- Provides a distinct color change from green to dark when processed.
- Writable strip for record keeping.
- ISO 11140-1 Type 1 Compliant.



Dry Heat Chemical Indicator Strips

| Ref. # | Description | Quantity |
|---------|-------------------------------------|----------|
| DIS-100 | 3.94 in x 0.62 in (1.6 cm x 0.2 cm) | 100/Box |

Instrument Packaging & Monitoring

AirView™ II Bowie-Dick Test Pack

- Specifically designed to easily detect presence of air in pre-vacuum steam sterilizers.
- Small, compact size reduces waste and storage space.
- Easy-to-read test sheet clearly identifies failures / can be filed for record keeping purposes.
- Lead-free ink reduces amount of hazardous waste / can be disposed in regular trash.
- Convenient exterior Chemical Process Indicators allow for easy verification of processed from unprocessed test packs.
- ISO 11140-5 Type 2 Compliant.




AirView™ II Bowie-Dick Test Pack

| Ref. # | Description | Quantity |
|--------|--|------------------------|
| MBD030 | 5 in x 3.75 in x 0.625 in (12.7 cm x 9.5 cm x 1.6 cm) | 30 Test Packs/ Case |

CSR Sterilization Wrap

- Highest bacterial filtration efficiency (BFE) available.*
- Highly permeable – efficient sterilant penetration.
- Petroleum- and polyolefin-free.
- Manufactured from environmentally-friendly cellulose material made from wood pulp, eliminating skin sensitivity issues associated with oily residues inherent in plastic-based material.
- Positive bacterial barrier to both airborne and waterborne contaminants.
- Exceptional tensile strength—wet or dry.
- For use with steam and EtO sterilization processes.
- EN ISO 11607-Part 1 and EN 868 – Part 2 Compliant.



CSR Sterilization Wrap

| Ref. # | Description | Quantity |
|--------|--|------------|
| SW12 | 12 in x 12 in (30.5 cm x 30.5 cm) | 1,000/Case |
| SW15 | 15 in x 15 in (38.1 cm x 38.1 cm) | 500/Case |
| SW18 | 18 in x 18 in (45.7 cm x 45.7 cm) | 500/Case |
| SW20 | 20 in x 20 in (50.8 cm x 50.8 cm) | 500/Case |
| SW24 | 24 in x 24 in (61 cm x 61 cm) | 500/Case |
| SW30 | 30 in x 30 in (76.2 cm x 76.2 cm) | 250/Case |
| SW36 | 36 in x 36 in (91.4 cm x 91.4 cm) | 250/Case |
| SW40 | 40 in x 40 in (101.6 cm x 101.6 cm) | 250/Case |
| SW45 | 45 in x 45 in (114.3 cm x 114.3 cm) | 250/Case |
| SW48 | 48 in x 48 in (121.9 cm x 121.9 cm) | 100/Case |
| SW54 | 54 in x 54 in (137.2 cm x 137.2 cm) | 100/Case |

*Data on File
Made in the USA from imported materials.



Instrument Packaging & Monitoring

Lead-Free Steam Indicator Tape

- Designed for use by health care providers to securely seal wrapped packs and to demonstrate that the unit has been exposed to a steam sterilization process.
- Made with adhesive and ink that is free of natural rubber latex, lead or solvents; safer for your patients, staff and the environment.
- Validated for use with all common dynamic-air-removal (pre-vacuum / SFPP) and gravity steam sterilization cycles, allowing simplified stocking of one tape for all steam sterilization needs.
- Made with a resilient pressure sensitive adhesive which securely seals sterilization packs, ensuring that sterile integrity is maintained during sterilization, transport and storage.
- Distinct ink color change from beige to dark brown/black when exposed to a steam sterilization process, easy to interpret indication of processed and unprocessed packs to avoid mix-ups.
- Saturated crepe paper provides the elongation needed for pack expansion during sterilization.
- Tape surface can be written on with indelible non-toxic ink to document appropriate information (e.g. date, sterilizer, contents, load/cycle).
- ISO 11140-1:2014 Type 1 Compliant.



Lead-Free Steam Indicator Tape

| Ref. # | Description | Quantity |
|----------|-----------------------------------|---------------|
| STLF24MM | 1 in x 60 yd (24 mm x 55 m) | 18 Rolls/Case |
| STLF18MM | 0.75 in x 60 yd (18 mm x 55 m) | 24 Rolls/Case |
| STLF12MM | 0.5 in x 60 yd (12 mm x 55 m) | 36 Rolls/Case |

Sterilization Pouches

Duo-Check™ Self-Seal Sterilization Pouches



- Internal and external indicators change color when processed, are easily visible, and are compatible with steam (indicators turn from pink to dark brown) and EtO (indicators turn from green to gold/brown).
- Blue-tinted transparent film provides easy instrument viewing and identification of punctures or tears in film.
- When properly sealed, sterility is maintained up the three (3) years provided the integrity of the pouch has not been compromised, i.e. pouch becomes wet, damaged or is opened during storage. Please refer to AAMI ST79 Section 11.1.3 Shelf Life, additionally, follow your facility policy on sterility maintenance.
- Meets ISO 11140 and 11607 standards for pouch performance.
- Meets ANSI/AAMI ST79 Comprehensive Guide to Steam Sterilization and Assurance in Healthcare Facilities.
- For use with steam and EtO sterilization processes only.
- Available in 15 sizes.
- Double pouch validation verified.

Duo-Check™ Self-Seal Sterilization Pouches

| Ref. # | Description | Quantity |
|---------|--|----------|
| SCZ | 2 in x 8 in (5.0 cm x 20.3 cm) | 200/Box |
| SCXX | 2.25 in x 4 in (5.7 cm x 10.2 cm) | 200/Box |
| SCX | 2.75 in x 9 in (7.0 cm x 22.9 cm) | 200/Box |
| SCQ | 3.25 in x 9 in (8.9 cm x 22.9 cm) | 200/Box |
| SCXS | 3.5 in x 5.25 in (8.9 cm x 13.3 cm) | 200/Box |
| SCS | 3.5 in x 9 in (8.9 cm x 22.9 cm) | 200/Box |
| SCS5 | 3.5 in x 9 in (8.9 cm x 22.9 cm) | 500/Box |
| SCL322 | 3.5 in x 22 in (8.9 cm x 55.9 cm) | 100/Box |
| SCB1X | 4.25 in x 11 in (10.8 cm x 27.9 cm) | 200/Box |
| SCM | 5.25 in x 10 in (13.3 cm x 25.4 cm) | 200/Box |
| SCL515 | 5.25 in x 15 in (13.3 cm x 38.1 cm) | 200/Box |
| SCL | 7.5 in x 13 in (19.0 cm x 33.0 cm) | 200/Box |
| SCL816 | 8 in x 16 in (20.3 cm x 40.6 cm) | 200/Box |
| SCL1015 | 10 in x 15 in (25.4 cm x 38.1 cm) | 100/Box |
| SCL1215 | 12 in x 15 in (30.5 cm x 38.1 cm) | 100/Box |
| SCL1218 | 12 in x 18 in (30.5 cm x 45.7 cm) | 100/Box |

Duo-Check™ Sterilization Pouches, Pink with a Purpose

| Ref. # | Description | Quantity |
|--------|--|----------|
| SCSPP | 3.5 in x 9 in (8.9 cm x 22.9 cm) | 200/Box |
| SCMPP | 5.25 in x 10 in (13.3 cm x 25.4 cm) | 200/Box |



Sterilization Pouches

Sure-Check™ Self-Seal Sterilization Pouches

- Be sure! Triple your confidence with the world's most advanced sterilization pouches with with Type 4 multi-parameter internal/external indicators for steam sterilization processes.
- The first pouch with built-in internal/external multi-parameter indicators that only change color when all 3 criteria for steam sterilization have been achieved: specific time, specific temperature, and specific exposure to sterilant.
- When used in steam processes pink indicators will turn dark brown. When used in EtO processes green indicator will turn gold/brown.
- Eliminates the extra step of adding separate Type 4 multi-variable internal indicators to each pouch.
- When properly sealed, sterility is maintained up the three (3) years provided the integrity of the pouch has not been compromised, i.e. pouch becomes wet, damaged or is opened during storage. Please refer to AAMI ST79 Section 11.1.3 Shelf Life, additionally, follow your facility policy on sterility maintenance.
- Meets ANSI/AAMI ST79 Comprehensive Guide to Steam Sterilization and Assurance in Healthcare Facilities and CDC Guidelines for Disinfection and Sterilization in Healthcare Facilities 2008.
- For use with steam and EtO sterilization processes only.
- Available in 14 sizes.
- Meets ISO 11140-1 Type 4 and 11607 standards for pouch performance.
- Double pouch validation verified.



Sure-Check™ Self-Seal Sterilization Pouches

| Ref. | Description | Quantity |
|----------|---|----------|
| SCXX2 | 2.25 in x 4 in (5.7 cm x 10.2 cm) | 200/Box |
| SCX2 | 2.75 in x 9 in (7.0 cm x 22.9 cm) | 200/Box |
| SCXS2 | 3.5 in x 5.25 in (8.9 cm x 13.3 cm) | 200/Box |
| SCS2 | 3.5 in x 9 in (8.9 cm x 22.9 cm) | 200/Box |
| SCS25 | 3.5 in x 9 in (8.9 cm x 22.9 cm) | 500/Box |
| SCL3222 | 3.5 in x 22 in (8.9 cm x 55.9 cm) | 100/Box |
| SCW2 | 5.25 in x 6.5 in (13.3 cm x 16.5 cm) | 200/Box |
| SCM2 | 5.25 in x 10 in (13.3 cm x 25.4 cm) | 200/Box |
| SCL5152 | 5.25 in x 15 in (13.3 cm x 38.1 cm) | 200/Box |
| SCL2 | 7.5 in x 13 in (19.0 cm x 33.0 cm) | 200/Box |
| SCL8162 | 8 in x 16 in (20.3 cm x 40.6 cm) | 200/Box |
| SCL10152 | 10 in x 15 in (25.4 cm x 38.1 cm) | 100/Box |
| SCL10172 | 10.5 in x 17 in (26.7 cm x 43.2 cm) | 100/Box |
| SCL12152 | 12 in x 15 in (30.5 cm x 38.1 cm) | 100/Box |
| SCL12182 | 12 in x 18 in (30.5 cm x 45.7 cm) | 100/Box |

Self-Sealing Autoclave Bags, Paper

- Unique hospital double-fold seal with pull-off adhesive tab.
- Gusset feature provides 50% more inner space.
- Color process indicators – steam.



Made outside the USA.

Self-Sealing Autoclave Bags, Paper

| Ref. # | Description | Quantity |
|------------|---|-----------------|
| OC | 2.5 in x 10.5 in +1.5 in Gusset (6.4 cm x 26.7 cm + 3.8 cm Gusset) | 1000/ Carton |
| OCS | 2.5 in x 8.75 in +1.5 in Gusset (6.4 cm x 22.2 cm + 3.8 cm Gusset) | 1000/ Carton |

Self-Seal Nylon Sterilization Pouches

- Manufactured for use with dry heat processes up to 420°F (215°C).
- Material is clear, strong and puncture-resistant with self-seal adhesive closure.
- Complimentary external chemical indicator labels are included with every box.
- Pouches are notched for ease of opening after processing.
- Available in 5 sizes.



Self-Seal Nylon Sterilization Pouches

| Ref. # | Description | Quantity |
|----------------|---------------------------------------|----------|
| NSP-400 | 2 in x 10 in (5.1 cm x 25.4 cm) | 100/Box |
| NSP-405 | 3 in x 10 in (7.6 cm x 25.4 in) | 100/Box |
| NSP-410 | 4 in x 10 in (10.2 cm x 25.4 cm) | 100/Box |
| NSP-420 | 7 in x 10.5 in (17.8 cm x 26.7 cm) | 100/Box |
| NSP-430 | 9.5 in x 13 in (24.1 cm x 33.0 cm) | 100/Box |

Steri-Coil™ Pouch Dividers

- Stainless steel dividers for placement of pouches in the sterilizer.
- Ensures proper loading of sterilizer and sterilization of contents.
- Available in 2 convenient sizes.



Steri-Coil™ Pouch Dividers

| Ref. # | Description | Quantity |
|----------------|--|----------|
| SCD-001 | 3.5 in x 3.5 in x 20 in (8.9 cm x 8.9 cm x 50.8 cm) | 1 Each |
| SCD-012 | 2.5 in x 6 in (6.4 cm x 15.2 cm) | 1 Each |



Sterilization Tubing

Sani-Roll™ Sterilization Tubing

Manufactured from medical-grade paper with external process indicators that change color when exposed to a steam and EtO sterilization processes.

Convenient, economical way to package instruments – use only as much as you need.

- Blue-tint laminate easily identifies if film has been compromised due to punctures and/or tears.
- For use with steam and EtO sterilization processes only.
- Available in multiple sizes and lengths.



Sani-Roll™ Sterilization Tubing

| Ref. # | Description | Quantity |
|----------|-----------------|---------------------|
| SCT2 | 2 in (5.1 cm) | 100 ft (30 m) Roll |
| SCT2200 | 2 in (5.1 cm) | 656 ft (200 m) Roll |
| SCT3 | 3 in (7.6 cm) | 100 ft (30 m) Roll |
| SCT3200 | 3 in (7.6 cm) | 656 ft (200 m) Roll |
| SCT4 | 4 in (10.2 cm) | 100 ft (30 m) Roll |
| SCT4200 | 4 in (10.2 cm) | 656 ft (200 m) Roll |
| SCT6 | 6 in (15.2 cm) | 100 ft (30 m) Roll |
| SCT6200 | 6 in (15.2 cm) | 656 ft (200 m) Roll |
| SCT8 | 8 in (20.3 cm) | 100 ft (30 m) Roll |
| SCT8200 | 8 in (20.3 cm) | 656 ft (200 m) Roll |
| SCT10 | 10 in (25.4 cm) | 100 ft (30 m) Roll |
| SCT12200 | 12 in (30.5 cm) | 656 ft (200 m) Roll |

Sani-Tube™ Nylon Sterilization Tubing

Manufactured from clear, puncture-resistant, 0.04 mm nylon.

- For use in dry heat sterilization processes only.
- Available in 4 sizes, with and without process indicators.



Sani-Tube™ Nylon Sterilization Tubing,
With Indicators

| Ref. # | Description | Quantity |
|--------|----------------|--------------------|
| SI2 | 2 in (5.1 cm) | 100 ft (30 m) Roll |
| SI3 | 3 in (7.6 cm) | 100 ft (30 m) Roll |
| SI4 | 4 in (10.2 cm) | 100 ft (30 m) Roll |
| SI6 | 6 in (15.2 cm) | 100 ft (30 m) Roll |

Sani-Tube™ Nylon Sterilization Tubing,
Without Process Indicator

| Ref. # | Description | Quantity |
|--------|----------------|--------------------|
| SN2 | 2 in (5.1 cm) | 100 ft (30 m) Roll |
| SN3 | 3 in (7.6 cm) | 100 ft (30 m) Roll |
| SN4 | 4 in (10.2 cm) | 100 ft (30 m) Roll |
| SN6 | 6 in (15.2 cm) | 100 ft (30 m) Roll |

Notes

Index

A

| | |
|--|----|
| Advantage Plus™ Cotton-Filled Sponges | 56 |
| Advantage Plus™ Non-Woven Sponges | 56 |
| Advantage Plus™ Towels | 37 |
| Advantage™ Non-Woven Sponges | 56 |
| Advantage™ Saliva Ejectors | 51 |
| Advantage™ Towels | 37 |
| AirView™ II Bowie-Dick Test Pack | 79 |
| Amalgam Squeeze Cloths | 55 |
| Ambi-ance™ Ambidextrous Nitrile Powder-Free Gloves | 38 |
| Aprons | 38 |
| Articulating Film | 70 |
| Articulating Foil | 70 |
| Articulating Paper | 70 |
| Axess™ Capnography Adapters | 25 |
| Axess™ Low Profile Nasal Mask | 21 |
| Axess™ Scavenging Circuits | 23 |

B

| | |
|--|----|
| Bard Parker® Blades & Scalpels | 65 |
| Barrier Film with Finger-Lift Edge | 62 |
| Barrier Gowns | 38 |
| Bibs | 36 |
| Biological Monitoring | 76 |
| BiteWing-Ease™ Barrier Sleeves | 58 |
| Bite Wing Tabs | 59 |
| Blades | 65 |
| Bottle Systems | 14 |
| Bracket Tray Covers | 60 |

C

| | |
|---|----|
| Catheters | 66 |
| Chair Sleeves | 62 |
| Chemical Monitoring | 78 |
| Chemistries | 4 |
| Cleaning Solutions & Detergents | 10 |
| ClearView™ Capnography Adapters | 25 |
| ClearView™ Classic Nasal Mask | 21 |
| ClearView™ Scavenging Circuits | 23 |
| Comfort Plus™ Premium Saliva Ejectors | 51 |
| ConFirm™ 10 In-Office Biological Monitoring System | 76 |
| ConFirm™ 24 In-Office Biological Monitoring System | 76 |
| ConFirm™ Mail-in Sterilizer Biological Monitoring Service | 77 |
| Cotton-Filled Sponges | 56 |
| Cotton & Non-Woven Sponge Chart | 57 |
| Cotton Roll Dispenser Deluxe | 54 |
| Cotton Rolls | 54 |

| | |
|--|----|
| Cotton Tipped Applicators | 54 |
| CrossZyme™ Concentrated Ultrasonic Cleaner | 10 |
| CSR Sterilization Wrap | 79 |
| Cup Dispenser | 73 |
| Cups | 73 |

D

| | |
|---|----|
| Dental Dam | 71 |
| Dental Water Management | 12 |
| Dental Water System Accessories | 17 |
| Dental Water Systems | 16 |
| DentaPure™ 60-Day Independent Water Bottle Cartridge | 14 |
| DentaPure™ 365-Day Independent Water Bottle Cartridge | 14 |
| DentaPure™ 365-Day Municipal Cartridge | 14 |
| Digi-Flo™ Automatic Switching Manifold | 20 |
| Digital Sensor Sleeves | 59 |
| Digital Ultra™ Flushmount Flowmeter | 20 |
| Disposable Scalpels | 65 |
| DrainHub™ Multi-Port Drain Adapter | 17 |
| Dry Heat Chemical Indicator Strips | 78 |
| Duo-Check™ Sterilization Pouches | 81 |

E

| | |
|---|----|
| Econoback™ Towels | 37 |
| Edge-Ease™ Cushions | 58 |
| EfferZyme™ Enzymatic Cleaning Tablets | 10 |
| Equipment | 18 |
| Evacuator Tips | 52 |
| E-Z Scrub™ Brushes | 69 |

F

| | |
|--------------------------|----|
| Face Shield | 33 |
| Flowmeters | 20 |
| Fluoride Gel Trays | 43 |
| Fluoride Varnish | 42 |

G

| | |
|--|----|
| Gauze | 55 |
| General Purpose Ultrasonic Cleaner, Non-ammoniated Concentrate | 11 |
| Gloves | 38 |
| Gowns | 38 |
| GumNumb™ Topical Anesthetic Gel | 46 |

H

| | |
|-------------------------------------|----|
| Handles for Scalpels & Blades | 65 |
| Hand Sanitizers | 39 |
| Headrest Covers | 61 |

| | |
|---|----|
| High-Level Disinfectants/Sterilants & Accessories | 9 |
| HVE Combo-Tip™ Evacuators..... | 52 |
| HVE Tips..... | 52 |
| Hygoformic™ Saliva Ejectors..... | 51 |

I

| | |
|--|----|
| iBarrier™ Sleeve for Tablets..... | 62 |
| Imaging Comfort | 58 |
| Iodine Test Strips | 15 |
| Isofluid™ Earloop Mask..... | 31 |
| Isofluid™ Earloop Mask with Secure Fit™ Mask Technology | 29 |
| Isofluid™ Plus Earloop Mask..... | 31 |
| Isofluid™ Plus Pure Earloop Mask with Secure Fit™ Mask Technology | 29 |
| Isolator™ Plus N95 Particulate Respirator | 32 |
| Isolite™ Earloop Mask..... | 32 |
| It's-A-Wrap™ Barriers..... | 58 |

L

| | |
|--------------------------------------|----|
| Lab Transport Bags..... | 64 |
| Lead-Free Steam Indicator Tape | 80 |
| Le Petit Masque..... | 32 |
| Light Sleeves, T-Handle | 63 |
| Liners | 36 |
| Liquid Ultra™ Solution..... | 15 |
| Luer-Lok® Needles..... | 66 |
| Luer-Lok® Syringes..... | 66 |

M

| | |
|---|----|
| Manifolds..... | 20 |
| Mask Feature Guide | 34 |
| Mask Ordering Information | 35 |
| Masks..... | 26 |
| Masks with Secure Fit™ Technology..... | 28 |
| Matrix Strips..... | 72 |
| MaxVac Plus™ Combo Tip™ Evacuators..... | 52 |
| Medicine/Mixing Cups | 73 |
| Municipally Plumbed Systems..... | 14 |

N

| | |
|--------------------------------|----|
| Nasal Masks..... | 21 |
| Needles | 66 |
| Non-Vented Evacuator Tips..... | 52 |
| Non-Woven Sponges | 56 |

O

| | |
|---------------------------|----|
| Operatory Protection..... | 60 |
|---------------------------|----|

P

| | |
|---|--------|
| Patient Barrier Aprons | 36, 38 |
| Patient Safety Mask with Shield..... | 33 |
| Patient Towels..... | 36 |
| Personal Inhaler Plus+™ Nasal Masks | 22 |
| Pillow Cases..... | 62 |
| PIP+™ Scavenging Circuits..... | 24 |
| Plastic Cups | 73 |
| Polyback™ Contour Bib Towels..... | 37 |
| Polyback™ Towels | 37 |
| Poly Coated Mixing Pads | 72 |
| Poly-Gard™ Deluxe Towels | 37 |
| Poly-Gard™ Towels..... | 37 |
| PPE | 26 |
| Premium™ All-Gauze Sponges..... | 55 |
| Premium™ Cotton-Filled Sponges | 56 |
| Premium™ Non-Woven Sponges..... | 56 |
| Premium™ Saliva Ejectors | 51 |
| Preventives | 40 |
| Proback™ Disposable Paper Towels..... | 37 |
| Procedural Earloop Mask..... | 30 |
| Procedural Earloop Mask with Secure Fit™ Mask Technology | 28 |
| Professional™ Regular Towels | 37 |
| Professional™ Thrift Towels | 37 |
| ProGuard™ 5% Sodium Fluoride Varnish | 42 |
| Prophy Angles | 44 |
| Prophy Paste..... | 45 |

R

| | |
|---|---|
| Rapicide™ OPA/28 High-Level Disinfectant..... | 9 |
| Rapicide™ OPA/28 Test Strips..... | 9 |

S

| | |
|--|----|
| Safe-Flo™ Saliva Ejectors and Valves | 50 |
| Safety IV Catheters..... | 66 |
| Saliva Ejectors..... | 50 |
| Sani-Cloth® AF3 Germicidal Disposable Wipe..... | 7 |
| Sani-Cloth® Bleach Germicidal Disposable Wipe..... | 7 |
| Sani-Cloth® HB Germicidal Disposable Wipe | 8 |
| Sani-Cloth® Plus Germicidal Disposable Cloth..... | 8 |
| Sani Glut™ Glutaraldehyde 3% | 9 |
| Sani-Hands® Instant Hand Sanitizing Wipes..... | 39 |
| Sani HLD Holder..... | 9 |
| Sani ProZyme™ Enzymatic Detergent | 10 |
| Sani-Roll™ Sterilization Tubing | 84 |
| Sani-Tab™ Chain-Free Patient Towel..... | 36 |

Index

| | |
|--|----|
| SaniTex Plus™ 2 Cleaner, Sanitizer & Disinfectant..... | 6 |
| SaniTex Plus™ Cleaner, Disinfectant & Deodorizer | 6 |
| SaniTex Plus™ Germicidal Disposable Wipe..... | 8 |
| Sani-Tube™ Nylon Sterilization Tubing | 84 |
| SaniTyz™ Hand Sanitizer..... | 39 |
| SaniTyz™ Wall-mounted Dispenser | 39 |
| Sani Vac™ Evacuation System Cleaner..... | 11 |
| Scalpels | 65 |
| Scavenging Circuits | 23 |
| See Clear® Lens Cleaning Wipes..... | 33 |
| Self-Sealing Autoclave Bags, Paper | 83 |
| Self-Seal Nylon Sterilization Pouches | 83 |
| Sensor Slippers™ Cushions | 58 |
| Sharps Containers | 68 |
| Shock/In-Line Treatment..... | 15 |
| Single-Use Disposables | 48 |
| Slip-Ease™ Cushions | 59 |
| Slotted Tip..... | 52 |
| Sparkle™ Air Water Syringe Tips..... | 53 |
| Sparkle™ Disposable Prophyl Angles | 44 |
| Sparkle™ EZ Contra Disposable Angles..... | 44 |
| Sparkle™ Free Prophyl Paste..... | 45 |
| Sparkle™ Prophyl Paste..... | 45 |
| Sparkle V™ 5% Sodium Fluoride Varnish | 42 |
| Standard Masks | 30 |
| Steam Indicator Tape..... | 80 |
| STEAMPlus™ Type 5 Integrators | 78 |
| Steri-Coil™ Pouch Dividers..... | 83 |
| Sterilization | 74 |
| Sterilization Pouches..... | 81 |
| Sterilization Tubing | 84 |
| Super Sani-Cloth® Germicidal Disposable Wipe..... | 7 |
| Sure-Check™ Sterilization Pouches | 82 |
| Surface Disinfectants | 6 |
| Surface Disinfection Microorganism Chart..... | 5 |
| Surgical Aspirator Tips..... | 53 |
| Surgical Molded Mask..... | 32 |
| Surgical Tie-on Mask | 31 |
| Surgical & Waste Accessories..... | 68 |
| Syringe-needle Combinations..... | 67 |
| Syringes | 66 |
| Syringe Sleeves | 63 |

T

| | |
|-----------------------------|----|
| Tape..... | 80 |
| Tartar & Stain Remover..... | 11 |

| | |
|---|----|
| Tongue Depressors..... | 72 |
| Topical Gel..... | 46 |
| Towels..... | 37 |
| Tray Sleeves..... | 64 |
| Treatment Accessories..... | 72 |
| Twist™ Disposable Prophyl Angles | 44 |
| Twist™ Plus Oscillating Disposable Prophyl Angles | 44 |
| Type 4 Steam Indicator Strips..... | 78 |

U

| | |
|---|----|
| Ultra-Gard™ Towels..... | 37 |
| Ultra-Gauze™ Non-Woven Sponges | 55 |
| Ultra-Pure™ Non-Woven Sponges | 57 |
| Ultra™ Earloop Mask | 30 |
| Ultra™ Earloop Mask with Secure Fit™ Mask Technology | 28 |
| Ultra™ Sensitive Earloop Mask..... | 30 |
| Ultra™ Sensitive Earloop Mask with Secure Fit™ Mask Technology | 28 |

V

| | |
|--|----|
| VistaCheck™ Dual Check Valve Backflow Preventer | 17 |
| VistaClear™ Centralized Water Filtration System..... | 16 |
| VistaCool™ Direct-to-Drain Cooling System for Autoclave Wastewater..... | 17 |
| VistaPure™ Water Purification System..... | 16 |

W

| | |
|---------------------------------|----|
| Wallmate™ Sharps Collector..... | 69 |
| Wrap-Ease™ Cushions..... | 59 |

X

| | |
|-------------------|----|
| X-Ray Sleeve..... | 63 |
|-------------------|----|

Z

| | |
|--|----|
| Zap™ Dual Arch Fluoride Gel Trays..... | 43 |
| Zap™ Fluoride Gel | 43 |
| Zap™ Neutral pH Fluoride Gel..... | 43 |

Notes

Notes



Dedicated to innovative, high-quality solutions, services,
and education that ensure maximum compliance and
improve outcomes for healthcare professionals and patients.



All product names are trademarks of Crosstex International, Inc., a Cantel Medical Company, its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions.

Autoguard™, EZ Scrub™, Instaflash™, Insyte™, Luer-Lok®, and Wallmate™ are trademarks and registered trademarks of Becton Dickinson and Company.

BRAVO™, HYDRIM® and STATIM® are trademarks and registered trademarks of SciCan Corporation Ltd.

Bard Parker® and Rib-Back® are registered trademarks of Aspen Surgical Products, Inc.

DrainHub™ is a trademark of Vista Water Group, LLC.

iBarrier® is a trademark of William Neumann.

Rinn®, XCP®, Uni-Grip® and Eezee-Grip® are registered trademarks of Dentsply Sirona Inc.

Sani-Cloth®, Sani-Hands®, and See Clear® are registered trademarks of Professional Disposables International, Inc.

Velcro® is a registered trademark of Velcro Industries B.V.

© 2019 Crosstex International, Inc. DLIT00866 Rev A 0419