


It's not worth the headache

Reduce post-lumbar puncture headaches (PLPH)

Introducing lumbar puncture and myelogram trays with atraumatic Whitacre needles



The atraumatic needle's pencil-point design is proven to reduce PLPH.²



The Safe-T PLUS adult lumbar puncture tray shown with the BD Whitacre spinal needle.

The American Academy of Neurology recommends atraumatic needles for lumbar puncture procedures.

CareFusion now offers standard and Safe-T™ PLUS lumbar puncture and myelogram trays with atraumatic BD™ Whitacre needles. The American Academy of Neurology recommends atraumatic needles for lumbar puncture procedures vs. traditional cutting needles.¹ Trusted by anesthesiologists performing spinal anesthesia for years, atraumatic needles are proven to reduce PLPH in lumbar puncture and myelogram patients.² The atraumatic BD Whitacre needle provides:

- Precise access and control.
- A pencil-point configuration designed to spread dura fibers, helping minimize cerebrospinal fluid (CSF) loss and incidence of PLPH.
- Clear visualization of CSF.

Convenient trays are designed to save time and money.

CareFusion lumbar puncture and myelogram trays contain the components you need to measure cerebrospinal fluid pressure and collect cerebrospinal fluid specimens. Our Safe-T PLUS and standard trays are designed to be convenient, easy to use and intuitive with your procedural flow.

Safe-T <i>PLUS</i> trays		
Cat. no.	Description	
4306CWSP	Safe-T <i>PLUS</i> adult lumbar puncture tray with 22 G x 3 ½" BD Whitacre spinal needle	
4324AWSP	Safe-T <i>PLUS</i> myelogram tray with 22 G x 3 ½" BD Whitacre spinal needle	
Safe-T <i>PLUS</i> tray components		
Safe-T <i>PLUS</i> tray components	4306CWSP	4324AWSP
BD SafetyGlide™ needle, 18 G x 1 ½" (3.8 cm)		•
BD SafetyGlide needle, 22 G x 1 ½" (3.8 cm)	•	•
BD SafetyGlide needle, 25 G x 1" (2.5 cm)	•	•
Point-Lok® needle retention safety device	2	2
Filter Straw® filter	•	•
Medication label sheet	•	•
Fine-tip marker	•	•
Time out label	•	•
Needle foam stop	•	•
Ear loop mask	•	•
PVP-I triple swab sticks	•	•
Standard components		
2 piece manometer, 550 mm	•	
3-way stopcock	•	
Extension tube, 5"	•	•
Extension tube, 20"		•
Prelabeled specimen vials with caps, 10 mL	4	3
Swab sticks	3	3
Lidocaine HCL, USP, 1%, 5 mL	•	•
Luer-lock syringe, 5 mL	•	•
Luer-lock syringe, 20 mL		•
Gauze pads	3	3
Bandage	•	•
Towel	•	•
Fenestrated drape with tape	•	•
Hospital wrap	•	•

Standard trays		
Cat. no.	Description	
4306CW	Adult lumbar puncture tray with 22 G x 3 ½" BD Whitacre spinal needle	
4324AW	Myelogram tray with 22 G x 3 ½" BD Whitacre spinal needle	
Standard tray components		
	4306CW	4324AW
Syringe, 3 mL with 25 G x 5/8" (1.6 cm) needle	•	
Needle, 25 G x 5/8" (1.6 cm)		•
Needle, 22 G x 1 ½" (3.8 cm)	•	•
Needle, 18 G x 1 ½" (3.8 cm)		•
3-way stopcock	•	
2 piece manometer, 550 mm	•	
Extension tube, 5"	•	•
Extension tube, 20"		•
Syringe, 5 mL		•
Plastic syringe, 20 mL		•
Prelabeled specimen vials, with caps, 10 mL	4	3
Gauze pads	3	3
Swab sticks	3	3
Lidocaine HCL, USP, 1%	2 mL	5 mL
Antiseptic prep well	•	•
Fenestrated drape	•	•
Towel	•	•
Bandage	•	•
Hospital wrap	•	•

References

1 Armon, C., Evans, RW. Addendum to assessment: Prevention of post-lumbar puncture headaches: Report of the Therapeutics and Technology Assessment Subcommittee of the American Academy of Neurology. *Neurology*. August 2005;65(4):510-512.
2 Arendt, K., Demaerschalk, D., Wingerchuk, M. et al. Atraumatic lumbar puncture needles: After all these years, are we still missing the point? *Neurologist* January 2009;15(1):17-20.

CareFusion
Vernon Hills, IL

carefusion.com

